IRMA ACTH

REF IM2030, B89463

TABLE OF CONTENTS

English	2	Slovenčina (SK)	34
Français (FR)	6	한국어 (KO)	38
Deutsch (DE)	10	Türkçe (TR)	42
Italiano (IT)	14	Русский (RU)	46
Español (ES)	18	Србија (SR)	50
Ελληνικά (EL)	22	APPENDIX	54
Polski (PL)	26		
Čeština (CZ)	30		

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IRMA ACTH

REF IM2030, B89463

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF ACTH IN HUMAN PLASMA For *in vitro* diagnostic use.

PRINCIPLE

The immunoradiometric assay of corticotropin hormone (ACTH) is a "sandwich" type assay where the ACTH is first modified chemically by succinylation, yielding succinyl ACTH (sACTH).

ACTH is highly conserved in different biological species and is unstable. For these two reasons it is particularly difficult to immunize animals against the human hormone for the production of antibodies with sufficiently high affinity for use in an immunoradiometric assay. The chemical modification of human ACTH by succinylation permits to differentiate it from the endogenous ACTH of the immunized animals and enhances its stability. Monoclonal antibodies with high affinity for succinylated human ACTH were thus obtained. The specificity of the antibodies utilized in the assay requires that the same chemical modification is applied in the case of the samples. This is achieved simply and highly reproducibly with the succinylation reagent furnished in the kit. Furthermore, the chemical modification of the samples stabilizes the ACTH and thus facilitates sample handling.

Three mouse monoclonal antibodies, directed against three different epitopes of the sACTH molecule and not competing with each other, are used in the kit. The modified samples and calibrators are incubated first in tubes coated with two monoclonal antibodies. Following the first incubation, the contents of the tubes are aspirated carefully. A solution containing a third monoclonal antibody, ¹²⁵I-labeled, is then added to each tube. Following a second incubation, the contents of the tubes are aspirated and, after washing twice, bound radioactivity is measured. Unknown values are determined by interpolation from a standard curve. The concentration of ACTH in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Bring all reagents to room temperature before pipeting.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Protection against ionizing radiation

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection.

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HbsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X)	DANGER		
	H360		May damage fertility or the unborn child.
	P201		Obtain special instructions before use.
	P280		Wear protective gloves, protective clothing and eye/face protection.
	P308+P313		IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - 0.3% Sodium Borate Decahydrate 0.1 - 0.3%
Succinylation Reagent	WARNING		
	H302		Harmful if swallowed.
	H319		Causes serious eye irritation.
	H335		May cause respiratory irritation.
	P261		Avoid breathing vapours.
	P280		Wear protective gloves, protective clothing and eye/face protection.
	P304+P340		IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.
	P312		IF eye irritation persists: Get medical advice/attention.
	P337+P313		Store in a well-ventilated place. Keep container tightly closed. Succinic Anhydride > 95%
	P403+P233		
Calibrators	WARNING		
	H315		Causes skin irritation.
	H319		Causes serious eye irritation.
	H335		May cause respiratory irritation.
	P261		Avoid breathing vapours.
	P271		Use only outdoors or in a well-ventilated area.
	P280		Wear protective gloves, protective clothing and eye/face protection.

P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P304+P340	IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P312	Call a POISON CENTER or doctor/physician if you feel unwell.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P337+P313	If eye irritation persists: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before use.
P403+P233	Store in a well-ventilated place. Keep container tightly closed.
P405	Store locked up.
P501	Dispose of contents/container in accordance with local/national regulations Dimethyl Sulfoxide Succinic Anhydride Potassium Hydroxide

In view of the stabilization of ACTH by succinylation, plasma samples may also be succinylated immediately after centrifugation following the protocol of paragraph Assay procedure. The succinylated samples can then be kept for 7 days at 2-8°C. It is recommended that for longer periods they be kept at below -18°C. It is not necessary to aliquot succinylated samples since succinylated ACTH is relatively insensitive to repeated freezing and thawing.

- If samples have concentrations greater than the highest calibrator, they must be succinylated and then diluted in the zero calibrator.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Kit for determination of ACTH, 100 tubes (Cat #IM2030)

Anti-sACTH monoclonal antibodies-coated tubes: 2 x 50 tubes (ready for use)

¹²⁵I-labeled monoclonal anti-sACTH antibody: one 11 mL vial (ready for use)

The vial contains 814 kBq, at the date of manufacture, of ¹²⁵I-labeled immunoglobulins in liquid form, proteins, sodium azide (<0.1%) and a dye.

Calibrators: five 2 mL vials + one 4 mL of "zero" calibrator (ready for use)

The calibrators contain human serum treated with succinic anhydride. They are thus ready for use and **must not be chemically modified further**.

The calibrator vials contain between 0 to approximately 1,500 pg/mL of synthetic sACTH in human serum and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against the international standard NIBSC 74/555.

Control serum: 2 vials (lyophilized)

Attention: Due to the instability of the ACTH molecule, two vials of the same control serum are supplied in each kit. After reconstitution, do not conserve the control serum untreated by succinic anhydride, even frozen. Use the second vial for another series.

The vials contain synthetic ACTH lyophilized in human serum. The expected values are in the concentration range indicated on a supplement.

Alkaline solution: one 15 mL vial (ready for use)

The buffer contains sodium azide (<0.1%).

Succinylation reagent: one vial (powder)

The vial contains 600 mg of succinic anhydride.

DMSO (Dimethylsulfoxide): one 3 mL vial (ready for use)

DMSO solution with phenol red. Let DMSO equilibrate at room temperature until completely melted. Mix gently before use.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

Kit for determination of ACTH, 50 tubes (Cat #B89463)

Anti-sACTH monoclonal antibodies-coated tubes: 1 x 50 tubes (ready for use)

¹²⁵I-labeled monoclonal anti-sACTH antibody: one 11 mL vial (ready for use)

Calibrators: five 2 mL vials + one 4 mL of "zero" calibrator (ready for use)

Control serum: 2 vials (lyophilized)

Alkaline solution: one 15 mL vial (ready for use)

Succinylation reagent: one vial (powder)

DMSO (Dimethylsulfoxide): one 3 mL vial (ready for use)

Wash solution (20x): one 50 mL vial

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- polypropylene tubes.
- precision micropipets (300 µL; 500 µL).
- semi-automatic pipets (50 µL; 100 µL; 250 µL, 2 mL).

DMSO

WARNING

H227	Combustible Liquid
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing vapours.
P280	Wear protective gloves, protective clothing and eye/face protection.
P304+P340	IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
P312	Call a POISON CENTER or doctor/physician if you feel unwell.
P337+P313	If eye irritation persists: Get medical advice/attention.
P403+P233	Store in a well-ventilated place. Keep container tightly closed. Dimethyl Sulfoxide > 95%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Sample blood into siliconized glass or plastic tubes containing EDTA. Do not use heparin as anti-coagulant. Do not use non-siliconized glass tubes; they adsorb ACTH.
- Separate plasma as soon as possible by centrifugation at 2-8°C.
- Since ACTH is very unstable, plasma samples must be stored at below -18°C unless the assay is done immediately. Samples should be aliquoted to avoid repeated freezing and thawing.

- vortex type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter set for 125 iodine.

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature.

Reconstitution of control sample

The content of one of the two vials is reconstituted with the volume of distilled water indicated on the label. Wait for 10 - 20 min (max) following reconstitution and mix gently to avoid foaming before dispensing. After reconstitution, do not conserve the control serum untreated by succinyl anhydrid, even frozen. Use the second vial for another series.

Reconstitution of the succinylation reagent

The content of the vial is reconstituted with the volume of DMSO indicated on the label. Mix thoroughly before use and visually check complete dissolution of succinyl anhydride. After reconstitution, the reagent is stable at 2-8°C until the expiry date of the kit. The slight tanning, sometimes noted during the storage, does not affect in anyway the performance of the succinylation reagent.

Preparation of wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution can be stored at 2-8°C until the expiry date of the kit.

Assay procedure

Succinylation of the control serum and samples

Attention: do not succinylate calibrators
To a polypropylene tube, add sequentially in this order
500 µL of sample or control serum 250 µL of alkaline solution 50 µL of succinylation reagent
Vortex immediately Incubate 5 minutes at 18-25°C

Remark 1: The volumes indicated above suffice for an assay in duplicate. For sample non-duplicate assays, sample, alkaline solution and succinylation reagent volumes may be reduced by one-half.

Remark 2: The color indicator of the succinylating reagent must turn yellow which indicates that the reaction has occurred normally. An abnormal reaction is indicated by persistence of the violet color.

Immunological step 1	Immunological step 2	Counting
To antibody-coated tubes, add:	To each tube, add:	Wash twice with 2 mL of wash solution.
300 µL of calibrator, or succinylated sample or succinylated control serum	100 µL of tracer	Count bound cpm (B) and total cpm (T) for 1 min.
Incubate 1 hour at 18-25°C with shaking (>280 rpm).	Incubate 2 hours at 18-25°C with shaking (>280 rpm).	
Aspirate the content of each tube carefully	Aspirate the content of each tube carefully (except the 2 tubes «total cpm»)	

*Add 100 µL of tracer to 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of ACTH concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{cal} - cpm_{cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (pg/mL).

Other data reduction methods may give slightly different results.

Total activity: 257,224 cpm				
Calibrators	ACTH (pg/mL)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	182	0.07	-
1	8.2	829	0.32	647
2	26.0	1,847	0.72	1,665
3	87.0	5,505	2.14	5,323
4	451	23,888	9.30	23,706
5	1,322	57,455	22.34	57,273

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding ACTH concentration on the horizontal axis.

Note: The concentration indicated on each calibrator vial already takes into account the sample dilution due to the succinylation step. This dilution must be neglected when calculating the concentration in sample (direct determination on the standard curve).

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

Number	Average	Median	2.5 th percentile	97.5 th percentile
42	18.5 pg/mL	15 pg/mL	7 pg/mL	41 pg/mL

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Sensitivity

Analytical sensitivity: 0.31 pg/mL

Functional sensitivity: 4.82 pg/mL

Specificity

The antibodies used in the immunoassay are highly specific for ACTH. Extremely low cross reactivities were obtained against several related molecules (ACTH 1-24; 1-10; 18-39; 11-24; αMSH and POMC) either in the absence (cross-reactivities) or the presence (interferences) of ACTH.

Precision

Intra-assay

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 10.0%.

Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 10.8%.

Accuracy

Dilution test

High-concentration samples were succinylated and then serially diluted with the zero calibrator. The recovery percentages obtained were between 82.3% and 119%.

Recovery test

Low-concentration samples were spiked with known quantities of ACTH. Samples were then succinylated and assayed. The recovery percentages obtained were between 97.2% and 116%.

Measurement range (from analytical sensitivity to highest calibrator): 0.31 to approximately 1,500 pg/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Heparin binds to ACTH and must therefore not be used. Samples should therefore be collected on EDTA.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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PERFORMANCE CHARACTERISTICS

Specificity

The specificity of the assay was determined by measuring the ACTH equivalents given by high concentrations of related molecules either in the absence (cross-reactivity) or the presence (interference) of ACTH.

These assays were done with plasma samples to which the related molecules, acylated according to the kit protocol, had been added.

Cross-reactivity

Related molecule	Concentration (pg/mL)	Concentration of ACTH measured (pg/mL)
ACTH 1-24	1,000	92
	500	54
	200	18
ACTH 1-10	1,000	0
	500	0
	200	0
ACTH 18-39	2,000	0
	1,000	0
	500	0
ACTH 11-24	100,000	0
	50,000	0
	1,000	0
α -MSH	1,000	0
	500	0
	200	0
POMC	1,000	4
	500	2
	200	<2

Interference

Spiking	Concentration of ACTH expected (pg/mL)	Concentration of ACTH measured (pg/mL)
ACTH + ACTH 1-24 (1,000 pg/mL)	263	359
ACTH + α -MSH (1,000 pg/mL)	263	257

Precision

Intra-assay

EDTA plasma	P1	P2	P3
Number of determinations	25	25	25
Mean value (pg/mL)	6.84	73.62	400.1
C.V., (%)	9.98	6.77	2.63

Inter-assays

EDTA plasma	P1	P2	P3
Number of determinations	10	10	10
Mean value (pg/mL)	16.73	373.7	579.9
C.V., (%)	10.8	6.53	6.86

Accuracy

Dilution test

Three samples containing ACTH were succinylated, diluted in the zero calibrator and then assayed according to the procedure of the kit.

EDTA plasma	Dilution factor	Measured conc. (pg/mL)	Expected conc. (pg/mL)	Ratio (%) Measured/Expected
P1	-	380.9	-	-
	1:2	226.5	190.4	118.9
	1:4	104.3	95.21	109.5
	1:8	51.78	47.61	108.8
	1:16	21.52	23.80	90.41
	1:32	11.87	11.90	99.73
P2	-	467.0	-	-
	1:2	273.1	233.5	117.0
	1:4	123.3	116.7	105.6
	1:8	54.42	58.37	93.23
	1:16	28.05	29.19	96.11
	1:32	12.01	14.59	82.30
P3	-	307.6	-	-
	1:2	159.9	153.8	104.0
	1:4	89.66	76.91	116.6
	1:8	38.22	38.45	99.39
	1:16	18.58	19.23	96.64

Recovery test

Known amount of ACTH was added to EDTA plasma samples succinylated and assayed according to the procedure of the kit.

EDTA plasma	Endogen. conc. (pg/mL)	Added conc. (pg/mL)	Expected conc. (pg/mL)	Measured conc. (pg/mL)	Ratio (%) Measured/Expected
P1	16.90	13.09	29.18	30.13	103.3
	16.02	18.62	34.64	33.67	97.20
	15.69	45.59	61.28	64.65	105.5
P2	59.10	18.62	77.72	87.96	113.2
	56.70	71.46	128.2	142.8	111.4
	54.49	120.2	174.7	203.3	116.4
P3	97.20	36.72	133.9	138.8	103.6
	94.58	71.46	166.0	170.4	102.6
	90.89	120.2	211.1	228.3	108.2

¹²⁵I Characteristics

$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$

¹²⁵ I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

M Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

Σ Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použiť pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

I Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

T Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo(i) di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(η) θερμοκρασίας / 溫度範圍 / Temperaturüros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Staklik aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

! Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

S Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kulanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

MW Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktivyvioji medžiaga / Radioaktiv / Radioaktivnyy / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZEŃSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLİKE / ОПАСНО / ОПАСНОСТ / 危險

Ag ^{125I}

Ab ^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracerlar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontroliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用說明

SOLN | WASH | 20x

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösingskoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalo koncentratas 20X / 20X mosóoldat-konzentrátum / Koncentrat 20X roztworu płuczacego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промычного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

SOLN | ALK

Alkaline Solution, Solution alcaline, Alkalische Lösung, Soluzione alcalina, Solución alcalina, Solução alcalina, Alkalisk lösning, Αλκαλικό διάλυμα, 碱性溶液, Alkalino tirpalas, Lúgos oldat, Roztwór zasadowy, Zásaditý roztok, Alkalický roztok, 알칼리성 용액, Alkalín Çözeltisi, Щелочной раствор, Алкален разтвор, 鹼性溶液

REAG | SUC

Succinylation Reagent, Réactif de succinylation, Succinylierungsreagenz, Reattivo di succinilazione, Reactivo de succinilación, Reagente de Succinilação, Succinylation Reagent, Αντιδραστήριο ακετυλίωσης, 琥珀酰化试剂, Succinylation Reagent, Succinylation Reagent, Odczynnik do sukcylnylacji, Sukcynylacni reagencie, Sukcynylačné činidlo, 석시닐화 시약, Süksinilleme Reaktifi, Реагент для сукцинирования, Реактив за сукциниране, 琥珀酰化試劑

DMSO

Dimethyl Sulfoxide, Diméthyle sulfoxyde, Dimethylsulfoxid, Dimetilsolfossido, Dimetilsulfóxido, Dimetilsulfóxido, Dimetilsulfoxid, Διμεθυλοσουλφοξίδιο, 二甲亚砜, Dimetilsulfoksidas, Dimetil-szulfoxid, Dimetylosulfoftlenek, Dimetilsulfoxid, Dimetyl sulfoxid, 디메틸 술폭시드, Dimetil Sülfoksit, Диметил сульфоксид, Диметилсулфоксид, 二甲基亞砜