Product information

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| Adiponectin human ELISA DEE009 | 96 Determinations |
|---------------------------------------|---|
| Principle of the test | Enzyme-linked Immunoassay |
| Duration (incubation period) | 1.75 h |
| Antibodies | monoclonal antibodies |
| Buffer | Ready for use and 20fold concentrate |
| Calibrator | 5 single calibrators: 2 - 100µg/L, native human Adiponectin |
| Assay Range | 0.27 - 31000µg/L |
| Control | 2 controls freeze-dried |
| Sample | human serum / plasma |
| Required sample volume | 10 µL |
| Sample dilution | 1:310 |
| Analytical sensitivity | $\emptyset \leq 0.27 \mu\text{g/L}$ |
| Intra- / Interassay Variance | $\emptyset < 10 \%$ |

1 INTENDED USE

The ELISA DEE009 is intended to be used for quantitative measurement of human Adiponectin in human serum and plasma samples.

2 INTRODUCTION

Adiponectin is a 30kDa protein which percentage in serum proteins is 0.01%. It is mainly synthesized by adipocytes, but also muscle cells and hepatocytes have the ability to synthesize Adiponectin. It consists of a Collagen-like N-terminal and a globular C-terminal domain [1]. In vivo Adiponectin appears with different oligomers. Beside the trimer and dimer also high molecular multimers exist [1-3]. Two different receptors are known, both receptors are ubiquitarily expressed, though the distribution in the tissues varies. The Adiponectin Receptor 1 (AdipoR1) is especially in muscle- and AdipoR2 in liver tissue synthesized [4].

Several studies show, that adiponectin correlates negatively with BMI and thus it could have relevance for the energy metabolism for example through the regulation of fatty acid oxidation. Beside the correlation with BMI, Adiponectin level is associated with the Insulin-Resistance [5-7] and so also linked with Type II Diabetes. Adiponectin is associated also with glucose- und lipometabolism [8, 9]. The formerly proposed diagnostic value of the high molecular weight form of adiponectin was not verified using a commercially available test system for the determination of HMW adiponectin [10]. Blueher et al. evaluated the Demeditec DEE009 regarding its diagnostic value in diagnosis of insulin resistance. Results of ROC analysis showed an area under curve of 0.92, which indicates a diagnostic value [10]. Furthermore adiponectin is involved in inflammatory processes [11-15] and therewith it is of importance for appearance of arteriosclerosis [4, 5, 16] and coronaritis [17, 18], thus the determination of Adiponectin level in plasma could serve to estimate the risk of coronary disease [19, 20]. Beside this Adiponectin influences further physiological processes as for example the angiogenesis [21, 22].

3 ASSAY PRINCIPLE

The Demeditec ELISA for Adiponectin DEE009 is a so-called Sandwich-Assay using two specific and highly affine antibodies. The Adiponectin in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-Adiponectin-Antibody binds in turn to the immobilised Adiponectin. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction, the turn of the colour will be catalysed quantitatively depending on the Adiponectin-level of the samples.

4 WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use only. For Professional use only.

The Demeditec kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Demeditec will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided. A Material Safety Data sheet is available on request.

Do not use obviously damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Controls** **CONTROL 1 LYO** and **CONTROL 2 LYO**, and **Calibrators** **CAL A – E LYO**.

Source human serum for the controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore, all components and patient's specimens should be treated as potentially infectious.

Reagents **CAL A – E LYO**, **ENZ CONJ**, **DIL**, **WASH SOLN 20x**

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one** and **2-methyl-4-isothiazolin-3-one** (<0.015%)

| | |
|-----------|---|
| H317 | May cause an allergic skin reaction. |
| P280 | Wear protective gloves/ protective clothing/ eye protection/ face protection. |
| P272 | Contaminated work clothing should not be allowed out of the workplace. |
| P261 | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. |
| P333+P313 | If skin irritation or rash occurs: Get medical advice/ attention. |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water. |
| P501 | Dispose of contents/ container in accordance with local/ regional/ national/ international regulations. |

Substrate **SUB TMB**

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine (<0.05%)

| | |
|------------|---|
| H315 | Causes skin irritation. |
| H319 | Causes serious eye irritation. |
| H335 | May cause respiratory irritation. |
| P261 | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes. |
| P338 | Remove contact lenses, if present and easy to do. Continue rinsing. |

Stop Solution **STOP SOLN**

The Stop solution contains 0.2 M acid sulphur acid (H₂SO₄)

| | |
|------------|---|
| H290 | May be corrosive to metals. |
| H314 | Causes severe skin burns and eye damage. |
| P280 | Wear protective gloves/ protective clothing/ eye protection/ face protection. |
| P301+P330+ | IF SWALLOWED: rinse mouth. |
| P331 | Do NOT induce vomiting. |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes. |
| P338 | Remove contact lenses, if present and easy to do. Continue rinsing. |
| P309+P310 | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |

4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

5.1 Sample type

Serum and Plasma

Serum and heparin plasma levels are comparable. In EDTA- and Citrate Plasma-samples levels were found approx. 18% lower, because of the relatively high amount of anticoagulant.

5.2 Specimen collection

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. Hemolytic reactions have to be avoided.

5.3 Required sample volume: 10 µL

5.4 Sample stability

In firmly closable sample vials

- Storage at 20-25°C: 2 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles max. 3

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized, 3 freeze-thaw cycles showed no effect on the measured adiponectin concentration.

5.5 Interference

Hemoglobin, triglyceride and bilirubin in the sample do not interfere to a concentration of **5 mg/mL**, **100 mg/mL** and **100 µg/mL** respectively. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

5.6 Sample dilution

- Dilution: **1:310** with Dilution Buffer **DIL**
- Dilute for example **300 µL** Dilution Buffer **DIL** in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add **10 µL** Serum- or Plasma (dilution: 1:31). Add **900 µL** Dilution Buffer **DIL** in another PE-/PP-tube and **100 µL** of the thoroughly mixed first dilution. After mixing, use 2×100 µL from this **1:310** diluted sample in the assay.
- Sample stability after dilution of the sample: maximum 1 hour at 20-25°C.

6 MATERIALS

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the calibration curve.

| | | |
|---------------------------------------|--|---------------------|
| SORB MT | Microtiter plate , ready for use, coated with mouse-anti-Adiponectin-antibody. Wells are separately breakable. | (8x12) wells |
| CAL A - E LYO | Calibrators , lyophilized, (native human Adiponectin), concentrations are given on vial labels and on quality certificate in ng/mL. | 5 x 750 µL |
| CONTROL 1 LYO | Control 1 , lyophilised, (human serum), concentration is given on quality certificate in ng/mL. | 1 x 500 µL |
| CONTROL 2 LYO | Control 2 , lyophilised, (human serum), concentration is given on quality certificate in ng/mL. | 1 x 500 µL |
| ENZ CONJ | Antibody-HRP-Conjugate , ready for use, mouse-anti-hAdiponectin-antibody biotinylated + streptavidin horseradish peroxidase conjugate | 1 x 12 mL |
| DIL | Dilution Buffer , ready for use | 1 x 125 mL |
| WASH SOLN 20x | Washing Buffer , 20-fold concentrated solution | 1 x 50 mL |
| SUB TMB | Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine. | 1 x 12 mL |
| STOP SOLN | Stop Solution , ready for use, 0.2 M sulphuric acid. | 1 x 12 mL |
| | Sealing Tape , for covering the microtiter plate . | 2 x |
| | Instructions for use | 1 x |
| | Quality Certificate | 1 x |

6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WASH** **SOLN** **20x** (A. dest.), 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** calibrators **CAL** **A - E** **LYO** and Controls **CONTROL** **1** **LYO** and **CONTROL** **2** **LYO** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WASH** **SOLN** **20x** is 4 weeks stable at 2-8°C.

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Calibrators **A – E** and Controls **CONTROL 1 LYO** and **CONTROL 2 LYO** are reconstituted with the Dilution Buffer **DIL**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Controls **CONTROL 1 LYO** and **CONTROL 2 LYO** with the Dilution Buffer **DIL** in the same ratio (1:310) as the sample. The required volume of Washing Buffer **WASH SOLN 20x** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

Assay Procedure

When performing the assay, Blank, Calibrators **CAL A – E LYO**, Controls **CONTROL 1 LYO**, **CONTROL 2 LYO** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **ENZ CONJ** as well as the succeeding Substrate **SUB TMB** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STOP SOLN** should be added to the plate in the same order as Substrate **SUB TMB**. All determinations (Blank, Calibrators **CAL A – E LYO**, Controls **CONTROL 1 LYO**, **CONTROL 2 LYO** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate **SUB TMB**, stabilised H₂O₂-Tetramethylbencidine, is photosensitive–store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WASH SOLN 20x** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistep device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8 ASSAY PROCEDURE

| Preparation of reagents | | Reconstitution: | Dilution |
|--|--|---|-----------------------------|
| CAL A - E LYO | Calibrators | in 750 µL Dilution Buffer DIL | - |
| CONTROL 1 LYO & CONTROL 2 LYO | Control 1 Control 2 | in 500 µL Dilution Buffer DIL | 1:310 with DIL |
| WASH SOLN 20x | Washing Buffer | - | 1:20 with Aqua dest. |
| Sample dilution: with Dilution Buffer DIL 1:310 | | | |
| Before assay procedure bring all reagents to room temperature 20-25°C . | | | |
| Assay Procedure in Double Determination: | | | |
| Pipette | Reagents | Position | |
| 100 µL | Dilution Buffer DIL (Blank) | A1/A2 | |
| 100 µL | Calibrator A (2 ng/mL) | B1/B2 | |
| 100 µL | Calibrator B (10 ng/mL) | C1/C2 | |
| 100 µL | Calibrator C (30 ng/mL) | D1/D2 | |
| 100 µL | Calibrator D (70 ng/mL) | E1/E2 | |
| 100 µL | Calibrator E (100 ng/mL) | F1/F2 | |
| 100 µL | Control CONTROL 1 LYO (1:310 diluted) | G1/G2 | |
| 100 µL | Control CONTROL 2 LYO (1:310 diluted) | H1/H2 | |
| 100 µL | Sample (1:310 diluted) | in the rest of the wells according the requirements | |
| Cover the wells with the sealing tape. | | | |
| Sample Incubation: 1 h at 20-25°C, 350 rpm | | | |
| 3 x 300 µL | Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer (diluted) / well | In each well | |
| 100 µL | Antibody-HRP-Conjugate ENZ CONJ | In each well | |
| Cover the wells with the sealing tape. | | | |
| Incubation: 30 Minutes at 20-25°C, 350 rpm | | | |
| 3 x 300 µL | Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer (diluted) / well | In each well | |
| 100 µL | Substrate SUB TMB | In each well | |
| Incubation: 15 Minutes in the Dark at 20-25°C | | | |
| 100 µL | Stop Solution STOP SOLN | In each well | |
| | Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength. | | |

9 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All calibrators and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

9.1 Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of calibrator E should be above 1.00. Samples, which yield higher absorbance values than calibrator E, should be re-tested with a higher dilution.

10 EVALUATION OF RESULTS

10.1 Establishing of the calibration curve

| Calibrator | A | B | C | D | E |
|------------|---|----|----|----|-----|
| ng/mL | 2 | 10 | 30 | 70 | 100 |

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The Adiponectin concentration in ng/mL of the samples and controls **CONTROL 1 LYO** and **CONTROL 2 LYO** can be calculated by **multiplication** with the respective **dilution factor**.

10.2 Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

| | Blank | A | B | C | D | E |
|----------------------------|-------|-------|-------|-------|-------|-------|
| ng/mL | 0.0 | 2 | 10 | 30 | 70 | 100 |
| OD _(450-620 nm) | 0.008 | 0.071 | 0.357 | 1.022 | 2.053 | 2.817 |


Figure 1 Exemplary calibration curve

The exemplary shown calibration curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

10.3 Exemplary calculation of Adiponectin concentrations

Sample dilution: 1:310

| | |
|------------------------------------|-------|
| Measured extinction of your sample | 0.408 |
| Measured extinction of the blank | 0.008 |

Your measurement program will calculate the Adiponectin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2nd degree).

In this exemplary case the following equation is solved by the program to calculate the adiponectin concentration in the sample:

$$\begin{aligned}0.400 &= -6 \times 10^{-5}x^2 + 0.0344x \\ 11.89 &= x\end{aligned}$$

If the dilution factor (1:310) is taken into account the adiponectin concentration of the undiluted sample is

$$11.89 \times 310 = 3685.9 \text{ ng/mL} = 3.69 \text{ } \mu\text{g/mL}$$

The exemplary shown calibration curve in Fig.1 cannot be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

10.4 Interpretation of results

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

10.5 Limitation of procedure

The Demeditec sensitive human Adiponectin ELISA, DEE009 is based on monoclonal antibodies. Generally, this technique could be sensible to heterophilic antibodies or rheumatic factors in the sample. Their influence is reduced by assay design, but cannot be excluded completely.

11 REFERENCE VALUES

The expected values for serum adiponectin, which were determined with the Demeditec ELISA DEE009 in healthy donors and analysed by Prof. Dr. J. Kratzsch, Department of Laboratory Medicine, University Hospital Leipzig, are given below (Tab. 1).

These data show significant correlation between Adiponectin-Serum values and age as well as gender of the probands, in turn the correlation between the respective BMI seems to be less significant. In the samples of neonatal cord blood very high values were found. Several different statistical analyses were performed to adapt for certain individual demands. The best suited data can be chosen respectively for interpretation of the own measurements.

Table 1a Expected values for adults, gender specific mean as well as median, 5. and 95. percentiles are given.

| Sex | Number [n] | Mean [µg/ml] | Median [µg/ml] | Standard Deviation | 5 th Percentile [µg/ml] | 95 th Percentile [µg/ml] |
|--------------|------------|--------------|----------------|--------------------|------------------------------------|-------------------------------------|
| Female | 101 | 10.2 | 9.1 | 4.6 | 4.0 | 19.4 |
| Male | 125 | 6.8 | 6.1 | 4.1 | 2.0 | 13.9 |
| total | 226 | 8.3 | 7.5 | 4.6 | 2.4 | 19.3 |

Table 1b Expected values for children, gender specific mean as well as median, 5. and 95. percentiles are given.

| Sex | Number [n] | Mean [µg/ml] | Median [µg/ml] | Standard Deviation | 5 th Percentile [µg/ml] | 95 th Percentile [µg/ml] |
|--------------|------------|--------------|----------------|--------------------|------------------------------------|-------------------------------------|
| Female | 131 | 8.71 | 8.18 | 4.32 | 3.05 | 15.6 |
| Male | 134 | 8.97 | 8.12 | 5.13 | 3.36 | 18.6 |
| total | 265 | 8.84 | 8.18 | 4.74 | 3.33 | 16.5 |

Table 1c Expected values for Adiponectin, age specific mean as well as median, 5. and 95. percentiles are given.

| Age group [a] | Number [n] | Mean [µg/ml] | Median [µg/ml] | 5. Percentile [µg/ml] | 95. Percentile [µg/ml] |
|---------------|------------|--------------|----------------|-----------------------|------------------------|
| < 7.99 | 46 | 12.82 | 11.7 | 2.33 | 26.5 |
| 8 – 9.99 | 40 | 8 | 8.09 | 3.96 | 14.9 |
| 10-11.99 | 55 | 8.02 | 7.14 | 3.36 | 13.8 |
| 12 – 13.99 | 26 | 8.21 | 7.54 | 4.5 | 13.2 |
| 14 – 15.99 | 59 | 8.12 | 8.09 | 3.67 | 13.7 |
| 16 – 19.99 | 41 | 7.97 | 7.79 | 2.74 | 13.3 |
| alle | 267 | 8.88 | 8.18 | 3.33 | 16.7 |
| <hr/> | | | | | |
| 20 – 29.99 | 47 | 6.72 | 6.38 | 2.50 | 12.25 |
| 30 – 39.99 | 38 | 7.38 | 6.69 | 1.98 | 19.29 |
| 40 – 49.99 | 55 | 8.42 | 8.20 | 2.41 | 17.85 |
| 50 – 59.99 | 47 | 9.61 | 8.55 | 2.15 | 19.85 |
| > 60 | 32 | 9.52 | 8.57 | 3.00 | 21.10 |
| alle | 226 | 8.33 | 7.5 | 2.41 | 19.29 |

Table 1d Expected values for Adiponectin, age as well as gender specific mean and median, BMI and 25. and 75. percentiles are given.

| Female | | | Adiponectin ($\mu\text{g/ml}$): | | | |
|---------------------------|----|---------------------|-------------------------------------|---------|-------------------------|--------------|
| Age (Years): | n: | BMI: AV \pm SD | AV \pm SD: | Median: | Percentile: 25.- 75. | Min. – Max.: |
| Newborn Cord blood | 19 | -- | 29.80 \pm 12.49 | 26.1 | 19.5 -35.2 | 16.9 - 61.4 |
| < 3.99 | 9 | 15.73 \pm 0.79 | 14.43 \pm 7.76 | 11.2 | 8.2 - 21.8 | 2.3 - 26.7 |
| 4.0 - 7.99 | 11 | 16.01 \pm 1.94 | 8.46 \pm 4.73 | 9.3 | 2.9 - 12.1 | 1.4 - 15.6 |
| 8.0 - 9.99 | 22 | 17.58 \pm 3.84 | 7.92 \pm 3.00 | 8.2 | 5.2 - 10.0 | 3.6 - 15.1 |
| 10.0 - 11.99 | 33 | 17.83 \pm 1.86 | 7.66 \pm 4.59 | 6.6 | 5.0 - 8.8 | 3.1 - 20.9 |
| 12.0 - 13.99 | 11 | 19.85 \pm 2.31 | 8.22 \pm 5.64 | 7.5 | 6.5 - 9.2 | 4.9 - 13.2 |
| 14.0 - 15.99 | 27 | 19.91 \pm 1.72 | 8.83 \pm 9.25 | 8.9 | 5.2 - 11.8 | 2.6 - 17.7 |
| 16.0 – 19.99 | 18 | 21.64 \pm 2.64 | 9.00 \pm 3.22 | 8.7 | 6.9 - 11.2 | 2.7 - 14.0 |
| 20.0 - 29.99 | 24 | 23.12 \pm 5.01 | 7.39 \pm 3.35 | 7.3 | 5.7 - 9.0 | 3.4 - 17.8 |
| 30.0 - 39.99 | 17 | 23.20 \pm 2.86 | 9.19 \pm 3.89 | 8.6 | 7.2 - 10.4 | 3.6 - 19.3 |
| 40.0 - 49.99 | 26 | 24.50 \pm 4.11 | 9.93 \pm 3.59 | 9.5 | 7.5 - 11.6 | 4.4 - 19.6 |
| 50.0 - 59.99 | 21 | 24.61 \pm 3.31 | 11.5 \pm 5.49 | 10.0 | 8.0 - 15.9 | 2.0 - 23.1 |
| >60.0 | 8 | 24.63 \pm 1.89 | 15.6 \pm 4.64 | 15.3 | 11.4 - 18.2 | 11.2 - 24.1 |

| Male | | | Adiponectin ($\mu\text{g/ml}$): | | | |
|---------------------------|----|---------------------|------------------------------------|---------|-------------------------|--------------|
| Age (Years): | n: | BMI: AV \pm SD | AV \pm SD: | Median: | Percentile: 25.- 75. | Min. – Max.: |
| Newborn Cord blood | 10 | -- | 27.80 \pm 7.68 | 26.7 | 22.2 -31.0 | 15.6 - 40.6 |
| < 3.99 | 14 | 16.17 \pm 1.81 | 16.57 \pm 6.55 | 14.3 | 11.6 - 21.2 | 5.8 - 40.3 |
| 4.0 - 7.99 | 12 | 15.69 \pm 1.05 | 11.24 \pm 5.43 | 9.7 | 8.9 - 15.9 | 3.5 - 18.6 |
| 8.0 - 9.99 | 18 | 16.45 \pm 1.76 | 8.11 \pm 2.93 | 7.6 | 6.2 - 9.1 | 5.00 - 15.4 |
| 10.0 - 11.99 | 21 | 18.34 \pm 2.18 | 8.43 \pm 3.91 | 7.8 | 5.2 - 10.9 | 3.4 - 20.2 |
| 12.0 - 13.99 | 14 | 18.61 \pm 2.11 | 7.59 \pm 2.86 | 7.1 | 6.0 - 10.3 | 2.4 - 12.2 |
| 14.0 - 15.99 | 32 | 19.86 \pm 2.00 | 7.53 \pm 2.52 | 7.4 | 5.1 - 9.3 | 3.8 - 15.4 |
| 16.0 -19.99 | 23 | 22.03 \pm 2.42 | 7.16 \pm 3.53 | 6.9 | 4.2 - 9.6 | 2.0 - 13.9 |
| 20.0 - 29.99 | 23 | 23.43 \pm 2.48 | 5.44 \pm 2.29 | 5.8 | 4.0 - 6.9 | 1.3 - 10.3 |
| 30.0 - 39.99 | 21 | 23.33 \pm 2.72 | 5.92 \pm 4.60 | 4.4 | 2.7 - 6.7 | 1.9 - 20.6 |
| 40.0 - 49.99 | 22 | 23.79 \pm 2.41 | 6.13 \pm 2.92 | 5.5 | 3.8 - 8.3 | 2.1 - 11.6 |
| 50.0 - 59.99 | 23 | 26.68 \pm 2.77 | 7.45 \pm 4.50 | 6.7 | 5.0 - 8.8 | 1.4 - 19.6 |
| >60.0 | 24 | 25.72 \pm 2.12 | 7.48 \pm 3.92 | 7.6 | 4.6 - 9.2 | 3.0 - 21.1 |

n= Number of Probands AV=Average Value BMI=Body Mass Index (kg/m^2) SD=Standard Deviation

12 PERFORMANCE CHARACTERISTICS

12.1 Sensitivity

The analytical sensitivity of the ELISA DEE009 was measured by the variability of the signal of the blank. Based on the twofold standard deviation of the blank the mean analytical sensitivity is < 0.27 ng/mL (Range 0.094 to 0.59 ng/mL).

12.2 Specificity

Adiponectin exists in different oligomeric forms: the high, medium and low molecular weight form. Different numbers of the adiponectin monomer aggregate specifically to form a complex. In Figure 2a the five different forms of human adiponectin are shown schematically. In parallel the results of a size-exclusion chromatography of human serum measured with the Demeditec DEE009 Adiponectin ELISA are shown.


Figure 2 Adiponectin Structure (A) proteinaceous structure of human adiponectin including posttranslational modifications and multimeric forms (taken from Nakano et al¹). (B) Results of a serum separation by size-exclusion chromatography. The sample was fractionated and the adiponectin content of each fraction was measured by Demeditec DEE009 Adiponectin, for comparison the corresponding seize in kDa is indicated by black arrows.

The results shown in Figure 2b clearly demonstrate that the Demeditec DEE009 Adiponectin ELISA detects all forms of Adiponectin present in human serum: the trimer at 65 kDa, the hexamer at 150 kDa and the high molecular weight forms of >280 kDa. The Demeditec DEE009 Adiponectin ELISA therefore measures total Adiponectin.

12.3 Precision

Intra-Assay Variance & Accuracy

Intra assay variance and accuracy is exemplarily shown with two samples. The adiponectin concentration of these samples was repeatedly measured in one assay.

Table 2 Intra-Assay Variation. Recombinant adiponectin was diluted in dilution buffer and the adiponectin concentration of the dilution was measured repeatedly within one assay.

| | Determinations [n] | Mean value [µg/L] | Standard deviation [µg/L] | VC [%] | Target Value [µg/L] |
|----------|--------------------|-------------------|---------------------------|--------|---------------------|
| Sample 1 | 8 | 7.108 | 0.22 | 3.14 | 6 |
| Sample 2 | 8 | 107.96 | 3.97 | 3.67 | 100 |

In both samples the variability is less than 5% and the deviation from the target value is <20%.

¹ Nakano Y, Tahima S, Yoshimi A, Akiyama H, Tsushima M, Tanioka T, Negoro T, Tomita M, Tobe T: A novel enzyme-linked immunosorbent assay specific for high-molecular-weight adiponectin. *J Lipid Res* 2006, 47(7):1572-1582

Inter-Assay Variance

Serum samples were repeatedly measured in independent assays of different lots. On average the coefficient of variation was 7.5% (SD 1.6). The results of 5 samples are shown in table 3.

Table 3 Inter-Assay Variation

| | Sample 1 | Sample 2 | Sample 3 | Sample 4 | Sample 5 |
|--------------|----------|----------|----------|----------|----------|
| Mean [µg/mL] | 4.72 | 5.25 | 8.36 | 5.45 | 22.29 |
| CV [%] | 8.16 | 8.14 | 6.93 | 8.05 | 7.30 |
| n | 99 | 68 | 62 | 174 | 62 |

12.4 Linearity

Linearity of sample dilution was tested by serial dilution (1:100 – 1:4000) of human serum samples and recalculation of the adiponectin content in comparison to the mean adiponectin concentration of all dilutions (Table 4), no diluted sample showed a deviation of >30%.

Table 4 Linearity. Human serum sample were diluted in **DIL** and adiponectin content was recalculated. Measurements results are shown in [mg/L]. No deviation of the mean >30% was detected.

| µg/mL | Mean | 1:100 | 1:200 | 1:400 | 1:800 | 1:1000 | 1:2000 | 1:4000 |
|----------|-------|-------|--------|--------|--------|--------|--------|--------|
| Sample 1 | 5.76 | 6.53 | 6.331 | 5.764 | 5.49 | 6.067 | 6.114 | 4.056 |
| Sample 2 | 11.53 | 10.93 | 12.107 | 11.395 | 11.454 | 11.567 | 12.884 | 10.362 |
| Sample 3 | 12.07 | 13.57 | 12.853 | 12.03 | 11.974 | 11.338 | 11.548 | 11.169 |
| Sample 4 | 4.89 | 4.659 | 4.886 | 4.384 | 4.425 | 5.851 | 5.13 | n/a |

Additionally, dilutions of 1:1500 to 1:96000 were evaluated with two samples. In Figure 3 the results are shown and demonstrate that in the tested samples no effect of dilution could be detected on measured adiponectin concentrations. The deviation of the target concentration of each dilution was less than 30% in all dilutions.


Figure 3 Linearity. Adiponectin concentration was measured in two human samples diluted 1:1500 to 1:96000.

12.5 Recovery

Trueness and traceability of the Demeditec Adiponectin ELISA DEE009 was evaluated by recovery of recombinant adiponectin in human serum. The recovery of recombinant Adiponectin yielded in a serum matrix on average 110%.

Table 5 Recovery of recombinant human Adiponectin in Serum. Recombinant Adiponectin was added in different amounts to human serum. The Adiponectin content of the so enriched samples was measured and recovery in comparison to enriched buffer calculated.

| R&D recombinant Adiponectin lot 1022911 | DIL | Serum | Recovery |
|---|-------|---------|------------|
| ng/mL | ng/mL | ng/mL | % |
| 0 | 0 | 0.00778 | --- |
| 75 | 87.92 | 95.71 | 109 |
| 37.5 | 51.84 | 59.98 | 116 |
| 18.75 | 26.55 | 26.7 | 101 |
| 9.375 | 13.35 | 15.15 | 113 |

12.6 Interference

Interference of bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing adiponectin. For comparison the same amount of buffer without any substance was also added to the serum. Table 6 demonstrates that neither bilirubin nor triglycerides or haemoglobin exert any influence on the measurement of adiponectin in human serum.

Table 6 Interference. Serum samples were enriched with different amount of triglycerides, bilirubin or haemoglobin. The relative amount of adiponectin measured in comparison with native serum is shown here [%].

| Triglyceride 100 mg/mL | Bilirubin 100 µg/mL | Hemoglobin 1 mg/mL | Hemoglobin 5 mg/mL |
|------------------------|---------------------|--------------------|--------------------|
| 94 | 96 | 90 | 109 |
| 90 | 93 | 97 | -- |
| 95 | 94 | 93 | -- |

13 COMPARISON STUDIES

Demeditec Adiponectin, DEE009 was compared with two different, commercially available test systems. In linear regression analysis both tests showed a good coefficient of determination ($R^2=0.896$ and $R^2=0.97$). Thus, the comparability of the results between the Demeditec DEE009 and other test systems is possible. In dependence on the respective system the absolute deviation of measurement results is different, but because of the excellent correlation the results are comparable after applying a factor. The results of both studies are shown in figure 4.


Figure 4 Comparison of Demeditec DEE009 with commercially available test systems (left) radioimmunoassay (n=14) and (right) ELISA (n=84)

Instructions for use for scientific application

14 SCIENTIFIC APPLICATION

In addition to serum and plasma samples Adiponectin can be determined in other human body fluids and in cell culture supernatants of human cell lines for research purposes.

14.1 Samples suitable for scientific application

Serum, plasma, saliva, urine, breast milk and cell culture supernatant of human cell lines

The recommended dilution for serum and plasma samples in Dilution Buffer **DIL** (1:310).

In the other samples, the Adiponectin levels can vary considerably, the optimal dilution must be found out by the customer.

Table 7 Results of sample matrix tests . Adiponectin was added to the respectively diluted samples. Enriched samples were measured without further dilution. Shown is the relative recovery of added Adiponectin of the concentration measured in respectively enriched assay buffer.

| Matrix Dilution | Urine | Saliva | Breast Milk | Cell Culture Medium 10% FCS | Cell Culture Medium |
|-----------------|-------|--------|-------------|-----------------------------|---------------------|
| 1:2 | 80 | 87 | 89 | 83 | 95 |
| 1:5 | 95 | 80 | 92 | 92 | 97 |
| 1:10 | 92 | 87 | - | 101 | 85 |
| 1:20 | 94 | 99 | - | 83 | 91 |

14.2 Species Cross-Reactivity

Serum of the cited species was diluted and used as sample in this assay system. No signal was detected in serum of the following species:



Horse, Cow, Chicken, Rabbit, Dog, Guinea pig, Sheep, Mouse, Goat, Donkey, Rat, Cat

Whether this obvious non-reactivity is species specific should be assessed individually by each customer.

15 LITERATUR / REFERENCES

1. Nakano, Y., et al., Isolation and characterization of GBP28, a novel gelatin-binding protein purified from human plasma. *J Biochem (Tokyo)*, 1996. 120(4): p. 803-12.
2. Pajvani, U.B., et al., Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. *J Biol Chem*, 2003. 278(11): p. 9073-85.
3. Tsao, T.S., et al., Role of disulfide bonds in Acrp30/adiponectin structure and signaling specificity. Different oligomers activate different signal transduction pathways. *J Biol Chem*, 2003. 278(50): p. 50810-7.
4. Shimada, K., T. Miyazaki, and H. Daida, Adiponectin and atherosclerotic disease. *Clin Chim Acta*, 2004. 344(1-2): p. 1-12.
5. Higashiura, K., et al., Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)*, 2004. 61(6): p. 753-9.
6. Spranger, J., et al., Adiponectin is independently associated with insulin sensitivity in women with polycystic ovary syndrome. *Clin Endocrinol (Oxf)*, 2004. 61(6): p. 738-46.
7. Zoico, E., et al., Adipocytokines, fat distribution, and insulin resistance in elderly men and women. *J Gerontol A Biol Sci Med Sci*, 2004. 59(9): p. M935-9.
8. Ye, J.M., et al., PPARalpha /gamma ragaglitazar eliminates fatty liver and enhances insulin action in fat-fed rats in the absence of hepatomegaly. *Am J Physiol Endocrinol Metab*, 2003. 284(3): p. E531-40.
9. Yamauchi, T., et al., Adiponectin stimulates glucose utilization and fatty-acid oxidation by activating AMP-activated protein kinase. *Nat Med*, 2002. 8(11): p. 1288-95.
10. Blüher (Blueher), M., et al., Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care*, 2007. 30(2): p. 280-5.
11. Delaigle, A.M., et al., Induction of adiponectin in skeletal muscle by inflammatory cytokines: in vivo and in vitro studies. *Endocrinology*, 2004. 145(12): p. 5589-97.
12. Winzer, C., et al., Plasma adiponectin, insulin sensitivity, and subclinical inflammation in women with prior gestational diabetes mellitus. *Diabetes Care*, 2004. 27(7): p. 1721-7.
13. Xydakis, A.M., et al., Adiponectin, inflammation, and the expression of the metabolic syndrome in obese individuals: the impact of rapid weight loss through caloric restriction. *J Clin Endocrinol Metab*, 2004. 89(6): p. 2697-703.
14. Motoshima, H., et al., Adiponectin suppresses proliferation and superoxide generation and enhances eNOS activity in endothelial cells treated with oxidized LDL. *Biochem Biophys Res Commun*, 2004. 315(2): p. 264-71.
15. Wolf, A.M., et al., Adiponectin induces the anti-inflammatory cytokines IL-10 and IL-1RA in human leukocytes. *Biochem Biophys Res Commun*, 2004. 323(2): p. 630-5.
16. Okamoto, Y., et al., Adiponectin reduces atherosclerosis in apolipoprotein E-deficient mice. *Circulation*, 2002. 106(22): p. 2767-70.
17. Schlegel, A., Adiponectin and risk of coronary heart disease. *Jama*, 2004. 292(1): p. 40; author reply 40.
18. Choi, K.M., et al., Inflammation, Insulin Resistance, and Glucose Intolerance in Acute Myocardial Infarction Patients without a Previous Diagnosis of Diabetes Mellitus. *J Clin Endocrinol Metab*, 2004.
19. Nakamura, Y., et al., Implications of plasma concentrations of adiponectin in patients with coronary artery disease. *Heart*, 2004. 90(5): p. 528-33.
20. Pischon, T., et al., Plasma adiponectin levels and risk of myocardial infarction in men. *Jama*, 2004. 291(14): p. 1730-7.
21. Shibata, R., et al., Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. *J Biol Chem*, 2004. 279(27): p. 28670-4.
22. Fernandez-Real, J.M., et al., Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care*, 2004. 27(3): p. 739-45.

SYMBOLS USED WITH DEMEDITEC ASSAYS

| Symbol | English | Deutsch | Française | Espanol | Italiano |
|---|------------------------------------|--|--|---|-------------------------------------|
|  | European Conformity | CE-Konformitäts-kennzeichnung | Conforme aux normes européennes | Conformidad europea | Conformità europea |
|  | Consult instructions for use | Gebrauchsanweisung beachten | Consulter les instructions d'utilisation | Consulte las Instrucciones | Consultare le istruzioni per l'uso |
|  | In vitro diagnostic device | In-vitro-Diagnostikum | utilisation Diagnostic in vitro | Diagnóstico in vitro | Per uso Diagnostica in vitro |
|  | For research use only | Nur für Forschungszwecke | Seulement dans le cadre de recherches | Sólo para uso en investigación | Solo a scopo di ricerca |
|  | Catalogue number | Katalog-Nr. | Référence | Número de catálogo | No. di catalogo |
|  | Lot. No. / Batch code | Chargen-Nr. | No. de lot | Número de lote | Lotto no |
|  | Contains sufficient for <n> tests/ | Ausreichend für "n" Ansätze | Contenu suffisant pour "n" tests | Contenido suficiente para <n> ensayos | Contenuto sufficiente per "n" saggi |
|  | Note warnings and precautions | Warnhinweise und Vorsichtsmaßnahmen beachten | Avertissements et mesures de précaution font attention | Tiene en cuenta advertencias y precauciones | Annoti avvisi e le precauzioni |
|  | Storage Temperature | Lagerungstemperatur | Température de conservation | Temperatura de conservacion | Temperatura di conservazione |
|  | Expiration Date | Mindesthaltbarkeitsdatum | Date limite d'utilisation | Fecha de caducidad | Data di scadenza |
| | Legal Manufacturer | Hersteller | Fabricant | Fabricante | Fabbricante |
| <i>Distributed by</i> | Distributed by | Vertrieb durch | Distribution par | Distribución por | Distribuzione da parte di |
|  | Version | Version | Version | Versión | Versione |
|  | Single-use | Einmalverwendung | À usage unique | Uso único | Uso una volta |