Angiotensin I RIA KIT

REF IM3518

TABLE OF CONTENTS

English	2	Magyar (HU)	35
Français (FR)	5	Polski (PL)	38
Deutsch (DE)	9	Čeština (CZ)	41
Italiano (IT)	13	Slovenčina (SK)	44
Español (ES)	16	한국어 (KO)	47
Português Portugal (PT-PT)	19	Türkçe (TR)	50
Dansk (DA)	22	Русский (RU)	53
Ελληνικά (EL)	25	Србија (SR)	57
中文 (ZH-CN)	29	APPENDIX	60
Lietuviškai (LT)	32		

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Angiotensin I RIA KIT

REF IM3518

RADIOIMMUNOASSAY OF ANGIOTENSIN I FOR THE IN VITRO DETERMINATION OF PLASMA RENIN ACTIVITY (PRA) IN HUMAN PLASMA For *in vitro* diagnostic use.

PRINCIPLE

The Angiotensin I RIA serves for the quantitative determination of plasma renin activity (PRA) by the radioimmunoassay of the product of the reaction, angiotensin I.

The generation of angiotensin I is the result of the enzymatic cleavage of the renin substrate, angiotensinogen, in plasma samples in the presence of ACE inhibitor (ACE - Angiotensin-Converting Enzyme), an enzymatic inhibitor that blocks the conversion of angiotensin I to angiotensin II.

The immunoassay of angiotensin I is a radioimmunological competition assay. Unknown samples, control, and calibrators are incubated in polyclonal antibody-coated tubes with ¹²⁵I-labeled angiotensin I as tracer. After incubation, the contents of the tubes are aspirated. The bound radioactivity is then determined in a gamma counter. A standard curve is established and unknown values are determined by interpolation from the standard curve.

WARNING AND PRECAUTIONS

General remarks:

- Enzymatic inhibitor solution, calibrators, control sample and analyzed samples must be cooled to 2-8°C before pipeting.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the immunoassay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

All plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Tracer	WARNING		
	H317		May cause an allergic skin reaction.
	H412		Harmful to aquatic life with long lasting effects.
	P273		Avoid release to the environment.
	P280		Wear protective gloves, protective clothing and eye/face protection.
	P333+P313		If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364		Take off contaminated clothing and wash it before use. reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.05%
Inhibitor	DANGER	 	
	H317		May cause an allergic skin reaction.
	H360		May damage fertility or the unborn child.
	P201		Obtain special instructions before use.
	P280		Wear protective gloves, protective clothing and eye/face protection.
	P308+P313		IF exposed or concerned: Get medical advice/attention.
	P333+P313		If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364		Take off contaminated clothing and wash it before use. 8-Hydroxyquinoline 0.1 - 0.2%
Wash Solution (20X)	DANGER		
	H360		May damage fertility or the unborn child.
	P201		Obtain special instructions before use.
	P280		Wear protective gloves, protective clothing and eye/face protection.
	P308+P313		IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - 0.3% Sodium Borate Decahydrate 0.1 - 0.3%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Plasma samples have to be collected into cold EDTA tubes.
- Separate plasma from cells by centrifugation at 2-8°C.
- Keep plasma samples frozen (<-20°C, 1 year maximum) if determination is not to be performed immediately, after aliquoting in order to avoid repeated freezing and thawing.

Note: The temperature of plasma samples must be kept at 2-8°C in the course of sampling. Avoid further manipulation to prevent both formation and decomposition of angiotensin I.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Anti-angiotensin I polyclonal antibody-coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled angiotensin I: one 11 mL vial (ready-to-use)

The vial contains 260 kBq, at the date of manufacture, of ¹²⁵I-labeled angiotensin I in buffer with bovine serum albumin and a dye.

Calibrators: six 1 mL vials (ready-to-use)

The calibrator vials contain from 0 to approximately 30 ng/mL of angiotensin I in buffer with bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against RP 86/536.

Control sample: one 1 mL vial (ready-to-use)

The vial contains angiotensin I in buffer with bovine serum albumin and sodium azide (<0.1%). The expected value is in the concentration range indicated in a supplement.

Enzymatic inhibitor: one vial (lyophilized)

It also contains sodium azide (<0.1%).

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (75 µL; 100 µL).
- adjustable dispensers (200 µL; 300 µL; 2 mL).
- water bath.
- ice bath.
- vortex-type mixer.
- horizontal or orbital shaker.
- aspiration system.
- Gamma counter set for ¹²⁵I

PROCEDURE

Preparation and storage of reagents

Preparation of enzymatic inhibitor solution

The content of the vials is reconstituted with the volume of cold distilled water (4°C) indicated on the label and mixed. The reconstituted enzymatic inhibitor may be stored at 2-8°C until the expiry date of the kit.

Preparation of wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Enzymatic step – generation of angiotensin I

Remarks and recommendations

- The enzymatic inhibitor has to be cooled to 4°C before addition to the sample.
- Both incubation temperatures (4°C and 37°C) must be adhered to strictly, even slight variations may cause severe errors in determination.

- The enzymatic incubation time at 37°C should be determined as precisely as possible and kept within narrow limits for the whole set of tubes.
- The promptness of the temperature increase from 4°C to 37°C and the following reverse drop are critical. A circulating water bath is convenient for warming and, the use of an iced-cooled water bath is advisable for cooling.
- The promptness of the temperature increase and drop may be improved by using tubes made of material with good thermal conductivity (glass).
- If low plasma renin activity of the sample is expected, the incubation time of the enzymatic step may be prolonged for up to 3 hours.

Enzymatic step – procedure

Attention: Do not treat the calibrators and the control sample.

- Add 200 µL of pre-cooled enzymatic inhibitor to 200 µL of each plasma sample and mix.
- Split each sample into two 200 µL aliquots.
- Place the first aliquot into an ice-cold water bath in a refrigerator (intended for the determination of background angiotensin I at 4°C).
- Place the second one into the water bath set for 37°C (intended for the determination of generated angiotensin I at 37°C).
- Incubate all aliquots for 1 hour.
- After incubation, cool samples from 37°C to 4°C rapidly using ice water bath.

Immunoassay procedure

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes, add successively: 75 µL of calibrator, control or sample after enzymatic incubation at 37°C and at 4°C respectively and 100 µL of tracer.** Mix.	Incubate 2 hours at 18-25°C with shaking (> 280 rpm).	Aspirate carefully the contents of 2 tubes (except the 2 tubes "total cpm"). Wash with 2 mL of wash solution. Aspirate twice. Determine activity (cpm) for 1 min.

*Calibrators, control sample and analyzed samples have to be cooled to 4°C before pipeting. Mix samples gently before they are added.

**Add 100 µL of tracer to 2 additional tubes to obtain «total cpm».

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of angiotensin I concentrations in samples assayed at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *weighted cubic regression* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

Total activity: 68,511 cpm				
Calibrators	Angiotensin I (ng/mL)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0.00	17 444	25.5	100
1	0.30	13 732	20.0	78.7
2	1.00	9 390	13.7	53.8
3	3.00	5 431	7.93	31.1
4	10.0	2 746	4.01	15.7
5	30.0	1 243	1.81	7.13

(Example of standard curve, do not use for calculation)

Samples

For the control and samples incubated at 4°C, or at 37°C, locate the B/T or the B/B₀ value on the vertical axis and read off the corresponding angiotensin I concentration in ng/mL on the horizontal axis.

Calculation of plasma renin activity

The determination of plasma renin activity is performed indirectly by the measurement of the in vitro generation of angiotensin I (A-I) per hour. Background A-I, determined on plasma samples incubated at 4°C, is subtracted from the A-I generated at 37°C for the calculation of PRA using the following equation:

$$\text{PRA ng of A-I /mL/hr} = \frac{[\text{A-I (37°C)} - \text{A-I (4°C)}] \times 2}{\text{Enzymatic incubation time (hrs)}}$$

Where

A-I (37°C): angiotensin concentration in ng/mL of sample incubated at 37°C

A-I (4°C) : angiotensin concentration in ng/mL of sample incubated at 4°C

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The PRA values presented below are indicative only.

N	Normal adult	2.5th-97.5th percentile (ng/mL/hr)	Median (ng/mL/hr)	Min-Max (ng/mL/hr)
38	Early Morning, Supine	0.32-1.84	0.79	0.30-1.90
41	Upright, 2 Hours	0.60-4.18	2.20	0.48-4.88

Detail information about expected values for children (sorted according to age) can be found in the data sheet "APPENDIX".

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.07 ng/mL

Functional sensitivity: 0.20 ng/mL

Specificity

The antibody used in the immunoassay is highly specific for angiotensin I. Extremely low cross reactivity was obtained against angiotensin II.

Moreover, the influence of possible interferences on PRA result is eliminated by subtraction of background.

Precision

Intra-assay

Samples were assayed in 25 replicates in the same series. The coefficients of variation were found below or equal to 11.3%.

Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 20.9%.

Accuracy

Dependence on time of enzymatic incubation

The samples were incubated with enzymatic inhibitor for 60, 120, and 180 minutes. No significant effect on PRA results was found.

Dilution test

Plasma samples were serially diluted in the zero calibrator. The recovery percentages were obtained between 78% and 99%.

Recovery test

Plasma samples were spiked with known quantities of angiotensin I. The recovery percentages were obtained between 104% and 123%.

Measurement range (from analytical sensitivity to highest calibrator): 0.07 to approximately 30 ng/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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No dependence of PRA on time was observed.

Time, minutes		60	120	180
PRA ng/mL/hour	Sample A	0.30	0.32	0.33
	Sample B	0.56	0.48	0.53
	Sample C	1.04	1.15	1.19
	Sample D	1.64	1.56	1.65
	Sample E	5.12	5.19	4.93

Dilution test

Three plasma samples were diluted with the zero calibrator and assayed according to the procedure of the kit.

Sample	Dilution	Theoretical conc. (ng/mL)	Observed conc. (ng/mL)	Recovery (%)
A	undiluted	-	2.45	-
	1/2	1.23	1.21	98.8
	1/4	0.61	0.56	91.4
	1/8	0.31	0.27	88.2
B	undiluted	-	3.05	-
	1/2	1.53	1.46	95.7
	1/4	0.76	0.64	83.9
	1/8	0.38	0.34	89.2
C	undiluted	-	9.21	-
	1/2	4.61	4.30	93.4
	1/4	2.30	1.82	79.0
	1/8	1.15	0.90	78.2
	1/16	0.58	0.45	78.2

Recovery test

Plasma samples were spiked with known quantities of angiotensin I and assayed according to the procedure of the kit.

Initial conc. (ng/mL)	Added angiotensin (ng/mL)	Observed conc. (ng/mL)	Observed addition (ng/mL)	Recovery (%)
0.88	1.01	1.89	2.09	111
	1.30	2.18	2.68	123
	1.93	2.81	3.20	114
1.39	1.01	2.40	2.64	110
	1.30	2.69	2.98	111
	1.93	3.32	3.59	108
1.54	1.01	2.54	2.84	112
	1.30	2.84	3.31	117
	1.93	3.46	4.03	116
2.60	1.01	2.54	2.84	112
	1.30	2.84	3.31	117
	1.93	3.46	4.03	116
2.62	1.01	3.62	3.75	104
	1.30	3.92	4.19	107
	1.93	4.55	4.82	106

Precision

Intra-assay

Samples	1	2
Number of determinations	25	25
PRA, ng/mL/hour	1.49	2.94
CV (%)	11.25	11.25

Inter-assay

After generation of Angiotensin I, samples were determined in duplicates in 10 different series according to the procedure of the kit. Plasma renin activity was obtained and used for the calculation of inter-assay precision.

Samples	A	B	C	D	E
Number of determinations	10	10	10	10	10
PRA, ng/mL/hour	0.69	3.37	7.08	15.1	24.2
CV (%)	20.9	9.57	8.72	9.74	11.8

Specificity

	Cross reactivity (%)
Angiotensin I	100
Angiotensin II	ND
Angiotensin III	ND
Tetradecapeptide	ND
Angiotensinogen	ND

Expected data for children

Results are sorted according to age.

Children Upright	N	Angiotensin (ng/mL.hr)				
		Min	Max	Median	2.5 th percentile	97.5 th percentile
2-9 years	27	0.60	7.38	3.20	0.76	6.64
10-15 years	16	0.58	3.98	1.52	0.64	3.93

125I Characteristics

T_{1/2} (125I) = 1443 h = 60.14 d

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / За ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

Σ Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

i Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expiração / 만료 날짜 / Son Kulanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

MW Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktivyvioj medžiaga / Radioaktiv / Radioaktivny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLÍKE / ОПАСНО / ОПАСНОСТ / 危險

Ag ^{125I}

Ab ^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikator / Indikator (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrator / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用說明

SOLN WASH 20x

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösingskoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalu koncentratas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztworu płuczacego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

INH

Inhibitor / Inhibiteur / Ihibitore / Inhibidor / Inibidor / Hämmare / Αναστολέας / 抑制剂 / Inhibitorius / Inhibitor / 억제인자 / Inhibitor / ингибитор / Инхибитор / 抑制劑