Free beta-hCG IRMA kit

REF IM2459

TABLE OF CONTENTS

English	2	Čeština	25
Français	4	Slovenčina	28
Deutsch	7	한국어	31
Italiano	10	Türkçe	33
Español	13	Русский	35
Ελληνικά	16	APPENDIX	38
Magyar	19		
Polski	22		

Free beta-hCG IRMA kit

REF IM2459

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF FREE β -hCG IN HUMAN SERUM OR PLASMA

For *in vitro* diagnostic use.

PRINCIPLE

The free β -hCG assay is a two-step "sandwich" type assay in which two monoclonal mouse antibodies, directed against two different epitopes of β -subunit of hCG molecule and hence not competing are employed.

The samples or calibrators are first incubated in tubes coated with the first monoclonal antibody. The contents of the tubes are then aspirated and the presence of free β -hCG in the sample is revealed by incubation with the second, ¹²⁵I labelled antibody. The contents of the tubes are aspirated and the bound radioactivity is determined in a gamma counter. The free β -hCG concentrations in the samples are obtained by interpolation from the standard curve. The concentration of free β -hCG in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash solution (20x) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in tubes without additives or in tubes containing EDTA.
- Separate serum or plasma from cells by centrifugation.
- Samples may be kept at 2-8 °C, if they are to be assayed within 24 hours. For longer storage keep frozen (at <-20 °C, 8 months maximum) in aliquots to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature. Do not leave the samples at room temperature, if not necessary.
- If samples have concentrations greater than the highest calibrator, they must be diluted with zero calibrator.

Serum and EDTA plasma values for 20 samples (serum values ranging from 3.94 to 18.69 ng/mL) were compared using the IM2459 Free β -hCG IRMA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.8674[\text{serum}] + 0.0827$$

$$R = 0.9961$$

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8 °C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Anti-free β -hCG monoclonal antibody-coated tubes: 2 x 50 tubes (ready-to-use)

Monoclonal ¹²⁵I-labelled anti-free β -hCG tracer antibody: one 22 mL vial (ready-to-use).

The vial contains 650 kBq at the date of manufacture of ¹²⁵I-labelled immunoglobulin in buffer containing bovine serum albumin, sodium azide (<0.1 %) and a red dye.

Free β -hCG calibrators: five 1.0 mL vials + one 5 mL vial with zero calibrator (ready-to-use).

The calibrator vials contain from 0 to approximately 200 ng/mL of free β -hCG in bovine serum and/or buffer containing bovine serum albumin with sodium azide (<0.1 %). The exact concentrations are indicated on each vial label. The calibrators were calibrated using the 1st IRP WHO 75/551 standard (1mIU = 1 ng).

Assay buffer (Free β -hCG - Buffer): one 11 mL vial; ready for use.

The vial contains buffer with bovine serum albumin and sodium azide (<0.1%) and a blue dye.

Control samples: two vials (lyophilised).

The vials contain free β -hCG lyophilised in human serum with sodium azide (<0.1 %). The expected values are in the concentration range indicated on a supplement.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (100 µL)
- semi-automatic pipets (100 µL, 200 µL and 2 mL)
- vortex type mixer
- orbital shaker
- aspiration system
- gamma counter set for 125 iodine

RESULTS**Standard curve**

The results in this package insert were calculated using a log-log curve fit ("spline" mode) with determined radioactivity (cpm_{cal} - $\text{cpm}_{\text{cal}0}$) on vertical axis and the free β -hCG concentration of the calibrators on the horizontal axis (ng/mL). Other data reduction methods may give slightly different results.

Total activity: 195,274 cpm				
Calibrators	Free β -hCG (ng/mL)	cpm (n=3)	B/Bmax	$\text{cpm}_{\text{cal}} - \text{cpm}_{\text{cal}0}$
0	0.0	40	-	-
1	0.32	959	0.53	919
2	1.60	4,507	2.57	4,467
3	8.00	21,401	12.26	21,361
4	40.0	87,281	50.08	87,241
5	200	174,243	100	174,203

(Example of standard curve, do not use for calculation)

Samples

Locate the ratio $\text{cpm}_{\text{sample}} - \text{cpm}_{\text{cal}0}$ for each sample on the vertical axis of the standard curve, read-off the free β -hCG concentration of the sample on the horizontal axis in ng/mL. The concentrations of diluted samples must be corrected by the dilution factor.

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

	N	Min, ng/mL	Max, ng/mL	Mean, ng/mL	Median, ng/mL	2.5th-97.5th percentile, ng/mL
Male and Female	100	0.00	0.23	0.03	0.02	0.00-0.13

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended to analyse their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE**Preparation of reagents**

Let all the reagents come to room temperature.

Reconstitution of control samples

The content of the vials is reconstituted with the volume of distilled water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8 °C for 1 week or at <-18 °C, until the expiry date of the kit.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8 °C until the expiry date of the kit.

Summary of assay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions	Step 2 1 st incubation
To coated tubes add successively: 100 µL of calibrator, control or sample and 100 µL of assay buffer. Mix.	Incubate for 1 hour at 18-25 °C with shaking at >280 rpm. Aspirate carefully the contents of each tube.

Step 3 2 nd incubation	Step 4 Washing and counting
Add 200 µL of tracer* to all tubes. Mix. Incubate for 30 min at 18-25 °C with shaking at >280 rpm.	Aspirate carefully the contents of each tube (except of 2 tubes for total activity T). Wash the tubes with 2 mL of wash solution and aspirate. Repeat washing and aspiration. Count bound cpm (B) and total cpm (T) for 1 min.

*Add 200 µL of tracer to 2 additional tubes to obtain total cpm (T).

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.03 ng/mL

Functional sensitivity: 0.12 ng/mL

Specificity

The antibodies used in this kit exhibit no cross-reaction with human TSH, FSH, LH and intact HCG.

Precision**Intra-assay**

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 4.0 %.

Inter-assay

Serum samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 9.5 %.

Accuracy**Dilution test**

High-concentration serum samples were serially diluted. The recovery percentages obtained were between 83.0 % and 102 %.

Recovery test

Low-concentration serum samples were spiked with known quantities of free β -hCG. The recovery percentages obtained were between 90.6 % and 104%.

Measurement range (from analytical sensitivity to highest calibrator): 0.03 to approximately 200 ng/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Accuracy

Dilution test

Serum samples were diluted in the zero calibrator and assayed according to the procedure of the kit.

Serum	Dilution factor	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Ratio (%) Measured/Expected
S1	Undil.	53.11	-	-
	1/2	22.64	26.56	85.26
	1/4	11.02	13.28	83.00
	1/8	5.65	6.64	85.11
	1/16	2.82	3.32	84.96
	1/32	1.41	1.66	84.96
S2	Undil.	38.87	-	-
	1/2	18.38	19.44	94.57
	1/4	9.04	9.72	93.03
	1/8	4.45	4.86	91.59
	1/16	2.20	2.43	90.56
	1/32	1.06	1.21	87.27
S3	Undil.	27.09	-	-
	1/2	13.46	13.55	99.37
	1/4	6.85	6.77	101.1
	1/8	3.46	3.39	102.2
	1/16	1.69	1.69	99.82
	1/32	0.83	0.85	98.04
S4	Undil.	33.05	-	-
	1/2	15.68	16.53	94.89
	1/4	7.56	8.26	91.50
	1/8	4.05	4.13	98.03
	1/16	1.95	2.07	94.40
	1/32	0.94	1.03	91.01
S5	Undil.	29.84	-	-
	1/2	14.63	14.92	98.06
	1/4	7.45	7.46	99.87
	1/8	3.56	3.73	95.44
	1/16	1.80	1.87	96.51
	1/32	0.85	0.93	91.15

Plasma samples were diluted in the zero calibrator and assayed according to the procedure of the kit.

EDTA Plasma	Dilution factor	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Ratio (%) Measured/Expected
P1	Undil.	42.97	-	-
	1/2	19.16	21.49	89.18
	1/4	9.66	10.74	89.92
	1/8	4.96	5.37	92.34
	1/16	2.37	2.69	88.25
	1/32	1.22	1.34	90.85
P2	Undil.	33.82	-	-
	1/2	16.21	16.91	95.86
	1/4	8.09	8.46	95.68
	1/8	3.91	4.23	92.49
	1/16	1.80	2.11	85.16
	1/32	0.94	1.06	88.94
P3	Undil.	27.67	-	-
	1/2	13.65	13.84	98.66
	1/4	6.89	6.92	99.60
	1/8	3.41	3.46	98.59
	1/16	1.67	1.73	96.57
	1/32	0.79	0.86	91.36
P4	Undil.	35.74	-	-
	1/2	16.86	17.87	94.35
	1/4	8.61	8.94	96.36
	1/8	4.09	4.47	91.55
	1/16	2.01	2.23	89.98
	1/32	0.98	1.12	87.74
P5	Undil.	32.41	-	-
	1/2	15.42	16.21	95.16
	1/4	7.63	8.10	94.17
	1/8	3.64	4.05	89.85
	1/16	1.90	2.03	93.80
	1/32	0.94	1.01	92.81

Recovery test

Free β -hCG was added to 5 serum samples and assayed according to the procedure of the kit.

Serum	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
S1	0.31	0.78	1.10	1.04	94.71
	0.31	1.09	1.40	1.45	103.5
	0.31	1.54	1.85	1.80	97.50
	0.73	0.78	1.51	1.51	100.0
S2	0.72	1.09	1.81	1.64	90.64
	0.71	1.54	2.25	2.27	100.9
	0.88	0.78	1.67	1.69	101.4
S3	0.87	1.54	2.40	2.30	95.68
	0.85	2.26	3.11	2.99	96.04
	1.16	0.78	1.94	1.85	95.30
S4	1.13	1.54	2.67	2.67	99.88
	1.11	2.26	3.38	3.36	99.49
	1.41	0.78	2.20	2.28	103.8
S5	1.38	1.54	2.92	2.93	100.2
	1.33	2.96	4.30	4.06	94.50

Free β -hCG was added to 5 EDTA-plasma samples and assayed according to the procedure of the kit.

EDTA plasma	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
P1	3.58	3.92	7.50	7.34	97.87
	3.55	5.45	9.00	9.26	102.9
	3.51	7.69	11.20	10.90	97.30
P2	0.33	0.78	1.12	1.01	90.37
	0.33	1.09	1.42	1.39	97.87
	0.33	1.54	1.87	1.75	93.81
P3	0.49	0.78	1.27	1.21	94.94
	0.48	1.54	2.02	1.93	95.58
	0.47	2.26	2.74	2.49	91.01
P4	0.65	0.78	1.43	1.41	98.51
	0.63	1.54	2.17	2.13	98.02
	0.62	2.96	2.89	2.76	95.61
P5	0.96	0.78	1.75	1.68	96.27
	0.94	1.54	2.48	2.44	98.36
	0.91	2.96	3.87	3.65	94.31

¹²⁵I Characteristics

$$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$$

¹²⁵ I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Уробца / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Rækker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 内容物足夠執行 <n> 次测试

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문외 / Kullanna Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sicaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Urozornění / Urozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / ТЕΗΛΙΚΕ / ОПАСНО / ОПАСНОСТ / 危險

Ag ^{125I}

Ab ^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Ateksamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL **0**

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibrátor / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controlla / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontroliné / Kontroll / Kontrola / Kontrola / Kontrola / 청도관리 / Kontrol / Контроль / Контролна / 质控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mėgintuvėliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Túpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

SOLN **WASH** **20x**

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösningskoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 濃縮清洗液 20X / Pļavimo tirpalo koncentrātas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztworu płuczącego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

BUF

Buffer, Tampon, Puffer, Tampone, Tampón, Tampão, Buffert, Ρυθμιστικό διάλυμα, 缓冲液, Buferinis tirpalas, Puffer, Bufor, Pufur, Tlmitivý roztok, 완충, Tampon, Буфер, буфер, 緩衝劑