C-peptide IRMA KIT

REF IM3639

TABLE OF CONTENTS

English	2	Slovenčina	29
Français	5	한국어	32
Deutsch	8	Türkçe	35
Italiano	11	Русский	38
Español	14	中文 ZH-TW	41
Ελληνικά	17	Србија	44
Lietuviškai	20	APPENDIX	47
Polski	23		
Čeština	26		

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C-peptide IRMA KIT

REF IM3639

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF C-PEPTIDE IN HUMAN SERUM, PLASMA AND URINE

For *in vitro* diagnostic use.

PRINCIPLE

The immunoradiometric assay of C-peptide is a "sandwich" type assay. The same kit may be employed for the measurement of C-peptide either directly in serum and plasma or, after dilution, in urine.

In the kit, mouse monoclonal antibodies directed against two different epitopes of C-peptide and hence not competing are used.

Serum, plasma or diluted urine samples, the controls and calibrators are incubated in tubes coated with the first monoclonal antibody in the presence of the second monoclonal antibody labeled with iodine 125. After incubation, the contents of the tubes are rinsed so as to remove unbound ¹²⁵I-labeled antibody. The bound radioactivity is then determined in a gamma counter. The C-peptide concentrations in the samples are obtained by interpolation from the standard curve. The concentration of C-peptide in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as short time as possible to avoid evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Material of human origin

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

Collection of serum and plasma samples

- Collect blood in dry tubes or in tubes containing EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -20°C for 2 months or preferably at < -70°C for 1 year maximum) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.

Human C-peptide is relatively unstable molecule in liquid serum or plasma (the considerable degradation may be observed already in a few day period). The recommendation concerning storage should be therefore strictly adhered.

- If samples have concentrations greater than the highest calibrator, they must be diluted in the dilution buffer.

Serum and EDTA plasma values for 20 samples (serum values ranging from 576.1 to 3,739 pM) were compared using the IM3639 C-Peptide IRMA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.9623[\text{serum}] + 49.99$$

$$R = 0.9960$$

Collection of urine samples

- Collect 24-hour urine in flask.
- Determine volume.
- Centrifuge or filter before assay, if necessary.
- Samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage, keep frozen at < -18°C for 1 year maximum, after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- The content of C-peptide in urine samples is often over the concentration of the higher calibrator. Therefore dilution 1:21 with dilution buffer is recommended, e.g. 25 µL of urine should be diluted with 500 µL of dilution buffer.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Anti-C-peptide monoclonal antibody coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal anti-C-peptide antibody: one 17 mL vial (ready-to-use)

The vial contains 640 kBq, at the date of manufacture, of ¹²⁵I-labeled immunoglobulins in buffer containing bovine serum albumin, sodium azide (<0.1%) and a dye.

C-peptide dilution buffer: one 25 mL vial (ready-to-use)

The solution contains bovine serum albumin and sodium azide (<0.1%).

Calibrators: six vials (lyophilized)

The calibrator vials contain from 0 to approximately 6,400 pM of C-peptide in buffer with bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against the international standard, WHO IRR C-Peptide 84/510.

Control sera: two vials (lyophilized)

The vials contain C-peptide lyophilized in bovine serum and sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (50 µL; 25 µL).
- semi-automatic pipet (150 µL).
- vortex-type mixer.
- horizontal or orbital shaker (> 280 rpm).
- aspiration system.
- gamma counter set for ¹²⁵I.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of C-peptide concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{cal} - cpm_{cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (pM).

Other data reduction methods may give slightly different results.

Total activity: 227,129 cpm				
Calibrators	C-peptide (pM)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	53	0.02	-
1	32	985	0.43	932
2	103	2,935	1.29	2,882
3	375	10,863	4.78	10,810
4	1,530	42,652	18.78	42,599
5	6,100	118,237	52.06	118,184

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

The concentrations of diluted samples must be corrected by the dilution factor.

To convert pmol/L (pM) into ng/mL, multiply results by **0.003**.

24 hours urine excretion

For diluted urines (dilution 1:21), the following formula must be used for the determination of 24 hours C-peptide excretion:

$$\text{C-Peptide (nmol/24h)} = \frac{\text{measured concentration (pM)} \times 21 \times V}{1000}$$

where V = 24 h urine volume in litres

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

Sera

- The serum samples of 156 fasting healthy individuals were analyzed.

	Mean	Median	Standard Deviation (SD)	Min-Max
pM	519	471	255	32.6-1,458
ng/mL	1.56	1.41	0.766	0.098-4.380

The 95 percentile range was **160-1,100 pM** (0.48-3.30 ng/mL).

- As an example, the average C-peptide concentrations in the OGTT of twelve young, healthy individuals are presented in the following table.

Time (min)	C-peptide (pM) Average ± 1 SD	C-peptide (ng/mL) Average ± 1 SD
0	520 ± 169	1.56 ± 0.507
35	2,274 ± 1,116	6.82 ± 3.348
65	3,067 ± 987	9.20 ± 2.961
95	3,488 ± 1,038	10.46 ± 3.114
125	3,008 ± 593	9.02 ± 1.779

Urines

- 48 samples of the first morning urine and 20 samples of a 24-hour urine collection were determined.

	Mean	Standard deviation	Min-Max
First morning urine (nmol/L)	20.12	16.70	1.22 - 66.39
24 hour urine collection (nmol/24 hrs)	23.13	12.77	7.76 - 49.38

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation and storage of reagents

Let all the reagents come to room temperature.

Reconstitution of calibrators

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one week or frozen at < -18°C for a longer time, until the expiry date of the kit.

Reconstitution of control samples

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one day or aliquoted at < -18°C for a longer time, until the expiry date of the kit.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Bring all reagents to room temperature before pipeting.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To coated tubes, add successively: 50 µL of calibrator, control or sample (serum, plasma, pre-diluted urine) and 150 µL of tracer. Mix.	Incubate 2 hours at 18-25°C with shaking (> 280 rpm).	Aspirate carefully the contents of tubes (except the 2 tubes «total cpm»). Wash twice with 2 mL of wash solution. Count activity (cpm) for 1 min.

*Add 150 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 3.79 pM (0.011 ng/mL)

Functional sensitivity: 13.9 pM (0.042 ng/mL)

Specificity

No cross-reactivity was observed with insulin, pancreatic polypeptide, somatostatin or glucagon. The cross-reaction with proinsulin is only 3% (w/w, equivalent to 10% on a molar basis).

Precision

Intra-assay

Samples were assayed in 25 replicates in the same series. The coefficients of variation were found below or equal to 3.1% for serum samples and to 3.2% for urine samples.

Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 5.2% for serum samples and to 4.8% for urine samples.

Accuracy

Dilution test

High-concentration serum and urine samples were serially diluted in the dilution buffer. The recovery percentages were obtained between 96% and 114% for serum samples and between 86% and 105% for urine samples.

Recovery test

Serum and urine samples were spiked with known quantities of C-peptide. The recovery percentages were obtained between 83% and 105% for serum samples and between 80% and 104% for urine samples.

Measurement range (from analytical sensitivity to highest calibrator): 3.79 to approximately 6,400 pM.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

"Hook effect" of this IRMA kit may appear in samples containing C-peptide concentration equal or higher than 50,000 pM. The samples with clinical suspicion to very high concentrations should be diluted with dilution buffer. Values obtained must be corrected by the dilution factor.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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