IRMA CA 125 antigen

REF IM2233, B64312

TABLE OF CONTENTS

English	2	Čeština	19
Français	4	Slovenčina	22
Deutsch	7	한국어	25
Italiano	10	Русский	27
Español	13	APPENDIX	30
Ελληνικά	16		

IRMA CA 125 antigen

REF IM2233, B64312

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF CA 125 ANTIGEN IN HUMAN SERUM AND PLASMA For *in vitro* diagnostic use.

PRINCIPLE

The CA 125 antigen assay is a one-step "sandwich" type assay in which two mouse monoclonal antibodies, directed against two different epitopes of the molecule, are employed.

Samples or calibrators are incubated in tubes, coated with the first monoclonal antibody, in presence of the second ¹²⁵I-labeled monoclonal antibody. Following incubation, the contents of the tubes are aspirated and unbound labeled antibody is eliminated by washing. The amount of bound reactivity measured in a gamma counter is proportional to the CA 125 concentration. The unknown values are determined by interpolation from a standard curve.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Material of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or heparinized tubes.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8 °C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at <-20 °C, 3 months maximum). Samples should be aliquoted to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Dilute samples showing a concentration exceeding that of the highest calibrator with zero calibrator.

Serum and heparin plasma values for 15 samples (serum values ranging from 31.1 to 117 U/mL) were compared using the IM2233 CA125 antigen IRMA kit. Results are as follows:

[Heparin-plasma] = 1.1666 [serum] - 2.764

R = 0.9868

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8 °C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Kit for determination of CA 125, 100 tubes (Cat. #IM2233)

Anti-CA 125 antibody-coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal antibody: one 33 mL vial (ready-to-use)

The vial contains 480 kBq at time of manufacture of ¹²⁵I- labeled monoclonal antibody with bovine serum albumine, sodium azide (<0.1%) and a dye.

Calibrators: four 1 mL vials + two zero calibrator 5 mL vials (ready-to-use)

The calibrator vials contain from 0 to approximately 500 U/mL of human CA 125 in buffer containing bovine serum albumine and sodium azide (<0.1%). The exact concentration is indicated on each vial label. After opening, the vials can be stored two weeks at 2-8 °C. For longer periods, store at <-18 °C after aliquoting. Calibrators are verified to an internal reference standard.

Control serum: 1 mL vial (ready-to-use)

The vial contains human CA 125 in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement. After opening, the vial can be stored two weeks at 2-8 °C. For longer periods, store at <-18 °C after aliquoting.

Wash solution (20 X): one 50 mL vial

Concentrated solution, has to be diluted with 950 mL of distilled water before use.

Kit for determination of CA 125, 50 tubes (Cat. #B64312)**Anti-CA 125 antibody-coated tubes: 1 x 50 tubes** (ready-to-use)**¹²⁵I-labeled monoclonal antibody: one 33 mL vial** (ready-to-use)**Calibrators: four 1 mL vials + two zero calibrator 5 mL vials** (ready-to-use)**Control serum: 1 mL vial** (ready-to-use)**Wash solution (20 X): one 50 mL vial****MATERIALS REQUIRED, BUT NOT PROVIDED**

In addition to standard laboratory equipment, the following items are required:

- precision micropipet (100 µL).
- repeating micropipets (300 µL, 2 mL).
- "vortex" type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of CA 125 concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the package insert were calculated using a log-log curve fit ("spline" mode) with B/T (%) or B/Bmax (%) on vertical axis and the CA 125 concentration of the calibrators on the horizontal axis (U/mL). Other data reduction methods may give slightly different results.

Total activity: 211,761 cpm				
Calibrators	CA 125 (U/mL)	cpm	B/Bmax (%)	B/T (%)
0	0	51	0.1	
1	15	2,296	4.0	1.1
2	50	7,008	12.5	3.3
3	200	26,787	47.9	12.6
4	500	55,823	100.0	26.4

*(Example of standard curve, do not use for calculation)***Samples**

Locate the B/T (%) or the B/Bmax (%) for each sample on the vertical axis and read off corresponding CA 125 concentration in U/mL on the horizontal axis.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The following values obtained from healthy subjects are indicative only.

Number of subjects	Average (x)	95th percentile	99th percentile
181	11 U/mL	31.2 U/mL	55.3 U/mL

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use one of the following e-mail addresses : In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE**Preparation of reagents**

Let all the reagents come to room temperature.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8 °C until the expiry date of the kit.

Step 1 Immunological step*	Step 2 Washing	Step 3 Counting
To antibody coated tubes add successively: 100 µL of calibrator, control or sample and 300 µL of tracer.** Vortex gently. Incubate 4 hours at 18-25 °C with shaking (>280 rpm)	Aspirate carefully the content of each tube (except "total cpm") Wash three times with 2 mL of wash solution. Immediately aspirate the contents of tubes.	Count for 1 minute bound cpm (B) and total cpm (T)

*An incubation time of 2 hours is sufficient if the test is performed automatically. 2 hour incubation values may be slightly lower (see appendix, § Correlation of 2-hour and 4-hour incubation procedure)
**Add 300 µL of tracer to 2 additional tubes to obtain "total cpm".

PERFORMANCE CHARACTERISTICS*(for more details, see the data sheet "APPENDIX")*

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity**Analytical sensitivity:** 0.69 U/mL**Functional sensitivity:** 3.09 U/mL**Specificity:** The kit is specific of CA 125.**Precision****Intra-assay**

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 3.11 %.

Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 4.57 %.

Accuracy**Dilution test**

High-concentration samples were serially diluted with the kit zero calibrator. The recovery percentages obtained were between 97.3 % and 116 %.

Recovery test

Low-concentration samples were spiked with know quantities of CA 125. The recovery percentages obtained were between 90.6 % and 101 %.

Measurement range (from analytical sensitivity to highest calibrator):

0.69 to approximately 500 U/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

Shortage of incubation time to 2 hours was tested on SR300 instrument. Performance characteristics of the assay are not guaranteed if different automate is used.

Interference

The assay is not affected by triglycerids (up to 7 mM). High concentrations of bilirubin (>100 µM) and hemoglobin (>0.5 g/L) may influence the assay.

Hook effect

No hook effect was observed until 20,000 U/mL.

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Accuracy

Dilution test

Three serum/heparinized plasma samples were diluted in zero calibrator and assayed according to the assay procedure of the kit.

Serum	Dilution	CA-125 (U/mL)		Ratio (%) Measured/ Expected
		Measured concentration	Expected	
S1	-	111.4	-	-
	1:2	56.37	55.68	101.2
	1:4	30.72	27.84	110.3
	1:8	15.45	13.92	111.0
	1:16	7.61	6.96	109.3
	1:32	3.43	3.48	98.6
S2	-	129.1	-	-
	1:2	70.60	64.53	109.4
	1:4	37.28	32.27	115.5
	1:8	18.58	16.13	115.2
	1:16	8.85	8.07	109.7
S3	-	299.3	-	-
	1:2	145.7	149.7	97.3
	1:4	73.53	74.83	98.3
	1:8	37.01	37.42	98.9
	1:16	19.05	18.71	101.8
	1:32	9.26	9.35	99.0

Heparin plasma	Dilution	CA-125 (U/mL)		Ratio (%) Measured/ Expected
		Measured	Expected	
		Concentration		
H1	-	226.8	-	-
	1:2	118.9	113.4	104.9
	1:4	62.20	56.70	109.7
	1:8	31.65	28.35	111.6
	1:16	15.82	14.18	111.6
H2	-	160.1	-	-
	1:2	82.09	80.05	102.5
	1:4	41.91	40.03	104.7
	1:8	21.57	20.01	107.8
	1:16	10.03	10.01	100.2
H3	-	274.0	-	-
	1:2	131.5	137.0	96.0
	1:4	70.89	68.51	103.5
	1:8	36.24	34.25	105.8
	1:16	18.17	17.13	106.1
H3	1:32	8.69	8.56	101.5

Recovery test

To three sera/heparinized plasma containing low CA 125 concentrations were added known concentrations of CA 125. Samples are assayed according to the kit procedure.

Serum	Endogen. conc. (U/mL)	Added conc. (U/mL)	Expected conc. (U/mL)	Measured conc. (U/mL)	Ratio (%) Measured/ Expected
S1	78.03	5.93	83.96	80.59	95.98
	75.93	19.23	95.16	92.09	96.77
	71.79	45.45	117.2	106.2	90.59
S2	49.99	5.93	55.92	56.60	101.2
	48.64	19.23	67.88	66.10	97.38
	45.99	45.45	91.45	85.57	93.57
S3	115.9	5.93	121.9	121.8	99.94
	112.8	19.23	132.0	132.0	99.93
	106.7	45.45	152.1	148.4	97.52

Heparin plasma	Endogen. conc. (U/mL)	Added conc. (U/mL)	Expected conc. (U/mL)	Measured conc. (U/mL)	Ratio (%) Measured/ Expected
H1	59.57	5.93	65.49	65.05	99.32
	57.96	19.23	77.19	75.52	97.83
	54.80	45.45	100.3	96.22	95.98
H2	124.4	5.93	130.3	129.7	99.57
	121.0	19.23	140.3	141.1	100.6
	114.4	45.45	159.9	162.9	101.9
H3	242.9	5.93	248.8	248.8	100.0
	236.3	19.23	255.5	253.5	99.20
	223.4	45.45	268.9	272.5	101.3

Correlation of 2-hour and 4-hour incubation procedure

Values of 105 samples (ranging from 19.6 to 483.9 U/mL) were determined using standard 4-hour and shortened 2-hour procedure on Stratec SR 300. Results were as follows:

$$[2\text{-hour procedure}] = 0.8918 \times [4\text{-hour procedure}] - 0.2198$$

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