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IMMUNOTECH s.r.o., Radiova 1, 102 27 Prague 10, Czech Republic

IRMA MUC-1 gene associated antigen (CA 15-3)

REF IM2397, B64314

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF MUC-1 GENE ASSOCIATED ANTIGEN (CA 15-3) IN HUMAN SERUM AND PLASMA For *in vitro* diagnostic use.

PRINCIPLE

The MUC-1 gene associated antigen (CA 15-3) assay is a two-step "sandwich" type assay in which two mouse monoclonal antibodies, directed against two different epitopes of the molecule, are employed.

Samples or calibrators are incubated in tubes coated with the first monoclonal antibody, the contents of the tubes are then aspirated and the presence of CA 15-3 in the sample is revealed by incubation with a second, ¹²⁵I-labeled monoclonal antibody. The contents of the tubes are aspirated after this second incubation and unbound labeled antibody is eliminated by washing. The amount of bound reactivity measured in a gamma counter is proportional to the CA 15-3 concentration. The unknown values are determined by interpolation from a standard curve.

This device is indicated for use in the measurement of CA 15-3 antigen to aid in the management of breast cancer patients. Elevated concentrations may be observed in the serum or plasma of patients with benign conditions or other non-cancer disorders, as well as in breast cancer and other malignant diseases.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Material of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes containing EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -20°C, 3 months maximum). Samples should be aliquoted to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Dilute samples showing a concentration exceeding that of the highest calibrator with kit diluent.

Serum and EDTA plasma values for 15 samples (serum values ranging from 11.86 to 165.2 U/mL) were compared using the IM2397 IRMA MUC-1 gene associated antigen kit. Results are as follows:

[EDTA-plasma] = 1.1233 [serum] -1.1026

R = 0.9942

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Kit for determination of CA 15-3, 100 tubes (Cat #IM2397)

Anti-CA 15-3 antibody-coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal antibody: one 22 mL vial (ready-to-use)

The vial contains 600 kBq at time of manufacture of ¹²⁵I- labeled monoclonal antibody with bovine serum albumine, sodium azide (<0.1%) and a dye.

Calibrators: four 2 mL vials + one zero 4 mL vial (ready-to-use)

The calibrator vials contain between 0 to approximately 250 U/mL of human CA 15-3 in buffer containing bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. After opening, the vials can be stored at 2-8°C until the expiration date of the kit. Calibrators are verified to an internal reference standard.

Control serum: one 1 mL vial (ready-to-use)

The vial contains human CA 15-3 in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement. When opened, the vial may be stored at 2-8°C for 2 weeks. For longer periods, freeze at < -18°C.

Diluent: one 50 mL vial (ready-to-use)

The vial contains bovine serum albumin in buffer. After opening, the vial can be stored at 2-8°C until the expiration date of the kit.

Kit for determination of CA 15-3, 50 tubes (Cat #B64314)

Anti-CA 15-3 antibody-coated tubes: 1 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal antibody: one 22 mL vial (ready-to-use)

Calibrators: four 2 mL vials + one zero 4 mL vial (ready-to-use)

Control serum: one 1 mL vial (ready-to-use)

Diluent: one 50 mL vial (ready-to-use)

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (10 µL, 200 µL, 500 µL)

- repeating micropipets (200 µL, 2 mL).
- "vortex" type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of CA 15-3 concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{cal} - cpm_{cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (U/mL).

Other data reduction methods may give slightly different results.

| Total activity: 241,627 cpm | | | | |
|-----------------------------|----------------|-----------|---------|--------------------------|
| Calibrators | CA 15-3 (U/mL) | cpm (n=3) | B/T (%) | $cpm_{cal} - cpm_{cal0}$ |
| 0 | 0 | 215 | 0.09 | - |
| 1 | 15 | 7,414 | 3.07 | 7,199 |
| 2 | 53 | 21,360 | 10.08 | 21,145 |
| 3 | 143 | 42,813 | 17.72 | 42,598 |
| 4 | 279 | 63,497 | 26.28 | 63,282 |

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

The concentrations of diluted samples must be corrected by the dilution factor.

Note: the concentration indicated on each calibrator vial label already takes into account the 1:51 predilution which therefore must be neglected when calculating the concentration in a non-diluted sample.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The following values obtained from healthy subjects are indicative only.

| Number of subjects | Average (x) | 95 th percentil | 99 th percentil |
|--------------------|-------------|----------------------------|----------------------------|
| 160 | 20.3 U/mL | 36 U/mL | 43 U/mL |

Diagnostic sensitivity: To evaluate the clinical utility, 170 samples of patients with non-metastatic breast cancer prior to any treatment were assayed. Within this population, values above 35 U/mL were found in 49 patients. That corresponds to 28.8% clinical sensitivity of the marker at the time of diagnosis of breast cancer.

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation of reagents

Let the reagents come to room temperature.

Dilution of samples and controls

Dilute samples and control to be assayed to 1/51, by adding 10 µL of sample or control into plastic or glass tubes and then adding to each tube 500 µL of diluent. Shake gently before the assay.

DO NOT DILUTE CALIBRATORS

Assay procedure

Bring all reagents to room temperature before pipeting.

| Immunological step 1 | Immunological step 2 | Counting |
|---|---|--|
| To coated tubes add 200 µL of calibrator, of diluted sample or of diluted control | Aspirate the content of each tube, then wash twice with 2 mL of distilled water | Aspirate the content of each tube, then wash twice with 2 mL of distilled water (except the 2 " total cpm ") |
| | Carefully aspirate, then add 200 µL of tracer in each tube | Carefully aspirate |
| Incubate 2 hours at 18-25°C with shaking (>280 rpm). | Incubate 1 hour at 18-25°C with shaking (>280 rpm). | Count for 1 minute bound cpm (B) and total cpm (T) |

*Add 200 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.15 U/mL

Functional sensitivity: 5.34 U/mL

Specificity: The kit is specific of the CA 15-3.

Precision

Intra-assay:

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 6.2%.

Inter-assay:

Serum samples were assayed in duplicate in 25 different series. Coefficients of variation were found below or equal to 8.4%.

Accuracy

Dilution test

High-concentration serum samples were serially diluted with the kit diluent. The recovery percentages obtained were between 80.9% and 101%.

Recovery test

Low-concentration serum samples were spiked with known quantities of CA 15-3. The recovery percentages obtained were between 86% and 108%.

Measurement range (from analytical sensitivity to highest calibrator):

0.15 to approximately 250 U/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

Interference

The assay is not affected by hemoglobin (5 mg/mL), triglycerids (30 mg/mL) and bilirubin (0.4 mg/mL).

Hook effect

There is no hook effect when the two-step procedure is used.

IRMA MUC-1 gene associated antigen (CA 15-3)

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Specificity

The immunoassay is specific of the MUC-1 gene associated antigen (CA15-3).

Precision

Intra-assay

| Serum | S1 | S2 | S3 |
|--------------------------|-------|-------|-------|
| Number of determinations | 25 | 25 | 25 |
| Mean (U/mL) | 24.95 | 47.43 | 197.9 |
| CV (%) | 4.14 | 6.15 | 5.00 |

| EDTA-Plasma | P1 | P2 | P3 |
|--------------------------|-------|-------|-------|
| Number of determinations | 25 | 25 | 25 |
| Mean (U/mL) | 20.33 | 24.04 | 102.3 |
| CV (%) | 3.26 | 6.57 | 4.81 |

Inter-assays

| Serum | S1 | S2 | S3 |
|--------------------------|-------|-------|-------|
| Number of determinations | 25 | 25 | 25 |
| Mean (U/mL) | 21.97 | 51.88 | 183.4 |
| CV (%) | 5.02 | 8.41 | 6.39 |

| EDTA-Plasma | P1 | P2 | P3 |
|--------------------------|-------|-------|-------|
| Number of determinations | 10 | 10 | 10 |
| Mean (U/mL) | 19.07 | 66.78 | 126.1 |
| CV (%) | 11.44 | 13.62 | 14.30 |

Accuracy

Dilution test

Three serum/EDTA plasma samples were diluted in kit diluent and assayed according to the assay procedure of the kit.

| Serum | Dilution | Measured conc. (U/mL) | Expected conc. (U/mL) | Ratio (%) Measured/Expected |
|-------|----------|-----------------------|-----------------------|-----------------------------|
| S1 | - | 268.4 | - | - |
| | 1:2 | 122.5 | 134.2 | 91.33 |
| | 1:4 | 57.85 | 67.09 | 86.23 |
| | 1:8 | 28.55 | 33.54 | 85.11 |
| | 1:16 | 15.46 | 16.77 | 92.18 |
| | 1:32 | 8.19 | 8.39 | 97.66 |
| S2 | - | 339.5 | - | - |
| | 1:2 | 171.7 | 169.8 | 101.2 |
| | 1:4 | 75.07 | 84.88 | 88.45 |
| | 1:8 | 36.45 | 42.44 | 85.89 |
| | 1:16 | 17.17 | 21.22 | 80.92 |
| | 1:32 | 9.29 | 10.61 | 87.56 |
| S3 | - | 228.3 | - | - |
| | 1:2 | 103.6 | 114.2 | 90.76 |
| | 1:4 | 51.42 | 57.08 | 90.09 |
| | 1:8 | 25.79 | 28.54 | 90.37 |
| | 1:16 | 14.22 | 14.27 | 99.65 |
| | 1:32 | 7.15 | 7.13 | 100.2 |

| EDTA Plasma | Dilution | Measured conc. (U/mL) | Expected conc. (U/mL) | Ratio (%) Measured/Expected |
|-------------|----------|-----------------------|-----------------------|-----------------------------|
| P1 | - | 214.4 | - | - |
| | 1:2 | 94.73 | 107.2 | 88.38 |
| | 1:4 | 45.03 | 53.59 | 84.03 |
| | 1:8 | 22.50 | 26.80 | 83.97 |
| P2 | - | 251.3 | - | - |
| | 1:2 | 133.8 | 125.7 | 106.4 |
| | 1:4 | 63.06 | 62.84 | 100.4 |
| | 1:8 | 30.58 | 31.42 | 97.33 |
| P3 | - | 215.8 | - | - |
| | 1:2 | 97.81 | 107.9 | 90.67 |
| | 1:4 | 44.96 | 53.94 | 83.36 |
| | 1:8 | 22.34 | 26.97 | 82.84 |
| | 1:16 | 11.08 | 13.48 | 82.17 |

Recovery test

To three sera/EDTA plasma samples containing low CA 15-3 concentrations were added known concentrations of CA 15-3. Samples are assayed according to the kit procedure.

| Serum | Endogen conc. (U/mL) | Added conc. (U/mL) | Expected conc. (U/mL) | Measured conc. (U/mL) | Ratio (%) Measured/Expected |
|-------|----------------------|--------------------|-----------------------|-----------------------|-----------------------------|
| S1 | 25.5 | 18.6 | 43.0 | 44.1 | 97.51 |
| | 25.5 | 48.0 | 63.0 | 73.5 | 85.71 |
| | 25.5 | 91.0 | 125 | 116.5 | 107.3 |
| S2 | 27.0 | 18.6 | 44.0 | 45.6 | 96.49 |
| | 27.0 | 48.0 | 74.0 | 75.0 | 98.67 |
| | 27.0 | 91.0 | 123 | 118 | 104.2 |
| S3 | 35.0 | 18.6 | 58.0 | 54.0 | 108.2 |
| | 35.0 | 48.0 | 76.0 | 83.0 | 91.57 |
| | 35.0 | 91.0 | 135 | 126 | 107.1 |

| EDTA Plasma | Endogen conc. (U/mL) | Added conc. (U/mL) | Expected conc. (U/mL) | Measured conc. (U/mL) | Ratio (%) Measured/Expected |
|-------------|----------------------|--------------------|-----------------------|-----------------------|-----------------------------|
| P1 | 9.43 | 4.32 | 13.75 | 14.61 | 106.3 |
| | 9.13 | 8.38 | 17.51 | 18.42 | 105.2 |
| | 8.85 | 12.18 | 21.04 | 21.92 | 104.2 |
| P2 | 13.07 | 9.10 | 22.17 | 22.90 | 103.3 |
| | 12.67 | 17.63 | 30.29 | 31.73 | 104.8 |
| | 12.28 | 25.64 | 37.92 | 40.25 | 106.2 |
| P3 | 20.47 | 9.10 | 29.56 | 29.29 | 99.07 |
| | 19.83 | 17.63 | 37.45 | 40.74 | 108.8 |
| | 19.23 | 25.64 | 44.86 | 43.94 | 97.94 |

¹²⁵I Characteristics

T_{1/2} (¹²⁵I) = 1443 h = 60.14 d

| ¹²⁵ I | E (MeV) | % |
|------------------|---------|-----|
| Y | 0.035 | |
| X | 0.027 | 114 |
| | 0.032 | 25 |

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Vyrobcce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženkla / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Стаклик аралıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kulanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Bīyolōjik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktiviojvi medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikator (tracer) / 트레이서 / Tracer 'lar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibrátor / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μέτρηση / 质控品 / Kontrollinè / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmvavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

Diluent, Diuuant, Diüent, Diluente, Diluyente, Diluente, Spädmedel, Αραιωτικό, 稀釋液, Skiediklis, Diluens, Rozcieńczalnik, Redidlo, Riediaci roztok, 희석액, Dilüent, Дилуент, Дилуент, 稀釋劑