## IRMA CA 19-9 antigen

**REF** IM3151, B64313

### TABLE OF CONTENTS

English . . . . .	2	Čeština . . . . .	19
Français . . . . .	4	Slovenčina . . . . .	21
Deutsch . . . . .	6	한국어 . . . . .	23
Italiano . . . . .	8	Türkçe . . . . .	25
Español . . . . .	10	Русский . . . . .	27
Ελληνικά . . . . .	12	Србија . . . . .	30
Lietuviškai . . . . .	15	APPENDIX . . . . .	32
Polski . . . . .	17		

---

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# IRMA CA 19-9 antigen

**REF** IM3151, B64313

## IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF CA 19-9 ANTIGEN IN HUMAN SERUM

For *in vitro* diagnostic use.

### PRINCIPLE

The CA 19-9 antigen assay is a two-step "sandwich" type assay in which the same mouse monoclonal antibody, directed against the antigenic determinant of the molecule, is employed.

Samples or calibrators are incubated in tubes coated with the CA 19-9 monoclonal antibody, the contents of the tubes are then aspirated and the presence of CA 19-9 in the sample is revealed by incubation with the <sup>125</sup>I-labeled anti CA 19-9 monoclonal antibody. The contents of the tubes are aspirated after this second incubation and unbound labeled antibody is eliminated by washing. The amount of bound reactivity measured in a gamma counter is proportional to the CA 19-9 concentration. The unknown values are determined by interpolation from a standard curve.

### WARNING AND PRECAUTIONS

#### General remarks

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

#### Material of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

PI-IM3151-B64313-03

## GHS HAZARD CLASSIFICATION

Not classified as hazardous

**SDS**

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

## SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes.
- Separate serum from cells by centrifugation.
- Serum samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -20°C, 3 months maximum). Samples should be aliquoted to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Dilute samples showing a concentration exceeding that of the highest calibrator with calibrator zero.

## MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

#### Kit for determination of CA 19-9, 100 tubes (Cat #IM3151)

**Anti-CA 19-9 antibody-coated tubes: 2 x 50 tubes** (ready-to-use)

**<sup>125</sup>I-labeled monoclonal antibody: one 11 mL vial** (ready-to-use)

The vial contains 310 kBq at time of manufacture of <sup>125</sup>I-labeled monoclonal antibody with bovine serum albumine, sodium azide (<0.1%) and a dye.

**Calibrators: five 0.5 mL vials + one zero 5 mL vial** (ready-to-use)

The calibrator vials contain between 0 to approximately 480 U/mL of human CA 19-9 in buffer containing bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. After opening, the vials can be stored at 2-8°C until the expiration date of the kit. Calibrators are verified to an internal reference standard.

**Control serum: one 1 mL vial** (ready-to-use)

The vial contains human CA 19-9 in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement. When opened, the vial may be stored at 2-8°C for 2 weeks. For longer periods, freeze at < -18°C.

**Diluent: one 25 mL vial** (ready-to-use)

The vial contains bovine serum albumin in buffer. After opening, the vial can be stored at 2-8°C until the expiration date of the kit.

#### Kit for determination of CA 19-9, 50 tubes (Cat #B64313)

**Anti-CA 19-9 antibody-coated tubes 1 x 50 tubes** (ready-to-use)

**<sup>125</sup>I-labeled monoclonal antibody: one 11 mL vial** (ready-to-use)

**Calibrators: five 0.5 mL vials + one zero 5 mL vial** (ready-to-use)

**Control serum: one 1 mL vial** (ready-to-use)

**Diluent: one 25 mL vial** (ready-to-use)

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (50 µL, 100 µL, 200 µL).
- repeating micropipets (100 µL, 200 µL, 2 mL).
- "vortex" type mixer.
- **rotative shaker.**
- aspiration system.
- gamma counter.

## RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of CA 19-9 concentrations in samples measured at the same time as the calibrators.

## Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ( $cpm_{cal} - cpm_{cal0}$ ) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (U/mL).

Other data reduction methods may give slightly different results.

Total activity: 116,470 cpm				
Calibrators	CA 19-9 (U/mL)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	1,250	1.07	-
1	14	2,302	1.98	1,052
2	27	3,343	2.87	2,093
3	63	5,379	4.62	4,129
4	122	9,156	7.86	7,906
5	449	25,685	22.1	24,435

(Example of standard curve, do not use for calculation)

## Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

The concentrations of diluted samples must be corrected by the dilution factor.

## EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The following values obtained from healthy subjects are indicative only.

Number of subjects	Average (x)	95 <sup>th</sup> percentil	99 <sup>th</sup> percentil
157	14.8 U/mL	36.1 U/mL	45.8 U/mL

## QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: [imunochem@beckman.com](mailto:imunochem@beckman.com)

## PROCEDURE

### Preparation of reagents

Let the reagents come to room temperature at least 30 minutes before using.

**Bring all reagents to room temperature before pipeting.**

Immunological step 1	Immunological step 2	Counting
To coated tubes add 200 $\mu$ L of diluent, 50 $\mu$ L of sample, calibrator or control and vortex gently. Incubate 2 hours at 18-25°C with shaking (>280 rpm).	Aspirate the content of each tube, then wash twice with 2 mL of distilled water. Carefully aspirate, then add 100 $\mu$ L of tracer in each tube. Incubate 1 hour at 18-25°C with shaking (>280 rpm).	Aspirate the content of each tube, then wash twice with 2 mL of distilled water (except the 2 "total cpm"). Carefully aspirate. Count for 1 minute bound cpm (B) and total cpm (T)

\*Add 100  $\mu$ L of tracer in 2 additional tubes to obtain "total cpm".

## PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

### Sensitivity

**Analytical sensitivity:** 2 U/mL

**Functional sensitivity:** 7.49 U/mL

**Specificity:** The kit is specific of the CA 19-9.

### Precision

#### Intra-assay

Serum samples were assayed 20 times in the same series. The coefficients of variation were found below or equal to 7.9%.

#### Inter-assay

Serum samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 9.6%.

### Accuracy

#### Dilution test

High-concentration serum samples were serially diluted with the kit zero calibrator. The recovery percentages obtained were between 88% and 113%.

#### Recovery test

Low-concentration serum samples were spiked with known quantities of CA 19-9. The recovery percentages obtained were between 85.6% and 113%.

**Measurement range** (from analytical sensitivity to highest calibrator): 2 to approximately 480 U/mL.

## LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use haemolyzed, icteric or lipemic samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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## PERFORMANCE CHARACTERISTICS

### Specificity

Antibodies are specific of Sialyl Lewis<sup>x</sup> and show no known cross-reactivity with related antigens like Sialyl Lewis<sup>x</sup>.

### Precision

#### Intra-assay

Samples	S1	S2	S3
Number of determinations	20	20	20
Mean (U/mL)	30	106	229
CV (%)	6.8	4.8	7.9

#### Inter-assay

Samples	S2	S2	S3
Number of determinations	10	10	10
Mean (U/mL)	36.15	176.6	349.9
CV (%)	9.64	7.22	6.27

### Accuracy

#### Dilution test

Three serum samples were diluted in zero calibrator and assayed according to the assay procedure of the kit.

Sample	Dilution factor	Measured conc. (U/mL)	Expected conc. (U/mL)	Recovery (%) Measured/Expected
S1	Undil.	220	-	-
	1/2	110	110	99.0
	1/4	54.0	55.0	98.0
	1/8	29.7	27.5	108
	1/16	15.6	13.7	113
S2	Undil.	229	-	-
	1/2	103	114	89.0
	1/4	52.0	57.0	90.0
	1/8	25.1	28.6	88.0
	1/16	12.6	14.3	88.0
S3	Undil.	172	-	-
	1/2	76.0	86.0	89.0
	1/4	40.0	43.0	93.0
	1/8	22.6	21.4	105
	1/16	11.0	10.7	102

#### Recovery test

To three sera containing low CA 19-9 concentrations were added known concentrations of CA 19-9. Samples are assayed according to the kit procedure.

Sample	Endogen. conc. (U/mL)	Added conc. (U/mL)	Expected conc. (U/mL)	Measured conc. (U/mL)	Ratio (%) Measured/Expected
S1	10.94	14.42	25.35	28.56	112.6
	10.59	27.94	38.53	39.06	101.4
	10.27	40.64	50.91	50.68	99.55
S2	52.30	14.42	66.72	57.08	85.56
	50.66	27.94	78.60	72.77	92.58
	49.13	40.64	89.76	85.76	95.54
S3	46.02	14.42	60.44	60.12	99.48
	44.58	27.94	72.52	69.99	96.52
	43.23	40.64	83.86	79.59	94.90

### <sup>125</sup>I Characteristics

$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$

Urine	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

## Symbols Key

**REF** Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

**IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

**CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

**Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó / Producent / Výrobce / Уробца / 제조 / Üretici / Изготовлено / Произведено от / 製造商**

**Σ** Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Rækker til "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

**CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

**SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

**Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明**

**Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sicaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍**

**Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / [İspējimas / Figyelem / Uwaga / Urozornění / Urozornenie / 주의 / Dikkat / Внимание / 注意**

**Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日**

**LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

**Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期**

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Ag <sup>125</sup>I

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Ab <sup>125</sup>I

CAL

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CAL 0

CTRL

Control / Contrôle / Kontrolle / Controlla / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 청도관리 / Kontrol / Контроль / Контролна / 质控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintuvéltai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Túpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

DIL

Diluent, Diluant, Diluent, Diluente, Diluyente, Diluente, Spädmittel, Αραιωτικό, 稀释液, Skiediklis, Diluens, Rozcieńczalnik, Redidlo, Riediaci roztok, 희석액, Dilüent, Дилуэнт, Дилуент, 稀釋劑

10 September 2018