# ACTIVE® Dihydrotestosterone RIA

**REF** DSL9600i

## TABLE OF CONTENTS

English . . . . .	2	Polski . . . . .	30
Français . . . . .	5	Čeština . . . . .	34
Deutsch . . . . .	9	Slovenčina . . . . .	37
Italiano . . . . .	13	Türkçe . . . . .	40
Español . . . . .	16	Русский . . . . .	43
Português (Portugal) . . . . .	20	APPENDIX . . . . .	47
Ελληνικά . . . . .	23		
Magyar . . . . .	27		

---

# ACTIVE® Dihydrotestosterone RIA

**REF** DSL9600i

## RADIOIMMUNOASSAY FOR THE IN VITRO QUANTITATIVE MEASUREMENT OF DIHYDROTESTOSTERONE (DHT) IN HUMAN SERUM OR PLASMA For *in vitro* diagnostic use.

### PRINCIPLE

The radioimmunoassay of dihydrotestosterone (5 $\alpha$ -Dihydrotestosterone; DHT; 17 $\beta$ -Hydroxy-5 $\alpha$ -androstane-3-one) is a competitive assay. The procedure follows the basic principle of radioimmunoassay where there is a competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of [125I]-labeled dihydrotestosterone bound to the antibody is inversely proportional to the concentration of unlabeled DHT present. The separation of the free and bound antigen is achieved by decanting or aspirating the antibody-coated tubes. A standard curve is constructed and unknown DHT values are obtained from the curve by interpolation.

For Summary and Explanation of the Test see APPENDIX.

### WARNING AND PRECAUTIONS

#### General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

#### Materials of human origin

Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination.

PI-DSL9600I-02

Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

### GHS HAZARD CLASSIFICATION

Oxidation Solution

H411 Toxic to aquatic life with long lasting effects.  
P273 Avoid release to the environment.  
P391 Collect spillage.  
Potassium Permanganate 1 - 5%

Sample Buffer

DANGER

H226 Flammable liquid and vapour.  
H314 Causes severe skin burns and eye damage.  
P210 Keep away from heat, hot surfaces, and sparks. No smoking.  
P280 Wear protective gloves, protective clothing and eye/face protection.  
P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.  
P303+P361+P353 IF ON SKIN (or hair): Rinse skin with water.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor/physician.  
Ethyl Alcohol 10 - 20%  
Sodium Hydroxide 20 - 30%

**SDS**

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

### SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum and plasma-EDTA are the recommended sample types.
- Allow serum samples to clot completely before centrifugation.
- Store separated serum or plasma at 2-8 °C for up to 24 hours. For longer storage keep frozen at <-20 °C, up to 1 year. It is recommended to prepare aliquots to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Samples MUST be extracted prior to assay. See Procedure.
- Any sample reading greater than the highest calibrator should be diluted appropriately with the 0 pg/mL calibrator and reassayed.
- Serum and EDTA plasma values for 19 samples (serum values ranging from 53.75 to 341.5 pg/mL) were compared using the Active Dihydrotestosterone RIA kit DSL9600(i). Results are as follows:

[EDTA-plasma] = 0.9592[serum] + 31.583 R = 0.9783

### MATERIALS PROVIDED

All reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8 °C. Expiry dates printed on vial labels apply to the

long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for opened reagents are indicated in appropriate paragraphs.

**Anti-dihydrotestosterone antibody coated tubes: 2 x 50 tubes** (ready-to-use)

Plastic tubes with polyclonal rabbit anti-DHT immunoglobulin immobilized to the inside wall of each tube.

**125I-labeled dihydrotestosterone tracer: one 55 mL vial** (lyophilized)

On the date of manufacture, the vial contains 185 kBq (<5 µCi) of 125I-labeled DHT in buffer with proteins (BSA) and sodium azide (<0.1 %). Reconstitute tracer with 55 mL of deionized water. After reconstitution, store at 2-8 °C for up to 14 days in order to achieve maximum binding (%B<sub>0</sub>/T) > 25%.

**Calibrators: one 50 mL vial labeled 0** (ready to use) **and seven vials labeled 1-7** (lyophilized).

The calibrator vials contain from 0 to approximately 2,500 pg/mL (0 to approximately 8,600 pmol/L) of dihydrotestosterone in buffer with proteins (BSA) and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrator values were established using a certified reference material (Cerilliant). Do NOT extract calibrators prior to assay.

**Controls: two vials labeled 1, 2** (lyophilized)

The vials contain dihydrotestosterone in human serum and sodium azide (<0.1%). The expected values are in a supplement found in the kit. Controls MUST be extracted prior to assay. (See Procedure).

**Oxidation Solution: two 25 mL bottles** (ready-to-use)

Bottles contain potassium permanganate solution (<5 %) in a buffer with sodium azide (<0.1 %).

**Sample buffer: one 5.5 mL vial** (ready-to-use)

The vial contains buffer with <23 % Ethyl Alcohol and <29 % Sodium Hydroxide.

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm or 13 x 100 mm or 16 x 150 mm **GLASS** tubes and safety caps.
- Test tube rack for 12 x 75 mm tubes.
- Deionized water.
- Precision micropipet (100 µL, 250 µL and 400 µL).
- Semi-automatic pipet (500 µL).
- Vortex type mixer.
- 5 mL graduated glass pipettes - **USED IN EXTRACTION.**
- Organic solvents: n-hexane and ethanol (HPLC grade).
- Centrifuge (1500 x g, preferably refrigerated).
- Nitrogen gas or a speed-vac assembly with heating (for extractions).
- Shaker capable of ≥ 180 rpm.
- A sponge rack or similar device for decantation.
- Aspiration system.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.
- Semi-log (log-linear) graph paper or computer RIA data analysis program.

## RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of dihydrotestosterone concentrations in samples measured at the same time as the calibrators.

### Standard curve

The results below were calculated using a log-linear curve fit ("spline" mode) with B/T (%) or B/B<sub>0</sub> (%) on vertical axis and the DHT concentration of the calibrators on the horizontal axis (pg/mL). Other data reduction methods may give slightly different results.

ED50 = 100.7 pg/mL

Total activity: 57,936 cpm				
Calibrators	DHT (pg/mL)	cpm (n=2)	B/T (%)	B/B <sub>0</sub> (%)
0	0	37.670	65.0	100
1	24.0	31.164	53.8	82.7
2	48.0	26.403	45.6	70.1
3	96.0	19.356	33.4	51.4
4	192	12.449	21.5	33.0
5	480	6.948	12.0	18.4
6	960	4.570	7.88	12.1
7	2.400	2.514	4.33	6.67

(Example of standard curve, do not use for calculation)

### Samples

For each sample locate the ratio B/T (%) or B/B<sub>0</sub> (%) on the vertical axis of the standard curve and read off the corresponding DHT concentration of the sample on the horizontal axis in pg/mL. To convert concentrations from pg/mL to pmol/L, multiply results by 3.44.

## EXPECTED VALUES

Each laboratory should establish its own reference ranges. An internal study of healthy adults gave the DHT concentrations shown in the table below and are provided for reference only.

	Age	n	Median (pg/mL)	2.5% Percentile (pg/mL)	97.5% Percentile (pg/mL)	Range (pg/mL)
Male	20 - 60	119	274.7	109.2	583.1	33.69 - 756.4
Female	20 - 60	122	82.81	33.28	196.5	17.67 - 245.8

(For more details, see APPENDIX)

## QUALITY CONTROL

Good laboratory practices require that control samples be used regularly to ensure the quality of the results obtained. These controls must be processed in exactly the same way as the patient samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling, or deterioration of reagents. In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: [imunochem@beckman.com](mailto:imunochem@beckman.com)

## PROCEDURE

### Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

### Reconstitution of tracer

Reconstitute tracer with 55 mL of deionized water. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store at 2-8 °C for up to 14 days.

### Reconstitution of calibrators and control samples

Reconstitute calibrators 1-7 and controls 1-2 with the volume of deionized water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store at <-20 °C until expiration date of kit.

### Extraction of samples

*Note: Extractions must be done using clean, preferably disposable, GLASS tubes and pipettes. Allow the Oxidation Solution and other extraction reagents to reach room temperature (18-25 °C) before use.*

*Samples and controls require extraction. Do NOT extract calibrators.*

- Number one tube for each sample or control.
- Add 400 µL of sample or control to a numbered glass tube and add 500 µL of Oxidation solution. Thoroughly vortex and incubate at room temperature (18-25 °C) for 15 minutes.
- Prepare the extraction mixture: 98% n-hexane and 2% ethanol.
- Using a glass pipette, extract the oxidized sample by adding 4.0 mL n-hexane-ethanol mixture (98% hexane: 2% ethanol). Vortex each sample immediately for 1 minute.

- Add 50 µL of DHT Sample Buffer, cap tubes and mix gently by inverting tubes 3-4 times.
- Centrifuge at 1500 X g for 15 minutes at 2-8 °C to separate the organic layer from the aqueous layer.
- Transfer 2.5 mL of the upper organic layer into appropriately labeled clean glass tubes and evaporate to dryness, using either nitrogen gas or a speed-vac assembly with heating.

Reconstitute the dried material with 250 µL of the DHT zero calibrator. Thoroughly vortex and keep at room temperature (18-25 °C) at least 1 hour.

**NOTE:**

- DO NOT EXTRACT KIT CALIBRATORS
- Use 12 x 75 mm or 13 x 100 mm glass tubes with safety caps. If safety caps are not available, 16 x 150 mm tubes can be used. After extraction, cover the 16 x 150 mm tubes with aluminum foil to avoid evaporation of organic solvent.
- Avoid using lipemic samples.

**Assay procedure**

Run Calibrators, Controls and patient samples in duplicate.

Step 1 Additions	Step 2 Incubation	Step 3 Counting
Extract samples and controls as directed in section Procedure.  To antibody coated tubes successively add:  100 µL of calibrator, control or sample and immediately add  500 µL of tracer.*  Vortex all tubes.	Cover and incubate 2 hours on a shaker (≥180 rpm) at room temperature (18-25°C).  Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.	Add 3.0 mL deionized water to all tubes, except «total cpm» tubes.  Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.  Strike the tubes sharply on absorbent material to facilitate complete drainage. Drain for >2 minutes. Blot the tubes.  Count bound cpm (B) and total cpm (T) for 1 minute.

\*Add 500 µL of tracer to 2 additional tubes to obtain total cpm.

**PERFORMANCE CHARACTERISTICS**

(For more details, see APPENDIX)

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

**Sensitivity**

**Analytical sensitivity:** 9.14 pg/mL

**Functional sensitivity:** 19.63 pg/mL

**Specificity**

The antibody used in the immunoassay is highly specific for DHT. Extremely low cross reactivities were obtained with several related molecules (androstenedione, estradiol, testosterone etc).

**Precision**

**Intra-assay**

Samples were assayed 25 times in the same run. The coefficients of variation were ≤ 7.9 %.

**Inter-assay**

Samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 7.1 %.

**Accuracy**

**Dilution test**

Five serum samples were serially diluted with zero calibrator. The recovery percentages ranged from 80.5 % to 114 %.

**Recovery test**

Five serum samples were spiked with known quantities of DHT. The recovery percentages ranged from 88.0 % to 107 %.

**Measurement range** (from analytical sensitivity to highest calibrator): 9.14 to approximately 2,500 pg/mL.

**LIMITATIONS**

- Failure to follow these instructions for use (IFU) may significantly affect results.
- Failure to blot the tubes adequately following decantation may result in poor replication and spurious results.
- Avoid repeated freezing and thawing of reagents or specimens.
- Do not use hemolytic, lipemic or icteric samples.
- If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial.

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#### Summary and explanation of the test

5 $\alpha$ -Dihydrotestosterone (DHT; 17 $\beta$ -Hydroxy-5 $\alpha$ -androstane-3-one), the most potent naturally-occurring androgen, is produced from testosterone through the action of cholesterone 5 $\alpha$ -reductase. The concentrations of 5 $\alpha$ -reductase are highest in certain peripheral tissues, including genital skin and hair follicles, and is localized intracellularly in apparent association with the nuclear membrane. DHT exerts its biological action by intracellular binding to the androgen receptor; this complex is then transferred to the nucleus where DNA-binding occurs with resultant effects on DNA transcription. Most of the residual DHT undergoes intracellular metabolism to 3 $\alpha$ -androstenediol and 3 $\alpha$ -androstenediol glucuronide. Only a small proportion of DHT escapes into the peripheral circulation, where it is present primarily complexed to sex-hormone binding globulin. Testosterone causes virilization of the Wolffian ducts during fetal life, while DHT is responsible for the development of the male external genitalia and prostate, and is primarily responsible for the physical changes which occur during male sexual maturation. An autosomal-recessive genetic deficiency of 5 $\alpha$ -reductase, sometimes called male pseudohermaphroditism or pseudovaginal perineoscrotal hypospadias, leads to inadequate differentiation of DHT-dependent peripheral tissues. Male infants with this disorder have ambiguous genitalia and are often raised as females, although significant virilization may occur later in life presumably due to the natural increase in testosterone levels. Measurement of DHT concentrations can be complicated by antibody cross-reactivity to testosterone. This DHT Radioimmunoassay utilizes a sample oxidation/extraction procedure to remove most of the testosterone, coupled with a relatively specific immunoassay for DHT.

#### Sensitivity

The analytical sensitivity, or minimum detection limit, calculated by the interpolation of the mean minus two standard deviations of 20 replicates of the 0 pg/mL DHT calibrator, is 9.14 pg/mL.

#### Specificity

The cross-reactivity of the DHT antiserum has been measured against the following compounds. The percent cross-reactivity is expressed as the ratio of the DHT concentration to the concentration of the reacting compound at 50% binding of the 0 pg/mL DHT calibrator.

COMPOUND	% CROSS-REACTIVITY
Dihydrotestosterone	100
Androstenedione	1.90
Estradiol	1.41
Testosterone	0.02*
Androstenediol	0.25
Androstenediol Glucuronide	0.19
Androsterone Glucuronide	ND
Dehydroepiandrosterone	ND
Cortisol	ND
11-Deoxycortisol	ND
17 $\alpha$ -OH-Progesterone	ND
Progesterone	ND

\*after extraction, ND - Not detectable (<0.01 %)

#### Precision

##### Intra-assay

The intra-assay precision (after extraction) was determined from the mean of 25 determinations.

Samples	S1	S2	S3
Number of determinations	25	25	25
Mean value, pg/mL	71.47	332.4	843.4
C.V., %	7.89	6.52	4.45

##### Inter-assay

The inter-assay precision (after extraction) was determined from the mean of average duplicates for 10 separate runs.

Samples	S1	S2	S3
Number of determinations	10	10	10
Mean value, pg/mL	54.65	315.7	813.4
C.V., %	5.17	7.08	2.47

#### Accuracy

##### Dilution test

Five human serum samples were diluted with 0 pg/mL DHT calibrator after extraction and assayed.

Sample	Dilution factor	Measured conc. (pg/mL)	Expected conc. (pg/mL)	Ratio (%) Measured/Expected
S1	-	301.2	-	-
	1:2	148.9	150.6	98.87
	1:4	75.46	75.30	100.2
	1:8	42.93	37.65	114.0
S2	-	620.9	-	-
	1:2	268.4	310.5	86.45
	1:4	132.1	155.2	85.12
	1:8	64.20	77.61	82.72
S3	-	409.9	-	-
	1:2	173.2	205.0	84.50
	1:4	86.97	102.5	84.87
	1:8	53.73	51.24	104.9
S4	-	481.9	-	-
	1:2	194.0	240.9	80.50
	1:4	107.8	120.5	89.48
	1:8	49.21	60.23	81.70
S5	-	1,126	-	-
	1:2	561.6	563.0	99.75
	1:4	239.1	281.5	84.94
	1:8	127.0	140.8	90.24
	1:16	64.92	70.38	92.25
	1:32	33.44	35.19	95.03

#### Recovery test

Five serum samples containing different levels of endogenous DHT were spiked with known amounts of DHT after extraction and assayed.

Sample	Endogen. conc. (pg/mL)	Added conc. (pg/mL)	Expected conc. (pg/mL)	Measured conc. (pg/mL)	Ratio (%) Measured/Expected
S1	270.7	119.0	389.8	370.3	95.01
	279.6	409.8	689.4	636.2	92.28
	275.1	806.5	1,082	1,049	97.03
S2	284.8	119.0	403.9	388.8	96.28
	294.2	409.8	704.0	684.1	97.18
	289.4	806.5	1,096	1,153	105.2
S3	371.7	156.3	528.0	501.6	95.00
	390.0	409.8	799.9	854.3	106.8
	383.7	806.5	1,190	1,201	100.9
S4	154.7	80.65	235.3	230.3	97.85
	147.5	192.3	339.8	299.1	88.02
	157.2	409.8	567.0	545.8	96.25
S5	316.8	156.3	473.1	439.1	92.81
	332.4	409.8	742.2	756.2	101.9
	327.0	806.5	1,134	1,158	102.1

#### Expected values

Male Age	n	Median (pg/mL)	2.5 % Percentile (pg/mL)	97.5 % Percentile (pg/mL)	Range (pg/mL)
20 - 30	26	343.0	168.6	691.7	150.7 - 756.0
31 - 40	32	283.2	162.8	596.5	155.0 - 655.0
41 - 50	29	255.6	76.86	440.3	70.62 - 454.4

Male					
Age	n	Median	2.5 % Percentile	97.5 % Percentile	Range
		(pg/mL)	(pg/mL)	(pg/mL)	(pg/mL)
51 - 60	32	228.1	93.45	412.9	33.69 - 515.2
20 - 60	119	274.7	109.2	583.1	33.69 - 756.4

Female					
Age	n	Median	2.5 % Percentile	97.5 % Percentile	Range
		(pg/mL)	(pg/mL)	(pg/mL)	(pg/mL)
20 - 30	30	103.4	46.39	217.0	33.16 - 245.8
31 - 40	31	95.80	46.02	201.1	45.24 - 244.5
41 - 50	32	76.40	37.27	122.3	21.24 - 127.3
51 - 60	29	54.70	18.52	78.55	17.67 - 84.36
20 - 60	122	82.81	33.28	196.5	17.67 - 245.8

### <sup>125</sup>I Characteristics

$$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$$

<sup>125</sup> I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

### Symbols Key

**REF** Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

**IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

**CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Уробца / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Rækker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 内容物足夠執行 <n> 次测试

**CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

**SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sicaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Urozornění / Urozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kulanma Tarihi / Срок годности / Срок на годност / 到期日

**LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Ag <sup>125</sup>I

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Ab <sup>125</sup>I

CAL

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CAL 0

CTRL

Control / Contrôle / Kontrolle / Controlla / Control / Controllo / Kontrolle / Μέτρησης / 质控品 / Kontrollinè / Kontroll / Kontrola / Kontrola / Kontrola / 청도관리 / Kontrol / Контроль / Контролна / 质控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνώρα / 试管 / Mégintuvéliei / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanna Talimati / Инструкции / Инструкции за употреба / 使用说明

BUF

Buffer, Tampon, Puffer, Tampone, Tampón, Tampão, Buffert, Ρυθμιστικό διάλυμα, 缓冲液, Buferinis tirpalas, Puffer, Bufor, Pufir, Tlmivý roztok, 완충, Tampon, Буфер, буфер, 緩衝劑

SOLN OX

Oxidation Solution, Solution oxydante, Oxidationslösung, Soluzione di ossidazione, Solución de oxidación, Solução de oxidação, Oxidationslösning, Διάλυμα οξειδωσης, 氧化溶液, Oksidacijos tirpalas, Oxidáló oldat, Roztwór do utleniania, Oxidační činidlo, Oxidačný roztok, 산화 용액, Oksidasyon Çözeltisi, Окисляющий раствор, Разтвор за оксидация, 氧化溶液

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZEŃSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / ТЕХЛІКЕ / ОПАСНО / ОПАСНОСТ / 危險

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda,  
Calçada Aldebarã, 39, Centro de Apoio 2 - Alphaville,  
CEP 06541-055 - Santana de Parnaíba, SP, Brasil  
Telefone: (11) 4154-8818 - CNPJ: 42.160.812/0001-44

IMMUNOTECH s.r.o., Radiova 1, 102 27 Prague 10, Czech Republic