ESTRONE SULFATE RIA

REF DSL5400

TABLE OF CONTENTS

English	2	Čeština (CZ)	20
Français (FR)	5	Русский (RU)	23
Deutsch (DE)	8	Србија (SR)	26
Italiano (IT)	11	APPENDIX	29
Español (ES)	14	REFERENCES	32
Magyar (HU)	17		

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ESTRONE SULFATE RIA

REF DSL5400

RADIOIMMUNOASSAY FOR THE QUANTITATIVE MEASUREMENT OF ESTRONE SULFATE IN HUMAN SERUM OR PLASMA THIS ASSAY IS INTENDED FOR *IN VITRO* DIAGNOSTIC USE

PRINCIPLE

The radioimmunoassay of estrone sulfate [1,3,5 (10)-estratrien-3-ol-17-one-3-sulfate] is a competitive assay [1]. The procedure follows the basic principle of radioimmunoassay where there is competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of [¹²⁵I]-labeled estrone sulfate bound to the antibody is inversely proportional to the concentration of unlabeled estrone sulfate present. The separation of the free and bound antigen is achieved by using a double antibody system. After centrifugation a standard curve is constructed and unknown estrone sulfate values are obtained from the curve by interpolation.

For Summary and Explanation of the Test see APPENDIX.

WARNING AND PRECAUTIONS

General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should take place in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided and in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up [2].

Materials of human origin

Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination.

PI-DSL5400-04

Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

The Materials Safety Data Sheet (MSDS) is available upon request.

GHS HAZARD CLASSIFICATION

Antiserum	DANGER	
		
	H360	May damage fertility or the unborn child.
	P201	Obtain special instructions before use.
	P280	Wear protective gloves, protective clothing and eye/face protection.
	P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - < 0.3% Sodium Borate Decahydrate 0.1 - < 0.3%
Tracer	WARNING	
	H316	Causes mild skin irritation.
	P332+P313	If skin irritation occurs: Get medical advice/attention. Acetic Acid 0.1 - < 1%
Precipitating Reagent	WARNING	
		
	H317	May cause an allergic skin reaction.
	H412	Harmful to aquatic life with long lasting effects.
	P273	Avoid release to the environment.
	P280	Wear protective gloves, protective clothing and eye/face protection.
	P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364	Take off contaminated clothing and wash it before use.
		reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) < 0.05%

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum or EDTA plasma are the recommended sample types.
- Allow serum samples to clot completely before centrifugation.
- Serum or plasma samples may be stored in glass tubes at 2-8°C, if the assay is to be performed within two days. For longer storage keep frozen at <-20°C for up to 2 months. It is recommended to prepare aliquots to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.
- Any sample reading greater than the highest calibrator should be diluted appropriately with the zero calibrator and reassayed.

Serum and EDTA plasma values for 15 samples (serum values ranging from 0.81 to 4.35 ng/mL) were compared using the DSL5400 Estrone Sulfate RIA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.9744[\text{serum}] + 0.0775$$

$$R = 0.9898$$

MATERIALS PROVIDED

All reagents in the kit are stable until the expiry date indicated on the kit label, when stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for opened reagents are indicated in appropriate paragraphs.

Estrone-S Antiserum (BLUE): one 11 mL vial (ready-to-use)

The vial contains rabbit anti-estrone sulfate (polyclonal) serum in a buffer with proteins (BSA), sodium azide (<0.1%) and a dye.

¹²⁵I-labeled Estrone-Sulfate tracer (YELLOW): one 55 mL bottle (ready-to-use)

At the time of manufacture, the vial contains 185 kBq (<5μCi) of [¹²⁵I]-labeled estrone sulfate in buffer with proteins (BSA), sodium azide (<0.1%) and a dye.

Calibrators: one 2.0 mL vial, labeled 0 and six 1.0 mL vials, labeled 1-6 (ready-to-use)

The calibrator vials contain from 0 to approximately 80 ng/mL of estrone sulfate in buffer with proteins (BSA) and sodium azide (<0.1%). The exact concentration is indicated on each vial label. Once opened, store at 2-8°C for up to 3 weeks, or at < -20°C until expiration date of kit. The calibrators values were established using an internal standard.

Controls: two 1.0 mL vials, labeled 1,2 (ready-to-use)

The vials contain estrone sulfate in buffer with proteins (BSA) and sodium azide (<0.1%). The expected values are indicated in a supplement found in the kit. Once opened, store at 2-8°C for up to 3 weeks, or at < -20°C until expiration date of kit.

Precipitating Reagent: one 100 mL bottle (ready-to-use)

The bottle contains goat anti-rabbit gamma globulin serum in a buffer with polyethylene glycol as a precipitating aid and sodium azide (<0.1%).

NOTE: A precipitate may be visible in the reagent. Mix the bottle thoroughly prior to use in the assay.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm plastic polypropylene test tubes.
- Test tube rack for 12 x 75 mm tubes.
- Deionized water.
- Precision micropipets (100 μL).
- Semi-automatic pipets (100 μL, 500 μL, 1.0 mL).
- Vortex type mixer.
- Centrifuge (1500 x g, preferably refrigerated).
- Shaker capable of 180 rpm.
- A sponge rack for decantation or similar device.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.
- Semi-log (log-linear) graph paper or computer RIA data analysis program.

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

Assay procedure

Run calibrators, controls and patients samples in duplicate.

Step 1 Additions	Step 2 Incubation & Centrifugation	Step 3 Counting
To the bottom of labeled tubes, successively add: 100 μL of calibrator, control or sample,*	Incubate 3 hours at room temperature (18 - 25°C) on a shaker (180 rpm). Add 1.0 mL of well mixed Precipitating Reagent to all tubes (except the «total cpm» tubes) and immediately vortex all tubes.	Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle. Drain on absorbent material for 1-2 min and gently blot the tubes.
500 μL of tracer, and**	Incubate all tubes for 10 - 15 min at room temperature (18 - 25°C). Centrifuge*** all tubes (except the «total cpm» tubes) for 15 min at 1500 x g.	Count bound cpm (B) and total cpm (T) for 1 minute.
100 μL of antiserum (except the 2 «total cpm» tubes and 2 tubes for NSB). Vortex gently for 1-2 seconds.		

*To two NSB tubes add 200 μL of the 0 ng/mL estrone sulfate calibrator.

**Add 500 μL of tracer to 2 additional tubes to obtain «total cpm».

*** preferably refrigerated; approx. 3000 rpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of estrone sulfate concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

$$ED_{50} = 1.39 \text{ ng/mL}$$

Total activity: 70,576 cpm				
Calibrators	Estrone sulfate (ng/mL)	cpm (n=2)	B/T (%)	B/B ₀ (%)
0	0.0	37,037	49.4	100
1	0.05	34,360	45.7	92.3
2	0.20	30,403	40.0	81.0
3	1.00	21,522	27.5	55.5
4	5.00	12,099	14.1	28.5
5	25.0	6,050	5.5	11.2
6	80.0	3,926	2.5	5.1

(Example of standard curve, do not use for calculation)

Samples

For each sample locate the ratio B/T or B/B_0 on the vertical axis of the standard curve and read off the corresponding estrone sulfate concentration of the sample on the horizontal axis in ng/mL.

EXPECTED VALUES

Each laboratory should establish its own reference ranges. The following values are indicative only.

	n	Mean (SD) (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
Male	20	1.73 (0.62)	1.66	0.97-3.93
Female				
Follicular phase	31	2.29 (0.88)	2.02	0.79-4.57
Luteal phase	30	3.88 (2.20)	3.10	0.83-8.81
Pregnancy				
1 st trimester	30	11.72 (6.15)	10.47	3.09-30.32
2 nd trimester	30	17.20 (9.40)	17.02	2.76-36.41
3 rd trimester	39	31.54 (23.29)	27.45	5.04-100.2
Oral Contraceptive	27	1.42 (1.26)	1.02	0.58-6.36
Postmenopausal (no HRT)	28	1.06 (0.40)	0.97	0.50-1.91
Postmenopausal (+HRT)	15	4.39 (4.84)	1.62	0.53-15.77

QUALITY CONTROL

Good laboratory practices require that control samples be used regularly to ensure the quality of the results obtained. These controls must be processed in exactly the same way as the patient samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Non-specific binding (NSB) greater than 7% may indicate reagent deterioration, contamination or technical error.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling or deterioration of reagents.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

In the US, contact the Beckman Coulter technical support at 1-800-854-3633; or by email at: immunoassay@beckman.com

PERFORMANCE CHARACTERISTICS

(For more details, see APPENDIX)

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

The analytical sensitivity is 0.01 ng/mL

Specificity

The antibody used in the immunoassay is highly specific for estrone sulfate. Low cross reactivities were obtained with following compounds (estrone and estrone glucuronide, 17 β -estradiol 3-sulfate, etc.).

Precision

Intra-assay

Human serum samples were assayed in 14 times in the same run. The coefficients of variation were \leq 9.2%.

Inter-assay

Human serum samples were assayed in duplicate in 8 different runs. The coefficients of variation were \leq 8.8%.

Accuracy

Dilution test

Five serum samples were serially diluted with zero calibrator. The recovery percentages ranged from 80.2% to 119%.

Recovery test

Low-concentration serum samples were spiked with known quantities of estrone sulfate. The recovery percentages ranged from 83% to 111%.

Measurement range (from analytical sensitivity to highest calibrator): 0.01 to approximately 80 ng/mL.

Method Comparison

n	Slope	Intercept	r
41	36.66	2.52	0.86

LIMITATIONS

- Failure to follow these instructions for use (IFU) may significantly affect results.
- Failure to blot tubes adequately following decantation may result in poor replication and spurious values.
- Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.
- The reagents supplied in this kit are optimized to measure estrone sulfate levels in serum or EDTA plasma.
- Repeated freezing and thawing of reagents supplied in the kit and of specimens must be avoided.
- Do not use hemolyzed, icteric or lipemic specimens.
- The possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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APPENDIX

PERFORMANCE CHARACTERISTICS

Summary and explanation of the test

Estrone sulfate [1,3,5(10)-estratrien-3-ol-17-one-3-sulfate] is the most abundant circulating estrogen in non-pregnant women as well as normal men [3,4]. It is found in peripheral circulation, due to sulfokinases, present in the peripheral tissues as well as gonad and adrenal cortex, as a major metabolite of estradiol and estrone, and also known to have a longer half-life in blood [4,5]. Circulating levels of estrone sulfate are high during fetal life in both sexes, decrease to very low levels within the first few days of life, remain relatively low during childhood and increase steadily during puberty [1]. Postnatal circulating levels in females are higher than males. In premenopausal women, estrone sulfate levels generally parallel those of estrone, rising gradually during follicular phase, peaking just prior to ovulation, with a secondary and smaller increase during the luteal phase [1,6]. In post-menopausal women, estrone sulfate could play a significant role, since the ovary secretes very little estrogen, and can be used as an important marker in monitoring post-menopausal women. Measurement of estrone sulfate may be complicated by assay interference due to the presence of other estrogens and related compounds in the sample.

Sensitivity

The analytical sensitivity, or minimum detection limit, calculated by the interpolation of the mean minus two standard deviations of 16 replicates of the 0 ng/mL estrone sulfate calibrator, is 0.01 ng/mL.

Specificity

The cross-reactivity of the Estrone Sulfate antiserum has been measured against the following compounds:

Compound	% Cross-reactivity
Estrone sulfate	100
Estrone	4.9
Estrone glucuronide	3.4
17 β -Estradiol 3-sulfate	1.0
17 β -Estradiol	0.3
17 β -Estradiol glucuronide	0.1
5-Androsten-3 β -ol-17-one (DHA)	ND
Δ -Androsten-3,17-dione (androstenedione)	ND
5 α -Androstane-3 β ,17 β -diol 3-glucuronide	ND
5 α -Androstan-3 α -ol-17-one sodium sulfate	ND
5 α -Androstan-3 α -ol-17-one glucuronide (androsterone glucuronide)	ND
5 α -Androstan-3 β -ol-17-one glucuronide	ND
5 β -Androstan-3 α -ol-17-one sodium sulfate	ND
Hydrocortisone	ND
Progesterone	ND
4-Pregnen-21-ol-3,20-dione	ND
5-Androsten-3 β -ol-17-one sulfate (dehydroisoandrosterone-3-sulfate)	ND
4-Estren-17 β -ol-3-one	ND
1,4-Pregnadien-11 β ,17 α ,21-triol-3,20-dione (prednisolone)	ND
5 α -Androstan-3 α -ol-17-one (androsterone)	ND
17 α -Hydroxy- progesterone	ND
Corticosterone	ND
5 β -Androstan-3 β -ol-17-one	ND
Norethindrone	ND
Medroxyprogesterone acetate	ND

Precision

Intra-assay

The intra-assay precision was determined from the mean of 14/25 replicates each with three human serum/EDTA plasma samples.

Serum	N	Mean (ng/mL)	C.V. (%)
S1	14	0.35	9.2
S2	14	8.98	4.6
S3	14	59.33	4.7

EDTA plasma	N	Mean (ng/mL)	C.V. (%)
P1	25	0.58	11.72
P2	25	3.16	5.47
P3	25	13.3	3.53

Inter-assay

The inter-assay precision was determined from the mean of average duplicates for 8/10 separate runs with human serum or EDTA plasma samples.

Serum	N	Mean (ng/mL)	C.V. (%)
S1	8	0.08	8.8
S2	8	0.49	5.1
S3	8	11.30	5.5

EDTA plasma	N	Mean (ng/mL)	C.V. (%)
P1	10	1.02	16.64
P2	10	2.54	8.17
P3	10	7.04	7.83

Accuracy

Dilution test

Five human serum/EDTA plasma samples were diluted with 0 ng/mL estrone sulfate calibrator and assayed.

Serum	Dilution factor	Measured (ng/mL)	Expected (ng/mL)	Ratio (%) Measured/Expected
S1	-	13.17	-	-
	1:2	5.28	6.59	80.18
	1:4	2.82	3.29	85.65
	1:8	1.46	1.65	88.69
S2	1:16	0.89	0.82	108.1
	-	14.89	-	-
	1:2	6.71	7.45	90.13
	1:4	3.61	3.72	96.98
S3	1:8	1.60	1.86	85.96
	1:16	0.91	0.93	97.78
	-	6.56	-	-
	1:2	3.54	3.28	107.9
S4	1:4	1.53	1.64	93.29
	1:8	0.82	0.82	100.0
	1:16	0.34	0.41	82.93
	-	9.00	-	-
S5	1:2	4.50	4.50	100.0
	1:4	2.65	2.25	117.8
	1:8	1.34	1.13	119.1
	1:16	0.50	0.56	88.89
S5	-	17.42	-	-
	1:2	8.60	8.71	98.74
	1:4	4.75	4.36	109.1
	1:8	2.47	2.18	113.4
1:16	1.27	1.09	116.6	

EDTA plasma	Dilution factor	Measured (ng/mL)	Expected (ng/mL)	Ratio (%) Measured/Expected
P1	-	10.97	-	-
	1:2	4.65	5.49	84.78
	1:4	2.86	2.74	104.3
	1:8	1.38	1.37	100.6
	1:16	0.67	0.69	97.72
P2	-	8.51	-	-
	1:2	3.47	4.26	81.55
	1:4	2.26	2.13	106.2
	1:8	1.21	1.06	113.7
	1:16	0.52	0.53	97.77
P3	-	9.16	-	-
	1:2	3.81	4.58	83.19
	1:4	2.20	2.29	96.07
	1:8	1.12	1.15	97.82
	1:16	0.58	0.57	101.3
P4	-	15.66	-	-
	1:2	6.63	7.83	84.67
	1:4	3.78	3.92	96.55
	1:8	2.16	1.96	110.3
	1:16	0.94	0.98	96.04
P5	-	15.65	-	-
	1:2	7.15	7.83	91.37
	1:4	3.36	3.91	85.88
	1:8	1.83	1.96	93.55
	1:16	0.87	0.98	88.95

EDTA plasma	Endogen. (ng/mL)	Added (ng/mL)	Expected (ng/mL)	Measured (ng/mL)	Ratio (%) Measured/Expected
R1	2.13	0.41	2.54	2.75	108.1
	2.10	0.81	2.91	3.14	108.0
	2.07	1.19	3.26	3.75	115.1
R2	1.75	0.41	2.16	2.30	106.4
	1.72	0.81	2.53	2.91	115.1
	1.70	1.19	2.89	3.44	119.2
R3	0.74	0.41	1.15	1.14	99.34
	0.73	0.81	1.53	1.76	114.9
	0.71	1.19	1.90	2.27	119.2
R4	1.62	0.41	2.03	2.06	101.3
	1.60	0.81	2.40	2.52	104.9
	1.57	1.19	2.76	3.00	108.6
R5	2.00	0.41	2.41	2.67	110.9
	1.96	0.81	2.77	3.15	113.7
	1.93	1.19	3.12	3.58	114.6

Method Comparison

The DSL5400 Estrone sulfate RIA was compared to a commercially available RIA kit [Method X]. Forty-one serum samples ranging in value from 0.5 to 77.57 ng/mL (estrone sulfate) were assayed and linear regression analysis of the results yielded the following:

n = 41

Regression: [DSL5400] = 36.66 [Method X] + 2.52

r = 0.86

¹²⁵I Characteristics

T_{1/2} (¹²⁵I) = 1443 h = 60.14 d

¹²⁵ I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Recovery test

Three human serum and five EDTA-plasma samples containing different levels of endogenous estrone sulfate were spiked with different amounts of estrone sulfate and assayed.

Serum	Endogen. (ng/mL)	Added (ng/mL)	Expected (ng/mL)	Measured (ng/mL)	Ratio (%) Measured/Expected
S1	0.30	13.05	13.35	13.07	98
		36.63	32.93	27.35	83
		65.25	65.55	66.99	102
S2	0.45	13.05	13.50	11.26	83
		36.63	33.07	28.43	86
		65.25	65.70	63.71	97
S3	1.95	13.05	15.0	16.58	111
		36.63	34.57	36.39	105
		65.25	67.20	70.45	105

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räkcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženkla / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postępujcie podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo sarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvijo medžiaga / Radioaktiv / Radioaktywny / Radioaktivni / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Ag ¹²⁵I

Ab ¹²⁵I

Tracer / Traceur / Marcato / Trazador / Marcador / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Kalibrator / Βοηθουομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / kalibrátor / 보정 물질 / Kalibratör / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Controllo / Kontrolle / Μέτροπος / 质控品 / Kontrollné / Kontroll / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

Ab

Antiserum / Antisérum / Antisiero / Antisero / Anti-soro / Αντισώμα / 抗血清 / Antiserumas / Antisérum / Antisurovica / antisérum / 항혈청 / Антисыворотка / Антисерум / 抗血清

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instruções de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用說明

REAG **PREC**

Precipitating Agent / Réactif précipitant / Ausfällungsreagenz / Reattivo precipitante / Reactivo de precipitación / Reagente de precipitação / Fällningsreagens / Αντιροσθηριο καθίζησης / 沉淀试剂 / Nusodinamasis reagentas / Precipitaló reagens / Odczynnik do precypitacji / Precipitační reagencie / Precipitačná reagentia / 침전 시약 / Çökeltme Reaktifi / Осаждающий реагент / Реактив за утаяване / 沉澱試劑

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLIKE / ОПАСНО / ОПАСНОСТ / 危險

23 March 2020

REFERENCES

1. Burtis CA, Ashwood ER: Tietz Textbook of Clinical Chemistry, 2nd ed. W.B. Saunders Company, Philadelphia, 1994, p.1863.
2. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
3. Stanczyk FZ: Steroid Hormones. In: Mishell DR, Davajan V, Lobo RA (eds.): Infertility Contraception and Reproductive Endocrinology, 3rd. ed. Cambridge, Massachusetts, Blackwell Scientific Publications, Inc. 1991, p.69.
4. Levant SG, Barnes RB: Pharmacology of estrogen. In: Lobo RA (ed.): Treatment of the Postmenopausal Woman. New York: Raven Press, 1993, p.57.
5. Buster JE: Estrogen metabolism. In: Sciarra JJ (ed.): Endocrinology, Infertility, Genetics, Volume 5, Philadelphia: Harper & Row, 1987.
6. Cassidenti DL, Vijad AF, Vijad MA, Stanczyk FZ, Lobo RA.: Short Term Effects of Smoking on the Pharmacokinetic Profiles of Micronized Estradiol in Postmenopausal Women. Am. J. Obstetrics and Gynecology. 163:1953. 1990.