Product information

Information about other products is available at: www.demeditec.com

Ferritin ELISA

REF DE4408

 **96**

Demeditec Diagnostics GmbH
Lise-Meitner-Strasse 2
24145 Kiel – Germany
www.demeditec.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.
Por favor, usar a versão válida das instruções de utilização fornecidas com o kit.**

Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas	
The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente : Foram efetuadas as seguintes alterações em comparação com a versão anterior:	
Detailed editorial revision. Ausführliche redaktionelle Überarbeitung. Revisione editoriale dettagliata. Revisión editorial detallada. Révision éditoriale détaillée. Revisão editorial detalhada.	
1 INTENDED PURPOSE:	Indication of clinical use
1.1 CLINICAL SIGNIFICANCE:	Added
4 WARNINGS:	Updated
3.1 Reagents and materials supplied in the kit:	Enzyme Conjugate volume increased to 21 mL (before: 12 mL)
8.1 Preparation of Calibrators and Controls:	Concentration for Standard 5 has changed to 800 ng/mL (before: 1000 ng/mL);
8.4 Procedure:	Pipetting volume for standards/samples decreased to 10 µL (before: 20 µL), Pipetting volume for conjugate increased to 200 µL (before: 100 µL), Incubation time for substrate increased to 15 minutes (before: 10 minutes);
9 QUALITY CONTROL:	Updated
11 METROLOGY AND TRACEABILITY	Standardised against internal reference standards (serum matrix) (before: calibrated against the WHO 1 st IS Ferritin 80/602).
13 EXPECTED VALUES:	Updated data determined during assay validation:
14 PERFORMANCE CHARACTERISTICS:	Updated and additional data.
18 BIBLIOGRAPHY:	Updated

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

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For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Ferritin ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of ferritin in human serum or plasma.

Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of diseases affecting iron metabolism.

1.1 CLINICAL SIGNIFICANCE

Ferritin is a protein mainly present in the cytoplasm but also found in serum and plasma as an iron-carrier protein. It contains, stores and releases iron in a controlled fashion. One molecule of ferritin is capable of binding between 4000 and 5000 atoms of iron, making ferritin the major iron storage protein for the body^{1,2}.

Iron is an essential micronutrient as it is required for an adequate erythropoietic function, oxidative mechanism and cellular immune response^{1,2,3}. Under physiologic conditions there is a balance between iron absorption, iron transport and iron storage in the human body. When this balance is interrupted, two main conditions may present: iron deficiency and iron overload^{1,2,4}. In humans, ferritin acts as a buffer against iron deficiency and overload.

Iron deficiency is the most common and widespread nutritional disorder in the world, affecting 2 billion people and in particular children, women of childbearing age with heavy menstrual flow and miscarriages, and premenopausal women. The WHO estimates that 30% of nonpregnant and 42% of pregnant women are affected by iron deficiency and anemia⁵⁻¹⁰.

Since the concentration of ferritin is directly proportional to the total iron stores in the body, serum ferritin concentrations are considered as a common diagnostic tool in the diagnosis of diseases affecting iron metabolism¹.

In addition to ferritin levels additional markers assessed as part of the routine diagnosis of iron deficiency are transferrin saturation and the assessment of haemoglobin levels. All these parameters should also be considered together with other parameters such as the familial history of the patients and the evaluation of clinical symptoms.

2 PRINCIPLE OF THE METHOD

The Ferritin ELISA is based on simultaneous binding of human Ferritin to two monoclonal antibodies, one immobilised on the microwell plate and the other conjugated with horseradish peroxidase (HRP). After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is proportional to the ferritin concentration of in the sample.

Ferritin concentration in the sample is calculated through a calibration curve.

3 REAGENTS, MATERIALS, AND INSTRUMENTATION

3.1 Reagents and materials supplied in the kit

- CAL 0** Zero Standard (1 vial) 3 mL and
CAL 1 - 5 Standard (Standard 1 - 5), 5x (1 vial = 1 mL)
- CONTROL** Ferritin Control (1 vial) 1 mL
Control concentration is indicated on the Certificate of Analysis
- ENZ CONJ** Enzyme Conjugate (1 vial) 21 mL
Anti-Ferritin antibody conjugated with Horseradish peroxidase (HRP)
- SORB MT** Microtiterwells (1 microplate breakable)
Anti-Ferritin antibody adsorbed on microplate
- SUB TMB** Substrate Solution (1 vial) 15 mL
H₂O₂-TMB 0.26 g/L (avoid any skin contact)
- STOP SOLN** Stop Solution (1 vial) 15 mL
Sulphuric acid 0.15 mol/L (avoid any skin contact)
- WASH SOLN 10x** Wash Solution 10X (1 vial) 50 mL
Phosphate buffer 0.2 M, pH 7.4

3.2 Materials required but not provided

Distilled water

3.3 Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

4 WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5 PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction for Use.
- All reagents should be stored refrigerated at 2 °C - 8 °C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 °C - 28 °C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6 REAGENT STORAGE AND STABILITY

Store the kit at 2 °C - 8 °C in the dark.

The kit is stable at 2 °C - 8 °C until the expiry date stated on the external kit label.

Once opened, the kit is stable at 2 °C - 8 °C for 5 months.

The diluted wash solution is stable for 30 days at 2 °C - 8 °C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

7 SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	72 hours
Freeze/thaw cycles	3 cycles

8 PROCEDURE

8.1 Preparation of Calibrators and Controls

The calibrators are ready to use and are provided at the following concentrations:

	S0	S1	S2	S3	S4	S5
ng/mL	0	5	20	100	400	800

The Control is ready to use.

8.2 Preparation of the Wash Solution

Dilute the content of the vial "Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

8.3 Preparation of Samples

Ferritin determination should be performed with human serum or plasma.

8.4 Procedure

- **Allow all reagents to reach room temperature (22 °C - 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 °C - 8 °C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 °C - 8 °C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate to improve accuracy of the test results, prepare two wells for each point of the calibration curve (S0 - S5), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		10 µL	
Standard S0 – S5	10 µL		
Conjugate	200 µL	200 µL	
Incubate for 1 hour at room temperature (22 °C - 28 °C). Remove the content from each well. Wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each manual washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate for 15 minutes at room temperature (22 °C - 28 °C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620 - 630 nm or against Blank within 5 minutes.			

9 QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of the control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DEMEDIATEC recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹².

10 CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A cubic spline or 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is recommended**.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

11 MEASURING RANGE

The assay measuring range (AMR) is 2.84 – 800 ng/mL.

Any value that reads below 2.84 ng/mL should be reported as “< 2.84 ng/mL”. Any value that reads above 800 ng/mL should be reported as “> 800 ng/mL”.

12 METROLOGY AND TRACEABILITY

The Ferritin ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to another commercially available test method.

13 EXPECTED VALUES

The following ranges were determined using the Ferritin ELISA and are provided for information only. The 95% reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

	n	Median (ng/mL)	Reference interval (ng/mL)
Women			
Pre-menopausal	42	33.5	8.4– 156.9
Post-menopausal	42	44.8	16.8 – 180.8
Men	84	122.5	21.7 – 276.4

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

14 PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

14.1 Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.35 ng/mL
Limit of Detection (LoD)	2.00 ng/mL
Limit of Quantitation (LoQ)	2.84 ng/mL

14.2 Trueness

Trueness has been demonstrated through method comparison of the Ferritin ELISA to a commercially available assay using native donor samples – refer to section 15.5

14.3 Precision

Precision of the Ferritin ELISA was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in five replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV
1	75	9.21	1.23	13.3%
2	75	68.30	2.47	3.6%
3	75	178.38	9.44	5.3%
4	75	380.97	26.54	7.0%
5	75	467.72	25.58	5.5%
6	75	619.17	51.33	8.3%

Reproducibility: A total of 6 serum samples were assayed in five replicates, once a day for 5 days by 3 operators.

Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV
1	150	8.96	1.21	13.5%
2	150	67.85	4.31	6.4%
3	150	177.92	14.14	7.9%
4	150	377.41	32.10	8.5%
5	150	460.17	46.83	10.2%
6	150	614.18	75.43	12.3%

14.4 Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For ferritin concentration by Ferritin ELISA, the measurement procedure shows linearity for the interval from 2.5 to 849 ng/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

14.5 Method comparison

The Ferritin ELISA was compared against a commercially available quantitative automated ferritin assay, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 100 samples, selected to represent a wide range of ferritin concentrations [15.65 to 782 ng/mL], were assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
100	0.96 [0.92 – 1.01]	0.24 [-5.46 – 3.07]	0.96

14.6 Analytical Specificity

The following substances do not interfere with a bias of $> \pm 15\%$ in the Ferritin ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	1000 mg/dL
Total Protein	10 g/dL
Triglyceride	3000 mg/dL

14.7 Serum-plasma study

The Ferritin ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples covering the assay range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.02 [0.99 – 1.25]	-1.02 [-4.98 – 0.28]	1.00
Lithium Heparin	1.04 [0.95 – 1.09]	-1.61 [-3.33 – 0.42]	1.00
Sodium Heparin	1.07 [0.98 – 1.21]	-1.92 [-4.51 – -0.22]	1.00
EDTA	0.97 [0.88 – 1.06]	-1.02 [-3.24 – 0.66]	1.00

14.8 Hook effect

No hook effect was observed with the Ferritin ELISA up to 10 000 ng/mL.

15 LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays¹³. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16 WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

17 PRODUCT COMPLAINTS

For a patient/user/third party in the European Union and in countries with similar regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

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
18 BIBLIOGRAPHY / BIBLIOGRAFIA / BIBLIOGRAFÍA / LITERATUR

1. Muñoz M, García-Erce JA, Remacha ÁF. Disorders of iron metabolism. Part II: iron deficiency and iron overload.
J Clin Pathol. 2011 Apr;64(4):287-96.
2. Moyer TP, Highsmith WE, Smyrk TC, Gross JB Jr. Hereditary hemochromatosis: laboratory evaluation.
Clin Chim Acta. 2011 Aug 17;412(17-18):1485-92.
3. Crownover BK, Covey CJ. Hereditary hemochromatosis.
Am Fam Physician. 2013 Feb 1;87(3):183-90.
4. Daru J, Colman K, Stanworth SJ, De La Salle B, Wood EM, Pasricha SR. Serum ferritin as an indicator of iron status: what do we need to know?
Am J Clin Nutr. 2017 Dec;106(Suppl 6):1634S-1639S.
5. Peyrin-Biroulet L, Williet N, Cacoub P. Guidelines on the diagnosis and treatment of iron deficiency across indications: a systematic review.
Am J Clin Nutr. 2015 Dec;102(6):1585-94.
6. De Franceschi L, Iolascon A, Taher A, Cappellini MD. Clinical management of iron deficiency anemia in adults: Systemic review on advances in diagnosis and treatment.
Eur J Intern Med. 2017 Jul;42:16-23.
7. Cappellini MD, Musallam KM, Taher AT. Iron deficiency anaemia revisited.
J Intern Med. 2020 Feb;287(2):153-170.
8. Friedman AJ, Shander A, Martin SR, Calabrese RK, Ashton ME, Lew I, Seid MH, Goodnough LT. Iron deficiency anemia in women: a practical guide to detection, diagnosis, and treatment.
Obstet Gynecol Surv. 2015 May;70(5):342-53.
9. Lopez A, Cacoub P, Macdougall IC, Peyrin-Biroulet L. Iron deficiency anaemia.
Lancet. 2016 Feb 27;387(10021):907-16.
10. Daru J, Allotey J, Peña-Rosas JP, Khan KS. Serum ferritin thresholds for the diagnosis of iron deficiency in pregnancy: a systematic review.
Transfus Med. 2017 Jun;27(3):167-174.
11. Basic QC Practices On-line Course; <http://www.Westgard.com>.
12. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'.
Clin Chem, 34, 1988, pp 27–33

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Française	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	utilisation Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di catalogo
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributed by	Vertrieb durch	Distribution par	Distribución por	Distribuzione da parte di
	Version	Version	Version	Versión	Versione
	Single-use	Einmalverwendung	À usage unique	Uso único	Uso una volta