Active® IGFBP-3 IRMA

REF DSL6600

TABLE OF CONTENTS

English	2	Slovenčina (SK)	26
Français (FR)	5	한국어 (KO)	29
Deutsch (DE)	8	Türkçe (TR)	32
Italiano (IT)	11	Русский (RU)	35
Español (ES)	14	中文 (ZH-TW)	38
Magyar (HU)	17	APPENDIX	40
Polski (PL)	20	REFERENCES	43
Čeština (CZ)	23		

IMMUNOTECH s.r.o., Radiova 1, 102 27 Prague 10, 捷克共和國

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Active® IGFBP-3 IRMA

REF DSL6600

IMMUNORADIOMETRIC KIT FOR THE QUANTITATIVE MEASUREMENT OF IGFBP-3 IN HUMAN SERUM THIS ASSAY IS INTENDED FOR IN VITRO DIAGNOSTIC USE

PRINCIPLE

The two-site immunoradiometric assay (IRMA) of Insulin-like Growth Factor Binding Protein-3 (IGFBP-3) is a sandwich-type assay [1]. In the kit, goat polyclonal antibody is used. Samples or calibrators are incubated in tubes pre-coated with polyclonal antibody, with 125I-labeled polyclonal antibody. After incubation, unbound reagents are removed by washing the tubes. The amount of 125I-labeled anti-IGFBP-3 bound to the tube is directly proportional to the concentration of IGFBP-3 present in the sample. A standard curve is constructed and unknown IGFBP-3 values are obtained from the curve by interpolation.

For Summary and Explanation of the Test see APPENDIX.

WARNING AND PRECAUTIONS

General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up [2].

Materials of human origin

All serum samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

The Materials Safety Data Sheet (MSDS) is available upon request.

GHS HAZARD CLASSIFICATION

Calibrators/Controls

WARNING

H317

May cause an allergic skin reaction.

H412

Harmful to aquatic life with long lasting effects.

P273

Avoid release to the environment.

P280

Wear protective gloves, protective clothing and eye/face protection.

P333+P313

If skin irritation or rash occurs: Get medical advice/attention.

P362+P364

Take off contaminated clothing and wash it before use.

reaction mass of:
5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)
<0.05%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum is the recommended sample type. Do not use heparinized or EDTA plasma.
- Allow serum samples to clot completely before centrifugation.
- Undiluted samples may be stored at 2-8°C for up to 48 hours or at <-20°C for up to 8 weeks, after aliquoting so as to avoid repeated freezing and thawing. Samples which have been diluted for the assay (see Procedure) may be stored at 2-8°C for up to 24 hours or at <-20°C for up to 4 weeks. Thawing of sample should be performed at room temperature.
- Typical serum samples should be diluted 1:100 with the IGFBP-3 zero calibrator prior to assay. Samples with very low known or expected values should be diluted 1:50 with IGFBP-3 zero calibrator prior to assay.
- Any sample reading greater than the highest calibrator should be further diluted appropriately with the IGFBP-3 zero calibrator and reassayed.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS PROVIDED

All reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reconstituted reagents are indicated in appropriate paragraphs.

Anti-IGFBP-3 antibody coated tubes: 2 x 50 tubes (ready-to-use)

Plastic tubes with goat-anti IGFBP-3 polyclonal immunoglobulin immobilized to the inside wall of each tube.

125I-labeled anti-IGFBP-3 tracer (YELLOW): one 22 mL vial (ready-to-use)

At the time of manufacture, the vial contains 370 kBq, (<10 µCi) of 125I-labeled goat-anti-IGFBP-3 polyclonal antibody in buffer with proteins (BSA), sodium azide (<0.1%) and a dye.

Calibrators: one 55 mL bottle labeled 0 (ready-to-use) and five vials labeled 1-5 (lyophilized)

The calibrator vials contain from 0 to approximately 100 ng/mL of IGFBP-3 in buffer with proteins (BSA) and ProClin 300. The exact concentration is indicated on each vial label. Calibrators are verified to an internal reference standard. The volume for reconstitution is indicated on the vial label. After reconstitution, calibrators may be stored at 2-8°C for up to 48 hours or at < -20°C for up to 2 weeks. Longer storage periods are not recommended.

Zero calibrator may be ordered separately, too (cat. #A99225)

Controls: two vials labeled 1, 2 (lyophilized)

The vials contain IGFBP-3 in buffer with proteins (BSA) and ProClin 300. The volume for reconstitution is indicated on the vial label, the concentration range after reconstitution is given on a supplement. After reconstitution, controls may be stored at 2-8°C for up to 48 hours or at < -20°C for up to 2 weeks. Longer storage periods are not recommended.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm polystyrene uncoated test tubes.
- Test tube rack for 12 x 75 mm tubes.
- Precision micropipets (10 µL, 50 µL and 1.0 mL).
- Semi-automatic pipets (200 µL, 3.0 mL).
- Shaker capable of ≥ 180 rpm.
- Vortex type mixer.
- A sponge rack for decantation or similar device or aspiration system.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.
- Log-log graph paper or computer with IRMA data analysis program.

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

Reconstitution of calibrators and controls

The content of the vials is reconstituted with the volume of distilled water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing.

Dilution of samples

Dilute human serum samples 1:100 with IGFBP-3 zero calibrator prior to assay. (example: 10 µL serum + 1.0 mL zero calibrator). Dilute human serum samples with known or expected very low concentrations of IGFBP-3 1:50 with zero calibrator prior to assay. (example: 10 µL serum + 500 µL zero calibrator).

Note: Do NOT dilute IGFBP-3 Controls or Calibrators. Do NOT multiply Control results by sample dilution factor.

Assay procedure

Bring all reagents to room temperature before pipeting.

Run Calibrators, Controls and patient samples in duplicate.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
Dilute human serum samples 1:100 with zero calibrator.	Incubate overnight (18-24 hours) at room temperature (18-25°C) on shaker set at ≥ 180 rpm.	Add 3.0 mL of distilled water to all tubes (except «total cpm» tubes), Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.
To antibody coated tubes successively add:	Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.	Repeat washing steps 2 more times (for a total of 3 wash steps).
50 µL of calibrator, control or sample, immediately add 200 µL of tracer.	Strike the tubes sharply on absorbent material to facilitate complete drainage and blot for 1-2 minutes.	Count bound cpm (B) and total cpm (T) for 1 minute.
Vortex gently.		

* Add 200 µL of tracer to 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of IGFBP-3 concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{cal} - cpm_{cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

Total activity: 146,587 cpm				
Calibrators	IGFBP-3 (ng/mL)	cpm (n=2)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	1,295	-	-
1	2.3	4,252	2.05	2,957
2	5.7	8,386	4.95	7,091
3	28.0	24,623	16.20	23,328
4	60.0	46,272	31.20	44,977
5	120.0	79,597	54.30	78,302

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

Multiply sample results by dilution factor (eg. 100).

Results may be converted to nmol/L based on the molecular weight of the non-glycosylated reference standard (28.75 kDa). To convert concentrations from ng/mL to nmol/L, multiply results by 0.035.

EXPECTED VALUES

Each laboratory should establish its own reference ranges. The pediatric data in the following table were obtained for children with normal height and growth patterns.

Age (years)	N	Mean (ng/mL)	SD (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
0 - 1	49	1,965	548	1,990	1,030 - 3,090
1 - 2	42	2,154	649	2,080	1,100 - 3,620
2 - 3	28	2,231	605	2,140	1,200 - 3,990
3 - 4	18	2,312	786	2,120	1,400 - 4,250
4 - 5	23	2,363	459	2,290	1,630 - 3,150
5 - 6	19	2,676	553	2,720	2,000 - 4,230
6 - 7	22	2,924	617	2,990	2,000 - 4,210

(For more details, see APPENDIX)

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling, or deterioration of reagents. In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.27 ng/mL

Functional sensitivity: 0.31 ng/mL

Specificity

1 µg/tube of the following peptide hormones gave concentrations of <3 ng/mL: hIGFBP-1 (purified, amniotic fluid), rhIGFBP-2, rhIGFBP-4, rhIGFBP-5, rhIGFBP-6, rhIGF-I, rhIGF-II.

Precision

Intra-assay

Samples were assayed 20 times in the same run. The coefficients of variation were ≤4.4%.

Inter-assay

Samples were assayed in duplicate in 15 different runs. The coefficients of variation were ≤13.5%.

Accuracy

Dilution test

Five serum samples were serially diluted with zero calibrator. The recovery percentages ranged from 82.3% to 120%.

Recovery test

Low-concentration serum samples were spiked with known quantities of IGFBP-3. The recovery percentages ranged from 80.1% to 116%.

Measurement range (from analytical sensitivity to highest calibrator): 0.27 to approximately 100.0 ng/mL.

LIMITATIONS

- Failure to follow these instructions for use (IFU) may significantly affect results.
- Failure to blot tubes adequately following decantation may result in poor replication and spurious values.
- Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.
- Use of EDTA or heparinized plasma is not recommended.
Serum containing IGFBP-3 protease activity may contain immunoreactive IGFBP-3 fragments. Measurement of intact IGFBP-3 in these samples may require preliminary size-exclusion chromatography [3].
- Avoid repeated freezing and thawing of reagents or specimens.
- Do not use haemolyzed, icteric or lipemic samples.
- The possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Summary and explanation of the test

The insulin-like growth factors, IGF-I and IGF-II, are potent mitogenic and anabolic peptides with structural homology to insulin [4]. The IGFs are produced by multiple body tissues under control of growth hormone (GH) and other regulatory factors, and act as autocrine, paracrine and endocrine regulators of cell growth. In serum and other biological fluids, the IGFs are complexed to six specific, structurally homologous binding proteins (IGFBPs) [5;6;7]. IGFBP-3 is the major IGFBP in postnatal serum, where its molar concentration in normal serum approximates the molar concentration of total IGF peptide (IGF-I + IGF-II) [8].

Human IGFBP-3 is a 264 amino acid peptide, with a 27 amino acid signal peptide and 18 cysteine residues [9]. Serum levels of IGFBP-3 rise gradually during childhood, increase more rapidly during puberty, reach peak levels at mid to late puberty, and fall during adult life [10]. Determination of serum IGFBP-3 levels promises to be of considerable value in assessment of disorders of the GH-IGF axis. Like IGF-I, IGFBP-3 is GH dependent; however, IGFBP-3 levels may be less age-dependent than IGF-I levels [10]. Total IGF levels show higher correlation with GH sufficiency than either IGF-I or IGF-II levels alone [11] and IGFBP-3 levels correspond with total IGF levels. In addition, IGF assays typically require an extraction procedure to remove interfering IGF-binding proteins; this step may be a source of assay error [12].

Assays for IGFBP-3 do not require prior sample extraction. These aspects of IGFBP-3 physiology and chemistry allow its measurement to reflect GH activity. Thus, in GH deficiency of either pituitary or hypothalamic origin, IGFBP-3 levels are uniformly low. Blum et al [10] have reported that 128/132 children with GH deficiency have serum IGFBP-3 concentrations below the 5th percentile for age. Similar data have been reported by Hasegawa et al. [13]. In situations of GH receptor deficiency, such as Laron syndrome, anti-GH antibodies, liver disease or malnutrition, IGFBP-3 levels are also reduced, even in the presence of elevated serum GH concentrations [3].

Sensitivity

The analytical sensitivity, or minimum detection limit, calculated as 2 standard deviations from the mean of 20 replicates of the 0 ng/mL IGFBP-3 calibrator, is 0.27 ng/mL.

Specificity

1 µg/tube of the following peptide hormones gave concentrations of <3 ng/mL:

PEPTIDE HORMONE
hIGFBP-1 (purified, amniotic fluid)
rhIGFBP-2
rhIGFBP-4
rhIGFBP-5
rhIGFBP-6
rhIGF-I
rhIGF-II

Interspecies Serum Comparisons:

NON-REACTIVE SERA	IMMUNOREACTIVE SERA
	Pig
	Cow
Sheep	Rabbit
Goat	Horse
Mouse	Dog
	Rat
	Guinea Pig

Since the actual concentrations of IGFBP-3 in the animal sera were not known and pure non-human IGFBP-3 was not available for testing, percent cross-reactivity could not be calculated. However, it was evident that none of the animal sera tested showed significant cross-reactivity or linearity in this assay.

Precision

Intra-assay

The intra-assay precision was determined from the mean of 20 replicates each.

Sample	N	Mean (ng/mL)	Standard deviation (ng/mL)	Coefficient of variation (%)
I	20	24.7	1.08	4.37
II	20	36.5	0.90	2.46
III	20	43.7	1.08	2.47

Inter-assay

The inter-assay precision was determined from the mean of 15 separate assays.

Sample	N	Mean (ng/mL)	Standard deviation (ng/mL)	Coefficient of variation (%)
I	15	17.6	1.79	10.2
II	15	44.2	5.97	13.5
III	15	61.5	3.83	6.23

Accuracy

Dilution test

Five serum samples were diluted with 0 ng/mL IGFBP-3 calibrator and assayed.

Sample	Dilution factor	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Ratio (%) Measured/Expected
S1	-	42.42	-	-
	1:2	24.23	21.21	114.2
	1:4	11.10	10.61	104.7
	1:8	4.66	5.30	87.88
	1:16	2.44	2.65	92.03
S2	1:32	1.33	1.33	100.3
	-	58.90	-	-
	1:2	32.85	29.45	111.5
	1:4	16.45	14.73	111.7
	1:8	6.94	7.36	94.26
S3	1:16	3.55	3.86	96.43
	1:32	1.91	1.84	103.8
	-	47.56	-	-
	1:2	25.31	23.78	106.4
	1:4	12.20	11.89	102.6
S4	1:8	5.21	5.95	87.64
	1:16	2.78	2.97	93.52
	1:32	1.65	1.49	111.0
	-	46.35	-	-
	1:2	27.60	23.18	119.1
S5	1:4	13.37	13.80	96.88
	1:8	5.68	6.90	82.32
	1:16	2.97	3.45	86.09
	1:32	1.56	1.73	90.43
	-	40.33	-	-
S5	1:2	24.16	20.17	119.8
	1:4	10.58	10.08	104.9
	1:8	4.62	5.04	91.64
	1:16	2.48	2.52	98.39
	1:32	1.50	1.26	119.0

Recovery test

IGFBP-3 was added to five serum samples and assayed according to the procedure of the kit.

Sample	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
S1	3.48	1.24	4.71	4.05	85.93
	3.42	3.05	6.47	5.25	81.10
	3.34	5.95	9.30	7.83	84.24
S2	4.94	1.85	6.78	5.60	82.56
	4.89	3.05	7.94	6.82	85.93
	4.73	7.08	11.80	12.18	103.2
S3	7.78	3.05	10.82	9.19	84.90
	7.59	5.95	13.54	13.86	102.3
	7.38	9.26	16.64	18.58	111.7
S4	4.80	3.05	7.85	6.97	88.80
	4.69	5.95	10.64	10.13	95.22
	4.51	10.32	14.83	17.27	116.4
S5	6.05	3.05	9.10	7.29	80.13
	5.90	5.95	11.86	10.39	87.63
	5.69	10.32	16.01	17.22	107.6

Expected values

Serum IGFBP-3 levels increase with age during childhood, peak in mid- to late-puberty and decline gradually during adult life. Levels differ between male and female within each age group and pubertal stage. In addition, other factors, such as nutritional status, may affect IGFBP-3 levels. Therefore, it is important that each laboratory carefully define the characteristics of the reference population prior to determining expected values. The data below are from multi-center studies.

The pediatric data in the following tables were obtained for children with normal heights and growth patterns.

IGFBP-3 levels by chronological age

Male

Age (years)	N	Mean (ng/mL)	SD (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
7 - 8	35	2,501	1,082	3,510	1,250 - 6,350
8 - 9	22	3,539	665	3,515	2,300 - 5,050
9 - 10	25	3,694	794	3,880	2,190 - 5,190
10 - 11	26	3,382	1,083	3,160	1,800 - 7,060
11 - 12	26	3,502	873	3,305	2,000 - 5,470
12 - 13	27	3,787	993	3,840	1,820 - 6,990
13 - 14	27	4,467	1,119	4,650	2,400 - 7,330
14 - 15	20	4,652	1,524	4,885	1,700 - 6,940
15 - 16	19	4,610	1,494	4,900	2,100 - 7,170
16 - 18	22	4,625	1,294	4,315	2,590 - 7,280
18 - 20	22	4,932	1,218	5,020	2,680 - 7,290
20 - 23	31	4,645	1,020	4,490	2,930 - 7,380
23 - 25	22	4,253	907	4,415	2,250 - 5,480
25 - 30	52	3,642	957	3,405	2,330 - 6,680
30 - 40	43	3,445	851	3,350	1,730 - 5,590

Female

Age (years)	N	Mean (ng/mL)	SD (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
7 - 8	21	4,067	1,137	4,020	2,060 - 6,530
8 - 9	17	3,923	830	4,070	2,560 - 5,530
9 - 10	25	4,440	1,020	4,350	2,860 - 7,740
10 - 11	24	4,116	1,115	4,105	2,690 - 7,200
11 - 12	20	4,771	1,447	4,200	2,300 - 7,740
12 - 13	23	4,710	1,699	4,400	1,800 - 8,410
13 - 14	32	4,310	1,225	4,170	2,000 - 7,110
14 - 15	24	4,371	1,089	4,220	2,600 - 7,320
15 - 16	23	4,440	944	4,630	2,400 - 5,980
16 - 18	24	4,319	1,084	4,190	2,000 - 6,470
18 - 20	21	4,430	1,530	4,840	2,310 - 7,480
20 - 23	28	4,668	1,446	4,340	2,760 - 7,350
23 - 25	20	4,567	1,319	4,185	2,920 - 7,000
25 - 30	48	4,220	1,191	3,985	2,050 - 7,600
30 - 40	29	3,789	1,114	3,780	2,300 - 7,260

Both Sexes

Age (years)	N	Mean (ng/mL)	SD (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
30 - 40	90	3,493	900	3,330	1,730 - 7,260
40 - 50	44	3,152	497	3,140	2,080 - 4,310
50 - 70	41	2,966	439	2,900	2,020 - 3,990

PATHOLOGIC SAMPLES

Diagnosis	Age (years, mean \pm SD)	N	Mean (ng/mL)	SD (ng/mL)	Median (ng/mL)	Absolute range (ng/mL)
Growth Hormone Deficiency	8.8 \pm 4.6	43	1,354	695	1,210	250 - 2,860
Constitutional Growth Retardation	10.8 \pm 3.2	18	2,607	719	2,530	1,520 - 4,320
Normal Variant Short Stature	7.3 \pm 4.2	20	1,994	601	1,830	1,270 - 2,990
Acromegaly	adult	57	5,777	1,374	5,540	1,550 - 9,370

¹²⁵I Characteristics

T_{1/2} (¹²⁵I) = 1443 h = 60.14 d

¹²⁵ I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirație / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktivyvioj medžiaga / Radioaktiv / Radioaktivny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLÍKE / ОПАСНО / ОПАСНОСТ / 危險

Ag^{125I}

Ab^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

SOLN|WASH|20x

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvätlösningskoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalu koncentratas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztvoru pľucšacego / Koncentrát mycího roztoku 20X / Концентрат промышленного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

REFERENCES

1. Miles LEM, Lipschitz DA, Bieber CP and Cook JD: Measurement of serum ferritin by a 2-site immunoradiometric assay. *Analyt Biochem* 61:209-224, 1974.
2. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
3. Guevara-Aguirre J, Rosenfeld R. et al.: Growth hormone receptor deficiency in Ecuador: clinical and biochemical phenotype in two populations. *J Clin Endocrinol Metab* 76: 417-423, 1993.
4. Daughaday W, Rotwein P: Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations. *Endocr Rev* 10:68-91, 1989.
5. Rosenfeld RG, Lamson L, Pham H, Oh Y, Conover C, De Leon DD, Donovan SH, Ocrant I, Giudice L.: Insulin-like growth factor-binding proteins. *Rec Progr Horm Res* 46:99-163, 1990.
6. Lamson G, Giudice L, Rosenfeld RG: Insulin-like growth factor binding proteins: structural and molecular relationships. *Growth Factors*. 5:19-28, 1991.
7. Cohen P, Fielder PJ, Hasegawa Y, Frisch H, Giudice LC, Rosenfeld RG: Clinical aspects of IGF binding proteins. *Acta Endocrinol (Copenh.)* 124 (suppl 2):74-85, 1991.
8. Baxter R, Martin J: Radioimmunoassay of growth hormone-dependent insulin like growth factor binding protein in human plasma. *J Clin Invest* 78:1504-1512, 1986.
9. Report on the nomenclature of the IGF binding proteins. *Endocrinology* 130: 1736-1737, 1992.
10. Blum WF, Ranke MB, Kietzman K, Gauggel E, Zeisel HJ, Bierich J: A specific radioimmunoassay for the growth hormone (GH)-dependent somatomedin binding protein: its use for diagnosis of GH deficiency. *J Clin Endocrinol Metab* 70: 1292-1298, 1990.
11. Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL: Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr* 109:428-33, 1986.
12. Powell DR, Rosenfeld RG, Baker BK, Liu F, Hintz RL: Serum somatomedin levels in adults with chronic renal failure: the importance of measuring insulin-like growth factor I (IGF-I) and IGF-II in acid chromatographed uremic serum. *J Clin Endocrinol Metab* 63:1186-1192, 1986.
13. Hasegawa Y, Hasegawa T, Aso T, Kotoh S, Tsuchiya Y, Nose O, Ohyama Y, Araki K, Tanaka T, Saisyo S, Yokoya S, Nishi Y, Miyamoto S, Sasaki N, Stene M: Usefulness and limitation of measurement of insulin-like growth factor binding protein-3 (IGFBP-3) for diagnosis of growth hormone deficiency. *Endocrinol Japon* 39(6): 585-591, 1992.