

# Thymidine Kinase REA KIT

Instruction for use in local language is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs).

## REVISION HISTORY

Previous version:	Current version:
PI-IM1948-05	IFU-IM1948-01
—	IVDR requirements incorporated
Chapter INTENDED USE removed	Chapter INTENDED PURPOSE added
—	Chapter APPENDIX: Interference data added

**REF** IM1948

## FOR PROFESSIONAL USE ONLY

### INTENDED PURPOSE

Thymidine Kinase REA KIT is an in vitro diagnostic manual medical device intended to be used by healthcare professionals for the quantitative measurement of thymidine kinase enzymatic activity in human serum and plasma. Measurement of thymidine kinase is intended to be used in the primary diagnosis and follow-up of patients with malignancies of the hematopoietic system and in the combination as a complementary tumour marker to the tumour cell proliferation assessment [1, 2, 3].

### PRINCIPLE

There are two thymidine kinase (TK) isoenzymes in eucaryotic cells, TK 1 and TK 2, differing biochemically and electrophoretically. This radioenzymatic assay has been optimised for TK 1 estimation (also known as fetal TK, dTK-F, or cytosol-TK).

The radioenzymatic assay of thymidine kinase is based on enzymatic phosphorylation of 5-[<sup>125</sup>I]-deoxyuridine (<sup>125</sup>I-Tracer), catalysed by thymidine kinase contained in the sample or calibrator. 5-[<sup>125</sup>I]-deoxyuridine is transformed to 5-iododeoxyuridine monophosphate (<sup>125</sup>I-d-UMP). <sup>125</sup>I-d-UMP is then separated from the reaction mixture by adsorption on ion exchange resin. The sorbent is rinsed to remove 5-[<sup>125</sup>I]-deoxyuridine. The radioactivity is then determined in a gamma counter. The TK enzymatic activity in the samples is obtained by interpolation from the standard curve. The enzymatic activity is directly proportional to the radioactivity.

### WARNING AND PRECAUTIONS

#### General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipetting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Sodium azide disposal must be in accordance with appropriate local regulations.

## Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All patient specimens should be handled as potentially infectious and waste should be discarded according to the country rules.

The summary of safety and performance for this in vitro diagnostic medical device is available to the public in the European database on medical device (EUDAMED) when this database is available, and the information has been uploaded by the Notified Body. The web address of the EUDAMED public web site is: <https://ec.europa.eu/tools/eudamed>.

To search the information about this product in EUDAMED, use BUDI-DI: 150995905IM19486P

## GHS HAZARD CLASSIFICATION

Not classified as hazardous



Safety Data Sheet is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs)

## SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum or EDTA plasma are the recommended sample types.
- Allow serum samples to clot completely before centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -20°C, 1 year maximum), after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Dilution of samples with concentration greater than the highest calibrator is not recommended.

Serum and EDTA plasma values for 24 samples (serum values ranging from 2.67 to 12.83 U/L) were compared using the IM1948 Thymidine Kinase REA KIT. Results are as follows:

[EDTA-plasma] = 0.8055 [serum] + 0.7026

R = 0.9271

## MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Once reconstituted, all the reagents should be processed. **Do not store reconstituted reagents.**

**<sup>125</sup>I-Tracer: one vial** (lyophilized)

The vial contains 190 kBq, at the date of manufacture, of <sup>125</sup>I-labeled deoxyuridine in buffer containing bovine serum albumin and sodium azide (<0.1%).

**Calibrator: one vial** (lyophilized)

The calibrator vial contains Eppendorf tube with thymidine kinase with activity 80 U/L, lyophilized in buffer with bovine serum albumin and sodium azide (<0.1%). The calibrator is traceable to an internal reference standard.

**Control sample: one vial** (lyophilized)

The control vial contains Eppendorf tube with thymidine kinase lyophilized in human serum with sodium azide (<0.1%). The activity range is indicated on a supplement. The control sample is traceable to an internal reference standard.

**Buffer: one vial** (lyophilized)

The vial contains lyophilized buffer with cofactors and sodium azide (<0.1%).

**BSA: one vial** (lyophilized)

The vial contains lyophilized BSA with sodium azide (<0.1%).

**Sorbent: one vial** (powder)

Sorbent suspension in distilled water must be prepared before the assay.

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Polypropylene tubes.

*Use only new polypropylene tubes. Polystyrene is not suitable because high acetone concentrations lead to its structural damage and, in consequence, to discrepant results.*

- Precision micropipette (25 µL).
- Semi-automatic pipette (300 µL).
- Dispensers (500 µL and 2 mL).
- Eppendorf tubes.

- Vortex type mixer.
- Aspiration system.
- Thermostat, 37°C.
- Magnetic stirrer.
- Acetone p.a., 250 mL.

Use only acetone of purity grade „puriss“ or „p.a.“. The use of lower grades (e.g. purum) may lead to improper sorbent decantation so that the background radioactivity reaches up to several thousands cpm and the standard curve cannot be established.

- Gamma counter set for <sup>125</sup>I.

## PROCEDURE

### Preparation of reagents

Let all the reagents come to room temperature.

### Reconstitution of BSA

The content of the vial with lyophilized BSA is reconstituted with 2 mL of distilled water. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing.

### Preparation of working buffer solution

The content of the vial with lyophilized buffer is reconstituted with 5 mL of distilled water. Wait for 10 minutes following reconstitution and mix gently to avoid foaming. Add to a beaker with 25 mL of distilled water and then add reconstituted BSA. The total volume of this working buffer solution is 32 mL.

### Reconstitution and preparation of calibrators

The content of the vial with lyophilized calibrator is reconstituted with the volume of working buffer solution indicated on the label. Wait for 10 minutes following reconstitution and mix gently to avoid foaming (see warning below this §). The enzyme activity of obtained thymidine kinase solution is 80 U/L. The other calibrators are prepared by further dilution in Eppendorf tubes. Six following levels 80, 40, 20, 10, 5, 2.5 U/L are recommended for the construction of standard curve.

Warning: Avoid foaming of the enzyme solution. It could cause enzyme denaturation. Repeating aspiration and re-addition of enzyme solution with pipette tip is satisfactory for proper mixing. Do not vortex.

### Preparation of reaction mixture

The content of the tracer vial with <sup>125</sup>I-deoxyuridine is reconstituted with 30 mL of working buffer solution. Wait for 10 min following reconstitution and mix gently to avoid foaming.

### Preparation of the sorbent suspension

Add the sorbent powder in 20 mL of distilled water. Prepare sorbent suspension at least 1 hour before use.

### Reconstitution of control sample

The content of the vial with lyophilized control sample is reconstituted with 100 µL of working buffer solution. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing.

### Preparation of 50% acetone

Pour 250 mL of distilled water into 250 mL of acetone and homogenize.

### Assay procedure

Step 1 Additions*	Step 2 Incubation
To polypropylene tubes add successively: 25 µL of working buffer solution (for calibrator 0), calibrator, control or sample and 500 µL of reaction mixture. Vortex gently 1-2 seconds.	Incubate 3.5 hours at 37°C in thermostat.

\*Add 500 µL of reaction mixture to 2 additional tubes to obtain total cpm.

Step 3 Adsorption	Step 4 Counting
Add 300 µL of the sorbent suspension to all tubes (except the 2 tubes «total cpm»). Vortex gently 1-2 seconds.” Incubate 30 minutes at 18-25°C without shaking. Aspirate the supernatant so that 2-3 mm of liquid remain above the sediment. Prevent aspiration of the sediment. The determination may be interrupted at this stage and the stoppered tubes may be stored at <-18°C until the next day.	Add 2 mL of 50% acetone to all tubes (except the 2 tubes «total cpm»). Vortex briefly. Let the sorbent sediment for 5 minutes and aspirate the supernatant. Repeat the sorbent washing by 50% acetone three more times.  Count bound cpm (B) and total cpm (T) for 1 minute.

\*\*As the sorbent settles quickly, keep it in suspension to maintain its concentration during pipetting. The use of magnetic stirrer is recommended. Add sorbent suspension into each tube separately, repeated addition with dispenser is not recommended.

## RESULTS

Results are obtained from the calibrator curve by interpolation. The curve serves for the determination of analyte activity in samples measured at the same time as the calibrators.

### Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ( $cpm_{cal} - cpm_{cal0}$ ) on the log vertical axis and analyte activity of the calibrators on the log horizontal axis.

Other calculation methods may give slightly different results.

Total activity: 118,187 cpm				
Calibrators	TK (U/L)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	220	0.19	-
1	2.5	1,501	1.27	1,281
2	5	2,426	2.05	2,206
3	10	4,081	3.45	3,861
4	20	7,399	6.26	7,179
5	40	13,840	11.71	13,620
6	80	25,105	21.24	24,885

(Example of standard curve, do not use for calculation)

### Samples

For each sample, locate the cpm ( $cpm_{sample} - cpm_{cal0}$ ) or B/T value on the vertical axis and read off the corresponding analyte activity on the horizontal axis.

The enzyme activity is given in units [U] per 1 litre of sample, where the unit is defined:

$$1 \text{ U} = 1.2 \times 10^{-12} \text{ mol/s}$$

## EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following serum values are indicative only.

Number of subjects	Average	95 <sup>th</sup> percentile	99 <sup>th</sup> percentile
	(U/L)		
45	3.26	6.29	10.13

## QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling or deterioration of reagents.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: [imunochem@beckman.com](mailto:imunochem@beckman.com)

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of EU Member State in which the user and/or patient is located.

## PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

### Sensitivity

**Analytical sensitivity:** 0.74 U/L

**Functional sensitivity:** 1.87 U/L

### Precision

#### Intra-assay

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 4.23%.

#### Inter-assay

Serum samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 6.49%.

### Accuracy

#### Dilution test

The dilution test is not applicable because no dilution matrix is recommended.

**Measurement range** (from analytical sensitivity to the highest calibrator): 0.74 to 80 U/L.

## **LIMITATIONS**

Failure to follow these instructions for use (IFU) may significantly affect results.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples. For more details, see Appendix, § Interference.

Increased values can also be found e.g. in infections or inflammatory diseases. If a temporary infection is suspected, it may be necessary to repeat the test at a later occasion.

---

## APPENDIX

### PERFORMANCE CHARACTERISTICS

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

#### Interference

Serum samples containing thymidine kinase enzymatic activities (low and high) were spiked with multiple concentrations of the substances listed below and assayed using Thymidine Kinase REA KIT. Values were calculated as described in CLSI EP07, 3<sup>rd</sup> ed. [4]. Interference was determined by testing controls (no interfering substance added) and matched test samples (with interfering substance added). No interference (defined as a shift in dose > 15 %) was found for addition of interferent up to concentration stated in the table below.

Interferent	Test concentration
Conjugated bilirubin	412.0 ug/mL
Hemoglobin	7,403 µg/mL
Triglycerides	22.70 mg/mL
Unconjugated bilirubin	351.3 µg/mL

In spite of hemoglobin, bilirubin (conjugated, unconjugated) and triglyceride interference data in the table, we advise to avoid using hemolyzed, lipemic or icteric samples.

#### Precision

##### Intra-assay

Serum	S1	S2	S3
Number of determinations	25	25	25
Mean value, (U/L)	6.88	18.96	53.78
C.V. (%)	4.23	4.04	3.68

EDTA plasma	P1	P2	P3
Number of determinations	25	25	25
Mean value, (U/L)	68.75	2.18	38.07
C.V. (%)	2.69	6.38	2.72

##### Inter-assay

Serum	S1	S2	S3
Number of determinations	10	10	10
Mean value, (U/L)	17.21	52.05	8.89
C.V. (%)	3.39	4.39	6.49




EDTA plasma	P1	P2	P3
Number of determinations	10	10	10
Mean value, (U/L)	11.28	54.16	6.55
C.V. (%)	9.01	5.76	6.49

#### <sup>125</sup>I Characteristics

$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$

<sup>125</sup> I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

## Symbols Key

	Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референца за производ / 產品參考
	In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználás / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Диагностика in vitro / За ин витро диагностика / 體外診斷
	Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄
	Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商
	Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Råcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 内容物足夠執行 <n> 次測試
	CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識
	Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表
	Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postępujcie podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanna Talimatna Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明
	Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo(i) di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度范围 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍 請參閱使用說明
	Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / İspijimas / Figelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意
	Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Dátum expirácie / 만료 날짜 / Son Kullanna Tarihi / Срок годности / Срок на годност / 到期日
	Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號
	Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

## REFERENCES

1. Topolcan O, Holubec L Jr. The role of thymidine kinase in cancer diseases. Expert Opin Med Diagn. 2008 Feb;2(2):129-41. doi: 10.1517/17530059.2.2.129. PMID: 23485133.
2. O'Neill KL, Buckwalter MR, Murray BK. Thymidine kinase: diagnostic and prognostic potential. Expert Rev Mol Diagn. 2001 Nov;1(4):428-33. doi: 10.1586/14737159.1.4.428. PMID: 11901857.
3. Alegria M, Robinson R, O'Neill K. Thymidine Kinase 1: A Universal Marker for Cancer. Cancer and Clinical Oncology, 2012; 2(1). doi:10.5539/cc.o.v2n1p159.
4. Approved Guideline - Interference Testing in Clinical Chemistry, EP07 3<sup>rd</sup> Edition. April 2018. Clinical and Laboratory Standards Institute.



IMMUNOTECH s.r.o., Radiova 1122/1, 102 00 Prague 10, Czech Republic

www.beckmancoulter.com