Product information

Information about other products is available at: www.demeditec.com

User's Manual

T3 free ELISA

Enzyme immunoassay for the quantitative measurement of free triiodothyronine (fT3) in serum

DE3801

96 wells

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Sommaire

1	INTRODUCTION	4
2	PRINCIPLE OF THE TEST	4
3	WARNINGS AND PRECAUTIONS	5
4	REAGENTS	6
5	SPECIMEN COLLECTION AND PREPARATION	7
6	ASSAY PROCEDURE	7
7	EXPECTED NORMAL VALUES	9
8	QUALITY CONTROL	9
9	PERFORMANCE CHARACTERISTICS	10
10	LIMITATIONS OF USE	11
11	LEGAL ASPECTS	11
1	EINLEITUNG	12
2	TESTPRINZIP	12
3	VORSICHTSMAßNAHMEN	12
4	BESTANDTEILE DES KITS	13
5	PROBENVORBEREITUNG	14
6	TESTDURCHFÜHRUNG	14
7	ERWARTETE WERTE	16
8	QUALITÄTSKONTROLLE	16
9	ASSAY-CHARAKTERISTIKA	16
10	GRENZEN DES TESTS	17
11	RECHTLICHE GRUNDLAGEN	17
1	INTRODUZIONE	18
2	PRINCIPIO DEL TEST	18
3	PRECAUZIONI	18
4	COMPONENTI DEL KIT	19
5	CAMPIONI	20
6	ATTUAZIONE DEL TEST	20
7	VALORI NORMALI	22
8	CONTROLLO QUALITÀ	22
9	CARATTERISTICHE DEL TEST	22
10	LIMITAZIONE DEL TEST	23
11	ASPETTI LEGALI	23

1	INTRODUCCIÓN	24
2	FUNDAMENTO DEL ENSAYO	24
3	PRECAUCIONES.....	24
4	COMPONENTES DEL KIT.....	25
5	MUESTRAS.....	26
6	PROCEDIMIENTO DE ENSAYO	26
7	VALORES ESPERADOS	28
8	CONTROL DE CALIDAD	28
9	CARACTERÍSTICAS DEL ENSAYO	28
10	LIMITACIONES DE USO	29
11	ASPECTOS LEGALES	29
1	INTRODUCTION.....	30
2	PRINCIPE DU TEST	30
3	PRECAUTIONS D'UTILISATION.....	30
4	COMPOSITION DU KIT	31
5	ECHANTILLON	32
6	RÉALISATION DU TEST	32
7	VALEURS ATTENDUES.....	34
8	CONTROLE DE QUALITE	34
9	CARACTERISTIQUES DU TEST	34
10	LIMITES D'UTILISATION.....	35
11	ASPECTS LEGAUX.....	35
12	REFERENCES / LITERATURE	36
	SYMBOLS USED WITH DEMEDITEC ASSAYS	36

1 INTRODUCTION

1.1 Intended Use

The **Demeditec T3 free ELISA** is an enzyme immunoassay for the quantitative **in vitro diagnostic** measurement of free triiodothyronine (fT3) in serum.

1.2 Summary and Explanation

Triiodothyronine is the most potent thyroid hormone of the human body. It effects various physiological functions including stimulation of metabolism and development of the foetal and neonatal brain (1,8). In blood, the vast majority of T3 (>99.5%) is bound to carrier proteins. The main transport protein is thyroxine-binding globulin (TBG). Only 0.2% - 0.5% of total T3 remains unbound (2). However, only the unbound (free) fraction of T3 exerts a biological activity. In many clinical conditions such as pregnancy, the concentration of carrier proteins changes, resulting in an alteration of total T3 levels while the concentration of free T3 remains constant (3). The biological activity of Triiodothyronine is regulated at the target tissue level by de-iodination and plasma membrane transport (4).

There is a gender-specific response to thyroid hormone with aging (5). The measurement of free T3 concentration more reliably correlates with the clinical status than total T3 levels.

Increased fT3 concentrations are observed with hyperthyroidism, particularly the isolated T3-hyperthyreosis. Decreased fT3 concentrations are observed with hypothyroidism, long-term therapy with thyreostatics, glucocorticoids, propranolol or amiodarone, chronically ill patients (low T3-T4 syndrome) or anorexia. Dialysis patients in general show marginally lower fT3 concentrations.

Free T3 determination provides an effective method for confirming manifest and isolated T3 hyperthyroidism (6). Moreover, free T3 can be used in subsequent treatment monitoring, as free T3 levels are generally elevated in patients with hyperthyroidism (7).

2 PRINCIPLE OF THE TEST

The Demeditec T3 free ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **principle of competitive binding**.

The microtiter wells are coated with a polyclonal (sheep) antibody directed towards a unique antigenic site of the T3 molecule.

Endogenous fT3 of a patient sample competes with a T3-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase conjugate is inversely proportional to the concentration of fT3 in the sample.

Having added the substrate solution, the intensity of colour developed is inversely proportional to the concentration of fT3 in the patient sample.

3 WARNINGS AND PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C - 8 °C in the sealed foil pouch and used in the frame provided.
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (21 °C - 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec.

4 REAGENTS

4.1 Reagents provided

1. **SORB MT Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-T3 antibody (polyclonal).
2. **CAL 0 – 5 Standard (Standard 0 - 5)**, 6 vials, 1 mL each, ready to use;
Concentrations: 0 – 1.0 – 2.5 – 5.0 – 10.0 – 20.0 pg/mL;
Conversion: 1 pg/mL × 1.536 = pmol/L;
Contain non-mercury preservative.
3. **CONTROL 1 & 2 Control Low & High**, 2 vials, 1 mL each, ready to use; For control values and ranges please refer to vial label or QC-Datasheet. Contain non-mercury preservative.
4. **ENZ CONJ Enzyme Conjugate**, 1 vial, 14 mL, ready to use; T3 conjugated with horseradish peroxidase;
Contains non-mercury preservative.
5. **SUB TMB Substrate Solution**, 1 vial, 25 mL, ready to use; Tetramethylbenzidine (TMB).
6. **STOP SOLN Stop Solution**, 1 vial, 14 mL, ready to use; Contains 0.5 M H₂SO₄, Avoid contact with the stop solution. It may cause skin irritations and burns.
7. **WASH SOLN 40x Wash Solution**, 1 vial, 30 mL (40X concentrated); See “Reagent Preparation”.

Note: Additional Standard 0 for sample dilution is available upon request.

4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 ± 10 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2-8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8 °C. Microtiter wells must be stored at 2-8 °C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Wash Solution

Add deionized water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated Wash Solution with 1170 mL deionized water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, Demeditec has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 4 days at 2-8 °C prior to assaying.

Specimens held for a longer time (up to 3 months) should be frozen only once at -20 °C prior to assay.

Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with Standard 0 and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

a) dilution 1:10: 10 µL sample + 90 µL Standard 0 (mix thoroughly)

b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL Standard 0 (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **50 µL** of each **Standard, Control** and **samples** with new disposable tips into appropriate wells.
3. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells.
Rinse the wells **4 times** with **300 µL - 400 µL** diluted Wash Solution per well. Strike the wells sharply on absorbent paper to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

6. Add **150 µL** of **Substrate Solution** to each well.
7. Incubate for **30 minutes** at room temperature.
8. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
9. Determine the absorbance (OD) of each well at **450 ± 10 nm** with a microtiter plate reader
It is recommended that the wells be read **within 10 minutes** after adding the Stop Solution.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 20 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard		Optical Units (450 nm)
Standard 0	0.0 pg/mL	1.92
Standard 1	1.0 pg/mL	1.69
Standard 2	2.5 pg/mL	1.43
Standard 3	5.0 pg/mL	1.12
Standard 4	10.0 pg/mL	0.81
Standard 5	20.0 pg/mL	0.55

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently healthy adults (18 to 81 years), using the Demeditec T3 free ELISA the following data were observed:

Population	n	Mean (pg/mL)	Median (pg/mL)	2.5 th - 97.5 th Percentile (pg/mL)	Range (min. - max.) (pg/mL)
Males	59	3.22	3.21	2.34 - 4.33	2.21 - 4.48
Females	58	3.09	2.99	2.19 - 4.88	1.99 - 5.58

The results correlate well to the ranges described in the literature (9, 10).

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels. The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results. Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or Demeditec directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.38 - 20 pg/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Sample	Concentration (ng/mL)	Cross-reactivity %
L-Thyroxine	1940	< 0.001
L-Tyrosine	100	< 0.001
3-Iodo-L-tyrosine	100	< 0.001
3,5-Diiodotyrosine	100	0.0012
Furosemide	100	< 0.001
Sodium Acetylsalicylate	100	< 0.001

9.3 Sensitivity

The analytical sensitivity of the Demeditec ELISA was calculated by subtracting 2 standard deviations to the mean of 20 replicate analyses of the Standard 0 and was found to be 0.381 pg/mL.

The Limit of Blank (LoB) is 0.327 pg/mL.

The Limit of Detection (LoD) is 0.536 pg/mL.

The Limit of Quantification (LoQ) is 0.990 pg/mL.

9.4 Reproducibility

9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	10	1.87	9.1
2	10	3.90	2.1
3	10	6.07	4.3
4	10	13.49	4.6

9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	30	1.95	9.5
2	30	4.09	12.0
3	30	6.04	9.9
4	30	13.34	7.4

9.4.3 Inter-Lot

The inter-assay (between-lots) variation was determined by repeated measurements of samples with 3 different kit lots.

Sample	n	Mean (pg/mL)	CV (%)
1	18	3.14	9.6
2	18	6.14	11.1
3	18	8.66	8.6
4	18	12.93	11.3

9.5 Recovery

Recovery of the Demeditec ELISA was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte.

	Sample 1	Sample 2	Sample 3	Sample 4
Concentration (pg/mL)	10.53	4.26	6.11	9.99
Average Recovery (%)	97.6	99.1	102.6	96.6
Range of Recovery (%)	from	96.5	96.4	100.9
	to	99.4	101.5	106.6

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of fT3 in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Demeditec.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement. Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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12 REFERENCES / LITERATURE

1. Bowen R (2010-07-24). "Physiologic Effects of Thyroid Hormones". Colorado State University. Retrieved: 2013-09-29.
2. Salter D R, Dyke C M, Wechsler A. Triiodothyronine (T3) and cardiovascular therapeutics: a review. *Journal of cardiac surgery*, 1992; 7(4), 363-374.
3. Osathanondh R, Tulchinsky D, Chopra, I J. Total and free thyroxine and triiodothyronine in normal and complicated pregnancy. *The Journal of Clinical Endocrinology & Metabolism*, 1976; 42(1), 98-104.
4. Visser TJ. Cellular uptake of thyroid hormones. In: De Groot LJ, et al. editors; 2016.
5. Suzuki S et al. Gender-specific regulation of response to thyroid hormone in aging. *Thyroid Res.* 2012; 26;5(1):1.
6. Saito O et al. (2012) Comparison between serum free triiodothyronine levels and body fluid distribution in hemodialysis patients. *Clin Exp Nephrol.*; 2012; 16(6), 952-8.
7. Bahn RS et al. Hyperthyroidism and other causes of thyrotoxicosis: management guidelines of the American thyroid association and the American association of clinical endocrinologists. *Thyroid*, 2011; 21, 593-641.
8. Sawin, C T et al. The free triiodothyronine (T3) index. *Annals of internal medicine*, 1978; 88(4), 474-477.
9. Kratzsch J et al. New reference intervals for thyrotropin and thyroid hormones based on National Academy of Clinical Biochemistry criteria and ultrasonography of the thyroid. *Clin. Chem.* 2005; 51:1480-6.
10. Kussmaul T et al. Thyroid analytes TSH, FT3 and FT4 in serum of healthy elderly subjects as measured by the Roche modular system: do we need age and gender dependent reference levels? *Clin Lab.* 2014; 60(9), 1551-9.

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Temperature de conservation	Temperatura de conservacion	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore