
Product information

Information about other products is available at: www.demeditec.com

User's Manual

Anti-R-TSH Ab RIA

Radioimmunoassay for the quantitative determination of thyrotropin receptor autoantibodies in human serum.

DETCT100

100 tubes

INTENDED USE

The Demeditec TSH receptor autoantibody (TRAb) Coated Tube (CT) kit is intended for use by professional persons only for the quantitative determination of thyrotropin receptor autoantibodies in human serum. Hyperthyroidism in Graves' disease is due to the presence of autoantibodies to the TSH receptor and measurement of these autoantibodies can be useful in disease diagnosis and management.

REFERENCES

J. Sanders et al. The Interaction of TSH Receptor Autoantibodies with 125 I-Labelled TSH Receptor
J Clin Endocrinol Metab 1999 84: 3797-3802

PATENTS

European patent EP 1021721 B1, US patent US 6,844,162 B1 and Japanese patent 4331403 apply.

ASSAY PRINCIPLE

In Demeditecs Anti-R-TSH RIA kit TSH receptor autoantibodies in patient sera, calibrators and controls are allowed to interact with TSH receptor coated onto plastic tubes. After a 2 hour incubation, the samples are discarded leaving TRAb bound to the immobilised TSH receptor. Porcine TSH labelled with 125I is added in a 2nd incubation step, where it interacts with TSH receptors which have not been blocked by bound TRAb. Any unbound 125I-labelled TSH is then removed from the tubes by a wash step prior to counting on a gamma counter. A lower level of radioactivity bound indicates the presence of TRAb in the test sample. The measuring range is 1 – 40 u/L (NIBSC 08/204).

STORAGE AND PREPARATION OF TEST SERUM SAMPLES

Sera to be analysed should be assayed soon after separation or stored, preferably in aliquots, at or below –20°C. 0.2 mL is sufficient for one assay (duplicate 100 µL determinations are recommended). Repeated freeze thawing or increases in storage temperature must be avoided. Incorrect storage of serum samples can lead to loss of TRAb activity. Do not use lipaemic or haemolysed serum samples. Do not use plasma in the assay. When required, thaw test sera at room temperature (20 – 25°C) and mix gently to ensure homogeneity. Centrifuge the serum prior to assay (preferably for 5 minutes at 10 – 15,000 rpm in a microfuge) to remove any particulate matter.

MATERIALS REQUIRED AND NOT SUPPLIED

- Tubes for total counts
- Suitable rack for assay tubes
- Pipettes capable of dispensing 50 µL, 100 µL and 1 mL
- Means of diluting wash buffer
- Pure water
- Orbital shaker
- Gamma counter

MATERIALS SUPPLIED IN 100 TUBE KITS

- 1x 11 mL ¹²⁵I-labelled TSH 180kBq/bottle (at manufacture)
- 5x 20 Tubes
- 1x10 mL Start Buffer
- 4x 0.7 mL Calibrators
- 1x 0.7 mL Positive Control
- 1x 0.7 mL Negative Control
- 1x 50 mL Concentrated Wash Solution

PREPARATION OF REAGENTS SUPPLIED

Store kits and all kit components at 2–8°C.

SORB CT 5 x 20 tubes	TSH Receptor Coated Tubes 20 tubes sealed in each foil bag. (A small amount of white residue may be present in the tubes but this does not affect assay performance.) After opening return any unused tubes to the original foil packet and seal. Then place foil bag in the self-seal plastic bag, with desiccant provided, and store at 2-8°C for up to 2 weeks.
START BUF 10 mL	Start Buffer. Coloured yellow. Ready for use.
CAL 1 - 4 4x 0.7 mL	Calibrators 1, 2, 8, and 40 units/L (units are NIBSC 90/672) Ready for use.
CONTROL - 1x 0.7 mL	Negative control. Ready for use.
CONTROL + 1x 0.7 mL	Positive control (See QC Datasheet for range) Ready for use.
TSH I-125 11 ml	¹²⁵I-Labelled TSH Coloured red. Ready for use.
WASH SOLN 10x 50 mL	Concentrated Wash Solution. Dilute to 500 mL with pure water before use. Store at 2–8°C up to kit expiry.

ASSAY PROCEDURE

Allow all reagents to stand at room temperature (20–25°C) for at least 30 minutes before use. A repeating type Eppendorf pipette is recommended for steps 1, 4, 5, 6 and 8.

1. Pipette **50 µL** of start buffer into each tube to be used.
2. Pipette **100 µL** of patient sera, calibrators (1-4) and controls into each tube.
3. Cover the tubes and incubate for **2 hours** at room temperature on an orbital shaker (200 shakes per min.).
4. Pipette **1 mL** of diluted wash solution into each tube and aspirate or decant each tube to waste. If decanting, allow the tubes to drain onto a clean, dry, absorbent surface for **2 minutes**.
5. Repeat wash step 4.
6. Pipette **100 µL** of ¹²⁵I labelled TSH into each tube and into two empty tubes for total counts. Cover the tubes and incubate for **1 hour** at room temperature on an orbital shaker (200 shakes per min)
7. Repeat wash step 4 twice (except tubes for total counts). If decanting, allow the tubes to drain as above for **5 minutes**.
8. Count each tube including total count tubes for ¹²⁵I for 1 minute.


RESULT ANALYSIS

Express the counts bound to each tube (B) as a percentage of the counts bound in the presence of the negative control (Bo) giving the % binding. A calibration curve can be established by plotting calibrator concentration on the x-axis (log scale) against the % binding on the y-axis (linear scale). Alternatively, % Inhibition (I) can be calculated using the formula;

$$\%I = 100 \times \left(1 - \frac{B}{B_0} \right)$$

Other data reduction methods can be used. Samples with high TRAb concentrations can be diluted in kit negative control. For example 20 µL of sample plus 180 µL of negative control to give a 10x dilution. Other dilutions (e.g. 100x) can be prepared from a 10x dilution or otherwise as appropriate. Some sera will not dilute in a linear way and we suggest that the dilution giving a value closest to 50% inhibition is used for calculation of TRAb concentration. The negative control has a concentration of 0 units/L, but can be assigned a value of 0.1 units/L to facilitate computer processing of data.

TYPICAL RESULTS (example only, not for calculation of actual results)

	% B	%I	units/L
Total Counts	83,978		
Control -	16.50	0	
CAL1	13.89	16	1
CAL2	12.24	26	2
CAL3	6.33	62	8
CAL4	2.17	87	40
Control +	9.53	42	3.9

ASSAY CUT OFF

	units/L
Negative	≤ 1
Borderline Positive	> 1-1.5
Positive	> 1.5

This cut off has been validated at Demeditec. However each laboratory should establish its own normal and pathological reference ranges for TRAb levels. Also it is recommended that each laboratory include its own panel of control samples in the assay.

CLINICAL EVALUATION

Clinical Specificity

242 samples from healthy blood donors were assayed in the TRAb CT kit. 242 (100%) were identified as being negative for TSH receptor autoantibodies.

Clinical Sensitivity

50 samples from patients diagnosed with Graves' disease were assayed using the TRAb CT kit. 46 (92%) were identified as being positive for TSH receptor autoantibodies.

Lower Detection Limit

The kit negative control was assayed 20 times and the mean and standard deviation calculated. The lower detection limit at 2 standard deviations was 0.33 units/L.

Inter Assay Precision (n=20)

Sample	units/L	CV (%)	% inhibition	CV (%)
1	0.9	17	12	13
2	2.2	13	27	7.9
3	5.2	7.1	57	5.2
4	13	15	74	1.9
5	24	16	84	2.1

Intra Assay Precision (n=25)

Sample	units/L	CV (%)	% inhibition	CV (%)
6	1.7	5.3	21	5.5
7	4.5	4.7	45	2.8

Clinical Accuracy

Analysis of sera from patients with autoimmune diseases other than Graves' disease indicated no interference from autoantibodies to thyroglobulin; thyroid peroxidase; glutamic acid decarboxylase; 21-hydroxylase; acetylcholine receptor; dsDNA or from rheumatoid factor.

Interference

No interference was observed when samples were spiked with the following materials; haemoglobin up to 0.5 mg/mL; bilirubin up to 0.20 mg/mL; Intralipid up to 1000 mg/dL; human LH up to 10 u/mL; hCG up to 160 u/mL; human FSH up to 15 u/mL and human TSH up to 0.3 u/L.

SAFETY CONSIDERATIONS**Start Buffer & Concentrated Wash Solution**

Signal word: Warning

Hazard statement(s)

H319: Causes serious eye irritation

Precautionary statement(s)

P280: Wear protective gloves/protective clothing/ eye protection/face protection

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

P337 + P313: If eye irritation persists: Get medical advice/attention

Follow the instructions carefully. Observe expiry dates stated on the labels and the specified shelf life for diluted reagents. Refer to Safety Data Sheet for more detailed safety information. The kit contains radioactive material. Users should make themselves aware of and observe any national and local legislation and codes of practice governing the use, storage, transportation and disposal of radioactive materials. Avoid all actions likely to lead to ingestion. Avoid contact with skin and clothing. Wear protective clothing and where appropriate personal dosimeters. Radioactive materials should only be used by authorised personnel and in designated areas. Wash hands thoroughly after handling. Monitor hands and clothing before leaving the designated area. Materials of human origin used in the preparation of the kit have been tested and found non reactive for HIV1 and 2, HCV antibodies and HbsAg but should none the less be handled as potentially infectious. Wash hands thoroughly if contamination has occurred and before leaving the laboratory. Sterilise all potentially contaminated waste, including test specimens before disposal. Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy. These materials should be handled as potentially infectious. Some components contain small quantities of sodium azide as preservative. With all kit components, avoid ingestion, inhalation, injection and contact with skin, eyes and clothing. Avoid formation of heavy metal azides in the drainage system by flushing any kit component away with copious amounts of water.

ASSAY PLAN

Allow all reagents and samples to reach room temperature (20–25 °C) before use	
Pipette:	50 µL Start buffer
Pipette:	100 µL Calibrators, Controls, Patient Sera
Incubate	2 hours at room temperature on an orbital shaker at 200 shakes/min
Pipette:	1 mL wash solution
Aspirate/Decant:	Tubes
Pipette:	1 mL wash solution
Aspirate/Decant:	Tubes (Drain for 2 minutes)
Pipette:	100 µL ¹²⁵ I TSH into each tube
Incubate:	1 hour at room temperature on an orbital shaker at 200 shakes/min
Pipette:	1 mL wash solution
Aspirate/Decant:	Tubes
Pipette:	1 mL wash solution
Aspirate/Decant:	Tubes (Drain for 5 minutes)
Count tubes in gamma counter for 1 minute	
For optimum results do not perform the assay at temperatures above 25 °C.	

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Temperature de conservation	Temperatura de conservacion	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore

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LITERATUR

J. Sanders et al. The Interaction of TSH Receptor Autoantibodies with ¹²⁵I-Labelled TSH Receptor
J Clin Endocrinol Metab 1999 84: 3797-3802

PATENTE

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European patent EP 1021721 B1, US patent US 6,844,162 B1 and Japanese patent 4331403.

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