# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

User's Manual

# Estrone-3-Sulfate equine ELISA

VET

REF DEV9933

Σ 96 wells

Demeditec Diagnostics GmbH  
Lise-Meitner-Strasse 2  
24145 Kiel – Germany  
[www.demeditec.com](http://www.demeditec.com)

**CONTENTS / INHALTSVERZEICHNIS**

1	INTRODUCTION.....	3
2	PRINCIPLE.....	3
3	WARNINGS AND PRECAUTIONS.....	3
4	REAGENTS.....	4
5	SAMPLES.....	5
6	ASSAY PROCEDURE.....	5
7	EXPECTED NORMAL VALUES.....	7
8	QUALITY CONTROL.....	8
9	PERFORMANCE CHARACTERISTICS.....	8
10	LIMITATIONS OF PROCEDURE.....	10
11	LEGAL ASPECTS.....	10
12	REFERENCES.....	11
13	SHORT INSTRUCTION.....	12
1	EINLEITUNG.....	13
2	METHODIK UND TESTPRINZIP.....	13
3	HINWEISE UND VORSICHTSMASSNAHMEN.....	13
4	KITBESTANDTEILE.....	14
5	PROBENENTNAHME UND -VORBEREITUNG.....	15
6	TESTDURCHFÜHRUNG.....	16
7	NORMALWERTE.....	18
8	QUALITÄTSKONTROLLE.....	19
9	TESTCHARAKTERISTIKA.....	19
10	LIMITATIONEN.....	21
11	RECHTLICHE GRUNDLAGEN.....	21
12	LITERATUR.....	22
13	KURZANLEITUNG.....	23
	Symbols used with Demeditec ELISAs.....	24

## 1 INTRODUCTION

### 1.1 Intended Use

The Demeditec Estrone-3-Sulfate equine ELISA is a competitive enzyme immunoassay for the quantitative determination of Estrone-3-Sulfate in mare serum.

### 1.2 Summary and Explanation

Estrone-3-Sulfate (E3S) is the predominant conjugated estrogen during pregnancy. It is produced by the fetus, possibly in association with the endometrium in the pregnant mare.

Different hormones are important for the complex events that occur during pregnancy in all mammals. In the mare these events include the maintenance of the corpus luteum function, formation of endometrial cups and development of secondary corpora lutea. Progesterone and PMSG (Pregnant Mare Serum Gonadotropine, eCG) and also free Estrogens, e.g. Estrone, are associated with these processes. It has been shown, that Estrone is rapidly conjugated after secretion and the ratio between conjugated and unconjugated estrogens is 100:1 in mare serum.

The conjugated estrogens, especially Estrone-3-Sulfate, provide the opportunity to improve the accuracy of pregnancy diagnosis, to monitor the pregnancy and to distinguish whether the fetal development is normal or impaired. The diagnosis of embryonic death is usually made by using techniques of palpation of the uterus per rectum or ultrasound echography. The determination of Estrone-3-Sulfate is an aid in the non-invasive diagnosis which allows monitoring of the fetoplacental unit during pregnancy. Only in mares with normal fetal development the values of Estrone-3-Sulfate show a tremendous increase between day 75 and 100 of gestation.

## 2 PRINCIPLE

The Demeditec Estrone-3-Sulfate equine ELISA test kit is a solid phase enzyme immunoassay (ELISA) in the microplate format, designed for the quantitative measurement of Estrone-3-Sulfate. The microplate is coated with a polyclonal antibody specific for Estrone-3-Sulfate.

Calibrators and samples as well as a biotin-labeled Estrone-3-Sulfate are pipetted into the antibody coated microplate. During one hour incubation the Estrone-3-Sulfate of the calibrator/sample competes with the biotin-labeled Estrone-3-Sulfate for the binding sites of the antibody fixed on the inner surface of the wells.

Afterwards, horseradish peroxidase-labeled streptavidine is added. During a 30 minutes incubation, this enzyme-labeled streptavidine binds to the biotinylated Estrone-3-Sulfate. All wells are washed to remove excess of non-bound enzyme-labeled streptavidine.

A chromogenic substrate, TMB (3,3',5,5'-Tetra-Methyl-Benzidine), is added to all wells. During a 30 minutes incubation, the substrate is converted to a colored end product (blue) by the fixed enzyme. Enzyme reaction is stopped by dispensing of hydrochloric acid as stop solution (change from blue to yellow). The color intensity is inversely related to the concentration of Estrone-3-Sulfate present in the sample. The optical density of the color solution is measured with a microplate reader at 450 nm.

## 3 WARNINGS AND PRECAUTIONS

1. This kit is strictly intended for veterinary use only. For use by staff, who is specially informed and trained in methods which are carried out by use of immunoassays only.
2. All blood components and biological materials should be handled as potentially hazardous in use and for disposal. Follow universal precautions when handling and disposing infectious agents.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2-8°C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (18-25°C) before starting the test. Temperature will affect the absorbance readings of the assay.

10. Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where samples or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may vary slightly.
17. Avoid contact with Stop Solution. It may cause skin irritations and burns.
18. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
19. For information please refer to Material Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH or on Demeditec homepage ([www.demeditec.com](http://www.demeditec.com)).

## 4 REAGENTS

### 4.1 Reagents Provided

1. **SORB MT** **Divisible Microplate**, 12 x 8 (break apart) strips with 96 wells, ready to use; coated with a polyclonal anti Estrone-3-Sulfate antibody
2. **CAL 0-6** **Calibrators**, 7 vials, lyophilized, for reconstitution see "Reagent preparation" (4.4). The Calibrators contain 0, 5, 10, 50, 100, 300 and 1000 ng/ml, respectively, in equine serum.
3. **BIOTIN CONJ** **Biotin-Labeled Estrone-3-Sulfate**, 1 vial, 11 ml, ready to use, clear; in a buffered solution.
4. **ENZ CONJ** **Enzyme-Labeled Streptavidine**, 1 vial, 3 ml, ready to use, blue; containing horseradish peroxidase-labeled streptavidine in a buffered matrix.
5. **SAM BUF** **Sample Buffer**, 1 vial, 11 ml, ready to use, yellow.
6. **SUB TMB** **TMB-Substrate Solution**, 1 vial, 22 ml, ready to use, clear. 3,3',5,5'-Tetramethyl-benzidine and hydrogen peroxide in a buffered solution.
7. **STOP SOLN** **Stop Solution**, 1 vial, 7 ml, ready to use; contains 2 M hydrochloric acid; clear. Avoid contact with the Stop Solution. It may cause skin irritations and burns.
8. **WASH SOLN 10x** **Wash Solution**, 1 vial, 50 ml, 10x concentrated; clear; see „Reagent preparation“.

### 4.2 Materials required but not provided

- Microplate reader capable for endpoint measurements at 450 nm
- Vortex mixer
- Microplate mixer operating at 900 rpm
- Distilled or deionized water
- Graduated cylinders for 500 ml
- Plastic container for storage of the wash solution
- Calibrated variable precision micropipettes
- Semi logarithmic graph paper or software for data reduction
- Absorbant paper and timer

### 4.3 Storage conditions

When stored at 2-8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. After first opening the reagents are stable for 30 days if used and stored properly. Keep away from heat and direct sunlight.

Microtiter wells must be stored at 2-8°C. Take care that the foil bag is sealed tightly.

Store Calibrators refrigerated, they will be stable at 2-8°C for 7 days after reconstitution. For longer storage aliquot and freeze at ≤-20°C.

Protect TMB-Substrate Solution from light.

#### 4.4 Reagent preparation

##### Calibrators:

Reconstitute lyophilized Calibrator 0 with 2.0 ml dist. water and Calibrator 1 through Calibrator 6 with 1.0 ml dist. water 30 min. before use. Mix gently.

##### Wash Buffer:

Dilute 50 ml of 10x concentrated Wash Solution with 450 ml dist. water to a final volume of 500 ml. The diluted Wash Solution is stable for at least 12 weeks at room temperature (18-25°C).

Precipitates may form when stored at 2-8°C, which should dissolve again by swirling at room temperature (18-25°C).

#### 4.5 Disposal of the kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

#### 4.6 Damaged test kit

In case of any severe damage of the test kit or components, Demeditec Diagnostics GmbH has to be informed in writing within one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

### 5 SAMPLES

For determination of Estrone-3-Sulfate equine serum is the preferred sample matrix. The procedure calls for 20 µl matrix per well. Blood collection should be as stress-free as possible. The samples may be stored refrigerated at 2-8°C for one week, or up to 6 months frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Samples expected to contain Estrone-3-Sulfate concentrations higher than the highest calibrator (1000 ng/ml) should be diluted with the Zero Calibrator before assay. The additional dilution step has to be taken into account for the calculation of the results.

Do not use hemolytic, icteric or lipemic specimens. Samples containing sodium azide should not be used in the assay.

### 6 ASSAY PROCEDURE

#### 6.1 General Remarks

- Do not interchange components of different lots.
- All components of the test kit supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent of foam formation.
- All reagents and specimens must be allowed to come to room temperature (18-25°C) before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross-contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Respect the incubation times as stated in this instructions for use.
- Calibrators, controls and samples should at least be assayed in duplicates.
- A calibrator curve must be established for every test run.
- For internal quality control we suggest to use **Estrone-3-Sulfate equine Control coded DEV9933C**. For more information please contact Demeditec.

## 6.2 Assay Procedure

### 1. Preparation of calibrators:

The calibrators are supplied lyophilized. At least 30 minutes before use, reconstitute the calibrator 0 with 2.0 ml of distilled water, and each of the remaining calibrators 1 through 6 with 1.0 ml of distilled water. Use volumetric pipettes and mix by gentle inversion.

	1	2	3	4	5	6	7	8	9	10	11	12
a	0	4	P2	P..								
b	0	4	P2	P..								
c	1	5	P3									
d	1	5	P3									
e	2	6	P4									
f	2	6	P4									
g	3	P1	P5									
h	3	P1	P5									

- Pipet **20 µl** of each **Calibrators, Controls** and **Samples** into the wells prepared.
- Add **100 µl** of **Biotin-Labeled Estrone-3-Sulfate** to every well.
- Add **100 µl** of **Sample Buffer** to every well.
- Incubate for **1 hour** at room temperature (18-25°C) on a plate shaker (900 rpm).  
**Afterwards do not wash or discard!**
- Add **25 µl** of **Enzyme-Labeled Streptavidine** to every well.
- Incubate for **30 minutes** at room temperature (18-25°C) on a plate shaker (900 rpm).
- Discard the content of the wells and wash each well **4 times** with **300 µl** buffered **Wash Solution**.
- Remove as much as possible Wash Solution by beating the microplate carefully.
- Add **200 µl** of **TMB-Substrate Solution** to all wells.
- Incubate without shaking for **30 minutes** at room temperature (18-25°C) in the dark.
- Add **50 µl** of **Stop Solution** to each well and mix carefully.
- Determine the absorbance of each well at 450 nm. It is recommended to read the wells within 15 minutes.

## 6.3 Calculation of Results

- Calculate the average absorbance values for each set of Calibrators, Controls and Samples.
- Using semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with concentration on the horizontal (X) axis and absorbance value on the vertical (Y) axis.
- Using the mean absorbance value for each sample, determine the corresponding concentration from the calibration curve.
- Automated method: The results in the package insert have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give different results.
- The concentration of the samples can be determined directly from this calibrator curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations, this dilution factor has to be taken into account.

## 6.4 Example of Typical Calibrator Curve

The figure below shows typical results for Demeditec Estrone-3-Sulfate equine ELISA tests. These data are intended for illustration only and should not be used to calculate results from another run.

Standard	Optical Density (450nm)
Calibrator 0 (0 ng/ml)	3.084
Calibrator 1 (5 ng/ml)	1.858
Calibrator 2 (10 ng/ml)	1.275
Calibrator 3 (50 ng/ml)	0.547
Calibrator 4 (100 ng/ml)	0.351
Calibrator 5 (300 ng/ml)	0.181
Calibrator 6 (1000 ng/ml)	0.100

## 7 EXPECTED NORMAL VALUES

In a reference range study equine serum samples were collected in the morning between 8 and 9 a.m. and in the evening between 5 and 6 p.m. Diurnal variations have not been observed. Analysis by the Demeditec Estrone-3-Sulfate equine test kit yielded the following results:

Group	Absolute Range (ng/ml)	n
Normal, non-pregnant mares	nondetectable - 10	22

During the course of pregnancy in 26 mares the Estrone-3-Sulfate serum concentrations were measured using the Demeditec Estrone-3-Sulfate equine ELISA.

Gestation Week	N	mean ng/ml	min ng/ml	max ng/ml
3	12	2.6	0.4	9.2
4	20	1.8	0.4	4.9
5	21	2.1	0.2	6.5
6	18	2.8	0.4	7.1
7	17	3.2	0.4	9.1
8	16	4.5	0.9	14.0
9	19	5.1	0.9	14.1
10	11	9.7	2.8	34.7
11	18	11	2.1	54.9
12	13	23	2.7	109
13	20	34	3.1	161
14	9	94	7.8	232
15	17	80	4.8	262
16	10	185	30.9	340
17	16	186	11.6	421
18	11	367	67.5	577
19	16	374	22.2	863
20	10	575	289	905
21	12	516	71.5	1233
22	9	497	247	889
23	11	571	106	861
24	12	634	263	1621
25	13	577	184	1229
26	9	528	255	949
27	12	523	213	1258
28	12	464	267	795
29	11	473	245	778
30	7	336	209	590
31	5	385	206	478
32	4	349	184	504
33	11	302	180	510
34	10	316	131	523
35	7	215	84	506
36	8	233	113	349
37	10	189	85	358
38	11	207	85	305
39	6	93	77	123
40	7	157	52	252
41	5	145	88	227
42	8	147	58	239
43	2	106	86	126
44	6	168	101	210
45	4	80	74	85
46	9	112	32	192
47	5	101	43	214
48	5	94	66	129

Because of differences which may exist between laboratories with respect of population, laboratory technique and selection of reference groups, it is recommended that each laboratory establishes its own normal and pathological ranges of equine Estrone-3-Sulfate. The reference ranges should be regarded as guidelines only.

## 8 QUALITY CONTROL

Good laboratory practice requires that controls are run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels. The controls and the corresponding results of the QC Laboratory are stated in the QC certificate included in the kit. The values and ranges stated on the QC certificate always refer to the current kit lot and should be used for direct comparison of the results.

Employ appropriate statistical methods for analyzing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices, microtiter plate reader, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or Demeditec directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Analytical Sensitivity

The analytical sensitivity was calculated by subtracting 2 standard deviations from the mean OD of 22 replicate analyses of Standard 0 and was found to be 0.424 ng/mL.

### 9.2 Specificity

The following materials have been evaluated for cross reactivity. The percentage indicates cross reactivity at 50% displacement compared to Estrone-3-Sulfate.

Compound	Concentration range (ng/ml)	Crossreactivity (%) at 50% binding
Estrone	1 – 10,000	10.44%
Estriol	1 – 10,000	<0.08%
Corticosterone	1 – 10,000	<0.08%
Equilin	1 – 10,000	<0.08%
Equilenin	1 – 10,000	0.70%
17-Hydroxyprogesterone	1 – 10,000	<0.08%
Androstendione	1 – 10,000	0.85%
5 $\alpha$ -Dihydrotestosterone	1 – 10,000	0.12%
Testosterone	1 – 10,000	0.08%
Progesterone	1 – 10,000	<0.08%
Androsterone	1 – 10,000	<0.08%
17 $\beta$ -Estradiol	1 – 10,000	<0.08%
Pregnenolone	1 – 10,000	<0.08%
$\beta$ -Estradiol-3-sulfate	1 – 10,000	0.68%
DHEA-S	1 – 10,000	0.10%
Cortisone	1 – 10,000	<0.08%

### 9.3 Reproducibility

#### 9.3.1 Intra-Assay

The intra-assay variation was determined by 20 replicate measurements of three serum samples within one run. The intra-assay variability is shown below:

	Serum 1	Serum 2	Serum 3
<b>Mean (ng/mL)</b>	30.1	176.4	472.6
<b>SD</b>	2.93	11.56	34.35
<b>CV (%)</b>	9.8	6.6	7.3
<b>n =</b>	20	20	20

#### 9.3.2 Inter-Assay

The inter-assay variation was determined by duplicate measurements of three serum samples in ten different runs.

	Serum 1	Serum 2	Serum 3
<b>Mean (ng/mL)</b>	28.2	161.4	424.7
<b>SD</b>	3.58	18.22	59.64
<b>CV (%)</b>	12.7	11.3	14.0
<b>n =</b>	10	10	10

### 9.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to three different samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was measured by the Demeditec Estrone-3-Sulfate equine ELISA. The percentage recoveries were determined by comparing observed and expected results of the samples.

Serum	Spiking (ng/ml)	Observed (ng/ml)	Expected (ng/ml)	Recovery (%)
1	-	0.3	-	-
	50	59.4	50.3	118
	250	265.2	250.3	106
	500	416.7	500.3	83
2	-	0.3	-	-
	50	40.7	50.3	81
	250	224.9	250.3	90
	500	364.7	500.3	73
3	-	43.1	-	-
	50	85.7	93.1	92
	250	287.8	293.1	98
	500	422.7	543.1	78

## 9.5 Linearity

Three serum samples were assayed both undiluted and diluted with Estrone-3-Sulfate Zero Calibrator. The observed and expected values are presented below in ng/ml.

Serum	Dilution	Observed (ng/ml)	Expected (ng/ml)	Linearity (%)
1	Native	201.6	-	-
	1 : 2	109.0	100.8	108
	1 : 4	56.4	50.4	112
	1 : 8	30.5	25.2	121
2	Native	401.0	-	-
	1 : 2	198.1	200.5	99
	1 : 4	101.4	100.3	101
	1 : 8	44.6	50.1	89
3	Native	476.2	-	-
	1 : 2	252.1	238.1	106
	1 : 4	147.5	119.0	124
	1 : 8	75.7	59.5	127

## 10 LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results

### 10.1 Interfering Substances

Do not use hemolytic, icteric or lipemic specimens. Samples containing sodium azide should not be used in the assay.

### 10.2 Drug Interferences

Any medication (cream, oil, pill, etc.) containing hormones could influence the measurement of this analyte. Clinically significant interactions cannot be ruled out with structurally similar drugs.

## 11 LEGAL ASPECTS

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include a sufficient number of controls within the test procedure for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Demeditec.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient therapeutic consequences should be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

## 12 REFERENCES

1. Bouhamidi R, Gaillard JL, Silberzahn P, Martin B  
Binding of Estrogen-3-Sulfates to stallion plasma and equine serum albumine.  
*J Steroid Biochem* 1992, 42: 345-349
2. Cox JE  
Oestrone and equilin in the plasma of pregnant mare  
*J Reprod Fert* 1975, Suppl 23: 463-468
3. Darenius K, Kindahl H, Knudsen O, Madej A, Edqvist LE  
PMSG, progesterone and oestrone sulfate during normal pregnancy and early fetal death.  
Departments of Obstetrics and Gynecology and Clinical Chemistry, Swedish University of  
Agricultural Sciences, Uppsala, Sweden
4. Hearp RB, Hamon M, Allen WR.  
Studies on oestrogen synthesis by the preimplantation equine conceptus  
*J Reprod Fert* 1982: Suppl 32: 343-352
5. Hughes JP, Stabenfeldt GH, Warren EJ  
Estrous cycle and ovulation in the mare  
*JAVMA* 1972 161:1367-1374
6. Hyland JH, Wright PJ, Manning SJ.  
An investigation of the use of plasma oestrone sulfate concentrations for the diagnosis of  
pregnancy in mares.  
*Australian Vet J* 1984, 61: 123
7. Jeffcott LB, Hyland JH, MacLean AA, Dyke T, Robertson-Smith G.  
Changes in maternal hormone concentrations associated with induction of fetal death at day  
45 of gestation in mares.  
*J Reprod Fert* 1987, Suppl 35: 461-467
8. Hyland JH, Langsford DA.  
Changes in urinary and plasma oestrone sulfate concentrations after induction foetal death  
in mares at 45 days of gestation  
*Australian Vet J* 1990, 67: 349-351
9. Kasman LH, Hughes JP, Stabenfeldt GH, Starr MD, Lasley BL  
Estrone sulfate concentrations as an indicator of fetal demise in horses  
*Am J Vet Res* 1988, 49: 184-187
10. Kindahl H, Knudsen O, Madej A, Edqvist LE.  
Progesterone, prostaglandin F-2 $\alpha$ , PMSG and oestrone sulfate during early pregnancy in the mare.  
*J Reprod Fert* 1982, Suppl 32: 353-359
11. Makawiti DW, Allen WE, Kilpatrick MJ  
Changes in oestrone sulfate concentrations in peripheral plasma of pony mares associated with follicular  
growth, ovulation and early pregnancy.  
*J Reprod Fert* 1983, 68: 481-487
12. Muyan M, Roser JF, Dybal N, Baldwin DM.  
Modulation of gonadotropin-releasing hormone-stimulated luteinizing hormone release in cultured male  
equine anterior pituitary cells by gonadal steroids  
*Biol of Reprod* 1993, 49: 340-345
13. Möstl E  
The horse feto-placental unit.  
*Exp Clin Endocrin* 1994, 102: 166-168
14. Ousey JC, Rossdale PD, Cash RSG, Worthy K.  
Plasma concentrations of progestagens, oestrone sulfate and prolactin in pregnant mares  
subjected to natural challenge with equid herpesvirus-1.  
*J Reprod Fert* 1987, Suppl 35: 519-528
15. Parkes RD, Blackmore DJ, Rance DA, Park BK, Dean PDG  
Plasma concentrations of equilin and oestrone in the assessment of fetoplacental function in the mare.  
*Vet Rec* 1977, 100: 511-512
16. Sanchi EM, LeBlanc MM, Weston PG.  
Progestagen, oestrone sulfate and cortisol concentrations in pregnant mares during medical  
and surgical disease  
*J Reprod Fert* 1991, Suppl 44: 627-634
17. Setchell BP, Cox JE.  
Secretion of free and conjugated steroids by the horse testis into lymph and venous blood  
*J Reprod Fert* 1982, Suppl 32: 123-127
18. Terqui M, Palmer E.  
Oestrogen pattern during early pregnancy in the mare  
*J Reprod Fert* 1979, Suppl 27: 441-446

**13 SHORT INSTRUCTION**(all sample sizes given in  $\mu\text{l}$ )

MP Well	ng/ml	0	1	2	3	4	5	6	Sample
		0	5	10	50	100	300	1000	
Steps	Solution								
Pipet	Calibrator	20	20	20	20	20	20	20	-
Pipet	Sample	-	-	-	-	-	-	-	20
Pipet	Biotin-Labeled Estrone-3-Sulfate	100	100	100	100	100	100	100	100
Pipet	Sample Buffer	100	100	100	100	100	100	100	100
Incubate for <b>1 hour</b> at RT (18-25°C) on a shaker (900 rpm) <b>Afterwards do <u>not</u> wash!</b>									
Pipet	Enzyme-labeled streptavidine	25	25	25	25	25	25	25	25
Incubate for <b>30 minutes</b> at RT (18-25°C) on a shaker (900 rpm)									
Decant Wash <b>4x with 300 <math>\mu\text{l}</math></b> of buffered wash solution									
Pipet	Substrate Solution	200	200	200	200	200	200	200	200
Incubate for <b>30 min</b> at RT (18-25°C) <b>in the dark</b>									
Pipet	Stop Solution	50	50	50	50	50	50	50	50
Read at $\lambda = 450 \text{ nm}$									

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
**12 LITERATUR**

1. Bouhamidi R, Gaillard JL, Silberzahn P, Martin B  
Binding of Estrogen-3-Sulfates to stallion plasma and equine serum albumine.  
J Steroid Biochem 1992, 42: 345-349
2. Cox JE  
Oestrone and equilin in the plasma of pregnant mare  
J Reprod Fert 1975, Suppl 23: 463-468
3. Darenius K, Kindahl H, Knudsen O, Madej A, Edqvist LE  
PMSG, progesterone and oestrone sulfate during normal pregnancy and early fetal death.  
Departments of Obstetrics and Gynecology and Clinical Chemistry, Swedish University of  
Agricultural Sciences, Uppsala, Sweden
4. Hearp RB, Hamon M, Allen WR.  
Studies on oestrogen synthesis by the preimplantation equine conceptus  
J Reprod Fert 1982: Suppl 32: 343-352
5. Hughes JP, Stabenfeldt GH, Warren EJ  
Estrous cycle and ovulation in the mare  
JAVMA 1972 161:1367-1374
6. Hyland JH, Wright PJ, Manning SJ.  
An investigation of the use of plasma oestrone sulfate concentrations for the diagnosis of  
pregnancy in mares.  
Australian Vet J 1984, 61: 123
7. Jeffcott LB, Hyland JH, MacLean AA, Dyke T, Robertson-Smith G.  
Changes in maternal hormone concentrations associated with induction of fetal death at day  
45 of gestation in mares.  
J Reprod Fert 1987, Suppl 35: 461-467
8. Hyland JH, Langsford DA.  
Changes in urinary and plasma oestrone sulfate concentrations after induction foetal death  
in mares at 45 days of gestation  
Australian Vet J 1990, 67: 349-351
9. Kasman LH, Hughes JP, Stabenfeldt GH, Starr MD, Lasley BL  
Estrone sulfate concentrations as an indicator of fetal demise in horses  
Am J Vet Res 1988, 49: 184-187
10. Kindahl H, Knudsen O, Madej A, Edqvist LE.  
Progesterone, prostaglandin F-2 $\alpha$ , PMSG and oestrone sulfate during early pregnancy in the mare.  
J Reprod Fert 1982, Suppl 32: 353-359
11. Makawiti DW, Allen WE, Kilpatrick MJ  
Changes in oestrone sulfate concentrations in peripheral plasma of pony mares associated with follicular  
growth, ovulation and early pregnancy.  
J Reprod Fert 1983, 68: 481-487
12. Muyan M, Roser JF, Dybal N, Baldwin DM.  
Modulation of gonadotropin-releasing hormone-stimulated luteinizing hormone release in cultured male  
equine anterior pituitary cells by gonadal steroids  
Biol of Reprod 1993, 49: 340-345
13. Möstl E  
The horse feto-placental unit.  
Exp Clin Endocrin 1994, 102: 166-168
14. Ousey JC, Rossdale PD, Cash RSG, Worthy K.  
Plasma concentrations of progestagens, oestrone sulfate and prolactin in preganant mares  
subjected to natural challenge with equid herpesvirus-1.  
J Reprod Fert 1987, Suppl 35: 519-528
15. Parkes RD, Blackmore DJ, Rance DA, Park BK, Dean PDG  
Plasma concentrations of equilin and oestrone in the assessment of fetoplacental function in the mare.  
Vet Rec 1977, 100: 511-512
16. Sanchi EM, LeBlanc MM, Weston PG.  
Progestagen, oestrone sulfate and cortisol concentrations in pregnant mares during medical  
and surgical disease  
J Reprod Fert 1991, Suppl 44: 627-634
17. Setchell BP, Cox JE.  
Secretion of free and conjugated steroids by the horse testis into lymph and venous blood  
J Reprod Fert 1982, Suppl 32: 123-127
18. Terqui M, Palmer E.  
Oestrogen pattern during early pregnancy in the mare  
J Reprod Fert 1979, Suppl 27: 441-446

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## SYMBOLS USED WITH DEMEDITEC ELISAS

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Temperature de conservation	Temperatura de conservacion	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore