Product information

Information about other products is available at: www.demeditec.com

Estrone-3-Sulfate equine ELISA

REF

DEV9933

96

Demeditec Diagnostics GmbH
Lise-Meitner-Strasse 2
24145 Kiel – Germany
www.demeditec.com

CONTENTS / INHALTSVERZEICHNIS

1	INTRODUCTION	3
2	PRINCIPLE	3
3	WARNINGS AND PRECAUTIONS	3
4	REAGENTS	4
5	SAMPLE COLLECTION AND PREPARATION.....	5
6	ASSAY PROCEDURE	5
7	EXPECTED NORMAL VALUES	7
8	QUALITY CONTROL	7
9	PERFORMANCE CHARACTERISTICS.....	8
10	LIMITATIONS OF PROCEDURE	9
11	LEGAL ASPECTS.....	10
12	REVISION HISTORY OF INSTRUCTION FOR USE.....	10
13	REFERENCES	11
14	SHORT INSTRUCTION.....	12
1	EINLEITUNG	13
2	METHODIK UND TESTPRINZIP.....	13
3	HINWEISE UND VORSICHTSMASSNAHMEN	13
4	KITBESTANDTEILE	14
5	PROBENENTNAHME UND -VORBEREITUNG	15
6	TESTDURCHFÜHRUNG.....	16
7	NORMALWERTE.....	18
8	QUALITÄTSKONTROLLE	18
9	LEISTUNGSSCHARAKTERISTIKA.....	19
10	LIMITATIONEN.....	21
11	RECHTLICHE GRUNDLAGEN	21
12	ÄNDERUNGSHISTORIE DER PACKUNGSBEILAGE	22
13	LITERATUR	22
14	KURZANLEITUNG.....	23

1 INTRODUCTION

1.1 Intended Use

The Demeditec Estrone-3-Sulfate equine ELISA is a competitive enzyme immunoassay for the quantitative determination of estrone-3-sulfate in horse serum and EDTA plasma.

1.2 Description of the analyte

Estrone-3-Sulfate (E3S) is the predominant conjugated estrogen during pregnancy. It is produced by the fetus, possibly in association with the endometrium in the pregnant mare.

Different hormones are important for the complex events that occur during pregnancy in all mammals. In the mare these events include the maintenance of the corpus luteum function, formation of endometrial cups and development of secondary corpora lutea. Progesterone and PMSG (Pregnant Mare Serum Gonadotropine, eCG) and also free estrogens, e.g. estrone, are associated with these processes. It has been shown that estrone is rapidly conjugated after secretion and the ratio between conjugated and unconjugated estrogens is 100:1 in mare serum.

The conjugated estrogens, especially estrone-3-sulfate, provide the opportunity to improve the accuracy of pregnancy diagnosis, monitor the pregnancy and distinguish whether the fetal development is normal or impaired. The diagnosis of embryonic death is usually made by using techniques of palpation of the uterus per rectum or ultrasound echography. The determination of estrone-3-sulfate is an aid in the diagnosis which allows monitoring of the feto-placental unit during pregnancy. Only in mares with normal fetal development the values of estrone-3-sulfate show a tremendous increase between day 75 and 100 of gestation.

2 PRINCIPLE

The Demeditec Estrone-3-Sulfate equine ELISA test kit is a solid phase enzyme immunoassay (ELISA) in the microtiterplate format, designed for the quantitative measurement of Estrone-3-Sulfate. The microplate is coated with a polyclonal antibody specific for Estrone-3-Sulfate.

Calibrators and samples as well as a Biotin-Labeled Estrone-3-Sulfate and Incubation Buffer are pipetted into the wells of an antibody-coated microtiterplate. During one hour incubation the estrone-3-sulfate of the Calibrator/sample competes with the Biotin-Labeled Estrone-3-Sulfate for the binding sites of the antibody fixed on the inner surface of the wells.

Afterwards, all wells are washed and the Enzyme-Labeled Conjugate - horseradish peroxidase conjugated to Streptavidine - is added. During a 30 minutes incubation, this Enzyme-Labeled Streptavidine binds to the biotinylated Estrone-3-Sulfate. All wells are washed again to remove excess of non-bound Enzyme-Labeled Conjugate.

A chromogenic substrate, TMB (3,3',5,5'-Tetramethylbenzidine), is added to all wells. During a 30 minutes incubation, the substrate is converted to a colored end product (blue) by the fixed enzyme. Enzyme reaction is stopped by dispensing of hydrochloric acid as stop solution (change from blue to yellow). The color intensity is inversely related to the concentration of estrone-3-sulfate present in the sample. The optical density of the color solution is measured with a microplate reader at 450 nm.

3 WARNINGS AND PRECAUTIONS

1. This kit is intended for laboratory use. Use by staff, who is specially informed and trained in methods which are carried out by use of immunoassays.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the instructions for use provided with the kit. Be sure that everything is understood.
3. The microtiterplate contains break apart strips. Unused wells must be stored at 2-8°C in the sealed foil pouch and used in the frame provided.
4. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
5. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
6. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
7. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
8. Allow the reagents to reach room temperature (18-25°C) before starting the test. Temperature will affect the absorbance readings of the assay.
9. Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
10. Do not smoke, eat, drink or apply cosmetics in areas where samples or kit reagents are handled.

11. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
12. Do not use reagents beyond expiry date as shown on the kit labels.
13. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
14. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
15. All blood components and biological materials should be handled as potentially hazardous in use and for disposal. Follow universal precautions when handling and disposing of infectious agents.
16. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
17. Avoid contact with Stop Solution. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, CMIT and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
20. For information, please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH or on Demeditec homepage (www.demeditec.com).
21. If product information, including labeling, is incorrect or inaccurate, please contact the kit manufacturer or supplier.

4 REAGENTS

4.1 Reagents Provided

1. **SORB MT Microtiterplate**, 12 x 8 (break-apart) strips with 96 wells, ready to use; coated with a polyclonal anti Estrone-3-Sulfate antibody.
2. **CAL 0-5 LYO Calibrators**, 6 vials, lyophilized, for reconstitution see 4.4 "Reagent preparation". The Calibrators contain 0, 5, 10, 50, 100 and 500 ng/ml, respectively, in equine serum.
3. **BIOTIN CONJ Biotin-Labeled Estrone-3-Sulfate**, 1 vial, 6 ml in a buffered matrix; ready to use, clear; contains a non-mercury preservative.
4. **INC BUF Incubation Buffer**, 1 vial, 6 ml, ready to use, yellow; contains a non-mercury preservative.
5. **ENZ CONJ Enzyme-Labeled Streptavidine**, 1 vial, 11 ml, ready to use, pale red; horseradish peroxidase conjugated to Streptavidine in a buffered matrix; contains a non-mercury preservative.
6. **SUB TMB TMB-Substrate Solution**, 1 vial, 22 ml, ready to use, clear. 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide in a buffered solution; contains 1-3% 2-Pyrrolidon as preservative.
7. **STOP SOLN Stop Solution**, 1 vial, 7 ml, ready to use; contains 2 M hydrochloric acid; clear. Avoid contact with the Stop Solution. It may cause skin irritations and burns.
8. **WASH SOLN 10x Wash Solution**, 1 vial, 50 ml, 10x concentrated; clear; contains a non-mercury preservative; see „Reagent preparation“.

4.2 Materials required but not provided

- Microplate reader capable for endpoint measurements at 450 nm
- Vortex mixer
- Microplate mixer operating at 900 rpm
- Distilled or deionized water
- Graduated cylinders for 500 ml
- Plastic container for storage of the wash solution
- Calibrated variable precision micropipettes, dispensers, multichannel pipettes
- Semi-logarithmic graph paper or software for data reduction
- Absorbant paper and timer

4.3 Storage conditions

When stored at 2-8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. After first opening the reagents are stable for 30 days if used and stored properly. Keep away from heat and direct sunlight.

Microtiter wells must be stored at 2-8°C. Take care that the foil bag is sealed tightly.

Store Calibrators refrigerated, they will be stable at 2-8°C for 7 days after reconstitution. For longer storage aliquot and freeze at <-20°C. Protect TMB-Substrate Solution from light.

4.4 Reagent preparation

Calibrators:

Reconstitute lyophilized Calibrator 0 **CAL 0 LYO** with 2.0 ml dist. water and Calibrator 1 through Calibrator 5 **CAL 1-5 LYO** with 0.5 ml dist. water 30 min. before use. Mix gently.

Wash Solution:

Dilute 50 ml of 10x concentrated Wash Solution **WASH SOLN 10x** with 450 ml dist. water to a final volume of 500 ml. The diluted Wash Solution **WASH SOLN 1x** is stable for at least 12 weeks at room temperature (18-25°C). Precipitates may form when stored at 2-8°C, which should dissolve again by swirling at room temperature (18-25°C).

4.5 Disposal of the kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged test kit

In case of any severe damage of the test kit or components, Demeditec Diagnostics GmbH has to be informed in writing within one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SAMPLE COLLECTION AND PREPARATION

For determination of Estrone-3-Sulfate equine serum and EDTA plasma are the preferred sample matrix. The procedure calls for 20 µl matrix per well. Blood collection should be as stress-free as possible. The samples may be stored refrigerated at 2-8°C for one week, or up to 6 months frozen <-20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Samples expected to contain Estrone-3-Sulfate concentrations higher than the highest Calibrator (500 ng/ml) should be diluted with the Zero Calibrator before assay. The additional dilution step has to be taken into account for the calculation of the results.

Do not use hemolytic, icteric or lipemic specimens. Samples containing sodium azide should not be used in the assay.

6 ASSAY PROCEDURE

6.1 General Remarks

- Do not interchange components of different lots.
- All components of the test kit supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent of foam formation.
- All reagents and specimens must be allowed to come to room temperature (18-25°C) before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each Calibrator, control or sample in order to avoid cross-contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Respect the incubation times as stated in these instructions for use.
- Calibrators, controls and samples should be assayed at least in duplicates.
- Microtiterplate washing is important. Improperly washed wells will give erroneous results. It is recommended to use a multichannel pipette or a multistepper, respectively, or an automatic microtiter plate washing system. Do not allow wells to dry between incubations. Do not scratch coated wells during rinsing and aspiration. Rinse and fill all reagents with care. While rinsing,

check that all wells are filled precisely with Wash Solution, and that there are no residues in the wells.

- A Calibrator curve must be established for every test run.
- For internal quality control we suggest to use **Estrone-3-Sulfate equine Control coded DEV9933C**. For more information, please contact Demeditec.

6.2 Assay Procedure

1. **Preparation of Calibrators:** The Calibrators are supplied lyophilized. At least 30 minutes before use, reconstitute the Calibrator 0 **CAL 0 LYO** with 2.0 ml of distilled water, and each of the remaining Calibrators 1 through 5 **CAL 1- 5 LYO** with 0.5 ml of distilled water. Use volumetric pipettes and mix by gentle inversion.
2. Prepare a sufficient number of microplate wells to accommodate Calibrators, controls and samples in duplicates.

	1	2	3	4	5	6	7	8	9	10	11	12
A	0	4	P3	P..								
B	0	4	P3	P..								
C	1	5	P4									
D	1	5	P4									
E	2	P1	P5									
F	2	P1	P5									
G	3	P2	P6									
H	3	P2	P6									

0-5: Calibrator 0-5, P1-P.. Sample

3. Pipet **20 µl** of **Calibrators** **CAL 0- 5**, **Controls** and **Samples** with new disposable tips into the appropriate wells of the microplate.
4. Add **50 µl** of **Biotin-Labeled Estrone-3-Sulfate** **BIOTIN CONJ** to each well.
5. Add **50 µl** of **Incubation Buffer** **INC BUF** to every well.
6. Incubate for **1 hour** at room temperature (18-25°C) on a plate shaker (900 rpm).
7. Discard the content of the wells and rinse the wells **4 times** with diluted **1x Wash Solution** **WASH SOLN 1x** (300 µl per well). Remove as much **Wash Solution** as possible by beating the microplate on absorbent paper.
8. Add **100 µl** of **Enzyme-Labeled Streptavidine** **ENZ CONJ** to each well.
9. Incubate for **30 minutes** at room temperature (18-25°C) on a plate shaker (900 rpm).
10. Discard the content of the wells and wash each well **4 times** with diluted **1x Wash Solution** **WASH SOLN 1x** (300 µl per well). Remove as much **Wash Solution** as possible by beating the microplate on absorbent paper.
11. Add **200 µl** of **TMB-Substrate Solution** **SUB TMB** to each well.
12. Incubate without shaking for **30 minutes** at room temperature (18-25°C) in the dark.
13. Add **50 µl** of **Stop Solution** **STOP SOLN** to each well and mix carefully.
14. Determine the optical density (OD) of each well at 450 nm. It is recommended to read the wells within 15 minutes.

6.3 Calculation of Results

1. Calculate the average OD values for each set of Calibrators, Controls and Samples.
2. The obtained OD values of the Calibrators (y-axis, linear) are plotted against their corresponding concentrations (x-axis, logarithmic) either on semilogarithmic paper or using an automated method.
3. Using the mean absorbance value for each sample, determine the corresponding concentration from the Calibrator curve.
4. Automated method: The results in the instructions for use have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give different results.
5. The concentration of the samples can be determined directly from this Calibrator curve. Samples with concentrations higher than that of the highest Calibrator have to be further diluted. For the calculation of the concentrations, this dilution factor has to be taken into account.

6.4 Example of Representative Calibrator Curve

OD results of a representative Calibrator curve for the Demeditec Estrone-3-Sulfate equine ELISA are given in the table below. These data are intended for illustration only and should not be used to calculate results from another run.

Calibrator	Optical Density (450nm)
Calibrator 0 (0 ng/ml)	2.963
Calibrator 1 (5 ng/ml)	2.121
Calibrator 2 (10 ng/ml)	1.661
Calibrator 3 (50 ng/ml)	0.677
Calibrator 4 (100 ng/ml)	0.425
Calibrator 5 (500 ng/ml)	0.163

7 EXPECTED NORMAL VALUES

In a reference range study twenty native mare serum samples were analyzed by Demeditec Estrone-3-Sulfate equine ELISA which were classified as non-gestating by LC-MS. The Upper Cut Off for non-gestating mares is found at 5.45 ng/ml Estrone-3-Sulfate.

Because of differences which may exist between laboratories with respect of population, laboratory technique, and selection of reference groups, it is recommended that each laboratory defines their own normal range with regard to their collective. The reference ranges should be regarded as guidelines only.

8 QUALITY CONTROL

Good laboratory practice requires that controls are run with each Calibrator curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels. Corresponding results of the external Estrone-3-Sulfate equine Control Set (product code. DEV9933C) are stated in the instructions for use respectively on the QC datasheet (Certificate of Analysis) included in the Estrone-3-Sulfate equine ELISA.

Employ appropriate statistical methods for analyzing control values and trends. If the results of the assay do not meet the established acceptable ranges of control materials, results should be considered invalid. In this case, check the following technical areas: Pipetting and timing devices, microtiter plate reader, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above-mentioned items without finding any error contact your distributor or Demeditec directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Limit of Blank (LoB)

The Limit of Blank (LoB) was found to be 0.36 ng/ml.

9.2 Specificity

The following materials have been evaluated for cross-reactivity used up to a concentration of 10,000 ng/ml. The percentage indicates cross-reactivity at 50% displacement compared to Estrone-3-Sulfate.

Compound	Cross-Reactivity (%) at 50% Binding
Estrone	0.45%
Estriol	< 0.13%
Corticosterone	< 0.13%
Equilin	< 0.13%
Equilenin	< 0.13%
17-OH-Progesterone	< 0.13%
Androstenedione	< 0.13%
5 α -Dihydrotestosterone	0.41%
Testosterone	< 0.13%
Progesterone	< 0.13%
Androsterone	< 0.13%
17 β -Estradiol	< 0.13%
β -Estradiol-3-Sulfate	0.29%
Pregnenolone	< 0.13%
DHEA-S	< 0.13%
Cortisone	< 0.13%

9.3 Reproducibility

9.3.1 Intra-Assay

The intra-assay variation was determined by 20 replicate measurements of three samples within one run.

	Sample 1	Sample 2	Sample 3
Mean (ng/ml)	224.1	69.2	24.1
SD	6.37	4.00	1.07
CV (%)	2.8	5.8	4.4
n =	20	20	20

9.3.2 Inter-Assay

The inter-assay variation was determined by duplicate measurements of three samples in eleven different runs.

	Sample 1	Sample 2	Sample 3
Mean (ng/ml)	224.4	75.0	24.6
SD	13.97	6.20	0.81
CV (%)	6.2	8.3	3.3
n =	11	11	11

9.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to three different samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was measured by the Demeditec Estrone-3-Sulfate equine ELISA. The percentage recoveries were determined by comparing observed and expected results of the samples.

Sample	Spiking (ng/ml)	Observed (ng/ml)	Expected (ng/ml)	Recovery (%)
1	native	26.9	-	-
	25	56.4	51.9	109
	100	139.7	126.9	110
	200	227.9	226.9	100
2	native	144.8	-	-
	25	193.3	169.8	114
	100	242.6	244.8	99
	200	310.3	344.8	90
3	native	30.8	-	-
	25	59.6	55.8	107
	100	161.1	130.8	123
	200	238.9	230.8	103

9.5 Linearity

Three samples were assayed both undiluted and diluted with Zero Calibrator. The percentage linearities were determined by comparing observed and expected results of the samples.

Sample	Dilution	Observed (ng/ml)	Expected (ng/ml)	Linearity (%)
1	native	263.2	-	-
	1 : 2	137.8	131.6	105
	1 : 4	66.7	65.8	101
	1 : 8	34.9	32.9	106
2	native	462.6	-	-
	1 : 2	254.4	231.3	110
	1 : 4	136.6	115.7	118
	1 : 8	68.8	57.8	119
3	native	260.5	-	-
	1 : 2	138.8	130.2	107
	1 : 4	64.5	65.1	99
	1 : 8	28.2	32.6	87

10 LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results

10.1 Interfering Substances

- Do not use any hemolytic, icteric, or lipemic specimens to avoid any interferences.
- Samples containing sodium azide should not be used.
- There is no significant interference of Biotin on the assay results up to 3000 ng/ml.
- Unspecific interferences with this immunoassay cannot be excluded (19-21). If unplausible results are suspected, they should be considered invalid and verified by further testing.

10.2 Drug Interferences

Any medication (cream, oil, pill, etc.) containing steroid hormones could influence the measurement of this analyte. Clinically significant interactions cannot be ruled out with structurally similar drugs.

10.3 High Dose Hook Effect

Up to a tested concentration of 25,000 ng/ml Estrone-3-Sulfate no High Dose Hook Effect could be observed for the Estrone-3-Sulfate equine ELISA.

11 LEGAL ASPECTS**11.1 Reliability of Results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include a sufficient number of controls within the test procedure for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern regarding results please contact Demeditec Diagnostics GmbH.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the therapeutic consequences should be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REVISION HISTORY OF INSTRUCTION FOR USE

Changes from the previous version 6-06/23 to actual version 7-11/25

Chap. 1	EDTA samples can also be used
Chap. 2	updated
Chap. 4	updated (Calibrators up to 500 ng/ml), updated volumes (4.1, 4.4)
Chap. 5	updated; EDTA samples can also be used. Information about limited use of Biotin due to interferences removed
Chap. 6.2	assay procedure updated
Chap. 7	expected normal values updated
Chap. 9	updated
Chap. 10.1	updated information about interferences with Biotin.
Chap. 14	updated
General	Editorial changes

13 REFERENCES

1. Bouhamidi R, Gaillard JL, Silberzahn P, Martin B: Binding of Estrogen-3-Sulfates to stallion plasma and equine serum albumine. *J Steroid Biochem* 1992, 42: 345-349
2. Cox JE: Oestrone and equilin in the plasma of pregnant mare. *J Reprod Fert* 1975, Suppl 23: 463-468
3. Darenius K, Kindahl H, Knudsen O, Madej A, Edqvist LE: PMSG, progesterone and oestrone sulfate during normal pregnancy and early fetal death. Departments of Obstetrics and Gynecology and Clinical Chemistry, Swedish University of Agricultural Sciences, Uppsala, Sweden
4. Hearp RB, Hamon M, Allen WR: Studies on oestrogen synthesis by the preimplantation equine conceptus. *J Reprod Fert* 1982: Suppl 32: 343-352
5. Hughes JP, Stabenfeldt GH, Warren EJ: Estrous cycle and ovulation in the mare. *JAVMA* 1972 161:1367-1374
6. Hyland JH, Wright PJ, Manning SJ: An investigation of the use of plasma oestrone sulfate concentrations for the diagnosis of pregnancy in mares. *Australian Vet J* 1984, 61: 123
7. Jeffcott LB, Hyland JH, MacLean AA, Dyke T, Robertson-Smith G: Changes in maternal hormone concentrations associated with induction of fetal death at day 45 of gestation in mares. *J Reprod Fert* 1987, Suppl 35: 461-467
8. Hyland JH, Langsford DA: Changes in urinary and plasma oestrone sulfate concentrations after induction foetal death in mares at 45 days of gestation. *Australian Vet J* 1990, 67: 349-351
9. Kasman LH, Hughes JP, Stabenfeldt GH, Starr MD, Lasley BL: Estrone sulfate concentrations as an indicator of fetal demise in horses. *Am J Vet Res* 1988, 49: 184-187
10. Kindahl H, Knudsen O, Madej A, Edqvist LE: Progesterone, prostaglandin F-2 α , PMSG and oestrone sulfate during early pregnancy in the mare. *J Reprod Fert* 1982, Suppl 32: 353-359
11. Makawiti DW, Allen WE, Kilpatrick MJ: Changes in oestrone sulfate concentrations in peripheral plasma of pony mares associated with follicular growth, ovulation and early pregnancy. *J Reprod Fert* 1983, 68: 481-487
12. Muyan M, Roser JF, Dybal N, Baldwin DM: Modulation of gonadotropin-releasing hormone-stimulated luteinizing hormone release in cultured male equine anterior pituitary cells by gonadal steroids. *Biol of Reprod* 1993, 49: 340-345
13. Möstl E: The horse feto-placental unit. *Exp Clin Endocrin* 1994, 102: 166-168
14. Ousey JC, Rossdale PD, Cash RSG, Worthy K: Plasma concentrations of progestagens, oestrone sulfate and prolactin in pregnant mares subjected to natural challenge with equid herpesvirus-1. *J Reprod Fert* 1987, Suppl 35: 519-528
15. Parkes RD, Blackmore DJ, Rance DA, Park BK, Dean PDG: Plasma concentrations of equilin and oestrone in the assessment of fetoplacental function in the mare. *Vet Rec* 1977, 100: 511-512
16. Sanchi EM, LeBlanc MM, Weston PG: Progestagen, oestrone sulfate and cortisol concentrations in pregnant mares during medical and surgical disease. *J Reprod Fert* 1991, Suppl 44: 627-634
17. Setchell BP, Cox JE: Secretion of free and conjugated steroids by the horse testis into lymph and venous blood. *J Reprod Fert* 1982, Suppl 32: 123-127
18. Terqui M, Palmer E: Oestrogen pattern during early pregnancy in the mare. *J Reprod Fert* 1979, Suppl 27: 441-446
19. Marks V.: False-Positive Immunoassay Results: A Multicenter Survey of Erroneous Immunoassay Results from Assays of 74 Analytes in 10 Donors from 66 Laboratories in Seven Countries. *Clinical Chemistry* 2002, 48:11: 2008-2016
20. Tate J, Ward G: Interferences in Immunoassays. *Clin. Biochem Rev* 2004, 25(2) : 105-120
21. Selby C: Interference in immunoassays. *Ann. Clin. Biochem* 1999, 36: 704-721

14 SHORT INSTRUCTION

all steps at room temperature (18-25°C)

all sample sizes given in µl

Steps	CAL 0-5 (0 - 500 ng/ml)	Sample
pipet CAL 0-5	20	-
pipet Samples	-	20
pipet BIOTIN CONJ	50	
pipet INC BUF	50	
incubate for 1 h on a shaker (900 rpm)		
decant wash 4x with WASH SOLN 1x	300	
pipet ENZ CONJ	100	
incubate for 30 min on a shaker (900 rpm)		
decant wash 4x with WASH SOLN 1x	300	
pipet SUB TMB	200	
incubate without shaking for 30 min in the dark		
pipet STOP SOLN	50	
read at $\lambda = 450 \text{ nm}$		

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13 LITERATUR

- 1 Bouhamidi R, Gaillard JL, Silberzahn P, Martin B: Binding of Estrogen-3-Sulfates to stallion plasma and equine serum albumine. J Steroid Biochem 1992, 42: 345-349
- 2 Cox JE: Oestrone and equilin in the plasma of pregnant mare. J Reprod Fert 1975, Suppl 23: 463-468
- 3 Darenius K, Kindahl H, Knudsen O, Madej A, Edqvist LE: PMSG, progesterone and oestrone sulfate during normal pregnancy and early fetal death. Departments of Obstetrics and Gynecology and Clinical Chemistry, Swedish University of Agricultural Sciences, Uppsala, Sweden
- 4 Hearp RB, Hamon M, Allen WR: Studies on oestrogen synthesis by the preimplantation equine conceptus. J Reprod Fert 1982: Suppl 32: 343-352
- 5 Hughes JP, Stabenfeldt GH, Warren EJ: Estrous cycle and ovulation in the mare. JAVMA 1972 161:1367-1374
- 6 Hyland JH, Wright PJ, Manning SJ: An investigation of the use of plasma oestrone sulfate concentrations for the diagnosis of pregnancy in mares. Australian Vet J 1984, 61: 123
- 7 Jeffcott LB, Hyland JH, MacLean AA, Dyke T, Robertson-Smith G: Changes in maternal hormone concentrations associated with induction of fetal death at day 45 of gestation in mares. J Reprod Fert 1987, Suppl 35: 461-467
- 8 Hyland JH, Langsford DA: Changes in urinary and plasma oestrone sulfate concentrations after induction foetal death in mares at 45 days of gestation. Australian Vet J 1990, 67: 349-351
- 9 Kasman LH, Hughes JP, Stabenfeldt GH, Starr MD, Lasley BL: Estrone sulfate concentrations as an indicator of fetal demise in horses. Am J Vet Res 1988, 49: 184-187
- 10 Kindahl H, Knudsen O, Madej A, Edqvist LE: Progesterone, prostaglandin F-2 α , PMSG and oestrone sulfate during early pregnancy in the mare. J Reprod Fert 1982, Suppl 32: 353-359
- 11 Makawiti DW, Allen WE, Kilpatrick MJ: Changes in oestrone sulfate concentrations in peripheral plasma of pony mares associated with follicular growth, ovulation and early pregnancy. J Reprod Fert 1983, 68: 481-487
- 12 Muyan M, Roser JF, Dybal N, Baldwin DM: Modulation of gonadotropin-releasing hormone-stimulated luteinizing hormone release in cultured male equine anterior pituitary cells by gonadal steroids. Biol of Reprod 1993, 49: 340-345
- 13 Möstl E: The horse feto-placental unit. Exp Clin Endocrin 1994, 102: 166-168
- 14 Ousey JC, Rosedale PD, Cash RSG, Worthy K: Plasma concentrations of progestagens, oestrone sulfate and prolactin in pregnant mares subjected to natural challenge with equid herpesvirus-1. J Reprod Fert 1987, Suppl 35: 519-528
- 15 Parkes RD, Blackmore DJ, Rance DA, Park BK, Dean PDG: Plasma concentrations of equilin and oestrone in the assessment of fetoplacental function in the mare. Vet Rec 1977, 100: 511-512
- 16 Sanchi EM, LeBlanc MM, Weston PG: Progestagen, oestrone sulfate and cortisol concentrations in pregnant mares during medical and surgical disease. J Reprod Fert 1991, Suppl 44: 627-634
- 17 Setchell BP, Cox JE: Secretion of free and conjugated steroids by the horse testis into lymph and venous blood. J Reprod Fert 1982, Suppl 32: 123-127
- 18 Terqui M, Palmer E: Oestrogen pattern during early pregnancy in the mare. J Reprod Fert 1979, Suppl 27: 441-446
- 19 Marks V.: False-Positive Immunoassay Results: A Multicenter Survey of Erroneous Immunoassay Results from Assays of 74 Analytes in 10 Donors from 66 Laboratories in Seven Countries. Clinical Chemistry 2002, 48:11: 2008-2016
- 20 Tate J, Ward G: Interferences in Immunoassays. Clin. Biochem Rev 2004, 25(2) : 105-120
- 21 Selby C: Interference in immunoassays. Ann. Clin. Biochem 1999, 36: 704-721

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SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Française	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	utilisation Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di catalogo
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta y advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributed by	Vertrieb durch	Distribution par	Distribución por	Distribuzione da parte di
	Version	Version	Version	Versión	Versione
	Single-use	Einmalverwendung	À usage unique	Uso único	Uso una volta