FSH IRMA KIT

REF IM2125, IM3301

TABLE OF CONTENTS

English	2	Čeština	17
Français	5	한국어	20
Deutsch	8	APPENDIX	22
Italiano	11		
Español	14		

 IMMUNOTECH s.r.o., Radiova 1, 102 27 Prague 10, Czech Republic

FSH IRMA KIT

REF IM2125, IM3301

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF FSH IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

PRINCIPLE

The immunoradiometric assay of FSH is a sandwich-type assay. In the kit, mouse monoclonal antibodies directed against two different epitopes of FSH and hence not competing are used. Samples or calibrators are incubated in tubes coated with the first monoclonal antibody in the presence of the second monoclonal antibody labeled with iodine 125. After incubation, the contents of the tubes are rinsed so as to remove unbound ¹²⁵I-labeled antibody. The bound radioactivity is then determined in a gamma counter. The FSH concentrations in the samples are obtained by interpolation from the standard curve. The concentration of FSH in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Tracer

DANGER

H360

May damage fertility or the unborn child.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned: Get medical advice/attention.

P405

Store locked up.

P501

Dispose of contents/container in accordance with local/national regulations
Boric Acid 1 - 2%

Wash Solution (20X)

DANGER

H360

May damage fertility or the unborn child.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned: Get medical advice/attention.
Boric Acid 0.1-0.3%
Sodium Borate Decahydrate 0.1-0.3%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in tubes containing no additives or EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -20°C, 1 year maximum), after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted in the zero calibrator.

Serum and EDTA plasma values for 18 samples (serum values ranging from 0.80 to 87.01 IU/L) were compared using the IM2125 FSH IRMA kit. Results are as follows:

[EDTA-plasma] = 1.0155[serum] - 0.1293

R = 0.9996

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Kit for the determination of FSH: 100 tubes (Cat. # IM2125)

Anti-FSH monoclonal antibody coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal anti-FSH antibody: one 5.5 mL vial (ready-to-use)

The vial contains 370 kBq, at the date of manufacture, of ¹²⁵I-labeled immunoglobulins in buffer with proteins, sodium azide (<0.1%) and a dye.

Calibrators: six 1 mL vials (ready-to-use)

The calibrator vials contain from 0 to approximately 180 IU/L of human FSH in human serum. The exact concentration is indicated on each vial. The calibrators were calibrated against the international standard WHO 94/632.

Control sera: two vials (lyophilised)

The vials contain human FSH lyophilized in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

Kit for determination of FSH: 400 tubes (Cat. # IM3301)**Anti-FSH monoclonal antibody coated tubes: 8 x 50 tubes** (ready-to-use)**¹²⁵I-labeled monoclonal antibody: four 5.5 mL vials** (ready-to-use)**Calibrators: six 1 mL vials** (ready-to-use)**Control sera: two vials** (lyophilised)**Wash solution (20X): two 50 mL vials****MATERIALS REQUIRED, BUT NOT PROVIDED**

In addition to standard laboratory equipment, the following items are required:

- Precision micropipet (100 µL).
- semi-automatic pipets (50 µL, 2 mL).
- Vortex type mixer.
- Horizontal or orbital shaker.
- Aspiration system.
- Gamma counter set for ¹²⁵I.

RESULTS

Results are obtained from the standard curve by interpolation. The standard curve serves for the determination of FSH concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{1cal} - cpm_{1cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (IU/L).

Other data reduction methods may give slightly different results.

Total activity: 127,770 cpm				
Calibrators	FSH (IU/L)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	130	0.10	-
1	0.53	324	0.25	194
2	2.60	1,179	0.92	1,049
3	15.7	6,330	4.95	6,200
4	61.0	22,118	17.3	21,988
5	180	52,660	41.2	52,530

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

The concentrations of diluted samples must be corrected by the dilution factor.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. FSH concentrations (IU/L) measured in serum samples from healthy subjects in the following table are indicative only.

Group	N	Min.	Max.	Median	2.5 th percentile	97.5 th percentile
Males	98	0.51	10.5	4.03	1.03	9.72
Females						
Follicular phase	351	2.27	29.7	6.07	2.89	14.5
Periovulatory phase	70	2.37	24.5	12.1	4.09	23.7
Luteal phase	309	1.00	21.6	4.15	1.58	14.2
Postmenopausal	48	24.0	141	68.14	28.13	130

FSH concentrations in sera from prepubertal children were found below 1IU/L.

(For more details, see APPENDIX)

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE**Preparation of reagents**

Let all the reagents come to room temperature.

Reconstitution of control sera

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 30 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one week or aliquoted at < -18°C for a longer time, until the expiry date of the kit.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution can be stored at 2-8°C until the expiry date of the kit.

Assay procedure

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes successively add: 100 µL of calibrator, control or sample 50 µL of tracer. Mix	Incubate 90 min at 18-25°C with shaking (>350 rpm) or overnight (18 hours) at 18-25°C without shaking.	Aspirate carefully the content of tubes (except of the 2 tubes «total cpm») Wash twice with 2 mL of wash solution Count bound cpm (B) and total cpm (T) for 1 minute.

*Add 50 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.17 IU/L

Functional sensitivity: 0.5 IU/L

Specificity

The antibody used in the immunoassay is highly specific for FSH. Extremely low cross reactivities were obtained with several related molecules (hLH, hTSH, hCG).

Precision**Intra-assay**

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 4.05%.

Inter-assay

Serum samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 8.20%.

Accuracy**Dilution test**

High-concentration serum samples were serially diluted with the zero calibrator. The recovery percentages were obtained between 81.6% and 117% for serum samples.

Recovery test

Low-concentration samples were spiked with known quantities of FSH. The recovery percentages were obtained between 90.75% and 103% for serum samples.

Measurement range (from analytical sensitivity to highest calibrator): 0.17 to approximately 180 IU/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use lipemic, icteric or haemolyzed sera in assay.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Specificity

The specificity of the assay is determined by measuring the "apparent FSH" response given by high concentrations of related hormones.

Cross-reactivity and interference

Hormone	Concentration	Apparent FSH value (IU/L)
hTSH	50 mIU/L	1.02 ± 0.82
hLH	50 IU/L	0.1 ± 0.84
hLH	1,260 IU/L	0.6 ± 0.03
hCG	900 IU/L	0 ± 0.031
hCG	4,500 IU/L	
hFSH + hCG	25 IU/L + 4,500 IU/L	26.9 ± 1.29

Accuracy

Dilution test

Serum samples were diluted in the zero calibrator and assayed in duplicate.

Serum	Dilution	IU/L		Ratio (%) Measured/ Expected
		Measured	Expected	
S1	-	19.01	-	-
	1:2	9.10	9.51	95.74
	1:4	4.10	4.75	86.27
	1:8	1.94	2.38	81.64
	1:16	0.98	1.19	82.48
	1:32	0.49	0.59	82.48
S2	-	22.79	-	-
	1:2	12.10	11.40	106.2
	1:4	6.20	5.70	108.8
	1:8	2.83	2.85	99.34
	1:16	1.40	1.42	98.29
	1:32	0.72	0.71	101.1
S3	-	24.20	-	-
	1:2	11.99	12.10	99.09
	1:4	5.59	6.05	92.40
	1:8	2.59	3.03	85.62
	1:16	1.25	1.51	82.64
	1:32	0.65	0.76	85.95
S4	-	49.80	-	-
	1:2	28.24	24.90	113.4
	1:4	14.53	12.45	116.7
	1:8	7.27	6.23	116.8
	1:16	3.42	3.11	109.9
	1:32	1.60	1.56	102.8
S5	-	91.19	-	-
	1:2	46.42	45.60	101.8
	1:4	23.64	22.80	103.7
	1:8	11.50	11.40	100.9
	1:16	5.39	5.70	94.57
	1:32	2.34	2.85	82.11

EDTA-plasma	Dilution	IU/L		Ratio (%) Measured/ Expected
		Measured	Expected	
P1	-	57.11	-	-
	1:2	30.26	28.56	106.0
	1:4	15.51	14.28	108.6
	1:8	7.56	7.14	105.9
	1:16	3.57	3.57	100.0
	1:32	1.69	1.78	94.69
P2	-	34.33	-	-
	1:2	18.40	17.17	107.2
	1:4	9.20	8.58	107.2
	1:8	4.27	4.29	99.50
	1:16	1.99	2.15	92.75
	1:32	0.97	1.07	90.42
P3	-	94.39	-	-
	1:2	46.05	47.20	97.57
	1:4	23.27	23.60	98.61
	1:8	11.62	11.80	98.49
	1:16	5.51	5.90	93.40
	1:32	2.46	2.95	83.40
P4	-	21.66	-	-
	1:2	11.08	10.83	102.3
	1:4	5.21	5.42	96.21
	1:8	2.41	2.71	89.01
	1:16	1.22	1.35	90.12
	1:32	0.56	0.68	82.73
P5	-	16.20	-	-
	1:2	7.88	8.10	97.28
	1:4	3.73	4.05	92.10
	1:8	1.74	2.03	85.93
	1:16	1.05	1.01	103.7
	1:32	0.42	0.51	82.96

Recovery test

Samples were spiked with known quantities of FSH and assayed.

Serum	Endogen. conc. IU/L	Added conc. IU/L	Expected conc. IU/L	Measured conc. IU/L	Ratio (%) Measured/Expected
S1	3.34	1.00	4.34	4.22	97.15
	3.30	2.97	6.28	5.84	93.02
	3.25	5.78	9.04	8.20	90.75
S2	10.24	3.60	13.84	14.15	102.2
	9.85	10.32	20.17	19.08	94.61
	9.50	16.36	25.86	23.98	92.72
S3	8.24	2.76	11.01	10.77	97.85
	8.00	7.98	15.98	15.05	94.18
	7.66	15.29	22.94	21.00	91.52
S4	8.24	1.71	6.60	6.75	102.2
	8.00	4.97	9.77	9.12	93.30
	7.66	9.67	14.35	13.21	92.06
S5	6.06	2.06	8.12	8.36	102.9
	5.92	6.05	11.98	11.24	93.85
	5.73	11.78	17.50	16.38	93.58

EDTA-plasma	Endogen. conc. (IU/L)	Added conc. (IU/L)	Expected conc. (IU/L)	Measured conc. (IU/L)	Ratio (%) Measured/Expected
P1	6.15	2.20	8.36	8.12	97.15
	6.01	6.39	12.40	11.71	94.45
	5.80	12.34	18.14	16.80	92.60
P2	5.44	1.71	7.15	7.12	99.59
	5.34	4.97	10.31	10.10	97.97
	5.19	9.67	14.87	13.82	92.94
P3	3.35	1.07	4.42	4.62	104.4
	3.31	3.18	6.49	6.33	97.49
	3.25	6.32	9.57	8.85	92.44
P4	5.72	1.85	7.57	7.97	105.2
	5.60	5.51	11.12	10.70	96.25
	5.43	10.76	16.20	15.38	94.95
P5	1.93	0.72	2.65	2.99	112.9
	1.92	2.13	4.05	4.05	99.97
	1.89	4.29	6.18	5.59	90.44

Precision

Intra-assay

Sample	Serum			EDTA-Plasma		
	1	2	3	1	2	3
Number of determinations	25	25	25	25	25	25
Mean value (IU/L)	2.26	10.67	48.79	6.62	46.44	3.58
C.V., %	4.05	2.62	1.12	2.97	1.52	4.53

Inter-assays

Sample	Serum			EDTA-Plasma		
	1	2	3	1	2	3
Number of determinations	10	10	10	10	10	10
Mean value (IU/L)	0.67	10.65	91.80	3.38	16.69	96.29
C.V., %	8.20	3.72	2.06	5.75	4.14	2.83

Expected values

Values of FSH in fertile females, detailed discrimination by menstrual cycle sub-phases.

Group	N	Min.	Max.	Median	2.5 th	97.5 th
					percentile	percentile
IU/L						
Follicular phase						
- Early	107	3.30	29.7	7.47	3.84	15.9
- Middle	141	2.64	20.6	6.11	3.22	13.9
- Late	103	2.27	15.6	4.85	2.40	8.71
Preovulatory peak	70	2.37	24.5	12.1	4.09	23.7
Luteal phase						
- Early	109	1.86	21.6	6.68	3.05	16.3
- Middle	163	1.02	15.1	3.57	1.57	9.28
- Late	37	1.00	13.5	2.95	1.29	12.3

¹²⁵I Characteristics

$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$

¹²⁵ I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferenz / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirație / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktiviojvi medžiaga / Radioaktiv / Radioaktivny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLÍKE / ОПАСНО / ОПАСНОСТ / 危險

Ag^{125I}

Ab^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmvavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

SOLN | **WASH** | **20x**

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösningkoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalo koncentratas 20X / 20X mosóoldat-konzentrátum / Koncentrat 20X roztworu płuczacego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X