# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

User's Manual

# IFN-gamma human ELISA

Enzyme immunoassay for the quantitative measurement of human interferon gamma (IFN-gamma) in serum and plasma

**DE4434**

**96 wells**

## 1. INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human interferon gamma (IFN-gamma) in serum and plasma.

## 2. CLINICAL BACKGROUND

### A. Biological activities

IFN- $\gamma$  (type 2, immune IFN) is structurally and functionally distinct from type 1 (alpha/beta) interferon and acts on a separate receptor. Only one IFN- $\gamma$  gene has been identified, coding for a 146 AA protein that is post-translationally processed into two glycosylated species of 20 and 25 Kd. Native IFN- $\gamma$  is pH2-labile, highly basic, and can aggregate to form dimers that are biologically active. IFN- $\gamma$  is a real lymphokine produced by activated T (and NK) cells. Despite its clear antiviral and cellular growth regulating activities, its immunomodulatory properties are believed to be the most important. IFN- $\gamma$  is the principal activator of macrophage function (Macrophage Activating Factor, MAF), and it also regulates the pathway of differentiation of myeloid cells. It plays an important role in the growth and differentiation of cytotoxic (and possibly suppressor) T cells, activates NK cells and acts as a B cell maturation factor. It regulates Ig isotype production and inhibits IgE responses. One of the modes of action of IFN- $\gamma$  is to induce the expression of membrane proteins, such as class 1 and class 2 MHC antigens and adhesion molecules on various cell types, high affinity Fc receptors for IgG on myelomonocytic cells, etc. Integrated in the cytokine network, IFN- $\gamma$  interacts with other cytokines, in either a synergistic (e.g. TNF) or antagonistic (e.g. IL-4) way.

### B. Clinical application

The precise role of IFN- $\gamma$  in human diseases and therapy is still poorly defined. Clearly, it is involved in the defense against parasites, intercellular pathogens and possibly tumour cells. Its therapeutic administration partially corrects the deficient immune response observed in lepromatous leprosy and the phagocyte defect of patients with X-Linked chronic granulomatous disease. A deficiency in IFN- $\gamma$  production has been related to persistent (e.g. EBV) viral infections, and a correlation could be established between the secretion of IFN- $\gamma$  by peripheral blood mononuclear cells during an herpetic infection and the time of a next recurrence. A defect in IFN- $\gamma$  production has also been recorded in several primary or secondary immunodeficiency states. IFN- $\gamma$  is seldom detected in the serum of healthy persons. Its production may be demonstrated "in situ" in several inflammatory disorders (Sarcoidosis, rheumatoid-arthritis, subacute thyroiditis, polymyositis, multiple sclerosis). Higher levels of serum IFN- $\gamma$  are measured during severe parasitic diseases (e.g. Plasmodium Falciparum malaria); during cytokine (IL-2) therapy; and after the first injections of OKT3.

## 3. PRINCIPLES OF THE METHOD

The IFN- $\gamma$ -ELISA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplate. The assay uses monoclonal antibodies (MAbs) directed against distinct epitopes of IFN- $\gamma$ . Calibrators and samples react with the capture monoclonal antibody (MAB 1) coated on microtiter well and with a monoclonal antibody (MAB 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAB 1 – human IFN- $\gamma$  – MAB 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the IFN- $\gamma$  concentration.

A calibration curve is plotted and IFN- $\gamma$  concentration in samples is determined by interpolation from the calibration curve. The use of the ELISA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

**4. REAGENTS PROVIDED**

Reagents	96 tests Kit	Reconstitution
<b>SORB MT</b> Microtiterplate with 96 anti IFN- $\gamma$ (monoclonal antibodies) coated wells	96 wells	<b>Ready</b> for use
<b>ENZ CONJ</b> Conjugate: HRP labelled anti-IFN- $\gamma$ (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin and thymol	1 vial 6 ml	<b>Ready</b> for use
<b>CAL 0</b> Zero calibrator in human serum, benzamidin and thymol	1 vial lyophil.	<b>Add</b> distilled water (see on the QC data sheet for the exact volume)
<b>CAL 1 - 5</b> Calibrator N = 1 to 5 (see exact values on QC data sheet) in human serum, benzamidin and thymol	5 vials lyophil.	<b>Add</b> 0.5 ml distilled water
<b>WASH SOLN 200x</b> Wash Solution (Tris-HCl)	1 vial 10 ml	<b>Dilute</b> 200x with distilled water (use a magnetic stirrer).
<b>CONTROL 1 &amp; 2</b> Controls - N = 1 or 2 in human serum, benzamidin and thymol	2 vials lyophil.	<b>Add</b> 0.5 ml distilled water
<b>SUB TMB</b> Chromogenic TMB (Tetramethylbenzidine)	1 vial 12 ml	<b>Ready</b> for use
<b>STOP SOLN</b> Stop Solution: HCl 1.0N	1 vial 12 ml	<b>Ready</b> for use

**Note:** 1. Use the zero calibrator for sample dilutions.  
2. 1 IU of the calibrator preparation is equivalent to 1 IU of the NIBSC Reference Reagent 87/586.

**5. SUPPLIES NOT PROVIDED**

The following material is required but not provided in the kit:

1. High quality distilled water
2. Pipettes for delivery of: 50  $\mu$ l, 200  $\mu$ l, 1 ml and 10 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Horizontal microtiterplate shaker capable of 700 rpm  $\pm$  100 rpm
6. Washer for Microtiterplates
7. Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)
8. Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph 10.A.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

**6. REAGENT PREPARATION**

- A. Calibrators:** Reconstitute the zero calibrator to the volume specified on the QC data sheet with distilled water and the other calibrators with 0.5 ml distilled water.
- B. Controls:** Reconstitute the controls with 0.5 ml distilled water.
- C. Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

## 7. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators and controls are stable for 4 days at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 2 months. Avoid successive freeze thaw cycles.
- The concentrated Wash Solution is stable at 18-28°C until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

## 8. SPECIMEN COLLECTION AND PREPARATION

- Serum must be removed as soon as possible from the clot of red cells after clotting and centrifugation, and kept at 4°C. If the samples are not used immediately, they must be kept at -20°C for maximum 2 months, and at -70°C for longer storage (maximum one year).
- Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at 18-28°C. It is recommended to vortex the samples before use.
- Sampling conditions can affect values, therefore, strict precautions have to be taken during sampling to avoid impurities contained in sampling materials that would stimulate IFN- $\gamma$  production by blood cells and thus falsely increase plasma IFN- $\gamma$  values.
- Collection tubes must be pyrogen-free. Plasma can be collected on sterile EDTA and rapidly separated after centrifugation. The use of heparin tubes is discouraged as batches of heparin are often contaminated with pyrogen.

## 9. PROCEDURE

### A. Handling notes

- Do not use the kit or components beyond expiry date.
- Do not mix materials from different kit lots.
- Bring all the reagents to 18-28°C prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
- Use a clean plastic container to prepare the Wash Solution.
- In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
- For the dispensing of the Revelation Solution and the Stop Solution avoid pipettes with metal parts.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section 12 paragraph E (Time delay).
- Prepare a calibration curve for each run, do not use data from previous runs.
- The Revelation Solution should be colourless. If a blue colour develops within a few minutes after preparation, this indicates that the reagent is unusable, and must be discarded.
- Dispense the Revelation Solution within 15 minutes following the washing of the microtiterplate.
- During incubation with Revelation Solution, avoid direct sunlight on the microtiterplate.

**B. Procedure**

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the strips into the holding frame.
3. Pipette 50 µl of each Calibrator, Control and Sample into the appropriate wells.
4. Pipette 50 µl of anti-IFN-γ-HRP conjugate into all the wells.
5. Incubate for 2 hours at 18-28°C on a horizontal shaker set at 700 rpm ± 100 rpm.
6. Aspirate the liquid from each well.
7. Wash the plate 3 times by:
  - Dispensing 0.4 ml of Wash Solution into each well
  - Aspirating the content of each well
8. Pipette 100 µl of the chromogenic Solution into each well within 15 minutes following the washing step.
9. Incubate the microtiterplate for 15 minutes at 18-28°C on a horizontal shaker set at 700 rpm ± 100 rpm, avoid direct sunlight.
10. Pipette 100 µl of Stop Solution into each well.
11. Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 30 minutes and calculate the results as described in section 10.

**10. CALCULATION OF RESULTS****A. Polychromatic Reading:**

1. In this case, the ELISA-AID™ software will do the data processing.
2. The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
3. A second reading is performed at 490 nm against the same reference filter.
4. The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
5. The principle of polychromatic data processing is as follows:
  - $X_i = \text{OD at 450 nm}$
  - $Y_i = \text{OD at 490 nm}$
  - Using a standard unweighted linear regression, the parameters A & B are calculated:  
 $Y = A \cdot X - B$
  - If  $X_i < 3$  OD units, then X calculated =  $X_i$
  - If  $X_i > 3$  OD units, then X calculated =  $(Y_i - B)/A$
  - A 4-parameter logistic curve fitting is used to build up the calibration curve.
  - The IFN-γ concentration in samples is determined by interpolation on the calibration curve.

**B. Bichromatic Reading**

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of IFN-γ (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

**11. TYPICAL DATA**

The following data are for illustration only and should never be used instead of the real time calibration curve.

	IFN-γ-ELISA	OD units Polychromatic model
Calibrator	0 IU/ml	0.03
	1 IU/ml	0.173
	2 IU/ml	0.339
	5 IU/ml	0.700
	10 IU/ml	1.353
	30 IU/ml	3.107

**12. PERFORMANCE AND LIMITATIONS****A. Detection Limit**

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.03 IU/ml.

**B. Specificity**

No significant cross-reaction was observed in presence of 50 ng of IL-1 $\alpha$ , IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-6, IL-7, IL-8, IL-10, TNF- $\alpha$ , TNF- $\beta$ , IFN- $\beta$ , TGF- $\beta$ , GM-CSF, OSM, MIP-1 $\alpha$ , MIP-1 $\beta$ , LIF, MCP-1, G-CSF and RANTES. This IFN- $\gamma$  assay is specific for human natural and recombinant IFN- $\gamma$ .

**C. Precision**

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> $\pm$ SD (IU/ml)	CV (%)	Serum	N	<X> $\pm$ SD (IU/ml)	CV (%)
A	10	1.26 $\pm$ 0.03	3.2	A	20	1.61 $\pm$ 0.03	5.8
B	10	12.28 $\pm$ 0.35	3.8	B	20	5.72 $\pm$ 0.51	8.8

SD: Standard Deviation; CV: Coefficient of variation

**D. Accuracy**

## RECOVERY TEST

Sample	Added IFN- $\gamma$ (IU/ml)	Recovered IFN- $\gamma$ (IU/ml)	Recovery (%)
Serum	20.5	20.3	99
	9.9	10.1	103
	4.7	4.8	101
	2.4	2.3	94
	1.0	1.0	99
Plasma	20.5	19.6	96
	9.9	9.5	97
	4.7	4.8	102
	2.4	2.4	101
	1.0	1.0	96

## DILUTION TEST

Sample	Dilution	Theoretical Conc. (IU/ml)	Measured Conc. (IU/ml)
Serum	1/1	-	18.0
	1/2	9.0	8.1
	1/4	4.5	4.4
	1/8	2.3	2.4
	1/16	1.1	1.3
	1/32	0.6	0.6
Plasma	1/1	-	11.7
	1/2	5.9	5.9
	1/4	2.9	3.2
	1/8	1.5	1.6
	1/16	0.7	0.8
	1/32	0.4	0.4

Samples were diluted with zero calibrator.

**E. Time delay between last calibrator and sample dispensing**

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

	TIME DELAY			
	T0	10 min	20 min	30 min
S1	2.7	2.7	2.7	2.7
S2	6.8	6.8	6.5	6.7
S3	4.2	4.1	4.0	4.0
N1	22.7	21.2	19.4	22.0
N2	24.9	23.7	21.7	20.1
N3	14.7	14.8	13.1	12.6
N4	17.9	15.3	13.8	13.9
N5	17.4	16.5	15.6	15.0

**F. Hook effect**

A sample spiked with IFN- $\gamma$  up to 500 000 IU/ml gives higher OD's than the last calibrator point.

**13. INTERNAL QUALITY CONTROL**

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the QC data sheet, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

**14. REFERENCE INTERVALS**

These values are given only for guidance; each laboratory should establish its own normal range of values.

For guidance, the mean of 30 normal serum samples was 0.28 IU/ml (SD = 0.15), ranging between 0 IU/ml and 0.77 IU/ml. This study was performed on samples from apparently healthy persons with low CRP levels.

The mean of 60 normal plasma, collected in strict sampling conditions, was 0.08 IU/ml (SD = 0.12), ranging between 0 IU/ml and 0.89 IU/ml).

**15. PRECAUTIONS AND WARNINGS****Safety**

For in vitro diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious. Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

**16. BIBLIOGRAPHY**

1. GUERRY D., ALEXANDER M.A., ELDER D.E. and HERLYN M.F. (1987) **IFN- $\gamma$  regulated the T cell response to precursor nevi and biologically early melanoma.** J. Immunol., 139 : 305-312.
2. LANDOLFO S., COFANO F., GIOVARELLI M., PRAT M., CAVALLO G. and FORNI G. (1985) **Inhibition of IFN- $\gamma$  may suppress allograft reactivity by T lymphocytes in vitro and in vivo.** Sci., 119 : 176-179.
3. MURRAY H.W. (1988) **Interferon gamma, the activated macrophage and host defense against microbial challenge.** Ann. Int. Med. , 108 : 595-608.
4. NATHAN C.F., KAPLAN G., LEVIS W.R., NUSRAT A., WITMER M.D., SHERWIN S.A., JOB C.K., HOROWITZ C.R., STEINMAN R.M. and COHN Z.A. (1986) **Local and systemic effects of intradermal rIFN gamma in patients with lepromatous leprosy.** N. Engl. J. Med., 315 : 6-12.
5. SNAPPER C.M. and PAUL W.E. (1987) **IFN- $\gamma$  and BSF-1 reciprocally regulate Ig Isotope production.** Sci., 236 : 944-946.
6. SUZUKI Y., ORELLANA M.A., SCHREIBER R.D. and REMINGTON J.S. (1988) **Interferon gamma : the major mediator of resistance against Toxoplasma gondii.** Sci., 240 : 516-518

**17. SUMMARY OF THE PROTOCOL**

	<b>CALIBRATORS (<math>\mu</math>l)</b>	<b>SAMPLE(S) CONTROLS (<math>\mu</math>l)</b>
Calibrators (0-5)	50	-
Samples, Controls	-	50
Anti-IFN- $\gamma$ -HRP conjugate	50	50
Incubate for 2 hours at 18-28°C with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 $\mu$ l of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 15 min at 18-28°C with continuous shaking at 700 rpm.		
Stop Solution	100	100
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

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