Product information

Information about other products is available at: www.demeditec.com

User's Manual

IGF-I 600 ELISA

Enzyme immunoassay for the quantitative measurement of Insulin-like Growth Factor 1 (IGF-1) in serum

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1 INTRODUCTION

1.1 Intended Use

The **DEMEDITEC IGF-1 600 ELISA** is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of Insulin-like Growth Factor 1 (IGF-1) in serum.

1.2 Summary and Explanation

Clinical Significance

Insulin-like growth factor 1 (IGF-1), also called somatomedin C, is a single-chain 70-amino acid polypeptide, that bears structural similarity to insulin. The IGFs (IGF-1 and IGF-2) are growth factors produced in the majority of the organs and body tissues mainly by the regulation of pituitary growth hormone. They have autocrine, paracrine and endocrine actions on metabolism and cell proliferation, growth and differentiation.

> 95% of serum IGF-1 circulates bound with high specificity and affinity to a family of 6 binding proteins, called IGFBPs (1 to 6) that modulate their bioactivity. IGFBP-3 is thought to be the major binding protein of IGF-1, forming a ternary complex of 140 000 Dalton with IGF-1 and an acid-labile subunit.

Most of the known IGF actions are mediated via IGF type 1 receptor (IGF1R). IGF-1 concentrations change with age, nutritional status, body composition and growth hormone secretion. Vitamin D has been shown to increase circulating IGF-1 and IGFBP-3.

Clinical Applications

The measurement of serum IGF-1 is of recognized value in the differential diagnosis of growth disturbances especially in children. IGF-1 is increased in acromegaly and decreased in impaired growth or dwarfism.

A single basal IGF-1 determination is useful in the assessment of short stature in children and in nutritional support studies of acutely ill patients. For the diagnosis of acromegaly, a single IGF-1 determination is considered more reliable than a random GH determination.

In chronic liver disease, IGF levels are decreased, and the circulating levels correlate to the extent of hepatocellular dysfunction. Patients with cirrhosis are characterized by a variety of metabolic disturbances, including nutritional and metabolic complications such as insulin resistance, malnutrition, osteopenia and hypogonadism, all related to IGF-1 deficiency.

2 PRINCIPLE OF THE TEST

The DEMEDITEC IGF-1 600 ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding.

Patient samples, standards and controls are acidified and neutralized prior to the assay procedure.

The microtiter wells are coated with a monoclonal antibody directed towards an antigenic site on the IGF-1 molecule.

The pre-treated sample is incubated at room temperature with Conjugate (biotinylated IGF-1). The wells are washed and then incubated with Enzyme Complex (Streptavidin-HRP-Complex).

After addition of the substrate solution, the intensity of colour developed is reverse proportional to the concentration of IGF-1 in the patient sample.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DEMEDITEC.

4 REAGENTS

4.1 Reagents provided

1. **HCl** **0.2 M HCl**, 2 vials, 3 mL each, ready to use; for sample acidification.
2. **BUF NEUTRA** **Neutralization Buffer**, 1 vial, 3 mL, ready to use; For neutralization of samples.
3. **SORB MT** **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells; Wells coated with anti-IGF-1 antibody (monoclonal).
4. **CAL 0 – 5** **Standard (Standard 0-5)**, 6 vials, 1 mL each, ready to use; Concentrations: 0 – 10 – 50 – 150 – 300 – 600 ng/mL; Conversion: 1 ng/mL x 0.13 = nmol/L; The standards are calibrated against the International Reference Reagent for IGF-1, NIBSC code: 02/254; Contain non-mercury preservative.
5. **CONTROL low & high** **Control Low & High**, 2 vials, 1 mL each, ready to use; For control values and ranges please refer to vial label or QC-Datasheet. Contain non-mercury preservative.
6. **ENZ CONJ** **Enzyme Conjugate**, 1 vial, 14 mL, ready to use, biotinylated IGF-1; Contains non-mercury preservative.
7. **ENZ COMP** **Enzyme Complex**, 1 vial, 20 mL, ready to use, Streptavidin-HRP Complex. Contains non-mercury preservative.
8. **SUB TMB** **Substrate Solution**, 1 vial, 14 mL, ready to use, Tetramethylbenzidine (TMB).
9. **STOP SOLN** **Stop Solution**, 1 vial, 14 mL, ready to use, contains 0.5 M H₂SO₄. Avoid contact with the stop solution. It may cause skin irritations and burns.
10. **WASH SOLN 40x** **Wash Solution**, 1 vial, 30 mL (40X concentrated), see „Preparation of Reagents“.

Note: Additional Standard 0 for sample dilution is available upon request

4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 nm, with reference wavelength at 620 nm - 630 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled water
- 1.5 mL Reaction Cups (e.g. from Eppendorf) for Sample Preparation (Acidification and Neutralization).
- Timer
- Semi logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Wash Solution

Add deionized water to the 40X concentrated Wash Solution. Dilute 30 mL of concentrated Wash Solution with 1170 mL distilled water to a final volume of 1200 mL. The diluted Wash Solution is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DEMEDITEC has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum can be used in this assay.

NOTE: In plasma significantly reduced values were observed.

Note: Samples containing sodium azide should not be used in the assay.

In general it should be avoided to use haemolytic, icteric or lipaemic specimens. For further information refer to chapter "Interfering Substances".

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 24 hours at 2 °C to 8 °C prior to assaying. Specimens held for a longer time (at least one year) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with Standard 0 and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:2: 50 µL serum + 50 µL Standard 0 (mix thoroughly)
b) dilution 1:10: 10 µL serum + 90 µL Standard 0 (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Acidification and Neutralization of Patient Samples, Standards and Controls

1. Pipette **50 µL Sample, Standard and Control** in 1.5 mL-Reaction caps (e.g. Eppendorf-Caps).
Please note: The standards should be acidified and neutralized too, according to the procedure described below
2. Add **50 µL 0.2 M HCl**.
3. Mix and incubate for 30 minutes.
4. For Neutralization add **10 µL Neutralization Buffer** to all caps and mix the solution.
A pH check and correction of pH is not necessary.
Immediately (within 10 minutes) continue with the test procedure in chapter 6.3.

6.3 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
 2. Dispense **20 µL** of each acidified and neutralized **Standard, Control** and **samples with new disposable tips** into appropriate wells.
 3. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
 4. Incubate for **120 minutes** at room temperature.
 5. Briskly shake out the contents of the wells.
Rinse the wells **3 times** with diluted Wash Solution (400 µL per well).
Strike the wells sharply on absorbent paper to remove residual droplets.
- Important note:**
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Dispense **150 µL Enzyme Complex** into each well.
 7. Incubate for **30 minutes** at room temperature.
 8. Briskly shake out the contents of the wells.
Rinse the wells **3 times** with diluted Wash Solution (400 µL per well).
Strike the wells sharply on absorbent paper to remove residual droplets.
 9. Add **100 µL of Substrate Solution** to each well.
 10. Incubate for **15 minutes** at room temperature.
 11. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
 12. Determine the absorbance (OD) of the solution in each well at **450 nm (reading) and at 620 - 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the Stop Solution.

6.4 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4 Parameter curve fit (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods). Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 600 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.4.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 ng/mL)	2.01
Standard 1 (10 ng/mL)	1.76
Standard 2 (50 ng/mL)	1.21
Standard 3 (150 ng/mL)	0.64
Standard 4 (300 ng/mL)	0.41
Standard 5 (600 ng/mL)	0.23

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy subjects, using the DEMEDITEC IGF-1 600 ELISA the following values are observed:

Age (years)	n	Mean (ng/mL)	2.5 th Percentile (ng/mL)	50 th Percentile (ng/mL)	97.5 th Percentile (ng/mL)	Min. Value (ng/mL)	Max. Value (ng/mL)
Newborn	55	101	65	94	172	65	184
1	10	75	45	77	115	44	119
2	10	91	67	92	120	64	126
3	8	143	103	133	204	103	206
4	11	121	71	119	195	65	208
5	10	165	130	160	232	126	248
6	11	162	80	170	231	62	236
7	10	171	102	169	236	95	242
8	12	200	156	199	249	150	251
9	13	191	61	207	270	44	275
10	11	227	132	237	273	108	278
11	13	227	139	246	308	136	315
12	9	284	174	279	372	158	375
13	3	295	203	265	411	200	419
14	6	269	153	255	416	151	424
15	11	269	119	303	375	111	386
16	9	229	106	228	379	88	407
17	6	251	175	243	349	174	355
18	4	244	186	222	338	184	346
19	5	256	168	246	387	163	402
20	3	317	246	346	363	240	363
21 - 25	23	191	97	175	304	92	304
26 - 30	7	211	137	218	278	133	284
31 - 35	8	184	120	193	229	115	229
36 - 40	5	232	190	220	294	188	300
41 - 45	14	154	107	156	216	105	228
46 - 50	8	150	86	164	196	84	200
51 - 55	17	148	102	140	212	100	214
56 - 60	18	151	96	141	217	91	218
61 - 65	11	130	81	124	216	79	217
66 - 70	10	119	80	117	176	79	183
71 - 75	15	124	57	124	199	46	212
76 - 80	14	122	74	115	192	73	205

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels. The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results. Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or DEMEDITEC directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 9.75 ng/mL - 600 ng/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Compound	% Cross reactivity
IGF-1	100
IGF-2	1.02
Insulin	3.3

9.3 Sensitivity

The analytical sensitivity of the DEMEDITEC ELISA was calculated by subtracting 2 standard deviations from the mean of 20 replicate analyses of the Standard 0 (S0) and was found to be 9.75 ng/mL.

9.4 Reproducibility

9.4.1 Intra-Assay

The within assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	20	89.34	7.39
2	20	227.39	6.93
3	20	390.82	6.39

9.4.2 Inter-Assay

The between assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	40	14.50	14.84
2	40	66.26	10.34
3	40	125.30	12.63

9.5 Recovery

Samples have been spiked by adding IGF-1 solutions with known concentrations in a 1:1 ratio. The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous IGF-1 + added IGF-1) / 2; because of a 1:2 dilution of serum with spike material).

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	
Concentration (ng/mL)	91.60	158.08	170.47	173.55	190.48	
Average Recovery (%)	103.5	102.0	95.8	109.2	95.9	
Range of Recovery (%)	from	96.8	96.6	86.1	98.8	87.5
	to	106.6	106.0	114.9	126.4	101.2

9.6 Linearity

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	
Concentration (ng/mL)	141.55	199.55	201.36	250.29	301.28	
Average Recovery (%)	94.4	94.3	97.7	95.1	97.4	
Range of Recovery (%)	from	88.2	85.4	87.5	89.8	91.4
	to	101.3	107.3	112.7	101.8	105.6

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.25 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

A biotin concentration of up to 1200 ng/mL in a sample has no influence on the assay results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of IGF-1 600 in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DEMEDITEC.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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
12 REFERENCES / LITERATURE

1. Daughaday E, Rotwein P: Insulin like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations. *Endocrin Rev* 10:68-91, 1989.
2. Baxter RC, Martin JL, Beniac VA: High molecular weight insulin-like growth factor binding protein complex. *J Biol Chem* 264:11843-11848, 1989.
3. Rechler M: Insulin-like growth factor binding proteins. *Vit Horm* 47:1-114, 1993.
4. Zapf J, Hauri C, Waldvogel M, Froesch ER: Acute metabolic effects and half-lives of intravenously administered insulin-like growth factors I and II in normal and hypophysectomized rats. *J Clin Invest* 77:1768-1775, 1986.
5. Guler HP, Zapf J, Froesch ER: Short-term metabolic effects of recombinant human insulin-like growth factor-I in healthy adults. *New Engl J Med* 317:1237-140, 1987.
6. Costigan DC, Guyda HJ, Posner BI: Free insulin-like growth factor I (IGF-1) and IGF-1I in human saliva. *J Clin Endocrinol Metab* 66:1014-1018, 1988.
7. Lewitt MS, Denyer GS, Cooney GJ, Baxter RC: Insulin-like growth factor-binding protein-1 modulates blood glucose levels. *Endocrinology* 129:2254-2256, 1991.
8. Lewitt MS, Saunders H, Baxter RC: Bioavailability of insulin-like growth factors (IGFs) in rats determined by the molecular distribution of human IGF-binding protein-3. *Endocrinology* 133:1797-1802, 1993
9. Lieberman SA et al.: Effects of recombinant human insulin-like growth factor-I (rhIGF-1) on total and free IGF-1 concentrations, IGF-binding proteins, and glycemic response in humans. *J Clin Endocrinol Metab* 75:30-36, 1992.
10. Schneiderman R, Maroudas A, Lee PDK: Concentrations of IGF-1 and its complexes in normal and osteoarthritic human cartilage: in situ values. *Orthopedic Res Soc*, submitted, 1994
11. Brabant G et al.: German KIMS Board. Serum insulin-like growth factor I reference values for an automated chemiluminescence immunoassay system: results from a multicenter study. *Horm Res.* 2003;60(2):53-60.
12. Elmlinger MW et al.: Reference ranges for two automated chemiluminescent assays for serum insulin-like growth factor I (IGF-1) and IGF-binding protein 3 (IGFBP-3). *Clin Chem Lab Med.* 2004;42(6):654-64.
11. Bonefeld K, Møller S: Insulin-like growth factor-I and the liver. *Liver Int.* 2011; 31(7):911-9.
12. Ameri P et al.: Interactions between vitamin D and IGF-I: from physiology to clinical practice. *Clin Endocrinol (Oxf).* 2013; 79(4):457-63.
13. Bidlingmaier M et al.: Reference Intervals for Insulin-like Growth Factor-1 (IGF-I) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-I Immunoassay Conforming to Recent International Recommendations *J Clin Endocrinol Metab*, 2014, 99(5):1712–1721.

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konfirmitäts-kennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Ussage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Temperature de conservation	Temperatura de conservacion	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore