Product information

Information about other products is available at: www.demeditec.com

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1 INTRODUCTION

1.1 INTENDED USE

The Demeditec PMN Elastase ELISA is an enzyme immunoassay for the manual quantitative determination of PMN Elastase in human EDTA or citrated plasma and seminal plasma. For in vitro diagnostic use only. For professional use in medical laboratories only. For single use.

1.2 DESCRIPTION OF THE ANALYTE

The human organism reacts with an inflammatory response to attacks of invading pathogens (micro-organisms and viruses) or damaged tissue (after accidents or surgery). Polymorphonuclear (PMN) granulocytes play an important role as primary defence cells in this inflammatory reaction. Different blood-stream mediators (cytokines, leukotrienes, complement factors, bacterial endotoxins, clotting and fibrinolysis factors) attract and stimulate these cells to phagocytize and destroy not naturally occurring agents.

PMN granulocytes use proteinases to digest these agents and tissue debris. One of these proteinases is PMN Elastase which is localised in the azurophilic granules of the polymorphonuclear granulocytes. During phagocytosis of foreign substances these enzymes are also partially excreted into the extracellular surrounding, where the activity of PMN Elastase is regulated by inhibitors (esp. the α_1 -proteinase inhibitor, α_1 -PI). An overwhelming release of PMN Elastase, however, can exceed the inhibitory potential of the α_1 -proteinase inhibitor. Thus, enzymatically active PMN Elastase, together with simultaneously produced oxidants (O_2 -radicals, H_2O_2 , OH-radicals), can cause local tissue injury.

Due to the bloodstream and lymphatic system, however, α_1 -PI is delivered subsequently and eventually able to form a complex with all excreted Elastase. Therefore, the concentration of the PMN Elastase/ α_1 -PI complex correlates with the released PMN Elastase and can be used as a measure for the activity of granulocytes during an inflammatory response.

Primarily, determinations of PMN Elastase find its application in observation of the course of trauma, shock and sepsis. Further indications are the areas of hemodialysis, infections by obstetrics, joint diseases, effusions of sport injuries, intestinal affection, pancreatitis and male adnex affections.

2 PRINCIPLE

The test kit is a solid phase enzyme immunometric assay (ELISA) in the microplate format, designed for the quantitative measurement of the complex of human PMN Elastase and α_1 proteinase inhibitor (α_1 -PI) in plasma. The microplate is coated with a first polyclonal antibody against human PMN Elastase (antigen).

Calibrators, controls and patient samples are pipetted into the antibody-coated microplate. During a 60 minutes incubation present antigens in the sample bind to the antibodies fixed on the inner surface of the wells. Non-reactive sample components are removed by a washing step.

Afterwards, a second polyclonal antibody against α_1 -PI, which is labeled with horseradish peroxidase, is added. During a 60 minutes incubation, the PMN Elastase/ α_1 -PI complex bound to the first antibody is specifically recognized by the enzyme labeled antibodies, and a sandwich complex is formed. An excess of enzyme conjugate is washed out.

A chromogenic substrate, TMB (3,3',5,5'-Tetramethylbenzidine), is added. During a 20 minutes incubation, the substrate is converted to a colored endproduct (blue) by the fixed enzyme. Enzyme reaction is stopped by dispensing of hydrochloric acid as stop solution (change from blue to yellow). The color intensity is direct proportional to the concentration of PMN Elastase present in the sample. The optical density of the color solution is measured with a microplate reader at 450 nm.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. All reagents containing human material should be handled as potentially infectious material.
4. The microplate contains break-apart strips. Unused wells must be stored at 2-8°C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing substrate solution that had previously been used for conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the washing steps.
9. Allow the reagents to reach room temperature (18-25°C) before starting the test. Temperature will affect the absorbance readings of the assay.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable protective gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may be slightly different.
17. Avoid contact with Stop Solution. It may cause skin irritation and burns.
18. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
19. For information please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH or on Demeditec homepage (www.demeditec.com).
20. All serious incidents occurring in relation to products made available on the EU market in accordance with Article 2(61) of Regulation (EU) 2017/746 shall be notified to the manufacturer and to the competent authority of the Member State where the user or patient is established in accordance with Article 82 of Regulation (EU) 2017/746.
21. If product information, including labeling, is incorrect or inaccurate, please contact the kit manufacturer or supplier.

4 REAGENTS

4.1 REAGENTS PROVIDED

1. **SORB MT Microtiterplate**, 12x8 (break-apart) strips, 96 wells.
Wells coated with polyclonal antibodies against PMN Elastase.
2. **CAL PMN Elastase Master Calibrator**, 1 vial (2 µg), lyophilized;
in plasma/buffer matrix containing PMN Elastase/ α1-PI complex
For reconstitution and concentrations see "Preparation of Reagents" (4.4).
3. **CONTROL 1-2 PMN Elastase Controls** in plasma/buffer matrix, 2 vials, lyophilized.
For control values and ranges please refer to QC-Datasheet.
For reconstitution see "Preparation of Reagents" (4.4).
4. **ENZ CONJ Enzyme Conjugate**, 1 vial, 16 ml, ready to use;
Enzyme-labeled anti-α₁-PI antibody, containing polyclonal antibodies labeled with horseradish peroxidase.
5. **SAM DIL Calibrator/Sample Diluent**, 1 vial, 50 ml, ready to use.
6. **SUB TMB Substrate Solution**, 1 vial, 22 ml, ready to use; Tetramethylbenzidine (TMB).
7. **STOP SOLN Stop Solution**, 1 vial, 7 ml, ready to use; contains 2 M hydrochloric acid.
Avoid contact with the stop solution. It may cause skin irritations and burns.
8. **WASH SOLN 10x Wash Solution**, 1 vial, 50 ml (**10X** concentrated);
see „Preparation of Reagents" (4.4).

4.2 MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable for endpoint measurements at 450 nm
- Vortex mixer
- Microplate mixer operating at 900 rpm
- Distilled water
- Graduated cylinders for 1000 ml
- Plastic containers for storage of the wash solution
- Calibrated variable precision micropipettes

4.3 STORAGE CONDITIONS

When stored at 2-8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2-8°C. After first opening the reagents are stable for 30 days if performed and stored properly. Keep reagents away from heat and direct sunlight.

The Wash Buffer is stable for 12 weeks after dilution.

Store Calibrators and Controls at ≤-20 °C (in aliquots), they will be stable for 30 days after reconstitution. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again. Protect TMB-Substrate Solution from light.

4.4 PREPARATION OF REAGENTS

PMN Elastase Master Calibrator:

For plasma: Reconstitute lyophilized Master Calibrator with **2 ml Calibrator/Sample Diluent** 30 min. before use (final concentration of 1000 ng/ml). Make a dilution serie with Calibrator/Sample Diluent to get calibrators with 1000; 500; 250; 125; 62.5; 31.3 and 15.6 ng/ml.

For seminal plasma: Reconstitute lyophilized Master Calibrator with **1 ml Calibrator/Sample Diluent** 30 min. before use (final concentration of 2000 ng/ml). Make a dilution serie with Calibrator/Sample Diluent to get calibrators with 2000; 1000; 500; 250; 125; 62.5 and 31.3 ng/ml.

PMN Elastase Controls:

For plasma and seminal plasma: Reconstitute with **1 ml Calibrator/Sample Diluent** 30 min before use.

Wash Solution:

Dilute 50 ml of 10X concentrated Wash Solution with 450 ml deionized water to a final volume of 500 ml. The diluted Wash Solution is stable for at least 12 weeks at room temperature (18-25°C).

Attention: For the determination of PMN Elastase in seminal plasma, please find detailed information about the reagent preparation in chapter 14 (see page 14).

4.5 DISPOSAL OF THE KIT

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 DAMAGED TEST KIT

In case of any severe damage of the test kit or components, DEMEDITEC DIAGNOSTICS GmbH has to be informed in writing within one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

For determination of PMN Elastase EDTA or citrated plasma are the preferred sample matrixes. For determination of PMN Elastase in seminal plasma please find a separate protocol in chapter 14 (see page 14).

Exudates, bronchial lavage and cerebrospinal fluid as sample material are not validated with this test. Therefore these sample material could be used for research purposes only, but not for use in in-vitro diagnostics.

Serum is not suitable, because during clotting PMN Elastase can be released *in vitro*. Culture supernatants are not suitable; because the assay detects only the PMN Elastase/ α_1 -PI complex and α_1 -PI is normally not present in culture medium.

All samples are prediluted 1:100 with Calibrator/Sample Diluent. Therefore 10 μ l of sample may be diluted with 990 μ L of Calibrator/Sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate plasma from the cells by centrifugation. Plasma samples can be stored at 2-8°C up to 5 days. For longer storage samples should be stored frozen at \leq -20°C. To avoid repeated thawing and freezing the samples should be aliquoted. Patient samples expected to contain higher PMN Elastase concentrations than the highest calibrator (1000 ng/ml) should be diluted in the Calibrator/Sample Diluent before further assaying. The additional dilution step has to be taken into account for the calculation of the results.

6 ASSAY PROCEDURE

6.1 GENERAL REMARKS

- Do not interchange components of different lots.
- All components of these test kits, supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent of foam formation.
- All reagents and specimens must be allowed to come to room temperature (18-25°C) before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross-contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Respect the incubation times as stated in this instructions for use.
- Calibrators, controls and samples should at least be assayed in double determinations.
- A calibrator curve must be established for every test run.

6.2 ASSAY PROCEDURE

1. Preparation of calibrators:

Label six tubes: G (500 ng/ml), F (250 ng/ml), E (125 ng/ml), D (62.5 ng/ml), C (31.3 ng/ml), and B (15.6 ng/ml). Pipet **0.5 ml** of the **Calibrator/Sample Diluent** into all tubes. Pipet **0.5 ml** of the **reconstituted PMN Elastase Master Calibrator** into tube G (500 ng/ml) and mix thoroughly. Transfer 0.5 ml from tube G (500 ng/ml) to tube F (250 ng/ml) and mix thoroughly. Repeat this process successively to complete the 2-fold dilution series. The reconstituted PMN Elastase Master Calibrator will serve as the highest calibrator H (1000 ng/ml). Use the PMN Elastase Calibrator/Sample Diluent as the zero calibrator A (0 ng/ml).

2. Dilute all patient samples 1:100 with **Calibrator/Sample Diluent** before assay. Therefore combine 10 µl of sample with 990 µl of Calibrator/Sample Diluent in a polystyrene tube. Mix well. Calibrators and controls are ready to use and must **not** to be diluted.
3. Prepare a sufficient number of microplate wells to accomodate calibrators, controls and prediluted patient samples in duplicates.

	1	2	3	4	5	6	7	8	9	10	11	12
a	A	E	C1	P..								
b	A	E	C1	P..								
c	B	F	C2									
d	B	F	C2									
e	C	G	P1									
f	C	G	P1									
g	D	H	P2									
h	D	H	P2									

4. For determination of PMN Elastase pipet **100 µl** of **Calibrators, Controls and prediluted patient samples** with new disposable tips into the wells according to the template.
5. Incubate for **60 minutes** at room temperature (18-25°C) on a plate mixer (900 rpm).
6. Discard the content of the wells and wash **4 times** with **300 µl** buffered wash solution. Remove as much wash solution as possible by beating the microplate carefully.
7. Pipet **150 µl** of **Enzyme Conjugate** into each well.
8. Incubate for **60 minutes** at room temperature (18-25°C) on a plate mixer (900 rpm).
9. Again discard the content of all wells and wash **4 times** with **300 µl** buffered wash solution. Remove as much wash solution as possible by beating the microplate carefully.
10. Dispense **200 µl** of **TMB Substrate Solution** into each well.
11. Incubate for **20 minutes** at room temperature (18-25°C) in the dark without shaking.
12. Add **50 µl** of **Stop Solution** to each well and mix carefully.
13. Determine the absorbance of each well at 450 nm. It is recommended to read the wells within 15 minutes.

For determination of PMN Elastase in seminal plasma please find a separate protocol in chapter 14 (see page 14).

6.3 CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of calibrators, controls and samples.
2. Using semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample, determine the corresponding concentration from the calibration curve.
4. Automated method: The results in the package insert have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be determined directly from this calibrator curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations, this dilution factor has to be taken into account.


6.4 EXAMPLE OF TYPICAL CALIBRATOR CURVE

The figure below shows typical results for Demeditec PMN Elastase test. These data are intended for illustration only and must not be used to calculate results from another run.

Calibrator		Optical Density (450 nm)
Calibrator A	0.0 ng/ml	0.098
Calibrator B	15.6 ng/ml	0.187
Calibrator C	31.3 ng/ml	0.292
Calibrator D	62.5 ng/ml	0.441
Calibrator E	125 ng/ml	0.699
Calibrator F	250 ng/ml	1.231
Calibrator G	500 ng/ml	1.916
Calibrator H	1000 ng/ml	2.751

7 EXPECTED NORMAL VALUES

In a normal range study with plasma samples from healthy blood donors (n = 57) the following ranges have been established with the Demeditec PMN Elastase test:

Frequency distribution of PMN elastase in citrated plasma of healthy blood donors (median = 35 ng/ml, 95. percentile = 64.9 ng/ml)

The results alone should not be the only reason for any diagnostic or therapeutical consequences. They have to be correlated to other clinical observations and diagnostic tests.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of plasma PMN Elastase. The reference ranges should be regarded as guidelines only.

8 QUALITY CONTROL

Good laboratory practice requires that controls are run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate included in the kit. The values and ranges stated on the QC certificate always refer to the current kit lot and should be used for direct comparison of the results. It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices, microtiter plate reader, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or Demeditec directly.

9 PERFORMANCE CHARACTERISTICS

9.1 ANALYTICAL SENSITIVITY

The lowest analytical detectable level of PMN Elastase that can be distinguished from the Zero Calibrator is 2.31 ng/ml at the 2SD confidence limit.

9.2 SPECIFICITY

The Demeditec PMN Elastase test is specific human PMN elastase only, respectively the PMN elastase/ α_1 -PI complex.

9.3 REPRODUCIBILITY

Statistics for Coefficients of Variation (CV) were calculated for each of three samples from the results of 10 determinations in a single run for Intra-Assay precision and the Inter-Assay precision was calculated from the results of 10 different runs of four samples:

Intra-Assay		
Sample No.	Mean (ng/ml)	CV (%)
1	80	5.2
2	241	4.7
3	358	4.6

Inter-Assay		
Sample No.	Mean (ng/ml)	CV (%)
1	128	5.7
2	216	6.4
3	346	4.4
4	681	5.7

9.4 RECOVERY

Three spiking solutions were prepared using the Calibrator/Sample Diluent (922, 615 and 478 ng/ml). A 50 µl aliquot of each solution (A, B, C) was spiked into 950 µl aliquots of three different patient plasma samples, for a spiking ratio of 1 to 20, leaving the plasma matrix of the spiked samples relatively intact. All samples were then assayed by the Demeditec PMN Elastase procedure.

Sample No.	Spiking Solution	measured Concentration [ng/ml]	expected Concentration [ng/ml]	Recovery [%]
1	-	23.2	-	-
	A	72.4	69.4	104
	B	59.3	54.1	109
	C	49.6	47.2	101
2	-	30.6	-	-
	A	73.4	76.7	96
	B	59.3	61.4	97
	C	56.8	54.4	104
3	-	61.7	-	-
	A	118.0	107.8	109
	B	100.8	92.5	109
	C	94.8	85.6	110

9.5 LINEARITY

In dilution experiments samples with high PMN Elastase concentrations were diluted with Calibrator/Sample Diluent and assayed in the Demeditec PMN Elastase test. The assay showed linearity over the full measuring range.

Sample No.	Dilution Factor	measured Concentration [ng/ml]	expected Concentration [ng/ml]	Recovery [%]
1	1:1	250.1	-	-
	1:2	145.2	125.1	116
	1:4	75.3	62.5	120
	1:8	37.1	31.3	119
2	1:1	465.5	-	-
	1:2	209.2	232.8	90
	1:4	111.1	116.4	95
	1:8	58.7	58.2	101

10 LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 INTERFERING SUBSTANCES

To simulate moderate and severe icterus, four samples were spiked with 100 and 200 mg bilirubin per liter. All samples were assayed, both spiked and unspiked, by the Demeditec PMN Elastase procedure, with the following results:

Sample No.	Unspiked	+100 mg/l Bilirubin	+200 mg/l Bilirubin
	[ng/ml]	[ng/ml]	[ng/ml]
1	100	106	104
2	249	245	261
3	572	575	534
4	903	964	910

The results show that severe icterus (bilirubin up to 200 mg/l) has no clinically significant effect on the Demeditec PMN Elastase procedure.

Icteric, lipaemic, and haemolytic samples and samples should not be used in this assay.

10.2 DRUG INTERFERENCES

Until today no substances (drugs) are known to us, which have an influence to the measurement of PMN Elastase in a sample.

11 LEGAL ASPECTS

11.1 RELIABILITY OF RESULTS

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Demeditec Diagnostics GmbH.

11.2 THERAPEUTIC CONSEQUENCES

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 LIABILITY

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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13 SHORT INSTRUCTION(all sample sizes given in μl)

MP Well	ng/ml	A	B	C	D	E	F	G	H	Control 1/2	Sample
		0	15.6	31.3	62.5	125	250	500	1000		
Steps	Solution										
Pipet	Calibrator	100	100	100	100	100	100	100	100	-	-
Pipet	Control	-	-	-	-	-	-	-	-	100	-
Pipet	Prediluted sample	-	-	-	-	-	-	-	-	-	100
Incubate for 60 min at RT (18-25°C) on a plate mixer (900 rpm)											
Decant Wash 4x with 300 μl of buffered wash solution											
Pipet	Enzyme conjugate	150	150	150	150	150	150	150	150	150	150
Incubate for 60 min at RT (18-25°C) on a plate mixer (900 rpm)											
Decant Wash 4x with 300 μl of buffered wash solution											
Pipet	Substrate Solution	200	200	200	200	200	200	200	200	200	200
Incubate for 20 min at RT (18-25°C) without shaking in the dark											
Pipet	Stop Solution	50	50	50	50	50	50	50	50	50	50
Read at $\lambda = 450 \text{ nm}$											

For a detailed description of the procedure see also page 7.

14 PROTOCOL FOR SEMINAL PLASMA

Separate the sperm by centrifugation (5 min).

Take the supernatant and freeze the seminal plasma at $\leq -20^{\circ}\text{C}$ for longer storage.

Preparation of the Calibrators and Controls

1. Reconstitute PMN Elastase Master Calibrator with **1 ml Calibrator/Sample Diluent** 30 min. before use (Calibrator H: final concentration of 2000 ng/ml). Label six tubes: G (1000 ng/ml), F (500 ng/ml), E (250 ng/ml), D (125 ng/ml), C (62.5 ng/ml), and B (31.3 ng/ml). Pipet **0.5 ml** of the **Calibrator/Sample Diluent** into all tubes. Pipet **0.5 ml** of the **reconstituted PMN Elastase Master Calibrator** into tube G (1000 ng/ml) and mix thoroughly. Transfer 0.5 ml from tube G (1000 ng/ml) to tube F (500 ng/ml) and mix thoroughly. Repeat this process successively to complete the 2-fold dilution series. The reconstituted PMN Elastase Master Calibrator will serve as the highest calibrator H (2000 ng/ml). Use the PMN Elastase Calibrator/Sample Diluent as the zero calibrator A (0 ng/ml).

G = 1000 ng/ml

F = 500 ng/ml

E = 250 ng/ml

D = 125 ng/ml

C = 62.5 ng/ml

B = 31.3 ng/ml

A = 0 ng/ml (Calibrator-/Sample Diluent)

2. Reconstitute controls with **1 ml** Calibrator-/Sample Diluent 30 min. before use and mix well. For storage keep the controls frozen in aliquots at $\leq -20^{\circ}\text{C}$.
3. Dilute seminal plasma **1:100** with Calibrator-/Sample Diluent (e.g. 10 μl seminal plasma + 990 μl Calibrator/Sample Diluent)
4. Pipet **100 μl Calibrators, Controls and diluted seminal plasma** into the wells according to the template.
5. Incubate for **60 minutes** at room temperature ($18-25^{\circ}\text{C}$) on a plate mixer (900 rpm).
6. Discard the content of the wells and wash **4 times** with **300 μl** buffered wash solution. Remove as much wash solution as possible by beating the microplate carefully.
7. Pipet **150 μl of Enzyme conjugate** into each well.
8. Incubate for **60 minutes** at room temperature ($18-25^{\circ}\text{C}$) on a plate mixer (900 rpm).
9. Again discard the content of all wells and wash **4 times** with **300 μl** buffered **Wash Solution**. Remove as much wash solution as possible by beating the microplate carefully.
10. Dispense **200 μl of TMB Substrate solution** into each well.
11. Incubate for **20 minutes** at room temperature ($18-25^{\circ}\text{C}$) in the dark without shaking.
12. Add **50 μl of Stop Solution** to each well and mix carefully.
13. Determine the absorbance of each well at 450 nm. It is recommended to read the wells within 15 minutes.

Interpretation of results:

Interpretation	Concentration
no inflammation	0 - 250 ng/ml
moderate inflammation	> 250 - 1,000 ng/ml
massive inflammation	> 1,000 ng/ml

The results alone should not be the only reason for any diagnostic or therapeutical consequences. They have to be correlated to other clinical observations and diagnostic tests.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of seminalplasma PMN Elastase. The given ranges should be regarded as guidelines only.

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