
# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

User's Manual

# Prolactin canine ELISA

**VET**

**REF** DEV9944

 96 Wells

**CONTENTS**

1 INTRODUCTION ..... 3  
2 PRINCIPLE ..... 3  
3 WARNINGS AND PRECAUTIONS ..... 4  
4 REAGENTS ..... 5  
5 SPECIMEN ..... 6  
6 ASSAY PROCEDURE ..... 6  
7 EXPECTED NORMAL VALUES ..... 8  
8 PERFORMANCE CHARACTERISTICS ..... 8  
9 LIMITATIONS OF PROCEDURE ..... 9  
10 REFERENCES ..... 10  
11 SHORT INSTRUCTION ..... 10

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## 1 INTRODUCTION

### 1.1 INTENDED USE

The Demeditec Prolactin canine ELISA is an enzyme immunoassay for the quantitative measurement of canine prolactin.

### 1.2 SUMMARY AND EXPLANATION

Canine prolactin (cPRL) is a single-chain polypeptide hormone of the canine anterior pituitary with a molecular mass of approx. 22,000. Prolactin from different species exhibits significant variations in the amino acid sequence. Canine prolactin differs from human prolactin at about 60 percent of all residues.

The secretion of cPRL from the pituitary is inhibited by hypothalamic prolactin-inhibitory factor (PIF). Although dopamine was long thought to be this PIF molecule, today it seems that there is a special peptide with prolactin-inhibiting activities. The release of prolactin is certainly stimulated by different peptides, particularly thyrotropin releasing hormone (TRH) and vasocative intestinal peptide (VIP). Estrogens and progesterone also seem to play a role in the secretion of prolactin, and neurogenic factors influence its release. Milking and suckling are immediately followed by an increase in serum cPRL.

The most important role of prolactin is stimulation of mammary gland growth and lactation. During pregnancy, prolactin levels in canine blood increase slightly; during lactation, significantly. Prolactin has a wide variety of other physiological actions. It affects water and electrolyte balance, metabolism and gonadal function; is an important stress hormone; and seems to play a role in the maintenance of the long interestrous interval in the bitch.

In dogs with pituitary-depent hyperadrenocorticism, prolactin levels in blood were higher than in healthy animals. Prolactin determinations can be used in the therapeutic control of hyperprolactinemia. During a pseudo pregnancy, prolactin is increased. Therapy with alkaloids like bromocriptine lowers PRL levels, and lactation and maternal behaviour are decreased.

The secretory capacity of the pituitary can be tested with the TRH stimulation test.

## 2 PRINCIPLE

The test kit is a solid phase enzyme immunometric assay (ELISA) in the microplate format, designed for the quantitative measurement of canine prolactin. The microplate is coated with a first monoclonal antibody specific for canine prolactin.

Calibrators and samples are pipetted into the antibody coated microplate. During a 2 hours incubation endogenous canine prolactin in the sample bind to the antibodies fixed on the inner surface of the wells. Non-reactive sample components are removed by a washing step.

Afterwards, a second polyclonal horseradish peroxidase-labeled antibody, directed against another epitope of the Prolactin molecule, is added. During an 1 hour incubation, a sandwich complex consisting of the two antibodies and the canine prolactin is formed. An excess of enzyme conjugate is washed out.

A chromogenic substrate, TMB (3,3',5,5'-Tetra-Methyl-Benzidine), is added to all wells. During a 30 minutes incubation, the substrate is converted to a colored end product (blue) by the fixed enzyme. Enzyme reaction is stopped by dispensing of hydrochloric acid as stop solution (change from blue to yellow). The color intensity is direct proportional to the concentration of canine prolactin present in the sample.

The optical density of the color solution is measured with a microplate reader at 450 nm. Bi-chromatic measurement with a 600 - 690 nm reference filter is recommended.

### 3 WARNINGS AND PRECAUTIONS

1. This kit is strictly intended for veterinary use only. Use by staff, who is specially informed and trained in methods which are carried out by use of immunoassays.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
4. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
5. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
6. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
7. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
8. Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
9. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
10. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
11. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
12. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
13. Do not use reagents beyond expiry date as shown on the kit labels.
14. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
15. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
16. Avoid contact with Stop Solution. It may cause skin irritation and burns.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
18. For information please refer to Material Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH.

## 4 REAGENTS

### 4.1 REAGENTS PROVIDED

1. **SORB MT** **Microtiter Plate**, 12 x 8 (break apart) strips with 96 wells, ready to use; coated with a monoclonal anti-canine prolactin antibody.
2. **CAL** **Master Calibrator**, 1 vial, 80 ng, lyophilized; in serum/buffer matrix containing highly purified canine prolactin; **For reconstitution see "Reagent preparation"**.
3. **SAM DIL** **Calibrator/Sample Diluent**, 1 vial, 6 ml, ready to use; canine prolactin free
4. **ENZ CONJ** **Enzyme Conjugate**, 1 vial, 22 ml, red, ready to use; containing horseradish peroxi-dase-labeled polyclonal anti canine prolactin antibody
5. **SAM BUF** **Sample Buffer**, 1 vial, 6 ml, yellow, ready to use;
6. **SUB TMB** **Substrate Solution**, 1 vial, 22 ml, ready to use; 3,3',5,5'-Tetra-Methyl-Benzidine in buffered peroxide solution
7. **STOP SOLN** **Stop Solution**, 1 vial, 7 ml, ready to use; contains 2 M hydrochloric acid
8. **WASH SOLN 10x** **Wash Solution**, 1 vial, 50 ml, 10 x concentrated; see "Reagent preparation".

### 4.2 MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable for endpoint measurements at 450 nm (optional reference filter in the range of 600 - 690 nm)
- Vortex mixer
- Microplate mixer operating more than 600 rpm
- Distilled or deionized water
- Graduated cylinders for 500 ml
- Plastic containers for storage of the wash solution
- Adjustable pipette for up to 1000 µl
- Dispenser or repeatable pipet for 25 µl, 50 µl and 200 µl.

### 4.3 REAGENT PREPARATION

#### Calibrators:

Reconstitute lyophilized Master Calibrator with **1 ml dest. water** 30 min. before use (end concentration of 80 ng/ml). Make a dilution serie with Calibrator/Sample Diluent to get calibrators with 80, 40, 20, 10, 5 and 2.5 ng/ml.

#### Wash Solution:

Dilute with 450 ml dist. water to a final volume of 500 ml.

*The diluted Wash Solution is stable for 12 weeks at room temperature.*

### 4.4 STORAGE CONDITIONS

When stored at 2°C to 8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2°-8°C. After first opening the reagents are stable for 30 days if used and stored properly.

Microtiter wells must be stored at 2°C to 8°C. Take care that the foil bag is sealed tightly. Store Calibrators refrigerated, they will be stable at 2°C to 8°C for 7 days after reconstitution. For longer storage aliquot and freeze at -20°C.

Protect Substrate Solution from light.

## 5 SPECIMEN

For determination of canine prolactin serum is the preferred sample matrix. The procedure calls for 25 µl matrix per well. Prolactin is one of the most sensitive stress hormones of the dog. Blood collection should therefore be as stress-free as possible.

The samples may be stored refrigerated at 2 - 8°C for one week, or up to 2 months frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Samples expected to contain canine prolactin concentrations higher than the highest calibrator (80 ng/ml) should be diluted with the Calibrator/Sample Diluent before assay. The additional dilution step has to be taken into account for the calculation of the results.

## 6 ASSAY PROCEDURE

### 6.1 GENERAL REMARKS

- Do not interchange components of different lots.
- All components should be at room temperature (21 – 26 °C) before use.
- All components of these test kits, supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent of foam formation.
- Use a disposable-tip micropipette to dispense plasma samples. Pipet directly to the bottom of the wells. Change the tip between samples, to avoid carryover contamination.
- For internal quality control we suggest to use **Canine Control coded DEV99CC**. For more information please contact DEMEDITEC.

### 6.2 ASSAY PROCEDURE

#### 1. Preparation of calibrators:

Label five tubes: F (40 ng/ml), E (20 ng/ml), D (10 ng/ml), C (5 ng/ml), and B (2.5 ng/ml). Pipet **0.1 ml** of the Calibrator/Sample Diluent into all tubes. Pipet 0.1 ml of the reconstituted Master Calibrator into tube F (40 ng/ml) and mix thoroughly. Transfer 0.1 ml from tube F (40 ng/ml) to tube E (20 ng/ml) and mix thoroughly. Repeat this process successively to complete the 2-fold dilution series. The reconstituted Calibrator will serve as the highest calibrator G (80 ng/ml). Use the Calibrator/Sample Diluent as the zero calibrator A (0 ng/ml).

	1	2	3	4	5	6	7	8	9	10	11	12
a	A	E	P2	P..								
b	A	E	P2	P..								
c	B	F	P3									
d	B	F	P3									
e	C	G	P4									
f	C	G	P4									
g	D	P1	P5									
h	D	P1	P5									

- Pipet 25 µl of each calibrator and patient sample into the wells prepared.
- Add 50 µl of Sample Buffer to every well.
- Rotate for 2 hours at room temperature (21 - 26 °C) on a plate mixer (600 – 900 rpm).
- Discard the content of the wells and wash 4 times with 300 µl buffered wash solution. Remove as much wash solution as possible by beating the microplate carefully.
- Add 200 µl of Enzyme Conjugate to all wells.
- Shake again for 1 hour at room temperature (21 - 26 °C) on a plate mixer (600 – 900 rpm).
- Discard the content of the wells and wash 4 times with 300 µl buffered wash solution. Remove as much wash solution as possible by beating the microplate carefully.
- Add 200 µl of liquid Substrate Solution to all wells.
- Incubate without shaking for 30 minutes in the dark.
- Add 50 µl of Stop Solution to each well and mix carefully.
- Read the optical density at 450 nm. Bi-chromatic measurement with a reference at 600-690 nm is recommended.

**The developed color is stable for at least 15 minutes. Read optical densities during this time.**

### 6.3 CALCULATION OF RESULTS

For evaluation of canine prolactin a 4-Parameter-Fit with lin-log coordinates for optical density (linear scale) and concentration (logarithmic scale) is recommended.  
Spline approximation with lin-log coordinates and log-log coordinates are also suitable.

### 6.4 EXAMPLE OF TYPICAL CALIBRATOR CURVE

The figure below shows typical results for canine prolactin test kits. These data are intended for illustration only and should not be used to calculate results from another run.

	Replicate (OD)	Mean (OD)	Binding (%)	canine prolactin (ng/ml)
Calibrators				
A	0.059 ----- 0.055	0.057	-	0
B	0.148 ----- 0.148	0.148	4.4	2.5
C	0.244 ----- 0.299	0.272	8.1	5
D	0.495 ----- 0.559	0.527	15.7	10
E	0.917 ----- 0.953	0.935	27.8	20
F	1.996 ----- 2.074	2.035	60.6	40
G (Bmax)	3.207 ----- 3.507	3.357	100	80
Unknown Samples				
X 001	0.795 ----- 0.772	0.784	23.3	15.8
X 002	1.703 ----- 1.730	1.717	51.1	33.9
X 003	2.384 ----- 2.446	2.415	71.9	50.0

## 7 EXPECTED NORMAL VALUES

In a reference range study canine serum samples were collected in the morning between 8 and 9 a.m. and in the evening between 5 and 6 p.m. Diurnal variations have not been observed. Analysis by the Demeditec Prolactin canine ELISA kit yielded the following results:

Group	Absolute Range (ng/ml)	n
Normal dogs	nd - 21	26
nd = non detectable		

Because of differences which may exist between laboratories with respect of population, laboratory technique and selection of reference groups, it is recommended that each laboratory establishes its own normal and pathological ranges of canine prolactin. The reference ranges should be regarded as guidelines only.

## 8 PERFORMANCE CHARACTERISTICS

### 8.1 ANALYTICAL SENSITIVITY

The lower detection limit for canine prolactin was 0.4 ng/ml.

### 8.2 SPECIFICITY

The antibodies in the Demeditec Prolactin canine ELISA procedure are highly specific for canine prolactin. Detectable crossreactivities to other hormones that may be present in serum samples are not known.

### 8.3 REPRODUCIBILITY

Statistics for Coefficients of variation (CV) were calculated for each of three samples from the results of 12 pairs of wells in a single run for Intra-Assay precision and the Inter-Assay precision was calculated from the results of 10 different runs of three samples:

Intra-Assay		
Sample No	Mean $\bar{x}$ (pg/ml)	CV (%)
1	6.8	7.1
2	29	6.0
3	50	7.4

Inter-Assay		
Sample No	Mean $\bar{x}$ (pg/ml)	CV (%)
1	8.8	9.2
2	15	6.9
3	32	5.5

#### 8.4 RECOVERY

Three spiking solutions were prepared using the Sample Diluent, to represent the 600, 800 and 1000 ng/ml, respectively. A 50 µl aliquot of each solution (A, B, C) was spiked into 950 µl aliquots of two different patient serum samples, for a spiking ratio of 1 to 20, leaving the serum matrix of the spiked samples relatively intact. All samples were then assayed by the Demeditec Prolactin canine ELISA procedure.

Sample	Diluted Solution	measured Concentration [ng/ml]	expected Concentration [ng/ml]	Recovery [%]
1	-	8.8	-	-
	A	35.3	38.4	92
	B	48.3	48.4	100
	C	56.2	58.4	96
2	-	6.1	-	-
	A	33.7	35.8	94
	B	49.7	45.8	109
	C	53.1	55.8	95

#### 8.5 LINEARITY

In dilution experiments sera with high antibody concentrations were diluted with sample diluent and assayed in the Demeditec Prolactin canine ELISA kit. The assay showed linearity over the full measuring range.

Sample	Dilution Factor	measured Concentration [ng/ml]	expected Concentration [ng/ml]	Recovery [%]
1	8 in 8	54.9	-	-
	4 in 8	24.4	27.5	89
	2 in 8	12.8	13.7	93
	1 in 8	6.6	6.9	96
2	8 in 8	54.2	-	-
	4 in 8	27.8	27.1	103
	2 in 8	15.4	13.6	113
	1 in 8	7	6.8	103

#### 9 LIMITATIONS OF PROCEDURE

The Demeditec Prolactin canine ELISA has no "high-dose hook" effect, even with samples containing more than 300 ng/ml of canine prolactin. However, this effect is characteristic of immunometric assays. Samples expected to contain canine prolactin concentrations greater than the highest calibrator (80 ng/ml) should be diluted with the Calibrator/Sample Diluent.

**10 REFERENCES**

1. Grünau, B., Nolte, I., Hoppen, H.-O.  
 Untersuchung zur Behandlung der Scheinträchtigkeit der Hündin mit den Prolaktinhemmern Metergolin und Bromocriptin.  
 Tierärztl Prax 1996; 24: 149-155.
2. Harvey, M.J.A., Cauvin, A., Dale, M., Lindley, S., Ballabio, R.  
 Effect and mechanisms of the anti-prolactin drug cabergoline on pseudopregnancy in the bitch.  
 Small Animal Practice 38 (1997): 336-339.
3. Cortese, L., Oliva, G., Verstegen, J., Ciaramella, P., Persechino, A.  
 Hyperprolactinaemia and galactorrhoea associated with primary hypothyroidism in a bitch.  
 Small Animal Practice 38 (1997): 572-575.
4. Hoppen, H.O., Grünau, B., Hayer, M., Günzel-Apel, A.R.  
 Prolactin in canine reproduction: normal values under various conditions.  
 Proceedings "Advances in Veterinary Endocrinology" Berlin, 1993 October 6; 18 - 19.

**11 SHORT INSTRUCTION**

(all sample sizes given in µl)

MP Well	ng/ml	0	1	2	3	4	5	6	Sample
		(A)	(B)	(C)	(D)	(E)	(F)	(G)	
Steps	Solution	0	2.5	5	10	20	40	80	
Pipet	Calibrator	25	25	25	25	25	25	25	-
Pipet	Sample	-	-	-	-	-	-	-	25
Pipet	Sample Buffer	50	50	50	50	50	50	50	50
Incubate for 2 hours at RT on a shaker									
Decant Wash 4x with 300 µl of buffered wash solution									
Pipet	Enzyme Conjugate	200	200	200	200	200	200	200	200
Incubate for 1 hour at RT on a shaker									
Decant Wash 4x with 300 µl of buffered wash solution									
Pipet	Substrate Solution	200	200	200	200	200	200	200	200
Incubate for 30 min at RT in the dark									
Pipet	Stop Solution	50	50	50	50	50	50	50	50
Read at λ = 450 nm									

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