Product information

Information about other products is available at: www.demeditec.com

User's Manual

Prolactin rat ELISA

Enzyme immunoassay for the quantitative measurement of prolactin in rat serum.

RUO

REF

DEV9966

96 Wells

CONTENTS

1. INTRODUCTION 3
2. PRINCIPLE 3
3. WARNINGS AND PRECAUTIONS..... 4
4. REAGENTS 5
5. SPECIMEN 6
6. ASSAY PROCEDURE 6
7. EXPECTED NORMAL VALUES..... 8
8. PERFORMANCE CHARACTERISTICS 8
9. LIMITATIONS OF PROCEDURE 9
10. REFERENCES 10
11. SHORT INSTRUCTION 11

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SYMBOLS USED WITH DEMEDITEC ELISAS.....24

1. INTRODUCTION

1.1. INTENDED USE

The Demeditec Prolactin rat ELISA is an enzyme immunoassay for the quantitative measurement of prolactin in rat serum.

1.2. SUMMARY AND EXPLANATION

Rat prolactin (rPRL) is a single-chain polypeptide hormone of the rat anterior pituitary with a molecule mass of approximately 23,000. Prolactin from different species exhibits significant variations in the amino acid sequence. Rat prolactin differs from human prolactin at about 50 percent of all residues. The secretion of rPRL from the pituitary is inhibited by hypothalamic prolactin-inhibitory factor (PIF). Although dopamine was long thought to be this PIF molecule, today it seems that there is a special peptide with prolactin-inhibiting activities.

The release of prolactin is certainly stimulated by different peptides, particularly thyrotropin releasing hormone (TRH) and vasoactive intestinal peptide (VIP). There is also evidence that rat posterior pituitary lobe contains a special prolactin releasing hormone.

The most important role of prolactin is stimulation of mammary gland growth and lactation. During pregnancy, blood prolactin levels climb, but the increases can differ enormously between rats. High prolactin levels are observed during lactation. Prolactin has a wide variety of other physiological actions, for example on the ovary. In the rat, prolactin has a luteotrophic effect which is not seen in many other species. Furthermore, prolactin is a stress hormone.

In rats, as in humans, prolactin exhibits a sleep-related diurnal variation. Peak values are seen in the late afternoon and nadir values in the morning.

Because of the variety of its actions, prolactin is one of the preferred hormones to monitor when testing the influence of new therapeutic agents and drugs on the endocrine system in the rat.

2. PRINCIPLE

The Demeditec Prolactin rat ELISA kit is a solid phase enzyme immunometric assay (ELISA) in the microplate format, designed for the quantitative measurement of rat prolactin. The microplate is coated with a first monoclonal antibody specific for rat prolactin.

Calibrators and samples are pipetted into the antibody coated microplate. During a 2 hours incubation endogenous rat prolactin in the sample bind to the antibodies fixed on the inner surface of the wells. Non-reactive sample components are removed by a washing step. Afterwards, a second polyclonal horseradish peroxidase-labeled antibody, directed against another epitope of the prolactin molecule, is added. During another 1 hour incubation, a sandwich complex consisting of the two antibodies and the rat prolactin is formed. An excess of enzyme conjugate is washed out.

A chromogenic substrate, TMB (3,3',5,5'-Tetra-Methyl-Benzidine), is added to all wells. During a 30 minutes incubation, the substrate is converted to a colored end product (blue) by the fixed enzyme. Enzyme reaction is stopped by dispensing of hydrochloric acid as stop solution (change from blue to yellow). The color intensity is direct proportional to the concentration of rat prolactin present in the sample. The optical density of the color solution is measured with a microplate reader at 450 nm. Bichromatic measurement with a 600 - 690 nm reference filter is recommended.

3. WARNINGS AND PRECAUTIONS

1. This kit is strictly intended for research use only. Use by staff, who is specially informed and trained in methods which are carried out by use of immunoassays.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
4. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
5. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
6. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
7. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
8. Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
9. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
10. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
11. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
12. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
13. Do not use reagents beyond expiry date as shown on the kit labels.
14. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
15. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
16. Avoid contact with Stop Solution. It may cause skin irritation and burns.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
18. For information please refer to Material Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH.

4. REAGENTS

4.1. REAGENTS PROVIDED

1. **SORB MT Microtiter Plate**, 12 x 8 (break apart) strips with 96 wells, ready to use
Removable wells coated with a monoclonal anti-rat prolactin antibody.
2. **CAL Master Calibrator**, 1 vial, 80 ng, lyophilized,
in serum/buffer matrix containing highly purified rat prolactin,
For reconstitution see "Reagent preparation".
3. **SAM DIL Calibrator/Sample Diluent**, 1 vial, 6 ml, ready to use
rat prolactin free
4. **ENZ CONJ Enzyme Conjugate**, 1 vial, 22 ml, red, ready to use
containing horseradish peroxidase-labeled polyclonal anti rat prolactin antibody in a buffered
solution with preservative
5. **SAM BUF Sample Buffer**, 1 vial, 6 ml, yellow, ready to use
6. **SUB TMB Substrate Solution**, 1 vial, 22 ml, ready to use
3,3',5,5'-Tetra-Methyl-Benzidine in buffered peroxide solution
7. **STOP SOLN Stop Solution**, 1 vial, 7 ml, ready to use
containing 2 M hydrochloric acid
8. **WASH SOLN 10x Wash Solution**, 1 vial, 50 ml, 10 x concentrated,
See "Reagent preparation".

4.2. MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable for endpoint measurements at 450 nm (optional reference filter in the range of 600 - 690 nm)
- Vortex mixer
- Microplate mixer operating more than 600 rpm
- Distilled or deionized water
- Graduated cylinder for 500 ml
- Plastic containers for storage of the wash solution
- Adjustable pipette for up to 1000 µl
- Dispenser or repeatable pipet for 25 µl, 50 µl, and 200 µl

4.3. REAGENT PREPARATION

Calibrators:

Reconstitute lyophilized Master Calibrator with **1 ml dest. water** 30 min. before use (end concentration of 80 ng/ml). Make a dilution series with Calibrator/Sample Diluent to get calibrators with 80, 40, 20, 10 and 5 ng/ml.

Wash Solution:

Dilute with 450 ml dist. water to a final volume of 500 ml.

The diluted Wash Solution is stable for 12 weeks at room temperature.

4.4. STORAGE CONDITIONS

When stored at 2°C to 8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2°-8°C. After first opening the reagents are stable for 30 days if used and stored properly.

Microtiter wells must be stored at 2°C to 8°C. Take care that the foil bag is sealed tightly. Store Calibrators refrigerated, they will be stable at 2°C to 8°C for 7 days after reconstitution. For longer storage aliquot and freeze at -20°C. Protect Substrate Solution from light.

5. SPECIMEN

For determination of rat prolactin serum is the preferred sample matrix. The procedure calls for 25 µl matrix per well.

Prolactin is one of the most sensitive stress hormones of the rat. Blood collection should therefore be as stress-free as possible. The samples may be stored refrigerated at 2 - 8 °C for one week, or up to 2 months frozen at -20 °C. To avoid repeated thawing and freezing the samples should be aliquoted.

Samples expected to contain rat prolactin concentrations higher than the highest calibrator (80 ng/ml) should be diluted with Calibrator/Sample Diluent before assay. The additional dilution step has to be taken into account for the calculation of the results.

6. ASSAY PROCEDURE

6.1. GENERAL REMARKS

- Do not interchange components of different lots.
- All components should be at room temperature (18 – 28 °C) before use.
- All components of this test kit, supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent of foam formation.
- Use a disposable-tip micropipette to dispense serum samples. Pipet directly to the bottom of the wells. Change the tip between samples, to avoid carryover contamination.
- For internal quality control we suggest to use **Rat Control Set coded DEV99RC. For more information please contact DEMEDITEC.**

6.2. ASSAY PROCEDURE

1. Preparation of calibrators:

Label four tubes: E (40 ng/ml), D (20 ng/ml), C (10 ng/ml) and B (5 ng/ml). Pipet **0.1 ml** of the Calibrator/Sample Diluent into all tubes. Pipet 0.1 ml of the reconstituted Master Calibrator into tube E (40 ng/ml), and mix thoroughly. Repeat this process successively to complete the 2-fold dilution series. The reconstituted Calibrator will serve as the highest calibrator F (80 ng/ml). Use the Calibrator/Sample Diluent as the zero Calibrator A (0 ng/ml).

	1	2	3	4	5	6	7	8	9	10	11	12
a	A	E	P3	P..								
b	A	E	P3	P..								
c	B	F	P4									
d	B	F	P4									
e	C	P1	P5									
f	C	P1	P5									
g	D	P2	P6									
h	D	P2	P6									

2. Pipet 25 µl of each calibrator, control and sample into the wells prepared.
3. Add 50 µl of Sample Buffer to every well.
4. Rotate for 2 hours at room temperature (18 - 28 °C) on a plate mixer (600 - 900 rpm).
5. Discard the content of the wells and wash 4 times with 300 µl buffered wash solution. remove as much wash solution as possible by beating the microplate carefully.
6. Add 200 µl of Enzyme Conjugate to all wells.
7. Shake again for 1 hour (600 – 900 rpm).
8. Discard the content of the wells and wash 4 times with 300 µl buffered Wash Solution. Remove as much Wash Solution as possible by beating the microplate carefully.
9. Add 200 µl of liquid Substrate Solution to all wells.
10. Incubate without shaking for 30 minutes in the dark.
11. Add 50 µl of Stop Solution to each well and mix carefully.
12. Read the optical density at 450 nm. Bi-chromatic measurement with a reference at 600-690 nm is recommended.

The developed color is stable for at least 30 minutes. Read optical densities during this time.

6.3. CALCULATION OF RESULTS

For evaluation of rat prolactin a 4-Parameter-Fit with lin-log coordinates for optical density (linear scale) and concentration (logarithmic scale) is recommended.

Spline approximation with lin-log coordinates and log-log coordinates are also suitable.

6.3.1. EXAMPLE OF TYPICAL CALIBRATOR CURVE

The figure below shows typical results for Demeditec Prolactin rat ELISA test kits. These data are intended for illustration only and should not be used to calculate results from another run.

	Replicate (OD)	Mean (OD)	Binding (%)	Rat Prolactin (ng/mL)
Calibrators				
A	0.079	0.073	-	0
	----- 0.067			
B	0.189	0.179	6.9	5
	----- 0.169			
C	0.328	0.342	12.7	10
	----- 0.357			
D	0.766	0.796	29.7	20
	----- 0.827			
E	1.591	1.635	61.0	40
	----- 1.679			
F	2.588	2.677	100	80
	----- 2.765			
Unknown Samples				
X 001	0.318	0.324	12,1	9.8
	----- 0.329			
X 002	0.577	0.590	22,0	16.4
	----- 0.603			
X 003	1,733	1,716	64,1	47.3
	----- 1,698			

7. EXPECTED NORMAL VALUES

In a reference range study rat serum samples were collected in the morning between 8 and 9 a.m. and in the evening between 5 and 6 p.m. Diurnal variations have not been observed. Analysis by the Demeditec Prolactin rat ELISA kit yielded the following results:

Group	Absolute Range (ng/ml)	n
Normal female rats	nd - 17.9	15
Normal male rats	nd- - 23.4	10
nd = nondetectable		

Because of differences which may exist between laboratories with respect of population, laboratory technique and selection of reference groups, it is recommended that each laboratory establishes its own normal and pathological ranges of rat prolactin. The reference ranges should be regarded as guidelines only.

8. PERFORMANCE CHARACTERISTICS

8.1. ANALYTICAL SENSITIVITY

The lower detection limit for rat prolactin is approximately 0.6 ng/ml.

8.2. SPECIFICITY

The antibodies in the Demeditec Prolactin rat ELISA procedure are highly specific for rat prolactin. Detectable crossreactivities to other hormones that may be present in serum samples are not known.

The following substances were tested:

	added quantity (ng/ml)	measured concentration (ng/ml)	cross-reactivity (%)
Rat TSH	5,000	1.5	0.03
Rat FSH	10,000	5.1	0.5
Rat LH	5,000	nd	-
Rat GH	4,000	nd	-

8.3. REPRODUCIBILITY

Statistics for Coefficients of variation (CV) were calculated for each of three samples from the results of 20 pairs of wells in a single run for Intra-Assay precision and the Inter-Assay precision was calculated from the results of 14 different runs of three samples:

Intra - Assay			
Serum No	Mean \bar{x} (ng/ml)	SD $\pm s$ (ng/ml)	CK (%)
1	6.6	0.20	3.0
2	19.6	0.77	3.9
3	33.0	1.80	5.5

Inter - Assay			
Serum No	Mean \bar{x} (ng/ml)	SD $\pm s$ (ng/ml)	CK (%)
1	6.5	0.28	4.3
2	20.5	0.71	3.5
3	33.4	1.48	4.4

8.4. RECOVERY

Three spiking solutions were prepared using the Sample Diluent, to represent the 600, 800 and 1,000 ng/ml, respectively. A 50 µl aliquot of each solution (A, B, C) was spiked into 950 µl aliquots of two different rat serum samples, for a spiking ratio of 1 to 20, leaving the serum matrix of the spiked samples relatively intact. All samples were then assayed by the Demeditec Prolactin rat ELISA procedure.

Demeditec Prolactin rat ELISA

Sample	Diluted Solution	measured Concentration	expected Concentration	Recovery
		[ng/ml]	[ng/ml]	[%]
1	-	26.1	-	-
	A	54.7	56.1	97
	B	66.6	66.1	100
	C	79.9	76.1	105
2	-	24.2	-	-
	A	52.6	54.2	103
	B	67.5	64.2	95
	C	81.1	74.2	91

8.5. LINEARITY

In dilution experiments sera with high rat prolactin concentrations were diluted with sample diluent and assayed in the Demeditec Prolactin rat ELISA kit. The assay showed linearity over the full measuring range.

Demeditec Prolactin rat ELISA

Sample	Dilution Factor	measured Concentration	expected Concentration	Recovery
		[ng/ml]	[ng/ml]	[%]
1	8 in 8	20.9	-	-
	4 in 8	10.6	10.5	101
	2 in 8	5.9	5.2	113
	1 in 8	3.0	2.6	115
2	8 in 8	34.7	-	-
	4 in 8	16.7	17.4	96
	2 in 8	9.0	8.7	103
	1 in 8	4.5	4.3	105

9. LIMITATIONS OF PROCEDURE

Effect of Anticoagulants

To determine whether anticoagulants interfere with the assay, blood was collected from 30 rats into plain and EDTA vacutainer tubes. All samples were assayed by the Demeditec Prolactin rat ELISA procedure, with the following results.

$$(\text{EDTA}) = 1.05 (\text{Serum}) - 9.3 \text{ ng/ml} \quad r = 0.978$$

Means: 3.68 ng/ml (Serum)
 3.78 ng/ml (EDTA)

A limited study with citrated and heparinized plasma show comparable results to EDTA plasma.

"High-Dose Hook"-Effect

Rat sera containing up to 300 ng/ml Prolactin were measured with the Demeditec Prolactin rat ELISA assay. A High-Dose Hook effect could not be observed.

10. REFERENCES

1. Leung, F. C., Russel, S. M.; Nicoll, C. S.
Relationship between bioassay and radioimmunoassay estimates of prolactin in rat serum.
Endocrinology 1978; 103: 1619 - 28.
2. Martinat, N., Hall, E., Ravault, J. P., Dubois, M. P.
Purification of rat prolactin: development of an homologous radioimmunological assay and comparison with the NIAMDD system.
Ann Biol Anim Bioch Biophys 1979; 19: 1771 - 48.
3. Beach, J. E., Miles, D. J., Lukes, Y. G., Vigersky, R. A.
Microplate solid-phase radioimmunoassay for rat prolactin.
J Lab Clin Med 1985; 105: 294 - 298.
4. Butcher, R. L., Collins, W. E., Fugo, N. W.
Plasma concentration of LH, FSH, prolactin, progesterone and estradiol-17 β throughout the 4-day estrous cycle of the rat.
Endocrinology 1974; 94: 1704 - 8.
5. Barbieri, R. L., Todd, R. B., Morishita, H., Ryan, K. J., Fishman, J, Naftolin, F.
Response of serum prolactin to catechol estrogen in the immature rat.
Fertil Steril 1980; 34: 391 - 3.
6. Wong, C. C.
Endogene und exogene Einflüsse auf die Variabilität der Hormonausschüttung bei der Ratte (Endogenous and exogenous influences on the variability of hormone release in the rat).
Thesis, Hannover (Germany): Univ. of Hannover, 1981.
7. Campbell, G. A., Kurcz, M., Marshall, S., Meites, J.
Effects of starvation in rats on serum levels of FSH, LH, TSH, growth hormone and prolactin: response to LH-releasing hormone and thyrotropin-releasing hormone.
Endocrinology 1977; 100: 580 - 7.
8. Haggi, E., Aoki, A.
Prolactin content in rat pituitary gland. RIA of prolactin after different extraction procedures.
Acta Endocrinol 1981; 97: 338 - 42.
9. Reichel, J.
Wirkungen gonadaler Steroide auf die adeno-hypophysären Thyreotropin-Releasing-Hormon-Rezeptoren, den Prolactinserumspiegel und die hormonelle Hypophysen-Schilddrüsen-Achse der Ratte (Effects of gonadal steroids on the pituitary TRH-receptors, the serum prolactin concentrations and the pituitary-thyroid axis in the rat)
Thesis.
Lübeck (Germany): University of Lübeck, 1990.
10. Moishige, W. K., Pepe, G. J., Rothschild, I.
Serum LH, prolactin and progesterone levels during pregnancy in the rat. Endocrinology 1973; 92: 1527 - 30.

11. SHORT INSTRUCTION

(all sample sizes given in µl)

MP Well Steps	ng/ml Solution	0 (A)	1 (B)	2 (C)	3 (D)	4 (E)	5 (F)	Sample
		0	5	10	20	40	80	
Pipet	Calibrator	25	25	25	25	25	25	-
Pipet	Sample	-	-	-	-	-	-	25
Pipet	Sample Buffer	50	50	50	50	50	50	50
Incubate for 2 hours at RT on a shaker								
Decant Wash 4x with 300 µl of buffered wash solution								
Pipet	Enzyme Conjugate	200	200	200	200	200	200	200
Incubate for 1 hour at RT on a shaker								
Decant Wash 4x with 300 µl of buffered wash solution								
Pipet	Substrate Solution	200	200	200	200	200	200	200
Incubate for 30 min at RT in the dark								
Pipet	Stop Solution	50	50	50	50	50	50	50
Read at $\lambda = 450$ nm								

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10. LITERATUR


1. Leung, F. C., Russel, S. M.; Nicoll, C. S.
Relationship between bioassay and radioimmunoassay estimates of prolactin in rat serum.
Endocrinology 1978; 103: 1619 - 28.
2. Martinat, N., Hall, E., Ravault, J. P., Dubois, M. P.
Purification of rat prolactin: development of an homologous radioimmunological assay and comparison with the NIAMDD system.
Ann Biol Anim Bioch Biophys 1979; 19: 1771 - 48.
3. Beach, J. E., Miles, D. J., Lukes, Y. G., Vigersky, R. A.
Microplate solid-phase radioimmunoassay for rat prolactin.
J Lab Clin Med 1985; 105: 294 - 298.
4. Butcher, R. L., Collins, W. E., Fugo, N. W.
Plasma concentration of LH, FSH, prolactin, progesterone and estradiol-17 β throughout the 4-day estrous cycle of the rat.
Endocrinology 1974; 94: 1704 - 8.
5. Barbieri, R. L., Todd, R. B., Morishita, H., Ryan, K. J., Fishman, J, Naftolin, F.
Response of serum prolactin to catechol estrogen in the immature rat.
Fertil Steril 1980; 34: 391 - 3.
6. Wong, C. C.
Endogene und exogene Einflüsse auf die Variabilität der Hormonausschüttung bei der Ratte (Endogenous and exogenous influences on the variability of hormone release in the rat).
Thesis, Hannover (Germany): Univ. of Hannover, 1981.
7. Campbell, G. A., Kurcz, M., Marshall, S., Meites, J.
Effects of starvation in rats on serumlevels of FSH, LH, TSH, growth hormone and prolactin: response to LH-releasing hormone and thyrotropin-releasing hormone.
Endocrinology 1977; 100: 580 - 7.
8. Haggi, E., Aoki, A.
Prolactin content in rat pituitary gland. RIA of prolactin after different extraction procedures.
Acta Endocrinol 1981; 97: 338 - 42.
9. Reichel, J.
Wirkungen gonadaler Steroide auf die adeno-hypophysären Thyreotropin-Releasing-Hormon-Rezeptoren, den Prolactinserumspiegel und die hormonelle Hypophysen-Schilddrüsen-Achse der Ratte (Effects of gonadal steroids on the pituitary TRH-receptors, the serum prolactin concentrations and the pituitary-thyroid axis in the rat)
Thesis.
Lübeck (Germany): University of Lübeck, 1990.
10. Moishige, W. K., Pepe, G. J., Rothschild, I.
Serum LH, prolactin and progesterone levels during pregnancy in the rat.
Endocrinology 1973; 92: 1527 - 30.

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SYMBOLS USED WITH DEMEDITEC ELISAS

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konformitäts-kennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Temperature de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore