PSA free IRMA KIT

REF IM2520**TABLE OF CONTENTS**

English	2	Čeština (CZ)	22
Français (FR)	4	Slovenčina (SK)	24
Deutsch (DE)	6	한국어 (KO)	26
Italiano (IT)	8	Türkçe (TR)	28
Español (ES)	10	Русский (RU)	30
Ελληνικά (EL)	13	Србија (SR)	33
Lietuviškai (LT)	16	APPENDIX	35
Magyar (HU)	18		
Polski (PL)	20		

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PSA free IRMA KIT

REF IM2520

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF FREE PROSTATE SPECIFIC ANTIGEN (PSA) IN HUMAN SERUM For *in vitro* diagnostic use.

PRINCIPLE

The determination of prostate specific antigen (PSA) is a "sandwich" type assay. In the kit, mouse monoclonal antibodies directed against two different epitopes of PSA and hence not competing are used. The samples the controls and calibrators are incubated in monoclonal antibody-coated tubes with the second ¹²⁵I-labelled antibody, specific for free PSA. After incubation, the content of tubes is aspirated and the tubes are rinsed so as to remove unbound ¹²⁵I-labeled antibody. The bound radioactivity is then determined in a gamma counter. The free PSA concentrations in the samples are obtained by interpolation from the standard curve. The concentration of free PSA in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- Do not mix the reagents from kits of different lots.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in tubes containing no additive.
- Separate serum from cells by centrifugation.
- Serum samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at <-70°C, 1 year maximum) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample must be performed at room temperature.
- The samples need not be diluted, as there is recommended to analyse samples containing concentrations of total PSA up to 10 ng/mL.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Anti-PSA monoclonal antibody-coated tubes: 50 tubes (ready-to-use)

Monoclonal ¹²⁵I-labelled anti-PSA antibody: one 5.5 mL vial; (ready-to-use)

The vial contains less than 275 kBq (at the date of manufacture) of ¹²⁵I-labelled antibody in buffer containing bovine serum albumin, sodium azide (<0.1%) and a dye.

Calibrators: five 1.0 mL vials and one 2.0 mL vial with «zero» calibrator (ready-to-use)

The vials contain from 0 to approximately 30 ng/mL of human free PSA in buffer with bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against the 1st IS NIBSC 96/668.

Control samples: two vials (lyophilized)

The vials contain human free PSA lyophilised in bovine serum and sodium azide (<0.1%). The volume for reconstitution is indicated on the vial label, the concentration range after reconstitution is given on a supplement.

Wash solution (500x): one 1 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipette (200 µL)
- Semi-automatic pipette (100 µL, 2 mL)
- Vortex-type mixer
- Horizontal or orbital shaker
- Aspiration system
- Gamma counter set for ¹²⁵I

PROCEDURE

Preparation and storage of reagents

Let all the reagents come to room temperature.

Reconstitution of control samples

The content of the vials is reconstituted with the volume of distilled water indicated on the vial label. Wait at least 10 minutes and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one day or frozen and aliquoted at <-18°C until the expiry date of the kit.

Preparation of the wash solution

- Pour the content of the vial into 500 mL of distilled water and homogenize.
- The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Assay procedure

Let all the reagents come to room temperature.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To coated tubes add successively: 200 µL of calibrator, control or sample and 100 µL of tracer. Mix.	Incubate for 2 hours at 18-25°C with shaking (> 280 rpm).	Aspirate carefully the content of tubes (except the 2 tubes «total cpm»). Wash twice with 2mL of wash solution and aspirate. Count bound cpm (B) and total cpm (T) for 1 min.

*Add 100 µL of tracer to 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of free PSA concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the package insert were calculated using a log-log curve fit ("spline" mode) with determined radioactivity ($cpm_{cal} - cpm_{cal.0}$) on vertical axis and the free PSA concentration of the calibrators on the horizontal axis (ng/mL). Other data reduction methods may give slightly different results.

Total activity: 239,264 cpm				
Calibrators	Free PSA (ng/mL)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal.0}$
0	0	300	0.13	—
1	0.3	1,598	0.67	1,298
2	1.0	4,437	1.85	4,137
3	3.0	12,350	5.16	12,050
4	10.0	39,754	16.6	39,454
5	30.0	102,910	43.0	102,610

(Example of standard curve, do not use for calculation)

Samples

For each sample or control, locate the ($cpm_{sample} - cpm_{cal.0}$) value on the vertical axis and read off the corresponding free PSA concentration in ng/mL on the horizontal axis.

EXPECTED VALUES

The clinical evaluation of results is based on the concentration ratio freePSA/totalPSA (in percents). It is suggested that each laboratory establishes its own range of values for the distinguishing benign and malign prostate disease, on the bases of the sufficient number of clinically proved samples. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

A retrospective clinical trial was conducted on 174 patient samples; specimen selection criterion was the concentration of total PSA ranging from 4 to 10 ng/mL at the time of diagnosis. 54 patients from this group were diagnosed as prostate cancer positive then. Cut-off values corresponding to 90% clinical sensitivity and specificity are shown at the table.

Cut-off	Sensitivity	Specificity
23%	90%	38%
16.3%	63%	90%

The probability of prostate cancer absence was significantly increased in patients with free/total PSA ratio higher than 23%. In such a case, and if other

clinical or laboratory findings and evidences do not suggest the presence of prostate cancer, it is possible to avoid unnecessary prostate biopsies in these patients.

The free/total PSA values between 16.3 and 23 % was assessed as equivocal; results below 16.3% indicated significantly increased probability of prostate cancer presence.

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended to analyse their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.01 ng/mL

Functional sensitivity: 0.13 ng/mL

Specificity

The antibodies used in this kit do not cross-react with human CEA, AFP, prolactin, hCG and PAP. The cross-reactivity with PSA-aa1-antichymotrypsin complex was determined to be less than 1%.

Precision

Intra-assay

Samples were assayed in 10 replicates in the same series. The coefficients of variation were found below or equal to 4.1% for serum samples.

Inter-assay

Samples were assayed in duplicate in 9 different series. Coefficients of variation were found below or equal to 5.8% for serum samples.

Accuracy

Dilution test

Five serum samples were diluted with the "zero" calibrator and assayed. The recovery percentages obtained were between 82% and 117%.

Measurement range (from analytical sensitivity to highest calibrator):

0.01 to approximately 30 ng/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Do not use hemolytic, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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