TOTAL T4 RIA KIT

REF IM1447, IM3286

TABLE OF CONTENTS

English	2	Slovenčina	22
Français	4	한국어	24
Deutsch	7	Türkçe	26
Italiano	9	Русский	28
Español	11	中文 ZH-TW	31
Português (Portugal)	13	Србија	33
Ελληνικά	15	APPENDIX	35
Polski	18		
Čeština	20		

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telephone: 0800-771-8818

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TOTAL T4 RIA KIT

REF IM1447, IM3286

RADIOIMMUNOASSAY FOR THE IN VITRO DETERMINATION OF TOTAL THYROXINE (TT4) IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

PRINCIPLE

The radioimmunoassay of total thyroxine (TT4) is a competition assay. Samples and calibrators are incubated with ¹²⁵I-labeled T4, as tracer, in antibody-coated tubes. After incubation, the liquid content of tubes is aspirated and the bound radioactivity is determined in a gamma counter. A standard curve is constructed and unknown values are obtained from the curve by interpolation.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes with EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage (up to 1 months) keep frozen at < -20°C after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.

Serum and EDTA plasma values for 15 samples (serum values ranging from 71.08 to 167.32 nmol/L) were compared using the IM1447 Total T4 RIA kit. Results are as follows:

[EDTA-plasma] = 0.6673 [serum] + 22.12

R = 0.9749

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution are indicated in paragraph Procedure.

Kit for determination of total T4, 100 tubes (Cat. #IM1447)

Anti-T4 monoclonal antibody-coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled T4 tracer: one 55 mL vial (ready-to-use)

The vial contains 110 kBq, at the date of manufacture, of ¹²⁵I-labeled T4 in buffer with proteins, sodium azide (<0.1%) and a dye.

Calibrators: six 0.5 mL vials (ready-to-use)

The calibrator vials contain from 0 to approximately 400 nmol/L of T4 in human serum with sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated using the standard IRMM-468.

Control serum: two vials (lyophilised)

The vials contain T4 lyophilised in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on the supplement.

Kit for determination of total T4, 400 tubes (Cat. #IM3286)

Anti-T4 monoclonal antibody-coated tubes: 8 x 50 tubes (ready-to-use)

¹²⁵I-labeled T4 tracer: four 55 mL vial (ready-to-use)

Calibrators: six 0.5 mL vials (ready-to-use)

Control serum: two vials (lyophilised)

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipets (20 µL).
- Semi-automatic pipets (500 µL).
- Vortex type mixer.
- Horizontal or orbital shaker.
- Aspiration system.
- Gamma counter set for 125 iodine.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of TT4 concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *weighted cubic regression* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (nmol/L).

Other data reduction methods may give slightly different results.

Total activity: 44,140 cpm				
Calibrators	TT4 (nmol/L)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0	34,959	77.2	100
1	26	31,590	71.6	90.4
2	53	26,213	59.4	75.0
3	105	17,048	38.6	48.8
4	210	9,727	22.0	27.8
5	420	5,187	11.8	14.8

(Example of standard curve, do not use for calculation).

Samples

Locate for each sample the ratio B/T or B/B_0 on the vertical axis of the standard curve and read-off the corresponding TT4 concentration of the sample on the horizontal axis in nmol/L.

To convert concentrations from nmol/L to ng/dL, multiply results by 77.7.

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. Normal concentration range of total T4 in serum for untreated euthyroid individuals was found:

	N	Mean, (nmol/L)	Median, (nmol/L)	Min-Max, (nmol/L)	2.5th-97.5th percentile (nmol/L)
Female	50	118.3	114.7	59.65 - 168.2	74.1 - 160.3
Male	50	96.42	96.76	59.13 - 141.4	68.91 - 122.3
Male and Female	100	107.4	101.9	59.13 - 168.2	69.32 - 159.7

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Reconstitution of control serum

The contents of the vials must be brought to room temperature before reconstitution with the volume of distilled water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C for one week or aliquoted at < -18°C until the expiry date of the kit.

Assay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes add successively: 20 µL of calibrator, control or sample and 500 µL of tracer Mix.	Incubate 1 hour at 18-25°C with shaking (>280 rpm).	Aspirate carefully the content of tubes (except the 2 tubes «total cpm»). Count bound cpm (B) and total cpm (T) for 1 min.

*Add 500 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 10.63 nmol/L

Functional sensitivity: 16.71 nmol/L

Specificity

The antibody used in the immunoassay is highly specific for T4. Extremely low cross reactivities were obtained with several related molecules (e.g. L-T3).

Precision

Intra-assay

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 3.29%.

Inter-assay

Serum samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 7.53%.

Accuracy

Dilution test

High-concentration samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 88.1% and 112% for serum.

Recovery test

Samples were spiked with known quantities of TT4. The recovery percentages were obtained between 81.0% and 107% for serum.

Measurement range (from analytical sensitivity to highest calibrator): 10.63 to approximately 400 nmol/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Specificity

Data on cross-reactivity with several hormones are presented in the following table:

Analogue	Cross reaction (%)
L-T4, (L-thyroxine), (3,3,5,5-tetraiodo-L-thyronine)	100
L-T3 - (3,3,5-triiodo-L-thyronine)	ND
3,5-diiodo-L-thyronine	ND
3,5-diiodo-L-tyrosine dihydrate	ND
3,3,5-triiodothyroacetic acid (TRIAc)	ND

Precision

Intra-assay

Serum	S1	S2	S3
Number of determinations	25	25	25
Mean (nmol/L)	89.35	122.0	265.2
CV (%)	3.29	2.64	2.55

Plasma-EDTA	P1	P2	P3
Number of determinations	25	25	25
Mean (nmol/L)	58.52	271.8	389.6
CV (%)	6.45	5.18	5.97

Inter-assay

Serum	S1	S2	S3
Number of determinations	10	10	10
Mean (nmol/L)	53.23	115.6	457.4
CV (%)	6.82	4.67	7.53

Plasma-EDTA	P1	P2	P3
Number of determinations	10	10	10
Mean (nmol/L)	40.78	117.3	353.6
CV (%)	10.70	8.68	4.74

Accuracy

Dilution test

Five serum and plasma samples were chosen for dilution test. Samples were serially diluted by calibrator 0 and assayed. Results are shown below:

Serum	Dilution	Serum		Ratio (%) Measured/ Expected
		Measured Concentration (nmol/L)	Expected	
S1	-	166.1	-	-
	1:2	86.84	83.03	104.6
	1:4	36.58	41.51	88.12
	1:8	21.17	20.76	102.0
S2	-	234.0	-	-
	1:2	116.6	117.00	99.69
	1:4	57.55	58.50	98.38
	1:8	28.47	29.25	97.34
	1:16	13.60	14.62	93.00
S3	-	257.1	-	-
	1:2	132.2	128.54	102.8
	1:4	65.56	64.27	102.0
	1:8	34.09	32.14	106.1
S4	-	15.87	16.07	101.5
	1:2	253.2	-	-
	1:4	131.0	126.59	103.5
	1:8	66.20	65.51	101.1
S5	-	30.94	32.76	94.46
	1:16	16.54	16.38	101.0
	1:2	140.9	138.88	101.5
	1:4	71.29	69.44	102.7
S5	1:8	37.29	34.72	107.4
	1:16	19.39	17.36	111.7

EDTA plasma	Dilution	EDTA-Plasma		Ratio (%) Measured/ Expected
		Measured Concentration (nmol/L)	Expected	
P1	-	166.1	-	-
	1:2	86.84	83.03	104.6
	1:4	26.01	22.54	115.4
P2	-	188.5	-	-
	1:2	88.03	94.25	93.40
	1:4	38.73	47.13	82.19
P3	-	26.68	23.56	113.2
	1:2	290.1	-	-
	1:4	158.3	145.1	109.1
P4	-	82.47	72.53	113.7
	1:8	33.54	36.26	92.49
	1:16	20.12	18.13	111.0
P5	-	328.0	-	-
	1:2	165.8	164.0	101.1
	1:4	78.29	82.00	95.48
P5	-	35.09	41.00	85.59
	1:8	20.27	20.50	98.88
	1:16	346.2	-	-
P5	-	185.4	173.1	107.1
	1:2	90.85	86.56	105.0
	1:4	46.53	43.28	107.5
P5	1:8	22.71	21.64	104.9

Recovery test

The recovery was performed on five serum/EDTA-plasma samples. Samples were spiked by calibrators with higher concentration. Volume of added calibrator was up to 10% of volume of samples. Samples were then assayed. Results are shown in the table below:

Serum	Endogen. conc. nmol/L	Added conc. nmol/L	Expected conc. nmol/L	Measured conc. nmol/L	Ratio (%) Measured/ Expected
S1	66.75	46.32	113.1	114.2	101.0
	65.46	90.87	156.3	142.1	90.91
	61.89	214.8	276.7	228.3	82.53
S2	67.76	46.32	114.1	122.2	107.1
	66.46	90.87	157.3	155.9	99.07
	62.84	214.8	277.6	256.9	92.55
S3	116.7	76.21	192.9	176.3	91.40
	113.0	147.7	260.7	219.9	84.37
	109.6	214.8	324.4	262.7	81.00
S4	108.2	76.21	184.4	172.8	93.71
	104.8	147.7	252.5	230.9	91.45
	101.6	214.8	316.4	268.8	84.95
S5	93.82	46.32	140.1	144.3	103.0
	92.02	90.87	182.9	170.3	93.14
	87.00	214.8	301.8	269.1	89.16

EDTA-plasma	Endogen. conc. nmol/L	Added conc. nmol/L	Expected conc. nmol/L	Measured conc. nmol/L	Ratio (%) Measured/ Expected
P1	61.57	9.27	70.84	71.75	101.3
	61.57	18.55	80.12	78.53	98.02
	61.57	37.09	98.66	92.74	94.00
P2	75.93	9.27	85.20	80.27	94.21
	75.93	18.55	94.48	90.96	96.27
	75.93	37.09	113.0	105.4	93.24
P3	73.95	9.27	83.22	82.00	98.53
	73.95	18.55	92.50	90.32	97.64
	73.95	37.09	111.0	103.4	93.08
P4	96.75	1.65	98.40	103.7	105.3
	76.75	3.24	99.99	114.0	114.0
	96.75	4.58	101.3	114.6	113.1
P5	69.49	9.27	78.76	83.01	105.4
	69.49	18.55	88.04	85.66	97.30
	69.49	37.09	106.6	110.5	103.6

Kinetics

Incubation time may be prolonged by one hour without effect on the results.

¹²⁵I Characteristics

$$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$$

¹²⁵ I	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25

Symbols Key

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktivtyoiji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 质控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provör / σωληνώρα / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Диагностика in vitro / Ин vitro диагностика / 體外診斷

Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó / Producent / Výrobce / 제조 / Üretici / Изготовлено / Произведено от / 製造商

Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Råcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida teste için yeterli / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Знаčka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문외 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижете Инструкциите за употреба / 請參閱使用說明

Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Eύρος(-η) θερμοκρασίας / 温度范围 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Ispėjimas / Figyelem / Uwaga / Urozozomení / Uprozozomenie / 주의 / Dikkat / Внимание / 注意

Expiration Date / Date d'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期