Product information

Information about other products is available at: www.demeditec.com

User's Manual

Anti-Ovarian Ab ELISA

Enzyme immunoassay for the quantitative measurement of antibodies directed against oocytes in human serum

DE2937

96 wells

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

| | | |
|---|--|----|
| 1 | INTENDED USE..... | 3 |
| 2 | PRINCIPLE OF THE TEST..... | 3 |
| 3 | WARNINGS AND PRECAUTIONS..... | 4 |
| 4 | REAGENTS..... | 5 |
| 5 | SPECIMEN COLLECTION AND PREPARATION..... | 6 |
| 6 | ASSAY PROCEDURE..... | 7 |
| 7 | EXPECTED NORMAL VALUES..... | 8 |
| 8 | QUALITY CONTROL..... | 8 |
| 9 | PERFORMANCE CHARACTERISTICS..... | 8 |
| 10 | LIMITATIONS OF USE..... | 9 |
| 11 | LEGAL ASPECTS..... | 9 |
| | | |
| 1 | ZWECKBESTIMMUNG..... | 10 |
| 2 | TESTPRINZIP..... | 10 |
| 3 | VORSICHTSMAßNAHMEN..... | 11 |
| 4 | BESTANDTEILE DES KITS..... | 12 |
| 5 | PROBENVORBEREITUNG..... | 13 |
| 6 | TESTDURCHFÜHRUNG..... | 13 |
| 7 | ERWARTETE WERTE..... | 15 |
| 8 | QUALITÄTSKONTROLLE..... | 15 |
| 9 | ASSAY-CHARAKTERISTIKA..... | 16 |
| 10 | GRENZEN DES TESTS..... | 16 |
| 11 | RECHTLICHE GRUNDLAGEN..... | 16 |
| | | |
| 1 | DESTINAZIONE D'USO..... | 17 |
| 2 | PRINCIPIO DEL TEST..... | 17 |
| 3 | PRECAUZIONI..... | 17 |
| 4 | COMPONENTI DEL KIT..... | 18 |
| 5 | CAMPIONI..... | 19 |
| 6 | ATTUAZIONE DEL TEST..... | 19 |
| 7 | VALORI NORMALI..... | 21 |
| 8 | CONTROLLO QUALITÀ..... | 21 |
| 9 | CARATTERISTICHE DEL TEST..... | 21 |
| 10 | LIMITAZIONE DEL TEST..... | 21 |
| 11 | ASPETTI LEGALI..... | 22 |
| | | |
| 1 | FINALIDAD PREVISTA..... | 23 |
| 2 | FUNDAMENTO DEL ENSAYO..... | 23 |
| 3 | PRECAUCIONES..... | 24 |
| 4 | COMPONENTES DEL KIT..... | 25 |
| 5 | MUESTRAS..... | 26 |
| 6 | PROCEDIMIENTO DE ENSAYO..... | 26 |
| 7 | VALORES ESPERADOS..... | 28 |
| 8 | CONTROL DE CALIDAD..... | 28 |
| 9 | CARACTERÍSTICAS DEL ENSAYO..... | 28 |
| 10 | LIMITACIONES DE USO..... | 28 |
| 11 | ASPECTOS LEGALES..... | 29 |
| | | |
| 12 | REFERENCES / LITERATURE..... | 30 |
| | | |
| SYMBOLS USED WITH DEMEDITEC ASSAYS..... | | 32 |

1 INTENDED USE

The **DEMEDIATEC Anti Ovarian Ab ELISA** is an enzyme immunoassay for the quantitative **in vitro diagnostic** measurement of antibodies directed against oocytes in human serum.

1.1 Summary and Explanation

Antibodies directed against oocyte antigens may cause infertility in women. The application of the Anti Ovarian Ab ELISA may serve as an aid for the diagnosis of immunologically caused disorders of fertility or premature ovarian failure (1). Auto-antibodies against oocyte antigens have been shown to be associated with poor ovarian response to controlled ovarian hyperstimulation (2). Furthermore, autoimmunity is responsible for approximately 4-30% of premature ovarian insufficiency (POI) (3) and is associated with polycystic ovary syndrome (PCOS) (4-6).

Unwanted childlessness is a growing problem with which up to 20% of all couples in the reproductive age are confronted temporarily or long-term. The definition of infertility according to the WHO (WHO Laboratory Manual for the Examination of Human Semen and Semen Cervical-Mucus Interaction, 1999) is the absence of a conception within 12 months of unprotected intercourse.

2 PRINCIPLE OF THE TEST

The DEMEDIATEC Anti Ovarian Ab ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle. The microtiter wells are coated with a mix of oocyte proteins.

During incubation, anti-oocyte antibodies in the samples (standards, controls, patient specimen) bind to the coated surface of the wells. A washing step removes unbound sample components.

Added enzyme conjugate binds to the immobilized antigen-antibody-complexes. The conjugate contains anti-human immunoglobulin antibodies, labelled with horseradish peroxidase (HRP).

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of color is proportional to the concentration of the analyte in the sample. A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay. However, values for the patient samples will not be affected.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DEMEDITEC.

4 REAGENTS

4.1 Reagents provided

1. **SORB MT** Microtiterwells, 12 x 8 (break apart) strips, 96 wells; Wells coated with a mix of oocyte antigen.(incl. 1 sealing film)
2. **DIL** Dilution Buffer / **Zero Standard**, 1 vial, 50 mL, ready to use; Concentration: 0 U/mL
Contains non-mercury preservative.
3. **CAL 1 - 4** **Standard (Standard 1 - 4)**, 4 vials, 0.5 mL each, ready to use;
Concentrations: 6 – 25 – 50 – 100 U/mL; Contain non-mercury preservative.
4. **CONTROL** **Quality Control**, 1 vials, 0.5 mL each, ready to use; For control values and ranges please refer to Certificate of Analysis (CoA). Contain non-mercury preservative.
5. **ENZ CONJ** **Enzyme Conjugate**, 1 vial, 8 mL, ready to use; Anti-human IgG antibody conjugated with horseradish peroxidase; Contains non-mercury preservative.
6. **SUB TMB** **Substrate Solution**, 1 vial, 14 mL, ready to use; Tetramethylbenzidine (TMB).
7. **STOP SOLN** **Stop Solution**, 1 vial, 14 mL, ready to use; Contains 0.5 M H₂SO₄.
Avoid contact with the stop solution. It may cause skin irritations and burns.
8. **WASH SOLN 40x** **Wash Solution**, 1 vial, 30 mL (40X concentrated); See "Reagent Preparation".

Note: Additional Dilution Buffer for sample dilution is available upon request.

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Incubator for 37 °C
- Tubes for sample dilution
- Calibrated variable precision micropipettes
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Wash Solution

Add distilled water to the 40X concentrated Wash Solution. Dilute 30 mL of concentrated Wash Solution with 1170 mL distilled water to a final volume of 1200 mL. The diluted Wash Solution is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DEMEDITEC must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum can be used in this assay.

Note: Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter "Interfering Substances".

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to assaying.

Specimens stored for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

Prior to assaying, dilute each patient specimen **1:100** with Dilution Buffer.

Example:

Dilution 1:100: 5 µL sample + 495 µL Dilution Buffer (mix thoroughly)

Note: The Quality Control is ready to use and must not be diluted!

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **50 µL** each of **Zero Standard, Standard, Quality Control** and **diluted sample with new disposable tips** into appropriate wells.
3. Cover with foil and incubate for **60 minutes** at **37 °C**.
4. Rinse the wells **3 times** with **400 µL** diluted Wash Solution per well, if a plate washer is used.
- OR -
Briskly shake out the contents of the wells.
Rinse the wells **3 times** with **300 µL** diluted Wash Solution per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
5. Dispense **50 µL Enzyme Conjugate** into each well.
6. Cover with foil and incubate for **60 minutes** at **37 °C**.
7. Rinse the wells **5 times** with **400 µL** diluted Wash Solution per well, if a plate washer is used.
- OR -
Briskly shake out the contents of the wells.
Rinse the wells **5 times** with **300 µL** diluted Wash Solution per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
8. Add **50 µL** of **Substrate Solution** to each well.
9. Incubate for **30 minutes** at room temperature.
10. Stop the enzymatic reaction by adding **50 µL** of **Stop Solution** to each well.
11. Determine the optical density of the solution in each well at **450 nm (reading)** and at **620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the Stop Solution.

6.3 Calculation of Results

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Using graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4-Parameter Rodbard or 4-Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve.

The standards are already pre-diluted, therefore the 1:100 dilution of the samples must not be taken into account for the final calculation of sample concentrations.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

| Standard | Optical Density (450 nm) |
|------------------------|--------------------------|
| Zero Standard (0 U/mL) | 0.03 |
| Standard 1 (6 U/mL) | 0.27 |
| Standard 2 (25 U/mL) | 0.94 |
| Standard 3 (50 U/mL) | 1.59 |
| Standard 4 (100 U/mL) | 2.50 |

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

Normal values 0 U/mL – 10 U/mL

Elevated values > 10 U/mL

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DEMEDITEC directly.

9 PERFORMANCE CHARACTERISTICS**9.1 Intra-assay**

Mean CV: 6.40 % (range from 5.50 % – 7.80 %)

For the determination of the intra-assay precision, 6 kits from 6 different batches (produced on different days) were used. One patient sample (OD > 1.0) was applied 96 times per testing procedure.

9.2 Inter-assay

Mean CV: 7.80 % (range from 5.50 % – 9.60 %)

For the determination of the inter-assay precision one strip each of 12 kits stemming from 6 different batches (produced on different days) were used. One patient sample (OD > 1.0) was applied 72 times per testing procedure.

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Severely haemolytic or lipaemic sera or sera from patients with liver diseases should not be used. Results may be adversely affected by certain pathologic conditions, such as poly- and monoclonal gammopathies, autoimmune diseases or by an altered immune status.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DEMEDITEC.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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12 REFERENCES / LITERATURE

1. Shamilova NN et al. The role of genetic and autoimmune factors in premature ovarian failure. *J Assist Reprod Genet.* 2013,30(5), 617-22.
2. Mohammadi Yeganeh L et al. The association of different auto-antibodies against ovarian tissues and gonadotropins and poor ovarian response in intracytoplasmic sperm injection cycles. *Hum Fertil (Camb).* 2017, 20(2), 126-131.
3. Deroux a et al. Female infertility and serum auto-antibodies: a systematic review. *Clin Rev Allergy Immunol.* 2017, 53(1), 78-86.
4. Kirshenbaum M, Orvieto R. Premature ovarian insufficiency (POI) and autoimmunity - an update appraisal. *J Assist Reprod Genet.* 2019; 36(11), 2207-2215.
5. Mobeen H, Afzal N, and Kashif M. Polycystic ovary syndrome may be an autoimmune disorder. *Scientifica (Cairo).* 2016; ID 4071735.
6. Petříková J, Iazúrová IL. Ovarian failure and polycystic ovary syndrome. *Autoimmun Rev.* 2012, 11(6-7), 471-8.

SYMBOLS USED WITH DEMEDITEC ASSAYS

| Symbol | English | Deutsch | Français | Espanol | Italiano |
|---|------------------------------------|--|--|---|-------------------------------------|
|  | European Conformity | CE-Konformitäts-kennzeichnung | Conforme aux normes européennes | Conformidad europea | Conformità europea |
|  | Consult instructions for use | Gebrauchsanweisung beachten | Consulter les instructions d'utilisation | Consulte las Instrucciones | Consultare le istruzioni per l'uso |
|  | In vitro diagnostic device | In-vitro-Diagnostikum | Usage Diagnostic in vitro | Diagnóstico in vitro | Per uso Diagnostica in vitro |
|  | For research use only | Nur für Forschungszwecke | Seulement dans le cadre de recherches | Sólo para uso en investigación | Solo a scopo di ricerca |
|  | Catalogue number | Katalog-Nr. | Référence | Número de catálogo | No. di Cat. |
|  | Lot. No. / Batch code | Chargen-Nr. | No. de lot | Número de lote | Lotto no |
|  | Contains sufficient for <n> tests/ | Ausreichend für "n" Ansätze | Contenu suffisant pour "n" tests | Contenido suficiente para <n> ensayos | Contenuto sufficiente per "n" saggi |
|  | Note warnings and precautions | Warnhinweise und Vorsichtsmaßnahmen beachten | Avertissements et mesures de précaution font attention | Tiene en cuenta advertencias y precauciones | Annoti avvisi e le precauzioni |
|  | Storage Temperature | Lagerungstemperatur | Temperature de conservation | Temperatura de conservacion | Temperatura di conservazione |
|  | Expiration Date | Mindesthaltbarkeitsdatum | Date limite d'utilisation | Fecha de caducidad | Data di scadenza |
|  | Legal Manufacturer | Hersteller | Fabricant | Fabricante | Fabbricante |
| <i>Distributed by</i> | Distributor | Vertreiber | Distributeur | Distribuidor | Distributore |